

REVEAL™ Distal Attachment Cap

Endoscope Reference Chart



Endoscope		Recommended REVEAL™ Distal Attachment Cap					
Model	Distal End Outer Diameter (mm)	BX00711770 (pink)	BX00711771 (blue)	BX00711772 (orange)	BX00711773 (light green)	BX00711774 (yellow)	BX00711775 (dark green)

OLYMPUS®

Gastroscope 1100 Series

GIF-1100	9.8		✓				
----------	-----	--	---	--	--	--	--

Gastrosopes 190 series

GIF-XP190N	5.4	Not Compatible					
GIF-H190	9.2	✓					
GIF-HQ190	9.9		✓				
GIF-1TH190	10			✓			

Gastrosopes 180 series

GIF-N180	4.9	Not Compatible					
GIF-XP180N	5.5	Not Compatible					
GIF-Q180	8.8	Not Compatible					
GIF-H180	9.8		✓				
GIF-H180J	9.9		✓				
GIF-2TH180	12.2	Not Compatible					

Gastrosopes 160 series

GIF-XP160	5.9	Not Compatible					
GIF-160	8.6	Not Compatible					
GIF-Q165	9.2	✓					
GIF-Q160	9.8		✓				
GIF-Q160Z	10.8	Not Compatible					
GIF-1TQ160	10.9				✓		
GIF-XTQ160	12.9					✓	
GIF-2T160	13.2					✓	

Gastrosopes 140 series

GIF-P140	8.7	Not Compatible					
GIF-XQ140	9.4	✓					
GIF-140	9.8		✓				
GIF-Q140	10.4			✓			
GIF-1T140	10.9				✓		

Gastrosopes 200 series

GIF-N260	4.9	Not Compatible					
GIF-XP260N	5	Not Compatible					
GIF-XP260	5	Not Compatible					
GIF-XQ260	9	✓					
GIF-Q260	9.2	✓					
GIF-XK240	9.7		✓				
GIF-H260	9.8		✓				
GIF-Q240Z	9.8		✓				
GIF-Q260J	9.9		✓				
GIF-H260Z	10.8	Not Compatible					
GIF-1T240	10.9				✓		
GIF-FQ260Z	11				✓		
GIF-2TQ260M	11.7				✓		
GIF-2T240	11.8				✓		

Please contact your local STERIS Sales Representative if you have any questions regarding the compatibility of your endoscope and the REVEAL Distal Attachment Cap.

Continued →

REVEAL™ Distal Attachment Cap Endoscope Reference Chart

Endoscope		Recommended REVEAL™ Distal Attachment Cap					
Model	Distal End Outer Diameter (mm)	BX00711770 (pink)	BX00711771 (blue)	BX00711772 (orange)	BX00711773 (light green)	BX00711774 (yellow)	BX00711775 (dark green)

OLYMPUS® (continued)

Gastrosopes 150 series

GIF-Q150	9.2	✓					
----------	-----	---	--	--	--	--	--

Colonoscope 1100 Series

CF-HQ1100	13.2					✓	
-----------	------	--	--	--	--	---	--

Colonoscopes 190 Series

PCF-PH190L/I	9.7		✓				
PCF-H190L/I	11.7				✓		
CF-HQ190L/I	13.2					✓	
CF-H190L/I	13.2					✓	

Colonoscopes 180 Series

PCF-Q180AL/I	11.3				✓		
PCF-H180AL/I	11.7				✓		
CF-Q180AL/I	13.2					✓	
CF-H180AL/I	13.9						✓
CF-H180DL/I	13.9						✓

Colonoscopes 160 Series

PCF-160AL/I	11.3				✓		
CF-Q160L/I	12.8					✓	
CF-Q160AL/I	12.8					✓	
CF-Q160S	12.8					✓	
CF-Q165L/I	12.8					✓	
CF-Q260DL/I	13.6						✓
CF-Q160DL/I	13.8						✓
CF-2T160L/L	13.8						✓
CF-FH260AZL/I	14.8	Not Compatible					
CF-Q160ZL/I	15	Not Compatible					

Colonoscopes 140 Series

PCF-140L/I	11.3				✓		
CF-140L/I	13.2					✓	
CF-140S	13.2					✓	
CF-Q140L/I	13.2					✓	
CF-1T140L/I	13.8						✓

Enteroscopes

SIF-Q260	9.2	✓					
SIF-Q180	9.2	✓					
SIF-Q240	9.8		✓				
SIF-Q140	10.5			✓			
SIF-100	10.9				✓		

PENTAX®

Gastrosopes i10 and i20 Series

EG16-K10	5.2	Not Compatible					
EG27-i10	9.2	✓					
EG27-i20c	9.2	✓					
EG29-i10	9.9		✓				
EG29-i20c	9.9		✓				

Please contact your local STERIS Sales Representative if you have any questions regarding the compatibility of your endoscope and the REVEAL Distal Attachment Cap.

Continued →

REVEAL™ Distal Attachment Cap Endoscope Reference Chart

Endoscope		Recommended REVEAL™ Distal Attachment Cap					
Model	Distal End Outer Diameter (mm)	BX00711770 (pink)	BX00711771 (blue)	BX00711772 (orange)	BX00711773 (light green)	BX00711774 (yellow)	BX00711775 (dark green)

PENTAX® (continued)

Gastrosopes 90i/90K Series

EG-1690K	5.4	Not Compatible					
EG-2490K	8	Not Compatible					
EG-2790K	9	✓					
EG-2790i	9	✓					
EG-2990Zi	9.8		✓				
EG-290Kp	9.8		✓				
EG-2990i	9.8		✓				
EG-2990K	9.8		✓				
EG-3490K	11.6				✓		
EG-3890TK	12.8					✓	

Gastrosopes 80i/80K Series

EG-1580K	5.1	Not Compatible					
EG-2985K	9.8		✓				
EG-2985	10.1			✓			
EG-3485	11.5				✓		
EG-3485K	11.6				✓		
EG-3885TK	12.8					✓	

Gastrosopes 70/72K Series

EG-1870K	6.1	Not Compatible					
EG-2470K	8	Not Compatible					
EG-2770K	9.2	✓					
EG-2970K	9.8		✓				
EG-3470K	11.5				✓		
EG-3470ZK	12.3	Not Compatible					
EG-3870CIK	12.8					✓	
EG-3870TK	13.4						✓

Gastrosopes 30/31/32 Series

EG-2530	8.4	Not Compatible					
EG-2731	9	✓					
EG-2930Q	9.8		✓				
EG-2931K	9.8		✓				
EG-2931	9.8		✓				
EG-3430Z	11.4				✓		
EG-3430K	11.4				✓		
EG-3431	11.4				✓		
EG-3831T	12.8					✓	

Gastrosopes 40 Series

EG-1540	5.1	Not Compatible					
EG-1840	6	Not Compatible					
EG-2540	8.4	Not Compatible					
EG-2940	9.8		✓				
EG-3440	11.4				✓		
EG-3840T	12.8					✓	

Please contact your local STERIS Sales Representative if you have any questions regarding the compatibility of your endoscope and the REVEAL Distal Attachment Cap.

Continued →

REVEAL™ Distal Attachment Cap Endoscope Reference Chart

Endoscope		Recommended REVEAL™ Distal Attachment Cap					
Model	Distal End Outer Diameter (mm)	BX00711770 (pink)	BX00711771 (blue)	BX00711772 (orange)	BX00711773 (light green)	BX00711774 (yellow)	BX00711775 (dark green)

PENTAX® (continued)

Colonoscopes i10 and i20 Series

EC34-i10L	11.5				✓		
EC34-i10F	11.5				✓		
EC34-i10M	11.5				✓		
EC34-i20c	11.5				✓		
EC38-i10L	13.2					✓	
EC38-i10F	13.2					✓	
EC38-i10F2	13.2					✓	
EC38-i10M	13.2					✓	
EC38-i10M2	13.2					✓	
EC38-i20c	13.2					✓	

Colonoscopes 90i/90K Series

EC-2990Li	10.2			✓			
EC-2990Fi	10.2			✓			
EC-2990Mi	10.2			✓			
EC-3490TLi	10.5			✓			
EC-3490TFi	10.5			✓			
EC-3490TMi	10.5			✓			
EC-3490Li	11.5				✓		
EC-3490Fi	11.5				✓		
EC-3490FK	11.5				✓		
EC-3490LK	11.5				✓		
EC-3890LZi	13					✓	
EC-3890FZi	13					✓	
EC-3890MZi	13					✓	
EC-3890Li	13.2					✓	
EC-3890LK	13.2					✓	
EC-3890Fi	13.2					✓	
EC-3890Fi2	13.2					✓	
EC-3890MK	13.2					✓	
EC-3890MK2	13.2					✓	
EC-3890FK	13.2					✓	
EC-3890FK2	13.2					✓	
EC-3890Mi	13.2					✓	
EC-3890Mi2	13.2					✓	
EC-3890TFK	13.2					✓	
EC-3890TLK	13.2					✓	

REVEAL™ Distal Attachment Cap Endoscope Reference Chart

Endoscope		Recommended REVEAL™ Distal Attachment Cap					
Model	Distal End Outer Diameter (mm)	BX00711770 (pink)	BX00711771 (blue)	BX00711772 (orange)	BX00711773 (light green)	BX00711774 (yellow)	BX00711775 (dark green)

PENTAX® (continued)

Colonoscopes 80/85 Series

EC-3485F	11.5				✓		
EC-3485K	11.6				✓		
EC-3485FK	11.6				✓		
EC-3485LK	11.6				✓		
EC-3885K	11.6				✓		
EC-3885K22	11.6				✓		
EC-3885MK2	12.8					✓	
EC-3885FK	12.8					✓	
EC-3885FK2	12.8					✓	
EC-3885LK	12.8					✓	
EC-3885TK	12.8					✓	
EC-3885TFK	12.8					✓	
EC-3885TLK	12.8					✓	
EC-3885F	12.8					✓	
EC-3885L	12.8					✓	
EC-380LKp	13					✓	
EC-380FKp	13					✓	
EC-380FK2p	13					✓	
EC-380MK2p	13					✓	
EC-3885M/F/L	14	Not Compatible					
EC-3885M2/F2	14	Not Compatible					

Colonoscopes 70/72 Series

EC-3470FK	11.5				✓		
EC-3470LK	11.5				✓		
EC-3470K	11.6				✓		
EC-3870K	11.6				✓		
EC-3870K2	11.6				✓		
EC-3870CILK	12.8					✓	
EC-3870CIFK	12.8					✓	
EC-3870CIFK2	12.8					✓	
EC-3870TK	12.8					✓	
EC-3870ZK	12.8					✓	
EC-3870LZK	13					✓	
EC-3872LK	13.4						✓
EC-3872LZK	13.4						✓
EC-3872TLK	13.4						✓

REVEAL™ Distal Attachment Cap Endoscope Reference Chart

Endoscope		Recommended REVEAL™ Distal Attachment Cap					
Model	Distal End Outer Diameter (mm)	BX00711770 (pink)	BX00711771 (blue)	BX00711772 (orange)	BX00711773 (light green)	BX00711774 (yellow)	BX00711775 (dark green)

PENTAX® (continued)

Colonoscopes 30/31/32 Series

EC-3430LK	11.6				✓		
EC-3430F	11.7				✓		
EC-3430L	11.7				✓		
EC-3430LZ	11.7				✓		
EC-3831F	12.8					✓	
EC-3830LZ	12.8					✓	
EC-3830TL	12.8					✓	
EC-3831L	12.8					✓	
EC-3832L	12.8					✓	
EC-3832TL	12.8					✓	
EC-3832LK	13.2					✓	

Colonoscopes 40 Series

EC-3440L	11.7				✓		
EC-3440F	11.7				✓		
EC-3840F	12.8					✓	
EC-3840TL	12.8					✓	
EC-3840L	12.8					✓	

Enteroscopes

VSF-2900	9.8		✓				
VSF-3430	11.7				✓		
VSF-3430K	11.6				✓		
VSF-3440	11.7				✓		

FUJIFILM®

Endoscopic Instrument

EI-740D/S	12.8					✓	
-----------	------	--	--	--	--	---	--

Gastrosopes 800 series

EG-840T	9.8		✓				
EG-840TP	7.9	Not Compatible					

Gastrosopes EVE 700 Series

EG-760R	9.2	✓					
EG-760CT	10.5			✓			
EG-740N	5.8	Not Compatible					

Gastrosopes EVE 700 Series

EG-760R	9.2	✓					
EG-760CT	10.5			✓			
EG-740N	5.8	Not Compatible					

Gastrosopes EVE 600 Series

EG-600WR	9.2	✓					
EG-600ZW	9.9		✓				

Please contact your local STERIS Sales Representative if you have any questions regarding the compatibility of your endoscope and the REVEAL Distal Attachment Cap.

Continued →

REVEAL™ Distal Attachment Cap Endoscope Reference Chart

Endoscope		Recommended REVEAL™ Distal Attachment Cap					
Model	Distal End Outer Diameter (mm)	BX00711770 (pink)	BX00711771 (blue)	BX00711772 (orange)	BX00711773 (light green)	BX00711774 (yellow)	BX00711775 (dark green)

FUJIFILM® (continued)

Colonoscopes EVE 700 Series

EC-740T/L	9.8			✓			
EC-760S	12.8					✓	
EC-760S-V/L	12.8					✓	
EC-760Z	12.8					✓	
EC-760P	11.1				✓		
EC-760ZP	11.7				✓		
EC-760P-V/L	11.8				✓		
EC-760R-V/L	12	Not Compatible					
EC-760R	12	Not Compatible					

Colonoscopes EVE 600 Series

EC-60WM/WI/WI	12	Not Compatible					
EC-600ZW/M, /L	12.8					✓	

Gastrosopes EVE 500 Series

EG-530NP	4.9	Not Compatible					
EG-530N	5.9	Not Compatible					
EG-530NW	5.9	Not Compatible					
EG-530FP	8.5	Not Compatible					
EG-530WR	9.4	✓					
EG-590WR	9.6		✓				
EG-530CT	10.8	Not Compatible					
EG-590ZW	10.8	Not Compatible					
EG-530D	11.5				✓		

Colonoscopes EVE 500 Series

EC-530MP	11				✓		
EC-530LP	11				✓		
EC-530LS	11.5				✓		
EC-530LS2	11.5				✓		
EC-530WM3/WI3/WL3	12.8					✓	
EC-590ZW/M	12.8					✓	
EC-590ZW/L	12.8					✓	
EC-590WI	12.8					✓	
EC-590W/L	12.8					✓	
EC-590WM	12.8					✓	
EC-530WM	12.8					✓	
EC-530WI	12.8					✓	
EC-530WL	12.8					✓	
EC-530MT	12.8					✓	
EC-530IT	12.8					✓	
EC-530FM	12.8					✓	
EC-530FI	12.8					✓	
EC-530FL	12.8					✓	
EC-530LT	12.8					✓	
EC-530HL	12.8					✓	
EC-530HL2	12.8					✓	
EC-530DM	12.8					✓	
EC-530DL	12.8					✓	

Please contact your local STERIS Sales Representative if you have any questions regarding the compatibility of your endoscope and the REVEAL Distal Attachment Cap.

Continued →

REVEAL™ Distal Attachment Cap Endoscope Reference Chart

Endoscope		Recommended REVEAL™ Distal Attachment Cap					
Model	Distal End Outer Diameter (mm)	BX00711770 (pink)	BX00711771 (blue)	BX00711772 (orange)	BX00711773 (light green)	BX00711774 (yellow)	BX00711775 (dark green)

FUJIFILM® (continued)

Gastrosopes EVE 400 Series

EG-470N5	5.9	Not Compatible					
EG-450PE5	8.2	Not Compatible					
EG-450WR5	9.4	✓					
EG-490WR5	9.8		✓				
EG-450ZW5	9.8		✓				
EG-490ZW5	10.8	Not Compatible					
EG-450CT5	10.8	Not Compatible					
EG-450D5	11.5				✓		

Gastrosopes EVE 200 Series

EG-270N5	5.9	Not Compatible					
EG-270NH5	5.9	Not Compatible					
EG-270N	5.9	Not Compatible					
EG-250PE5	8.2	Not Compatible					
EG-250WR5	9.4	✓					
EG-201FP	9.8		✓				
EG-250CT5	10.8	Not Compatible					
EG-250D5	11.5				✓		

FUJIFILM® (continued)

Colonoscopes EVE 400 Series

EC-490WR5	9.3	✓					
EC-450BI5	9.4	✓					
EC-450MP5	11				✓		
EC-450LP5	11				✓		
EC-450LS5	11.5				✓		
EC-490ZW5L	12.8					✓	
EC-490ZW5M	12.8					✓	
EC-450ZW5L	12.8					✓	
EC-450ZW5M	12.8					✓	
EC-450DL5	12.8					✓	
EC-450DM5	12.8					✓	
EC-450WI5	12.8					✓	
EC-450WL5	12.8					✓	
EC-450WM5	12.8					✓	
EC-450HL5	12.8					✓	

REVEAL™ Distal Attachment Cap Endoscope Reference Chart

Endoscope		Recommended REVEAL™ Distal Attachment Cap					
Model	Distal End Outer Diameter (mm)	BX00711770 (pink)	BX00711771 (blue)	BX00711772 (orange)	BX00711773 (light green)	BX00711774 (yellow)	BX00711775 (dark green)

FUJIFILM® (continued)

Colonoscopes EVE 200 Series

EC-250LP5	11				✓		
EC-250MP5	11				✓		
EC-250LS5	11.5				✓		
EC-250DL5	12.8					✓	
EC-250DM5	12.8					✓	
EC-250HL5	12.8					✓	
EC-250WI5	12.8					✓	
EC-250WL5	12.8					✓	
EC-250WM5	12.8					✓	
EC-201WI	12.8					✓	
EC-201WL	12.8					✓	
EC-201WM	12.8					✓	

Enteroscopes

EN-450P5/20	8.5	Not Compatible					
EN-450T5	9.4	✓					
EN-450T5/W	9.4	✓					
EN-580T	9.4	✓					
EN-200WM23	10.5			✓			
EN-410WM5/23	10.7	Not Compatible					

AMBU®

Gastrosopes

aScope™ Gastro	9.9			✓			
aScope™ Gastro Large	11.5				✓		

Please contact your local STERIS Sales Representative if you have any questions regarding the compatibility of your endoscope and the REVEAL Distal Attachment Cap.

The compatibility information contained in this Endoscope Reference Chart is based on endoscope manufacturers' published specifications. The Reference Chart is to be used for information purposes only. Because the Reference Chart may not include all available endoscopes, contact STERIS with any compatibility questions. Due to variability in manufacturing and repairs, the REVEAL™ Distal Attachment Cap may not be compatible with a listed endoscope.

REVEAL™ is a registered trademark of STERIS.

All other product and company names referenced are trademarks of their respective owners.