

July 14, 2025

Dear Valued Customer,

This letter is to inform you of a component change to the Rapicide™ PA High-Level Disinfectant Test Strips (Item ML02-0118). The new item is fully equivalent in function to the previous product, with the following changes to the dimensions:

- The test strip bottle is larger, and is intended to facilitate ease of holding and one-handed opening
- The test strip is longer, and is intended to allow for easier holding and use
- The test strip reaction pad is larger, and is intended to allow for an easier reading


There are no changes to the chemical formulation, instructions for use, packaging configuration, part number or other aspects of the product. Please continue to use Rapicide™ PA High-Level Disinfectant Test Strips in the same manner moving forward. A visual comparison is below for your reference.



Left: Previous test strip bottle. Right: New test strip bottle.

Please don't hesitate to contact your STERIS representative if you have any questions.

Sincerely,

A handwritten signature in black ink, appearing to read "Michael McGrath".

Michael McGrath  
Vice President, Canada