



STERIS

Surgical
Instruments

Surgical Instrument Sales

STERIS Instrument Management Services have over 20 years of experience in the management of surgical instruments and devices, working in partnership with the NHS and other clinical providers to maximise patient safety and drive improvements in quality outcomes for our Customers.

STERIS Instrument Management Services over a range of over 4000 instruments.

We are constantly expanding our range and utilise precision engineering and lean principles in all our processes. We strive towards continuous improvement enabling us to successfully support both the public and private healthcare sector in delivering safe and efficient patient care.



Surgical Instrument Repair



Surgical Device Repair



Surgical Instrument sales & health consumables




Surgical instrument Sterilisation



Endoscope Decontamination

Our skilled team of technicians offer one of the most comprehensive surgical instrument portfolios and repair services in the UK. Based at our repair facility in Hertfordshire, we offer reliable and high-quality services and products. We are also able to advise on aspects of maintenance and longevity to ensure you get the maximum life from your instruments.



Thank you for picking up our catalogue and choosing STERIS Instrument Management Services as your surgical instrument supplier.

You can contact us

STERIS IMS – Instrument Sales And Repairs

14 Pindar Road
Hoddesdon
Hertfordshire
EN11 0BZ

You can order today

Telephone orders

+44 (0)3452 414818

Email orders

instrumentinfo@steris.com



Insulation

Diathermy instruments are tested for integrity of diathermy insulation.

Accuracy

We ensure all moving parts fit easily together, that jaws are square, parallel and mate. Additionally, all serrations must mesh accurately.

Tungsten Carbide inserts

All tungsten carbide inserts are inspected under the microscope where we ensure all soldering is completed without any fractures or blow holes.



Quality First



Workmanship

We thoroughly control the workmanship of the instruments. We ensure that all screw joints & Rivets are immovable, that there are no movements or crevices in joints.

Under the Microscope

We use a microscope to inspect the jaws for all types of instruments. Magnification allows us to identify even the most minor imperfections.

Reliability

In accordance with BS5194 and ISO 13402 we carry out corrosion resistance testing. This ensures our instruments are resistant to corrosion during re-processing.



Inspections for Scissors, Forceps and Needle Holders

Elasticity

Artery forceps and needle holders are tested for elasticity in accordance with BS5194-2.

Functionality

We complete a clamp test for all needle holders as per BS5194-2 in order to ensure optimal functionality.

Sharpness

A cutting test is carried out to ensure all scissors are fit for purpose.


Material Composition

We complete a hardness test using the Rockwell hardness test method. This test confirms the instruments are the correct hardness for their intended purpose as specified within BS 5194.

We ensure that the hardness between mating surfaces never differs more than 4 units on the Rockwell scale. Mating surfaces need to be compatible in order to prevent damage.

Consistency in Processes

The size and shape of finished instruments are checked against the original drawing or master sample to ensure our processes remain consistent.



A Perfect Finish

We ensure the surface finish conforms to the specific requirements of each instrument and conforms to BS5194

Attention to Detail

We ensure that the surface condition of the instruments are free from pores, crevices and grinding marks. Before dispatch we check that all polishing material is removed.

Why Choose STERIS?

We know you have plenty of choice when it comes to selecting your surgical instrument provider, let us tell you why we believe STERIS can provide your organisation with highly reliable and safe products.

Quality First


We combine our quality management system with our corporate culture to create an environment that favours and promotes quality standards set out for our industry. This enables us to ensure that the safety of our surgical instruments is set at the heart of all of our operational efforts.

At STERIS Instrument Management Services, are certified to ISO 13485. We know that the services and products we supply are crucial to patient safety and therefore we never stop seeking excellence in quality through continuous improvement initiatives throughout our business activities.

This ensures that all instruments comply to the highest standards in terms of reliability and materials. All stainless steels used in our manufacturing meet the requirements set out in BS EN ISO 7153-1 as appropriate to the instrument and our inspection protocols ensures that the specifications within the appropriate parts of BS 5194 are met.

Our commitment to compliance and improvements is integral to our business and the sector in which we operate.

Order Today

 +44 (0)3452 414818

Email instrumentinfo@steris.com

Contents

General Scissors 11

| | |
|---------------------------|----|
| General Scissors | 12 |
| Scissors Tungsten Carbide | 20 |
| Scissors SuperCut | 23 |

Pliers and Wire Cutters 27

| | |
|-------------------------------|----|
| Pliers and Wire Cutters | 28 |
| Wire Cutters Tungsten Carbide | 30 |
| Scissors SuperCut | 23 |

Chiropody 33

Forceps 37

| | |
|-------------------------|----|
| Tissue Forceps | 38 |
| Artery Forceps | 41 |
| Dissecting Forceps | 52 |
| Suture Forceps | 61 |
| Aural Forceps | 62 |
| Nasal Forceps | 65 |
| Tonsil Forceps | 66 |
| Cholecystectomy Forceps | 67 |
| Lung Forceps | 67 |
| Bone Cutting Forceps | 68 |
| Rectal Forceps | 69 |

| | |
|-------------------|----|
| Endoscope Forceps | 69 |
|-------------------|----|

| | |
|--------------------------|----|
| Purse String Applicators | 72 |
|--------------------------|----|

| | |
|-------------------|----|
| Vulsellum Forceps | 72 |
|-------------------|----|

| | |
|-------------------|----|
| Tenaculum Forceps | 73 |
|-------------------|----|

Clamps, Clips & Pins 75

| | |
|---------------|----|
| Tubing Clamps | 76 |
|---------------|----|

| | |
|---------------------|----|
| Clips & Safety Pins | 77 |
|---------------------|----|

| | |
|-------------|----|
| Towel Clips | 78 |
|-------------|----|

Needle Holders 81

| | |
|----------------|----|
| Needle Holders | 82 |
|----------------|----|

Aneurysm Needles, Skin Hooks, Probes & Fistula Directors 91

| | |
|------------------|----|
| Aneurysm Needles | 92 |
|------------------|----|

| | |
|--------|----|
| Probes | 92 |
|--------|----|

| | |
|-------------------|----|
| Fistula Directors | 94 |
|-------------------|----|

| | |
|------------|----|
| Skin Hooks | 95 |
|------------|----|

Retractors 97

| | |
|------------|----|
| Retractors | 98 |
|------------|----|

Contents

Oral Instruments 115

| | |
|-------------------|-----|
| Tongue Depressors | 116 |
| Mouth Gags | 116 |

Tracheal Instruments 118

Aural Instruments 121

| | |
|---------------------|-----|
| Aural Curettes | 122 |
| Aural Wool Carriers | 125 |
| Aural Speculums | 126 |
| ENT Knives | 127 |

Nasal, Tonsil & Adenoid Instruments 129

| | |
|-----------------------|-----|
| Adenoid Curettes | 130 |
| Nasal Speculums | 130 |
| Tonsil Snare | 131 |
| Nasal Cartilage Knife | 132 |
| Nasal Chisels | 133 |

Dissectors & Elevators 135

Intestinal Clamps 139

| | |
|-------------------|-----|
| Intestinal Clamps | 140 |
|-------------------|-----|

Bone Instruments 143

| | |
|------------------------------|-----|
| Bone Rongeurs | 144 |
| Bone Elevators | 149 |
| Bone Rugine/Rasps | 152 |
| Bone Curettes/Hooks | 153 |
| Hammers and Mallets | 154 |
| Bone Chisels | 156 |
| Gouges | 157 |
| Osteotomes | 158 |
| Rib Shears and Rib Spreaders | 161 |

Plaster Cast 165

Rongeurs 169

| | |
|--------------------------|-----|
| Rongeurs | 170 |
| Semi-detachable Rongeurs | 172 |

Diathermy 177

| | |
|-------------------|-----|
| Diathermy | 178 |
| Diathermy Quivers | 179 |

Miscellaneous Instruments 181


| | |
|-----------------------|-----|
| Miscellaneous | 182 |
| Rasps | 186 |
| Speculums | 186 |
| Rectal Dilator | 187 |
| Urethral Forceps | 187 |
| Penile Clamp | 187 |
| Urethral Bougies | 188 |
| Catheter Introducers | 189 |
| Gall Stone | 190 |
| Dilators | 190 |
| Uterine Dilators | 190 |
| Hysterectomy Clamp | 192 |
| Myoma Screw | 193 |
| Uterine Curettes | 194 |
| Vaginal Speculum | 195 |
| Biopsy Punches | 201 |
| Gigli Saw and Handles | 204 |

Sterilisation Baskets & Trays 205

| | |
|-------------------------------|-----|
| Sterilisation Baskets & Trays | 206 |
| Trinket Boxes | 210 |
| Bowls | 211 |

Full product listings 212


| | |
|---------------------|-----|
| Product Codes Index | 212 |
| Product Names A-Z | 268 |



General Scissors

Dressing Stitch Scissors Sharp/Sharp Straight

| | |
|---------|-------|
| PH70200 | 100mm |
| PH70201 | 115mm |
| PH70202 | 130mm |
| PH70203 | 140mm |
| PH70204 | 150mm |
| PH70205 | 180mm |




Dressing Stitch Scissors Sharp/Sharp Curved

| | |
|---------|-------|
| PH70208 | 100mm |
| PH70209 | 115mm |
| PH70210 | 130mm |
| PH70211 | 140mm |
| PH70212 | 150mm |
| PH70213 | 180mm |



Dressing Stitch Scissors Sharp/Blunt Straight

| | |
|---------|-------|
| PH70216 | 100mm |
| PH70217 | 115mm |
| PH70218 | 130mm |
| PH70219 | 140mm |
| PH70220 | 150mm |
| PH70221 | 180mm |




Dressing Scissors Sharp/Blunt Curved

| | |
|---------|-------|
| PH70224 | 100mm |
| PH70225 | 115mm |
| PH70226 | 130mm |
| PH70227 | 140mm |
| PH70228 | 150mm |
| PH70229 | 180mm |



Dressing Scissors Blunt/Blunt Straight

| | |
|---------|-------|
| PH70232 | 100mm |
| PH70233 | 115mm |
| PH70234 | 130mm |
| PH70235 | 140mm |
| PH70236 | 150mm |
| PH70237 | 180mm |



Dressing Scissors Blunt/Blunt Curved

| | |
|---------|-------|
| PH70240 | 100mm |
| PH70241 | 115mm |
| PH70242 | 130mm |
| PH70243 | 140mm |
| PH70244 | 150mm |
| PH70245 | 180mm |



Dissecting Scissors with Open Shanks Straight

| | |
|---------|-------|
| PH70260 | 115mm |
| PH70261 | 130mm |
| PH70262 | 150mm |

Aufricht Scissors Straight

| | |
|---------|-------|
| PH70275 | 150mm |
|---------|-------|

Aufricht Scissors Curved

| | |
|---------|-------|
| PH70276 | 150mm |
|---------|-------|

Iris Scissors Fine Pointed Straight

| | |
|---------|-------|
| PH70281 | 100mm |
| PH70282 | 115mm |
| PH70283 | 130mm |

Iris Scissors Fine Pointed Curved


| | |
|---------|-------|
| PH70284 | 100mm |
| PH70285 | 115mm |
| PH70286 | 130mm |

Lister Bandage Scissors Angled

| | |
|---------|-------|
| PH70289 | 115mm |
| PH70290 | 140mm |
| PH70291 | 185mm |
| PH70292 | 200mm |
| PH70293 | 230mm |

Mayo Scissors Straight

| | |
|---------|-------|
| PH70300 | 140mm |
| PH70301 | 165mm |
| PH70302 | 190mm |
| PH70303 | 215mm |
| PH70304 | 230mm |




Mayo Scissors Curved

| | |
|---------|-------|
| PH70305 | 140mm |
| PH70306 | 165mm |
| PH70307 | 190mm |
| PH70308 | 215mm |
| PH70309 | 230mm |




Metzenbaum Scissors Straight

| | |
|---------|-------|
| PH70350 | 140mm |
| PH70351 | 165mm |
| PH70352 | 190mm |
| PH70353 | 215mm |
| PH70354 | 240mm |




Metzenbaum Scissors Curved

| | |
|---------|-------|
| PH70360 | 140mm |
| PH70361 | 165mm |
| PH70362 | 190mm |
| PH70363 | 215mm |
| PH70364 | 240mm |



McIndoe Scissors Curve

| | |
|---------|-------|
| PH70370 | 115mm |
| PH70371 | 190mm |



Kilner Scissors Straight

| | |
|---------|-------|
| PH70380 | 115mm |
|---------|-------|




Kilner Scissors Curved

PH70381 115mm



Strabismus Scissors Straight

PH70382 115mm



Strabismus Scissors Curved

PH70383 115mm



Universal Scissors for Ligature and Wire Cutting Angled

PH70402 115mm



Nelson Scissors Straight (Thoracic)

PH70420 230mm

PH70421 280mm



Nelson Scissors Curved (Thoracic)

PH70422 230mm



Wilson Tonsil Scissors Straight

PH70430 200mm



Wilson Tonsil Scissors Curved

PH70431 200mm



Carless Light Weight Scissors Straight

PH70454 165mm

PH70455 215mm



Carless Light Weight Scissors Curved


PH70457 165MM

PH70458 215MM




Pott Smith Scissors Angled to Side at 25°

PH70465 190mm




Pott Smith Scissors Angled to Side at 45°

PH70466 190mm



Pott Smith Scissors Angled on Flat



PH70467 190mm

Lawson Tait Scissors (Barnes) (Episiotomy) Angled



PH70470 130mm

Braun Stadler Scissors (Episiotomy) Angled




PH70474 150mm

PH70475 190mm

PH70476 200mm

Doyen Scissors (Uterine) Straight



PH70478 170mm

PH70480 200mm

Doyen Scissors (Uterine) Curved



PH70479 170mm

PH70481 200mm

Fergusson Scissors (Uterine) Straight



PH70483 210mm

Fergusson Scissors (Uterine) Angled to Flat


PH70484 210mm



Sims Scissors (Uterine) Straight

PH70488 200mm

PH70489 230mm



Sims Scissors (Uterine) Curved

PH70490 200mm

PH70491 230mm



Sims Scissors (Uterine) Angled to One Side

PH70492 200mm

PH70493 230mm



Did you know we
offer instrument
repair services?

We offer flexible and comprehensive services including a rapid repair service. All our services are ISO 13485 accredited with warranties provided for all repairs.

Call us today to find out more
(0)3452 414818

 **STERIS**

American Pattern Umbilical Cord Cutter

PH70500 105mm



STERIS IMS 03452 414818

Universal Tough Cut Scissors (Tubing Scissors)



| | |
|---------|-------|
| PH70520 | 140MM |
| PH70521 | 190MM |

Capsular Scissors Angled On Flat



| | |
|---------|-------|
| PH70522 | 250mm |
|---------|-------|

American Pattern Fistula Scissors Curved



| | |
|---------|-------|
| PH70530 | 180mm |
|---------|-------|

McKissock Dural Scissors Straight



| | |
|---------|-------|
| PH70538 | 150mm |
|---------|-------|

McKissock Dural Scissors Curved



| | |
|---------|-------|
| PH70539 | 150mm |
|---------|-------|


Lahey Scissors Curved



| | |
|---------|-------|
| PH70570 | 140mm |
|---------|-------|

Dressing/Stitch Scissors Straight Tungsten Carbide Jaws Sharp/Sharp

| | |
|---------|-------|
| PH70600 | 150mm |
| PH70601 | 180mm |




Dressing/Stitch Scissors Curved Tungsten Carbide Jaws Sharp/Sharp

| | |
|---------|-------|
| PH70605 | 150mm |
| PH70606 | 180mm |



Dressing/Stitch Scissors Straight Tungsten Carbide Jaws Sharp/Blunt

| | |
|---------|-------|
| PH70610 | 150mm |
| PH70612 | 180mm |



Dressing/Stitch Scissors Curved Tungsten Carbide Jaws Sharp/Blunt

| | |
|---------|-------|
| PH70615 | 150mm |
| PH70616 | 180mm |




Dressing/Stitch Scissors Straight Tungsten Carbide Jaws Blunt/Blunt

| | |
|---------|-------|
| PH70620 | 150mm |
| PH70621 | 180mm |



Iris Scissors Straight Tungsten Carbide Jaws Fine Pointed

| | |
|---------|-------|
| PH70625 | 115mm |
|---------|-------|




Iris Scissors Curved Tungsten Carbide Jaws Fine Pointed



PH70626 115mm

Kilner Scissors Straight Tungsten Carbide Jaws



PH70630 115mm

Kilner Scissors Curved Tungsten Carbide Jaws



PH70631 115mm

Strabismus Scissors Straight Tungsten Carbide Jaws



PH70635 115mm

Strabismus Scissors Curved Tungsten Carbide Jaws



PH70636 115mm

Mayo Scissors Straight Tungsten Carbide Jaws



| | |
|---------|-------|
| PH70640 | 140mm |
| PH70641 | 150mm |
| PH70642 | 180mm |
| PH70643 | 200mm |
| PH70644 | 230mm |

Mayo Scissors Curved Tungsten Carbide Jaws

| | |
|---------|-------|
| PH70650 | 140mm |
| PH70651 | 150mm |
| PH70652 | 180mm |
| PH70653 | 200mm |
| PH70654 | 230mm |



Lister Bandage Scissors Angled Tungsten Carbide Jaws

| | |
|---------|-------|
| PH70655 | 115mm |
| PH70656 | 140mm |
| PH70657 | 185mm |
| PH70658 | 200mm |
| PH70659 | 230mm |



Did you know we offer instrument repair services?

We offer flexible and comprehensive services including a rapid repair service. All our services are ISO 13485 accredited with warranties provided for all repairs.

Call us today to find out more
(0)3452 414818

 **STERIS**


Metzenbaum Scissors Straight Tungsten Carbide Jaws

| | |
|---------|-------|
| PH70660 | 140mm |
| PH70661 | 180mm |
| PH70662 | 200mm |
| PH70663 | 230mm |



Metzenbaum Scissors Curved Tungsten Carbide Jaws

| | |
|---------|-------|
| PH70670 | 140mm |
| PH70671 | 180mm |
| PH70672 | 200mm |
| PH70673 | 230mm |



Universal Cutters with Tungsten Carbide Jaws For Ligature and Wire Cutting

| | |
|---------|-------|
| PH70680 | 115mm |
|---------|-------|



STERIS IMS 03452 414818

Iris SuperCut Scissors Curved Fine Pointed Black Handles Straight



PH70701 115mm

Iris SuperCut Scissors Curved Fine Pointed Black Handles



PH70702 115mm

Kilner SuperCut Scissors Straight Black Handles



PH70705 115mm

Kilner SuperCut Scissors Curved Black Handles




PH70706 115mm

Stevens SuperCut Tenotomy Scissors Straight Black Handles



PH70710 115mm

Stevens SuperCut Tenotomy Scissors Curved Black Handles




PH70711 115mm

PH70713 135mm

Mayo Scissors Straight SuperCut Black Handles

| | |
|---------|-------|
| PH70715 | 140mm |
| PH70716 | 150mm |
| PH70717 | 180mm |
| PH70718 | 200mm |




Mayo Scissors Curved SuperCut Black Handles

| | |
|---------|-------|
| PH70720 | 140mm |
| PH70721 | 150mm |
| PH70722 | 180mm |
| PH70723 | 200mm |



Metzenbaum Scissors Straight SuperCut Black Handles

| | |
|---------|-------|
| PH70730 | 150mm |
| PH70731 | 180mm |
| PH70732 | 200mm |




Metzenbaum Scissors Curved SuperCut Black Handles

| | |
|---------|-------|
| PH70735 | 150mm |
| PH70736 | 180mm |
| PH70737 | 200mm |



Jameson Reynolds Scissors Curved SuperCut Black Handles

| | |
|---------|-------|
| PH70738 | 160mm |
|---------|-------|



McIndoe Scissors Curved SuperCut Black Handles

| | |
|---------|-------|
| PH70739 | 180mm |
|---------|-------|



Strabismus Scissors Straight SuperCut Black Handles

PH70745 230mm

Strabismus Scissors Curved SuperCut Black Handles

PH70746 115mm

Gorney/Freeman Face Lift Scissors Curved Tungsten Carbide SuperCut

PH70755 190mm

PH70757 190mm

Gorney/Freeman Face Lift Scissors Straight Tungsten Carbide SuperCut

PH70756 230mm

PH70758 230mm

McIndoe Scissors with Tungsten Carbide Jaws Straight

PH70800 140mm

PH70809 165mm

PH70801 180mm

PH70802 200mm

PH70803 230mm

McIndoe Scissors with Tungsten Carbide Jaws Curved

PH70804 140mm

PH70808 165mm

PH70805 180mm

PH70811 190mm

PH70806 200mm

PH70807 230mm

PH70810 280mm

Pliers and Wire Cutters

Dental Cutter/Plier

PH70679 XXmm



Small Wire Cutters Tungsten Carbide Jaws Straight

PH70685 125mm



Compound Action Wire Top Cutter Tungsten Carbide Jaws

PH70686 150mm



Compound Action Wire Cutter Tungsten Carbide Jaws Angled

PH70690 180mm



Compound Action Side Cutters Tungsten Carbide Jaws Straight

PH70691 200mm



Compound Action Top Wire Cutter Tungsten Carbide Jaws Straight

PH70692 200mm




Compound Action Side Wire Cutter Tungsten Carbide Jaws Straight

PH70693 230mm




Compound Action Large Wire Cutter To Cut 5.5mm diameter max with Tungsten Carbide Jaws

PH70694 500mm



Maun Plier Serrated Jaws and Side Cutter

PH70695 180mm



Power Cutter Tungsten Carbide Jaws Cutting up to 3.5mm Wire

PH70696 240mm



Small Wire Plier

PH70698 140mm



Large Wire Plier

PH70700 175mm



Chiroprody

Nail Nipper Straight Double Spring with Lock Straight

PH80526 140mm



Nail Nipper Concave Double Spring with Lock Curved

PH80527 140mm

PH80528 160mm



Nail Nipper Compound Action Straight

PH80534 150mm



Dovo German Pattern Nail Nipper Straight

PH80536 140mm




Thwaites Nail Nipper Pistol Grip Single Spring

PH80541 140mm



Double Ended Spatula and Packer with Solid Handle Straight

PH80546 155mm




Fickling Nail Elevator Straight

PH80549 155mm



Francis Chisel

PH80550 115mm



Blacks File Double Ended Fine Cut

PH80551 150mm



Blacks File Double Ended Course Cut

PH80553 140mm



Forceps

Collingwood Stewart Hernia Forceps

PH483470 130mm




Allis Tissue Forceps 3 into 4 Teeth

PH503550 150mm



Allis Tissue Forceps 4 into 5 Teeth

PH503551 150mm



Allis Tissue Forceps 5 into 6 Teeth

PH503552 150mm

PH503553 180mm

PH503554 200mm

PH503555 250mm

PH503556 280mm



Babcock Tissue Forceps

PH503558 160mm

PH503559 180mm

PH503560 200mm


PH503561 230mm



Collin Duval Tissue Forceps

PH503567 180mm

PH503568 200mm



Judd Allis Tissue Forceps 3 into 4 Teeth



| | |
|----------|-------|
| PH503575 | 150mm |
| PH503576 | 190mm |

Lahey Tissue Forceps 3 into 3 Teeth



| | |
|----------|-------|
| PH503577 | 150mm |
| PH503578 | 200mm |

Lane Tissue Forceps 1 into 2 Teeth



| | |
|----------|-------|
| PH503580 | 130mm |
| PH503581 | 150mm |
| PH503582 | 195mm |
| PH503583 | 200mm |

Littlewood Tissue Forceps 2 into 3 Teeth



| | |
|----------|-------|
| PH503590 | 190mm |
|----------|-------|

Poirier Tissue Forceps 2 into 3 Teeth



| | |
|----------|-------|
| PH503600 | 130mm |
|----------|-------|


Rutherford Morison Tissue Forceps 4 into 5 Teeth



| | |
|----------|-------|
| PH503605 | 160mm |
|----------|-------|

Stiles Tissue Forceps 6 into 7 Teeth

PH503607 200mm



Craft Vasectomy Forceps

PH503611 150mm




Ring Type Vasectomy Clamps

PH503630 3.0mm Ring 140mm

PH503631 3.5mm Ring 140mm

PH503632 4.0mm Ring 140mm




Scalpelless Vasectomy Forceps Curved Sharp

PH503635 140mm



Adson Artery Forceps Straight




PH90555 180mm

Adson Artery Forceps Curved




PH90556 180mm

Birkett Artery Forceps Straight



PH90558 190mm

Birkett Artery Forceps Curved




PH90559 190mm

Bland Sutton Artery Forceps Curved



PH90561 190mm

Cairn Artery Forceps Straight



PH90570 145mm

Cairn Artery Forceps Curved

PH90571 145mm



Cairn Artery Forceps Curved To One Side

PH90572 145mm



Crile Artery Forceps Straight

PH90582 100mm

PH90583 115mm

PH90584 130mm



Crile Artery Forceps Curved

PH90585 140mm

PH90586 160mm

PH90587 180mm



Crile Rankin Artery Forceps Straight

PH90588 140mm

PH90589 160mm




Crile Rankin Artery Forceps Curved

PH90590 140mm

PH90591 160mm



Cushing Artery Forceps Straight



PH90592 145mm

Cushing Artery Forceps Curved



PH90593 145mm

Cushing Artery Forceps Curved to One Side



PH90594 145mm

Denis Browne Tonsil Forceps Angled to One Side




PH90595 200mm

Dunhill Artery Forceps Straight



PH90605 130mm

Dunhill Artery Forceps Curved



PH90607 130mm

PH90608 180mm

Fickling Artery Forceps Angled Plain

PH90612 180mm



Fickling Artery Forceps Angled 1 into 2 Teeth


PH90613 180mm



Grey Turner Artery Forceps Straight

PH90618 145mm

PH90620 180mm



Grey Turner Artery Forceps Curved

PH90619 145mm

PH90621 180mm



Hartmann Artery Forceps Extra Fine With Box Joint Straight

PH90622 90mm




Hartmann Artery Forceps Extra Fine With Box Joint Curved

PH90623 90mm



Hartmann Artery Forceps Extra Fine With Box Joint Straight



PH90624 100mm

Halsted Mosquito Artery Forceps With Box Joint Straight



PH90626 130mm

Did you know we offer instrument repair services?

We offer flexible and comprehensive services including a rapid repair service. All our services are ISO 13485 accredited with warranties provided for all repairs.



Contact us today to find out more

+44 (0)3452 414818



Halstead Mosquito Artery Forceps Box Joint



PH90628 125mm

Halstead Mosquito Artery Forceps Straight 1 into 2 Teeth



PH90629 130mm


Halstead Mosquito Artery Forceps Curved 1 into 2 Teeth

PH90630 130mm



Micro Halstead Mosquito Artery Forceps Straight

PH90631 130mm




Micro Halstead Mosquito Artery Forceps Curved

PH90632 130mm



Heiss Half Curve Artery Forceps

PH90636 200mm



Heiss Full Curve Artery Forceps

PH90637 200mm




Halsted Mosquito Artery Forceps Curved Atraumatic Jaws

PH90639 130mm



Jolls Artery Forceps Straight




PH90645 145mm

Jolls Artery Forceps Curved



PH90646 145mm

Kelly Artery Forceps Straight



PH90650 140mm

PH90651 160mm

Kelly Artery Forceps Curved



PH90652 140mm


PH90653 160mm

Kelly Rankin Artery Forceps With Box Joint Curved



PH90655 160mm

Kelly Frazier Artery Forceps Curved



PH90657 180mm

Kocher Artery Forceps Straight 1 into 2 Teeth

| | |
|---------|-------|
| PH90682 | 140mm |
| PH90683 | 150mm |
| PH90684 | 180mm |
| PH90685 | 200mm |
| PH90686 | 230mm |



Kocher Artery Forceps Curved 1 into 2 Teeth

| | |
|---------|-------|
| PH90691 | 130mm |
| PH90692 | 140mm |
| PH90693 | 150mm |
| PH90694 | 180mm |
| PH90695 | 200mm |
| PH90696 | 230mm |



Baby Mixer Artery Forceps Curved

| | |
|---------|-------|
| PH90708 | 140mm |
| PH90709 | 180mm |



Mixer Artery Forceps Curved

| | |
|---------|-------|
| PH90712 | 200mm |
| PH90713 | 230mm |
| PH90714 | 280mm |



Mixer O'Shaughnessy Artery Forceps Curved

| | |
|---------|-------|
| PH90716 | 150mm |
| PH90717 | 190mm |



Moynihan Artery Forceps Curved

| | |
|---------|-------|
| PH90722 | 145mm |
| PH90723 | 180mm |




Negus Short Curve Artery Forceps



PH90730 190mm

Negus Long Curve Artery Forceps



PH90731 190mm

Rochester Pean Artery Forceps Straight




| | |
|---------|-------|
| PH90760 | 130mm |
| PH90761 | 140mm |
| PH90762 | 160mm |
| PH90763 | 180mm |
| PH90764 | 190mm |
| PH90765 | 200mm |
| PH90766 | 230mm |

Rochester Pean Artery Forceps Curved




| | |
|---------|-------|
| PH90770 | 130mm |
| PH90771 | 140mm |
| PH90772 | 160mm |
| PH90773 | 180mm |
| PH90774 | 190mm |
| PH90775 | 200mm |
| PH90776 | 230mm |

Robert Artery Forceps Curved



PH90804 230mm

Robert Artery Forceps Curved with Screw Joint



| | |
|---------|-------|
| PH90805 | 230mm |
| PH90815 | 300mm |

Sawtell Artery Forceps Curved

PH90807 180mm



Wilson (Tonsil) Artery Forceps Curved with Angled Shanks

PH90809 200mm



Did you know we offer instrument repair services?

We offer flexible and comprehensive services including a rapid repair service. All our services are ISO 13485 accredited with warranties provided for all repairs.

Call us today to find out more
(0)3452 414818



Spencer Wells Artery Forceps Straight

PH90836 130mm
PH90837 150mm
PH90838 180mm
PH90839 200mm
PH90840 230mm
PH90841 300mm



Spencer Wells Artery Forceps Curved

PH90847 150mm
PH90848 180mm
PH90849 200mm
PH90850 230mm
PH90851 300mm




Schmidt Artery Forceps Curved

PH90860 190mm




Sarrot Artery Forceps Straight




PH90861 240mm

Sarrot Artery Forceps Curved




PH90862 240mm

Belcher Gigli Saw Passing Forceps Curved



PH90870 300mm

Lloyd Davies Artery Forceps Straight



PH90880 250mm

Kilner Artery Forceps Straight



PH90890 145mm

Kilner Artery Forceps Curved



PH90891 145mm

Dissecting Forceps Bayonet Shanks and Serrated Jaws

| | |
|----------|-------|
| PH110940 | 180mm |
| PH110941 | 200mm |




Dissecting Forceps Turned Over End and Serrated Jaws

| | |
|----------|-------|
| PH110950 | 130mm |
| PH110951 | 150mm |
| PH110952 | 180mm |
| PH110953 | 200mm |
| PH110954 | 230mm |
| PH110955 | 250mm |
| PH110956 | 300mm |



Dissecting Forceps Block End and Serrated Jaws Medium Point

| | |
|----------|-------|
| PH110960 | 130mm |
| PH110961 | 150mm |
| PH110962 | 180mm |
| PH110963 | 200mm |
| PH110964 | 230mm |
| PH110965 | 250mm |
| PH110966 | 300mm |



Continental Pattern Dissecting Forceps Serrated Jaws and Narrow End

| | |
|----------|-------|
| PH110967 | 130mm |
| PH110968 | 140mm |
| PH110969 | 160mm |
| PH110970 | 180mm |
| PH110971 | 200mm |
| PH110972 | 250mm |



Continental Pattern Dissecting Forceps 1 into 2 Teeth and Broad End

| | |
|----------|-------|
| PH110982 | 130mm |
| PH110983 | 140mm |
| PH110984 | 160mm |
| PH110985 | 180mm |
| PH110986 | 200mm |
| PH110987 | 250mm |
| PH110988 | 300mm |



Adson Micro Dissecting Forceps Serrated Jaws

| | |
|----------|-------|
| PH111016 | 130mm |
|----------|-------|



Adson Micro Toothed Dissecting Forceps 1 into 2 Teeth

PH111017 130mm



Adson Browne Dissecting Forceps Serrated Jaws

PH111019 120mm



Adson Dissecting Forceps Serrated Jaws

PH111020 120mm

PH111022 180mm



Adson Dissecting Forceps 1 into 2 Teeth

PH111021 130mm

PH111023 180mm



Double Action Adson Dissecting Forceps 1-2 Teeth

PH111021DA 130mm



Bickford Dissecting Forceps With Serrated Jaws

PH111032 230mm



Bickford Dissecting Forceps 1 into 2 Teeth

PH111033 230mm



Bonney Dissecting Forceps Serrated Jaws

PH111035 180mm



Bonney Dissecting Forceps 1 into 2 Teeth

PH111036 180mm



Cairn Dissecting Forceps 2 into 3 Teeth

PH111045 165mm

PH111046 205mm



Charnley Suture Dissecting Forceps with Shutter and Slot

PH111051 145mm

PH111052 180mm



College Pattern Forceps For Dental Tweezer Curved

PH111057 150mm



We also offer decontamination and sterilisation services.

With over 20 years' experience and more than 100 million instruments processed each year, we can offer a fully outsourced service or contingency for your sterile services department.



Contact us today to find out more
+44 (0)3452 414818



Debaquey Tissue Dissecting Forceps Atraumatic Jaws



| | |
|----------|-------|
| PH111059 | 150mm |
| PH111060 | 180mm |
| PH111061 | 200mm |
| PH111062 | 230mm |

Dott Dissecting Forceps Serrated Jaws Straight



| | |
|----------|-------|
| PH111064 | 180mm |
|----------|-------|

Emmet Dissecting Forceps with Serrated Jaws Straight



| | |
|----------|-------|
| PH111175 | 200mm |
| PH111176 | 240mm |

Gillies Dissecting Forceps Serrated Jaws

PH111179 150mm



Gillies Dissecting Forceps 1 into 2 Teeth

PH111180 150mm



Iris Dissecting Forceps Serrated Jaws Straight

PH111185 100mm

PH111189 115mm



Iris Dissecting Forceps 1 into 2 Teeth Straight

PH111186 100mm

PH111190 115mm



Iris Dissecting Forceps Serrated Jaws Curved

PH111187 100mm

PH111191 115mm



Jackson Burrows Dissecting Forceps Serrated 1 into 2 Teeth Straight

PH111195 140mm



Jean Dissecting Forceps 1 into 2 Teeth Straight



| | |
|----------|-------|
| PH111202 | 140mm |
| PH111203 | 180mm |

Jefferson Dissecting Forceps 1 into 2 Teeth Straight



| | |
|----------|-------|
| PH111211 | 180mm |
|----------|-------|

Kurihara Endo Cervical Dissecting Forceps Satin Finish Straight



| | |
|----------|-------|
| PH111220 | 250mm |
|----------|-------|

Lane Dissecting Forceps 1 into 2 Teeth Straight



| | |
|----------|-------|
| PH111240 | 130mm |
| PH111241 | 140mm |
| PH111242 | 150mm |
| PH111243 | 180mm |
| PH111244 | 200mm |

Lane Dissecting Forceps 2 into 3 Teeth Straight



| | |
|----------|--------|
| PH111245 | 140mm |
| PH111246 | 150 mm |
| PH111247 | 180mm |
| PH111248 | 200mm |

Lane Dissecting Forceps 3 into 4 Teeth Straight



| | |
|----------|-------|
| PH111249 | 150mm |
| PH111250 | 180mm |
| PH111251 | 200mm |


McIndoe Dissecting Forceps Serrated Jaws

PH111264 150mm



McIndoe Dissecting Forceps 1 into 2 Teeth

PH111265 150mm



Moorfield Conjunctival Dissecting Forceps Straight

PH111268 115mm



Officer Dissecting Forceps 3 into 4 Teeth Straight

PH111280 200mm

PH111281 250mm



Ramsay Dissecting Forceps Serrated Jaws Straight

PH111290 180mm



Ramsay Dissecting Forceps 1 into 2 Teeth Straight

PH111291 180mm



Ramsay Dissecting Forceps 2 into 3 Teeth Straight



PH111292 180mm

Russian Pattern Dissecting Forceps Oval Toothed End Straight



PH111295 150mm

PH111296 200mm

PH111297 250mm


PH111298 180mm

St Marks Hospital Dissecting Forceps Serrated Jaws Straight



PH111301 290mm

Treves Dissecting Forceps 1 into 2 Teeth



PH111310 130mm

PH111311 150mm

PH111312 180mm

PH111313 200mm

Waugh's Dissecting Forceps Serrated Jaws Straight



PH111315 150mm

PH111316 180mm

PH111317 200mm

Waugh's Dissecting Forceps 1 into 2 Teeth Straight



PH111318 150mm

PH111319 180mm

PH111320 200mm

Adson Dissecting Forceps Tungsten Carbide Jaws Straight

| | |
|----------|-------|
| PH111330 | 120mm |
|----------|-------|



Adson Dissecting Forceps TC Jaws 1 into 2 Teeth

| | |
|----------|-----------|
| PH111331 | 120mm STR |
|----------|-----------|



Debaquey Dissecting Forceps With Tungsten Carbide Jaws Straight

| | |
|----------|-------|
| PH111350 | 150mm |
| PH111351 | 200mm |
| PH111352 | 230mm |



Martin Splinter Forceps Straight

| | |
|----------|-------|
| PH131405 | 75mm |
| PH131406 | 115mm |
| PH131407 | 130mm |



Lister Sinus Forceps Straight

| | |
|----------|-------|
| PH141455 | 130mm |
| PH141456 | 150mm |
| PH141457 | 180mm |
| PH141458 | 200mm |



Ramplsey Dressing And Sponge Holding Forceps Box Joint Straight

| | |
|----------|-------|
| PH161482 | 180mm |
| PH161483 | 240mm |

Denis Browne Suture Forceps Collar Bead Crusher

| | |
|----------|-------|
| PH161866 | 145mm |
|----------|-------|

Guys Galabin/Suture Forceps

| | |
|----------|-------|
| PH161867 | 140mm |
|----------|-------|

Tudor Edwards Scapula Retractor With 45mm Wide Blade x 50mm Deep

| | |
|----------|-------|
| PH252499 | 160mm |
|----------|-------|



Did you know we offer instrument repair services?

We offer flexible and comprehensive services including a rapid repair service. All our services are ISO 13485 accredited with warranties provided for all repairs.



Contact us today to find out more
+44 (0)3452 414818

Beale Aural Forceps Crocodile Action with Fine Serrated Jaw

PH342850 70mm Shoulder Length



Cawthorne Aural Forceps Crocodile Action with Fine Smooth Jaw

PH342851 70mm Shoulder Length



Hartmann Aural Forceps Crocodile Action with Extra Fine Serrated Jaw

PH342852 70mm Shoulder Length



Ormerod Aural Forceps Crocodile Action with Extra Fine Elongated Cup

PH342860 70mm Shoulder Length



Ormerod Aural Forceps Crocodile Action with Extra Fine Cup Jaw

PH342861 70mm Shoulder Length



Shea Aural Forceps Crocodile Action Straight with Oval Cup Jaw

PH342862 70mm Shoulder Length





Shea Aural Forceps Crocodile Action with Oval Cup Jaw Curved Upwards

PH342863 70mm Shoulder Length



Shea Aural Forceps Crocodile Action with Oval Cup Jaw Curved to the Right


PH342864 70mm Shoulder Length

Did you know we offer instrument repair services?

We offer flexible and comprehensive services including a rapid repair service. All our services are ISO 13485 accredited with warranties provided for all repairs.

Call us today to find out more
(0)3452 414818

 **STERIS**




Tilley Aural Dressing Forceps

PH342879 130mm



Wilde Aural Dressing Forceps

PH342880 130mm




Hartmann Crocodile Action Aural Forceps

PH342883 85mm

Hartmann Aural Forceps Crocodile Action 1 x 2 Teeth

PH342884 85mm Shoulder Length



Henckel Crocodile Action Aural Forceps

PH342887 75mm Shoulder Length




Dental Cartridge Syringe 2.2cc

PH383083



Tilley Aural Dressing Forceps

PH413176



Weil Blakley Nasal Turbinate Forceps Straight



PH413179 4mm Bite

Weil Blakley Nasal Turbinate Forceps 40 Degree Up



PH413180 4mm Bite

Weil Blakley Nasal Turbinate Forceps 90 Degree Up



PH413181 4mm Bite

Irwin Moore Nasal Turbinate Forceps



PH423211 190mm

Tilley Henckel Nasal Turbinate Forceps with Ethmoid Punches



| | |
|----------|-------|
| PH423215 | 4.5mm |
| PH423216 | 5.5mm |
| PH423217 | 6.5mm |
| PH423218 | 8mm |

Jansen Middleton Nasal Turbinate Septum Forceps with Double Action

PH423220 190mm



Luc Nasal Turbinate Forceps

| | | |
|----------|---------------------|-------|
| PH423221 | Oval Small (6mm) | 180mm |
| PH423222 | Oval Medium (7.5mm) | 180mm |
| PH423223 | Oval Large (9mm) | 180mm |



Walsham Nasal Septum Forceps

| | |
|----------|----------------|
| PH433255 | 220mm |
| PH433256 | Left 220mm |
| PH433257 | Right 220mm |



Did you know we offer instrument repair services?

We offer flexible and comprehensive services including a rapid repair service. All our services are ISO 13485 accredited with warranties provided for all repairs.

Call us today to find out more
(0)3452 414818

 **STERIS**

Denis Browne Tonsil Seizing Forceps

PH473401 215mm




Lahey Cholecystectomy Forceps Curved



| | |
|----------|-------|
| PH523745 | 190mm |
| PH523746 | 230mm |

Mixer Cholecystectomy Forceps Curved



| | |
|----------|-------|
| PH523759 | 160mm |
| PH523760 | 190mm |

Moynihan Cholecystectomy Forceps



| | |
|----------|-------|
| PH523763 | 230mm |
|----------|-------|

O'Shaughnessy Cholecystectomy Forceps



| | |
|----------|-------|
| PH523765 | 200mm |
|----------|-------|

Duval Lung Forceps

| | | |
|----------|----------------|-------|
| PH533801 | 13mm Wide Jaws | 180mm |
| PH533802 | 13mm Wide Jaws | 200mm |
| PH533803 | 19mm Wide Jaws | 200mm |
| PH533804 | 19mm Wide Jaws | 230mm |

| | | |
|----------|----------------|-------|
| PH533805 | 25mm Wide Jaws | 200mm |
| PH533806 | 25mm Wide Jaws | 230mm |
| PH533807 | 14mm Wide Jaws | 230mm |



Liston Bone Cutting Forceps Straight

| | |
|----------|-------|
| PH574070 | 130mm |
| PH574071 | 150mm |
| PH574072 | 180mm |
| PH574073 | 200mm |
| PH574074 | 230mm |



McIndoe Bone Cutting Forceps Angled on Flat with Compound Action

| | |
|----------|-------|
| PH574112 | 170mm |
| PH574113 | 190mm |



Ruskin Liston Bone Cutting Forceps Straight with Compound Action

| | |
|----------|-------|
| PH574115 | 150mm |
| PH574116 | 180mm |
| PH574118 | 200mm |



Ruskin Liston Bone Cutting Forceps Curved with Compound Action

| | |
|----------|-------|
| PH574119 | 150mm |
| PH574120 | 180mm |
| PH574121 | 200mm |



Stamm Bone Cutting Forceps

| | |
|----------|-------|
| PH574126 | 150mm |
|----------|-------|



Horsley Bone Cutting Forceps Curved with Compound Action

| | |
|----------|-------|
| PH574136 | 270mm |
|----------|-------|



Officer Rectal Biopsy Forceps with 8mm Cup Crocodile Action

PH675436 430mm Effective Shaft Length



Chevalier Jackson Endoscope Forceps with Crocodile Action and Serrated Jaw

| | |
|----------|------------------------------|
| PH685440 | 300mm Effective Shaft Length |
| PH685441 | 400mm Effective Shaft Length |
| PH685442 | 500mm Effective Shaft Length |
| PH685443 | 600mm Effective Shaft Length |




Chevalier Jackson Vocal Nodule Endoscope Forceps with Crocodile Action and Curved Up

| | |
|----------|------------------------------|
| PH685470 | 300mm Effective Shaft Length |
| PH685471 | 400mm Effective Shaft Length |
| PH685472 | 500mm Effective Shaft Length |
| PH685473 | 600mm Effective Shaft Length |




Chevalier Jackson Vocal Nodule Endoscope Forceps with Crocodile Action and Straight

| | |
|----------|------------------------------|
| PH685475 | 300mm Effective Shaft Length |
| PH685476 | 400mm Effective Shaft Length |
| PH685477 | 500mm Effective Shaft Length |
| PH685478 | 600mm Effective Shaft Length |



Purse String Applicator

PH252508 250mm



Purse String Applicator


PH252509 225mm



We also offer repairs of both flexible and rigid endoscopes

At STERIS we also offer repairs of both flexible and rigid scopes from our UK based repair facilities. All our fully accredited repairs are carried out in-house by experienced technicians. With us, your devices are in skilled and safe hands.

Contact us today to find out how we can help you reduce cost and improve turnarounds.

 **+44 (0)3452 414818**

 **steris-ims-instruments.com**

 **STERIS**

Teale Uterine Vulsellum Forceps CVD with 3 into 4 Teeth and Box Joint

PH725903 180mm

PH725904 230mm



Schroeder Uterine Vulsellum Forceps with 2 into 2 Teeth and Box Joint

PH725906 240mm STR



Luer Uterine Vulsellum Forceps with 1 into 1 Teeth and Box Joint

PH725911 240mm STR




Pratt Vulsellum Forceps with Teeth and Screw Joint Straight

PH725925 280mm



Lewer Vulsellum Forceps with 2 into 2 Teeth and Box Joint Curved

PH725932 240mm




Rotunda Hospital Pattern Vulsellum Forceps with 1 into 1 Teeth and Screw Joint Straight

PH725934 250mm




Schroeder Tenaculum Forceps with 1 into 1 Teeth and Box Joint Straight




PH725964 250mm

Green Armytage Uterine Dressing Forceps with Wide Serrated Jaws



PH725984 210mm

Green Armytage Uterine Dressing Forceps Angled Ends with Wide Serrated Jaws



PH725986 210mm

Bonney Uterine Polypus Forceps with Fenestrated Serrated Jaws and a Box Joint



PH726000 240mm

Clamps, Clips and Pins

Presbyterian Tubing Clamp Forceps with Serrated Jaws Straight

| | |
|----------|-------|
| PH100901 | 165mm |
| PH100902 | 195mm |



St Bartholomews Forceps for Tubing Clamp with Square Jaws Straight

| | |
|----------|-------|
| PH100903 | 150mm |
|----------|-------|



Square Jaw Tubing Clamp Straight

| | |
|----------|-------|
| PH100904 | 130mm |
| PH100906 | 180mm |
| PH100908 | 200mm |



Square Jaw Tubing Clamp Angled

| | |
|----------|-------|
| PH100905 | 115mm |
| PH100907 | 180mm |





X-Ray Bag Clip

PH100915



Zimmer Clip

PH100917



Instrument Pin

| | |
|----------|-------|
| PH100925 | 115mm |
| PH100926 | 130mm |
| PH100927 | 140mm |
| PH100928 | 150mm |

Backhaus Towel Clip Curved

| | |
|----------|-------|
| PH171550 | 90mm |
| PH171551 | 115mm |
| PH171552 | 135mm |



Mayo (Backhaus) Towel Clip Curved

| | |
|----------|-------|
| PH171553 | 135mm |
|----------|-------|



Towel Clip Non-Tearing Ball and Socket Curved

| | |
|----------|-------|
| PH171558 | 90mm |
| PH171559 | 100mm |
| PH171560 | 115mm |
| PH171561 | 130mm |



Towel Clip Non-Tearing Serrated Jaws Curved

| | |
|----------|-------|
| PH171562 | 100mm |
| PH171563 | 115mm |
| PH171564 | 140mm |



Towel Clip English Pattern Cross Action

| | |
|----------|------|
| PH171565 | 90mm |
|----------|------|



Schaedel Towel Clip Cross Action

| | |
|----------|------|
| PH171567 | 90mm |
|----------|------|





Lorna (Edna) Towel Clip Curved

PH171574 140mm



Robin Anchoring Towel Clip
130mm Curved

PH171585 5mm

PH171587 13mm

PH171589 25mmr

Did you know we offer instrument repair services?

We offer flexible and comprehensive services including a rapid repair service. All our services are ISO 13485 accredited with warranties provided for all repairs.



Contact us today to find out more
+44 (0)3452 414818

 **STERIS**

Needle Holders

Baumgartner Needle Holder

PH181650 140mm



Crile Wood Needle Holder

PH181684 150mm



Derf Needle Holder

PH181691 120mm



Foster Needle Holder And Scissors

PH181706 135mm



Fry Needle Holder

PH181707 115mm



Gillies Needle Holder with Right Hand Scissor

PH181710 160mm





Halsey Needle Holder
with Plain Jaw

PH181712 130mm



Halsey Needle Holder
with Serrated Jaw

PH181713 130mm

Did you know we
offer instrument
repair services?

We offer flexible and comprehensive services including a rapid repair service. All our services are ISO 13485 accredited with warranties provided for all repairs.

Call us today to find out more
(0)3452 414818

 **STERIS**



Kilner Needle Holder

PH181716 135mm



Lawrence Needle Holder

PH181722 150mm



Lawson Tait Needle Holder

PH181724 130mm


Mayo Hegar Needle Holder

| | | | |
|----------|-------|----------|-------|
| PH181770 | 130mm | PH181775 | 190mm |
| PH181771 | 140mm | PH181776 | 200mm |
| PH181772 | 150mm | PH181777 | 240mm |
| PH181773 | 165mm | PH181778 | 300mm |
| PH181774 | 180mm | | |



Naunton-Morgan Needle Holder with Multiple Action

| | |
|----------|-------|
| PH181785 | 180mm |
| PH181786 | 250mm |



Nievert Needle Holder with Offset Ring For Thumb

| | |
|----------|-------|
| PH181787 | 130mm |
|----------|-------|



Olsen-Hegar Needle Holder with Scissors

| | |
|----------|-------|
| PH181790 | 135mm |
| PH181791 | 160mm |



Thomson Walker Needle Holder Curved Shanks

| | |
|----------|-------|
| PH181821 | 200mm |
|----------|-------|





Wright Needle Holder

PH181828 130mm



Baumgartner Needle Holder

PH191550 135mm



Bozemann Needle Holder with Tungsten Carbide Jaws

PH191555 200mm



Crile Wood Needle Holder with Tungsten Carbide Jaws

| | | | |
|----------|-------|----------|-------|
| PH191560 | 140mm | PH191564 | 230mm |
| PH191561 | 150mm | PH191565 | 250mm |
| PH191562 | 180mm | PH191566 | 300mm |
| PH191563 | 200mm | | |



DeBakey Needle Holder with Tungsten Carbide Jaws

| | |
|----------|-------|
| PH191570 | 150mm |
| PH191571 | 180mm |
| PH191572 | 200mm |
| PH191573 | 230mm |
| PH191574 | 250mm |
| PH191575 | 300mm |



Derf Needle Holder with Tungsten Carbide Jaws

PH191585 130mm

Kilner Needle Holder with Tungsten Carbide Jaws

PH191590 135mm



Foster Gillies Needle Holder with Tungsten Carbide Jaws and Box Joint Angled

PH191595 135mm



Fry Needle Holder with Tungsten Carbide Jaws

PH191598 130mm



Gillies Needle Holder Tungsten Carbide Jaws and Box Joint with Right Hand Scissor

PH191600 130mm



Halsey Needle Holder with Tungsten Carbide Jaws


PH191605 130mm



Heaney Needle Holder with Tungsten Carbide Jaws

PH191610 200mm

PH191611 250mm



Lawrence Needle Holder with Tungsten Carbide Jaws



PH191613 150mm

Mason Needle Holder with Tungsten Carbide Jaws



PH191615 265mm

Mayo Hegar Needle Holder with Tungsten Carbide Jaws Box joint 205mm Straight



| | | | |
|----------|-------|----------|-------|
| PH191620 | 130mm | PH191624 | 200mm |
| PH191621 | 140mm | PH191625 | 230mm |
| PH191622 | 150mm | PH191626 | 250mm |
| PH191623 | 180mm | PH191627 | 300mm |

Naunton Morgan Needle Holder with Tungsten Carbide Jaws



PH191640 250mm

Nievert Needle Holder with Tungsten Carbide Jaws and Offset Ring



PH191645 130mm

Olsen Hegar Needle Holder with Tungsten Carbide Jaws and Scissors



| | |
|----------|-------|
| PH191650 | 140mm |
| PH191651 | 160mm |
| PH191652 | 180mm |

Ryder Needle Holder French Eye with Tungsten Carbide Jaws

| | |
|----------|-------|
| PH191660 | 130mm |
| PH191661 | 150mm |
| PH191662 | 180mm |
| PH191663 | 200mm |
| PH191664 | 230mm |



Ryder Needle Holder Micro French Eye with Tungsten Carbide Jaws and Box Joint Straight

| | |
|----------|-------|
| PH191675 | 130mm |
| PH191676 | 150mm |
| PH191677 | 180mm |
| PH191678 | 200mm |
| PH191679 | 230mm |



Sternum Needle Holder/Wire Twister with Tungsten Carbide Heavy Jaws

| | |
|----------|-------|
| PH191705 | 130mm |
| PH191706 | 180mm |
| PH191707 | 200mm |



Stratte Needle Holder with Tungsten Carbide Jaws and Curved Shanks

| | |
|----------|-------|
| PH191715 | 230mm |
|----------|-------|



Webster Needle Holder with Tungsten Carbide

| | |
|----------|-------|
| PH191750 | 130mm |
|----------|-------|



Castroviejo Needle Holder with Tungsten Carbide Straight Smooth Jaws and Centre Lock

PH665257 140mm



Castroviejo Needle Holder with Tungsten Carbide Curved Smooth Jaws and Centre Lock

PH665258 140mm



Castroviejo Needle Holder with Tungsten Carbide Cross Serrated Straight Jaws and Centre Lock

PH665259 140mm



Castroviejo Needle Holder with Tungsten Carbide Cross Serrated Curved Jaws and Centre Lock

PH665260 140mm



Did you know we offer instrument repair services?

We offer flexible and comprehensive services including a rapid repair service. All our services are ISO 13485 accredited with warranties provided for all repairs.



Contact us today to find out more
+44 (0)3452 414818



Anuerysm Needles, Skin Hooks, Probes and Fistula Directors

Aneurysm Needle with Flat Handle and Curved Laterally

PH201885 140mm



Aneurysm Needle with Oval Handle and Curved Laterally

PH201886 140mm



Negus Knot Tier or Pusher

PH201910 200mm



Symes Aneurysm Needle Flat Handle and Curved Laterally

PH201915 160mm



Probe with Eye, Silver Plated

PH212000 130mm

PH212001 150mm

PH212002 180mm

PH212003 200mm

PH212004 250mm

PH212005 300mm



Probe with Eye, Stainless Steel

PH212010 130mm

PH212011 150mm

PH212012 180mm


PH212013 200mm

PH212014 250mm

PH212015 300mm



Double Ended Probe, Stainless Steel




| | |
|----------|-------|
| PH212020 | 130mm |
| PH212021 | 150mm |
| PH212022 | 180mm |
| PH212023 | 200mm |
| PH212024 | 250mm |

St Bartholomews Wax Hook 180mm Angled



| | |
|----------|-------|
| PH212025 | 180mm |
|----------|-------|

Jobson Horne Probe with Ring



| | |
|----------|-------|
| PH212029 | 140mm |
| PH212030 | 180mm |

Leibrichs Probe



| | | | |
|----------|-----------------|----------|------------|
| PH212031 | Size 000 - 0000 | PH212035 | Size 1 - 3 |
| PH212032 | Size 00 - 0 | PH212036 | Size 4 - 5 |
| PH212033 | Size 1 - 2 | PH212037 | Size 5 - 6 |
| PH212034 | Size 2 - 3 | PH212038 | Size 3 - 4 |

Bowman Lachrymal Probe



| | |
|----------|-----------------|
| PH232081 | 0.6mm and 0.7mm |
| PH232082 | 0.7mm and 0.8mm |
| PH232083 | 0.8mm and 0.9mm |
| PH232084 | 0.9mm and 1.1mm |
| PH232085 | 1.1mm and 1.3mm |
| PH232086 | 1.3mm and 1.4mm |
| PH232087 | 1.4mm and 1.5mm |
| PH232088 | 1.5mm and 1.6mm |
| PH232089 | 1.6mm and 1.8mm |
| PH232090 | 1.8mm and 1.9mm |
| PH232091 | 0.9mm and 1.0mm |
| PH232092 | 1.1mm and 1.2mm |
| PH232093 | 1.7mm and 1.8mm |

Brodie Fistula Director Grooved

| | |
|----------|-------|
| PH222050 | 130mm |
| PH222051 | 140mm |
| PH222052 | 160mm |



Brodie Fistula Director Grooved with Probe End

| | |
|----------|-------|
| PH222053 | 130mm |
| PH222054 | 140mm |
| PH222055 | 60mm |



Lockhart Mummery Fistula Director Straight (No.1)

| | |
|----------|-------|
| PH222060 | 155mm |
|----------|-------|



Lockhart Mummery Half Curved Fistula Director (No.2, 45° Curve)

| | |
|----------|-------|
| PH222061 | 150mm |
|----------|-------|



Lockhart Mummery Curved Fistula Director (No.3, 90° Curve)

| | |
|----------|-------|
| PH222062 | 140mm |
|----------|-------|



Lockhart Mummery Full Curved Fistula Director (No.4, Retrograde)

| | |
|----------|-------|
| PH222063 | 160mm |
|----------|-------|



Gillies Skin Hook Fine Point 165mm

GILLIES SKIN HOOK FINE POINT 165MM

| | |
|----------|--------|
| PH262500 | Small |
| PH262501 | Large |
| PH262549 | Medium |



Kilner Detachable Skin Hook and Handle

KILNER DETACHABLE SKIN HOOK AND HANDLE

| | |
|----------|-------|
| PH262502 | 155mm |
|----------|-------|



McIndoe Double Prong Skin Hooks Wide End 150mm

MCINDOE DOUBLE PRONG SKIN HOOKS WIDE END 150MM


| | |
|----------|------|
| PH262505 | 9mm |
| PH262506 | 10mm |
| PH262507 | 11mm |



Frazier Small Skin Hook

FRAZIER SMALL SKIN HOOK

| | |
|----------|-------|
| PH262510 | 180mm |
|----------|-------|



Retractors


Charnley Pin Retractor T-Handle Complete with Pins and Chain

PH252007



Teales Gorget Small Tapered with Probe from 11mm to 4.5mm

PH252008



Teales Gorget Large Tapered with Probe from 17mm to 7mm

PH252009



Did you know we
offer instrument
repair services?

We offer flexible and comprehensive services including a rapid repair service. All our services are ISO 13485 accredited with warranties provided for all repairs.

Call us today to find out more
(0)3452 414818

 **STERIS**

Awty Channel Retractors


PH252018 12mm Blade Width

PH252019 17mm Blade Width



Baby Hohman Bone Levers 6mm


PH252022



STERIS IMS 03452 414818

Ronald Edwards Retractor with a Rack Measuring 450mm and Bladed Measuring 120mm Deep and 75mm Wide

PH252025



Hohmann Muller Spike 240mm



| | |
|----------|------------------|
| PH252100 | 8mm Blade Width |
| PH252101 | 18mm Blade Width |
| PH252102 | 43mm Blade Width |
| PH252103 | 70mm Blade Width |

Retractor Single Hook Straight



| | | |
|----------|-------|-------|
| PH252125 | 160mm | Sharp |
| PH252126 | 160mm | Blunt |

Retractor Double Hook Straight



| | | |
|----------|-------|-------|
| PH252127 | 160mm | Sharp |
| PH252128 | 160mm | Blunt |

Retractor Hook and 1 Prong Straight



| | | |
|----------|-------|-------|
| PH252129 | 160mm | Sharp |
| PH252133 | 160mm | Blunt |

Retractor Hook and 2 Prongs Straight

| | | |
|----------|-------|-------|
| PH252130 | 160mm | Sharp |
| PH252134 | 160mm | Blunt |



Retractor Hook and 3 Prongs Straight

| | | |
|----------|-------|-------|
| PH252131 | 160mm | Sharp |
| PH252135 | 160mm | Blunt |



Retractor Hook and 4 Prongs Straight

| | |
|----------|-------|
| PH252132 | 160mm |
|----------|-------|



Retractor Blunt with Hook And 4 Prongs Straight

| | |
|----------|-------|
| PH252136 | 160mm |
|----------|-------|



Austin Kay Right Angled Retractor with Lip Angled

| | |
|----------|-------|
| PH252138 | 150mm |
|----------|-------|



Allison Retractor

| | | |
|----------|-------|-------------|
| PH252141 | 280mm | Infant size |
| PH252142 | 280mm | Child size |
| PH252143 | 280mm | Adult size |



Aufricht Nasal Retractor with 50mm Blade and Chrome Plated Handle Angled



PH252144 160mm

Cairn (Rake) Scalp Retractor 3 Prongs 16mm Sharp Straight



PH252145 210mm

Cairn (Rake) Scalp Retractor 4 Prongs 19mm Sharp Straight



PH252146 210mm

Czerny Retractor Double Ended



PH252160 155mm

Deaver Retractor with Hollow Handle



| | |
|----------|--------------------|
| PH252170 | 18mm Blade x 180mm |
| PH252171 | 25mm Blade x 230mm |
| PH252172 | 25mm Blade x 250mm |
| PH252173 | 38mm Blade x 300mm |
| PH252174 | 50mm Blade x 300mm |
| PH252175 | 75mm Blade x 300mm |

Deaver Retractor (Without Handle) Plain Strip



| | |
|----------|--------------------|
| PH252176 | 25mm Blade x 230mm |
| PH252177 | 25mm Blade x 250mm |
| PH252178 | 38mm Blade x 300mm |
| PH252179 | 50mm Blade x 300mm |
| PH252180 | 75mm Blade x 300mm |

Doyen Abdominal Retractor 215mm

| | |
|----------|--------------------------|
| PH252185 | Small 90mm Blade Width |
| PH252186 | Medium 100mm Blade Width |
| PH252187 | Large 115mm Blade Width |




Doyen Retractor Ring Handles 215mm

| | |
|----------|--------------------------|
| PH252188 | Small 95mm width blade |
| PH252189 | Medium 105mm width blade |
| PH252190 | Large 115mm width blade |



Desmarres Retractor (Eye Lid) 160mm

| | |
|----------|-------------------|
| PH252193 | Small 9mm x 9mm |
| PH252194 | Medium 9mm x 12mm |
| PH252195 | Large 11mm x 14mm |



Bowdler Henry Rake Retractor for Lower Wisdom Teeth

| | |
|----------|-------|
| PH252240 | 190mm |
|----------|-------|



O'Connell Double Ended Retractor for Nerve Root

| | |
|----------|-------|
| PH252251 | 215mm |
|----------|-------|



Ward Double Ended Cheek Retractor used in the Removal of Wisdom Teeth

| | |
|----------|-------|
| PH252256 | 150mm |
|----------|-------|



Kelly Retractor



| | |
|----------|-----------------------------|
| PH252259 | Small Size 38mm Blade 190mm |
| PH252260 | Large Size 57mm Blade 190mm |

Kilner Cat Paw Retractor



| | |
|----------|-------|
| PH252261 | 150mm |
|----------|-------|

Kilner Cheek Retractor



| | |
|----------|-----------------------------------|
| PH252263 | Small 25mm Tapering to 20mm 125mm |
| PH252264 | Large 30mm Tapering to 28mm 125mm |

Kilner ALAE Retractor



| | |
|----------|----------------|
| PH252265 | 8mm Wide 90mm |
| PH252266 | 10mm Wide 90mm |
| PH252267 | 13mm Wide 90mm |

Kocher Retractor



| | | |
|----------|-----------|-------|
| PH252269 | 10mm Wide | 215mm |
| PH252270 | 15mm Wide | 215mm |
| PH252271 | 20mm Wide | 215mm |

Laster Cheek Retractor



| | |
|----------|-------|
| PH252278 | 100mm |
|----------|-------|

Langenbeck Retractor 215mm

| | |
|----------|-------------------------|
| PH252285 | Small 23mm x 7mm |
| PH252286 | Medium 44mm x 13mm |
| PH252287 | Large 65mm x 25mm |
| PH252288 | Extra Large 90mm x 20mm |



Morris Retractor Single Ended 230mm

| | |
|----------|-----------------------|
| PH252300 | 38mm Wide x 50mm Deep |
| PH252301 | 50mm Wide x 52mm Deep |
| PH252302 | 65mm Wide x 50mm Deep |
| PH252303 | 50mm Wide x 50mm Deep |
| PH252304 | 50mm Wide x 60mm Deep |
| PH252305 | 50mm Wide x 65mm Deep |
| PH252306 | 50mm Wide x 70mm Deep |
| PH252307 | 65mm Wide x 70mm Deep |



Morris Retractor Double Ended 230mm

| | |
|-----------|-----------------------------------|
| PH252308 | 32mm Wide x 38mm Wide x 57mm Deep |
| PH252309 | 4mm Wide x 50mm Wide x 57mm Deep |
| PH252310 | 57mm Wide x 64mm Wide x 57mm Deep |
| PH252310A | 30mm Wide x 40mm Wide x 70mm Deep |
| PH252310B | 45mm Wide x 50mm Wide x 70mm Deep |
| PH252310C | 60mm Wide x 65mm Wide x 70mm Deep |



We also offer decontamination and sterilisation services.

With over 20 years' experience and more than 100 million instruments processed each year, we can offer a fully outsourced service or contingency for your sterile services department.



Contact us today to find out more
+44 (0)3452 414818

Park Anal and Sphincterotomy Combination Retractor

PH252333 Adult Size 93mm Blade Length



Ragnel Davis Retractor

PH252338 150mm



St Marks and Lloyd Davis Retractor with Lip

PH252348 Blade Length 130mm x 60mm Width

PH252349 Blade Length 170mm/64mm x 51mm Width



St Marks and Lloyd Davis Retractor

PH252350 Blade Length 130mm x 60mm Width



St Marks and Lloyd Davis Retractor

PH252351 Blade Length 170mm/64mm x 51mm Width



Volkman Retractor Sharp 1 Prong

| | | |
|----------|-------|-------|
| PH252352 | 215mm | Sharp |
| PH252360 | 215mm | Blunt |



Volkman Retractor Sharp 2 Prong

| | | |
|----------|-------|-------|
| PH252353 | 215mm | Sharp |
| PH252361 | 215mm | Blunt |



Volkman Retractor 3 Prongs

| | | |
|----------|-------|-------|
| PH252354 | 215mm | Sharp |
| PH252362 | 215mm | Blunt |



Volkman Retractor 4 Prongs

| | | |
|----------|-------|-------|
| PH252355 | 215mm | Sharp |
| PH252363 | 215mm | Blunt |



Volkman Retractor 5 Prongs

| | | |
|----------|-------|-------|
| PH252356 | 215mm | Sharp |
| PH252364 | 215mm | Blunt |



Volkman Retractor 6 Prongs

| | | |
|----------|-------|-------|
| PH252357 | 215mm | Sharp |
| PH252365 | 215mm | Blunt |



Baby Senn Miller Retractor Sharp Double Ended

PH252368



Did you know we offer instrument repair services?

We offer flexible and comprehensive services including a rapid repair service. All our services are ISO 13485 accredited with warranties provided for all repairs.


Call us today to find out more
(0)3452 414818

 **STERIS**

Alms Skin Retractor Sharp

PH252370 70mm

PH252371 100mm



Alms Skin Retractor Blunt

PH252372 70mm

PH252373 100mm




Charnley Initial Incision Retractor

PH252374




Balfour Abdominal Retractor

PH252375 90mm Deep x 35mm Wide




Balfour Baby Abdominal Retractor

PH252376 27mm Deep



Goligher Abdominal Retractor

PH252383



Joll Thyroid Retractor

PH252404



Cone Laminectomy Retractor 3 into 4 Teeth with Hinged Arms

| | | |
|----------|-------------|-------|
| PH252407 | Sharp Teeth | 125mm |
| PH252408 | Blunt Teeth | 125mm |
| PH252409 | Sharp Teeth | 160mm |
| PH252410 | Blunt Teeth | 160mm |
| PH252411 | Sharp Teeth | 200mm |
| PH252412 | Blunt Teeth | 200mm |
| PH252413 | Sharp Teeth | 260mm |
| PH252414 | Blunt Teeth | 260mm |



Jackson Burrows Retractor

| | |
|----------|-------|
| PH252420 | 175mm |
|----------|-------|



Norfolk and Norwich Retractor 4 into 5 Teeth

| |
|----------|
| PH252424 |
|----------|



Norfolk And Norwich Self Retaining Retractor Long Teeth With 4 Into 5 Teeth

| |
|----------|
| PH252425 |
|----------|




Gelpi Vulva Retractor

| | |
|----------|------------------|
| PH252437 | Large Size 200mm |
| PH252438 | Small Size 170mm |



Irwin Retractor 4 into 4 Teeth

| | |
|----------|-------|
| PH252440 | 180mm |
|----------|-------|





Mollison Retractor 4 into 4 Teeth Adult

PH252443



Mollison Retractor 2 into 2 Teeth Child

PH252444



Travers Retractor 4 into 5 Teeth

PH252446



West Retractor 2 into 3 Teeth

PH252450 Blunt 115mm

PH252451 Sharp 115mm



West Retractor 3 into 4 Teeth

PH252452 Blunt 140mm

PH252453 Sharp 140mm



Weitlander Retractor 2 into 3 Teeth

PH252455 Sharp 115mm

PH252456 Sharp 140mm

PH252457 Sharp 160mm

PH252458 Blunt 115mm

PH252459 Blunt 140mm

PH252460 Blunt 160mm

Doyen Rectangular Blade Retractor

| | |
|----------|--------------------|
| PH252465 | Small 48mm x 60mm |
| PH252466 | Medium 48mm x 90mm |
| PH252467 | Large 48mm x 125mm |



Landon Vaginal Retractor

| | |
|----------|-------------------|
| PH252470 | 90mm x 25mm Blade |
|----------|-------------------|



Tudor Edwards Scapula Retractor

| | |
|----------|-------|
| PH252499 | 160mm |
|----------|-------|



Kocher Retractor

| | |
|----------|------------|
| PH252426 | 45mm Blade |
| PH252427 | 54mm Blade |
| PH252428 | 62mm Blade |



We also offer decontamination and sterilisation services.

With over 20 years' experience and more than 100 million instruments processed each year, we can offer a fully outsourced service or contingency for your sterile services department.

Contact us today to find out how we can help you reduce cost and improve turnarounds.

 **+44 (0)3452 414818**

 **steris-ims-instruments.com**

 **STERIS**

Oral Instruments

Lack Tongue Depressor

| | |
|----------|-----------------------|
| PH282535 | Infant Size 12mm Wide |
| PH282536 | Child Size 19mm Wide |
| PH282537 | Adult Size 22mm Wide |



Mayo Tongue Depressor

| | |
|----------|------------------|
| PH282538 | Child Size 140mm |
| PH282539 | Adult Size 200mm |



Did you know we offer instrument repair services?

We offer flexible and comprehensive services including a rapid repair service. All our services are ISO 13485 accredited with warranties provided for all repairs.

Call us today to find out more
(0)3452 414818



Davis Boyle Mouth Gag Frame

| | |
|----------|--------------------|
| PH292546 | Child 135mm x 70mm |
| PH292547 | Adult 150mm x 90mm |



Davis Boyle Mouth Gag Tongue Plate

| | |
|----------|-------|
| PH292552 | 65mm |
| PH292553 | 75mm |
| PH292554 | 90mm |
| PH292555 | 100mm |
| PH292556 | 115mm |




Daughty Mouth Gag Tongue Plate

| | |
|----------|-------|
| PH292562 | 64mm |
| PH292563 | 76mm |
| PH292564 | 89mm |
| PH292565 | 102mm |
| PH292566 | 114mm |



Draffin Bipod



| | |
|------------|---|
| PH292567Ad | Adult Fitted with 4 Rings For Varying Heights 480mm |
| PH292567Ch | Child Fitted with 4 Ring For Varying Heights 375mm |

Maguaran Plate For Holding Draffin Bipod



| | |
|----------|-------|
| PH292568 | 407mm |
|----------|-------|

Laryngeal Handle with Mirror



| | |
|-------------|------|
| PH322800-8 | 8mm |
| PH322800-10 | 10mm |
| PH322800-12 | 12mm |
| PH322800-14 | 14mm |
| PH322800-16 | 16mm |
| PH322800-18 | 18mm |
| PH322800-20 | 20mm |
| PH322800-22 | 22mm |
| PH322800-24 | 24mm |
| PH322800-26 | 26mm |
| PH322800-28 | 28mm |
| PH322800-30 | 30mm |

LARYNGEAL HANDLE ONLY

| | |
|----------|-------|
| PH322801 | 120mm |
|----------|-------|

Dental Mirror with Handle



| | |
|-------------|------|
| PH322810 -8 | 8mm |
| PH322810-10 | 10mm |
| PH322810-12 | 12mm |
| PH322810-14 | 14mm |
| PH322810-16 | 16mm |
| PH322810-18 | 18mm |
| PH322810-20 | 20mm |
| PH322810-22 | 22mm |
| PH322810-24 | 24mm |
| PH322810-26 | 26mm |
| PH322810-28 | 28mm |
| PH322810-30 | 30mm |

Tracheal Instruments



Bowlby Tracheal Dilating Forceps Curved

PH302595 Adult 140mm

PH302596 Child 115mm



Magil Intra Tracheal Dilating Forceps

| | |
|----------|-------------|
| PH302601 | Adult 240mm |
| PH302602 | Child 200mm |




Aural Instruments

House Aural Curette Double Ended Straight

PH362920 1.0mm and 1.2mm Spoon 180mm


PH362921 1.5mm and 1.8mm Spoon 180mm



House Aural Curette Double Ended Angled

PH362922 1.0mm and 1.2mm Spoon 180mm

PH362923 1.5mm and 1.8mm Spoon 180mm



MacEwen Double Ended Round Curette and Seeker (McEwan)

PH362930 185mm



MacEwen Double Ended Oval Curette and Seeker (McEwan)

PH362931 185mm



Did you know we offer instrument repair services?

We offer flexible and comprehensive services including a rapid repair service. All our services are ISO 13485 accredited with warranties provided for all repairs.




Contact us today to find out more

+44 (0)3452 414818


 **STERIS**

Cawthorne Hook Single End No.1 Slight Curve


PH362940 165mm

Cawthorne Hook Single End No.2 Medium Curve


PH362941 165mm

Cawthorne Hook Single End No.3 Full Curve

PH362942 165mm

Cawthorne Hook Single End No.4 Double Curve


PH362943 165mm

Cawthorne Hook Single End No.5 90° Angle 0.5mm Point

PH362944 165mm


Cawthorne Hook Single End No.6 90° Angle 1.0mm Point

PH362945 165mm




Cawthorne Hook Single Ended No.7 Retrograde

PH362946 165mm




Cawthorne Hook Single End No.8 45° Angle 0.5mm Point

PH362947 165mm



Cawthorne Hook Single End No.9 Straight 0.5mm Point

PH362948 165mm



Volkmann Aural Curette Double Ended Scoops Straight

| | | |
|----------|-------|--------|
| PH362960 | 215mm | Small |
| PH362961 | 215mm | Medium |
| PH362962 | 215mm | Large |




Volkmann Aural Curette Single End Straight

| | | |
|----------|---------|---------------------|
| PH362963 | Size 00 | Hollow Handle 180mm |
| PH362964 | Size 0 | Hollow Handle 180mm |
| PH362965 | Size 1 | Hollow Handle 180mm |
| PH362966 | Size 2 | Hollow Handle 180mm |
| PH362967 | Size 3 | Hollow Handle 180mm |
| PH362968 | Size 4 | Hollow Handle 180mm |
| PH362969 | Size 5 | Hollow Handle 180mm |
| PH362970 | Size 6 | Hollow Handle 180mm |



Musgrove Aural Double Ended Wool Carrier and Cerumen Pick

PH362981 170mm



St Barts Blunt Hook Wool Carrier

PH362982 170mm



Formby Cerumen Scoop and Hook

PH362983 180mm



Gruber Aural Speculum Round End

| | |
|----------|-----|
| PH362986 | 4mm |
| PH362987 | 5mm |
| PH362988 | 6mm |
| PH362989 | 7mm |



Gruber Aural Speculum Oval End

| | |
|----------|-----|
| PH362991 | 4mm |
| PH362992 | 5mm |
| PH362993 | 6mm |
| PH362994 | 7mm |



Shea Speculum

| | |
|----------|--------------|
| PH362995 | Set of 5 |
| PH362996 | No.1 - 3.7mm |
| PH362997 | No.2 - 4.0mm |
| PH362998 | No.3 - 5.0mm |
| PH362999 | No.4 - 6.0mm |
| PH363000 | No.5 - 7.0mm |



Rosen Aural Speculum Round End

| | |
|----------|-----|
| PH363021 | 4mm |
| PH363022 | 5mm |
| PH363023 | 6mm |
| PH363024 | 7mm |



We also offer decontamination and sterilisation services.

With over 20 years' experience and more than 100 million instruments processed each year, we can offer a fully outsourced service or contingency for your sterile services department.



Contact us today to find out more
+44 (0)3452 414818

Rosen Aural Speculum Oval End



| | |
|----------|-----|
| PH363026 | 4mm |
| PH363027 | 5mm |
| PH363028 | 6mm |
| PH363029 | 7mm |

Toynbee Aural Speculum



| | |
|----------|-----|
| PH363031 | 4mm |
| PH363032 | 5mm |
| PH363033 | 6mm |
| PH363034 | 7mm |

Tumarkin Nasal Speculum



| | | | |
|----------|------------------------|----------|------------|
| PH453335 | Full Set of 8 Sizes | PH453339 | 3 = 7.5mm |
| PH453336 | 0 = 4.5mm | PH453340 | 4 = 8.0mm |
| PH453337 | 1 = 6.5mm | PH453341 | 5 = 8.5mm |
| PH453338 | 2 = 7.0mm | PH453342 | 6 = 9.0mm |
| | | PH453343 | 7 = 11.0mm |

Agnew Myringotome Ear Knife Angled

| | |
|----------|-------|
| PH373051 | 190mm |
|----------|-------|



Nasal, Tonsil and Adenoid Instruments

St Clair Thomson Adenoid Curette

| | | |
|----------|----------------|-------|
| PH443272 | Size 0 - 8mm | 220mm |
| PH443273 | Size 1 - 10mm | 220mm |
| PH443274 | Size 1½ - 12mm | 220mm |
| PH443275 | Size 2 - 14mm | 220mm |
| PH443276 | Size 3 - 16mm | 220mm |
| PH443277 | Size 4 - 18mm | 220mm |



Beckmann Adenoid Curette

| | | |
|----------|---------------|-------|
| PH443279 | Size 0 - 8mm | 210mm |
| PH443280 | Size 1 - 10mm | 210mm |
| PH443282 | Size 2 - 14mm | 210mm |
| PH443283 | Size 3 - 16mm | 210mm |
| PH443284 | Size 4 - 18mm | 210mm |




Cottle Nasal Speculum

| | |
|----------|-------------|
| PH453310 | No.1 - 20mm |
| PH453311 | No.2 - 32mm |
| PH453312 | No.3 - 41mm |
| PH453313 | No.4 - 60mm |



Killian Nasal Speculum

| | |
|----------|-------------------|
| PH453314 | 38mm Blade Length |
| PH453315 | 64mm Blade Length |
| PH453316 | 76mm Blade Length |
| PH453317 | 89mm Blade Length |



St Clair Thomson Nasal Speculum

| | |
|----------|-------------------------------|
| PH453320 | Small 38mm Blade Length |
| PH453321 | Medium 51mm Blade Length |
| PH453322 | Large 64mm Blade Length |
| PH453323 | Extra Large 89mm Blade Length |




Thudichum Nasal Speculum

| | |
|----------|--------|
| PH453325 | Size 1 |
| PH453326 | Size 2 |
| PH453327 | Size 3 |
| PH453328 | Size 4 |
| PH453329 | Size 5 |
| PH453330 | Size 6 |
| PH453331 | Size 7 |



Eve Tonsil Snare Sliding Type

| | |
|----------|-------|
| PH463377 | 280mm |
|----------|-------|



Ballenger Nasal Cartilage Knife

PH463389 3mm Wide 200mm

PH463390 5mm Wide 200mm

PH463391 8mm Wide 200mm



We also offer decontamination and sterilisation services.

With over 20 years' experience and more than 100 million instruments processed each year, we can offer a fully outsourced service or contingency for your sterile services department.



Contact us today to find out more
+44 (0)3452 414818

 **STERIS**

Silver Nasal Chisel Straight

PH483450 180mm




Silver Nasal Chisel Curved to the Right

PH483451 180mm



Silver Nasal Chisel Curved to the Left


PH483452 180mm



Joseph Nasal Saw Vertical Shaft and Bayonet Crank

PH483458 Right 185mm

PH483460 Left 185mm




Freer Double Ended Separator and Elevator

PH493482 200mm



Hurd Double Ended Pillar Retractor and Tonsil Dissector

PH493486 200mm



Dissectors and Elevators

Gwynne Evans Double Ended 6mm Wide Dissector

PH493495 190mm



Hill Double Ended Elevator Curved to Side with 5mm Wide Rugine End

PH493497 220mm




Howarth Double Ended Elevator Curved on Flat with Rugine End

PH493498 220mm




MacDonald Double Ended Dissector

PH493505 190mm



Syme Double Ended Dissector One Blunt End and One Semi Sharp End

PH493508 190mm



Watson Cheyne Double Ended Elevator Fine with Probe

| | |
|----------|-------|
| PH493510 | 130mm |
| PH493511 | 170mm |
| PH493512 | 230mm |



Warwick James Single End Elevator Straight

| | |
|----------|-------|
| PH493515 | 170mm |
|----------|-------|



Warwick James Single End Elevator Curved Right

| | |
|----------|-------|
| PH493516 | 170mm |
|----------|-------|




Warwick James Single End Elevator Curved Left

| | |
|----------|-------|
| PH493517 | 170mm |
|----------|-------|



Cryers Elevator Left

| | |
|----------|-------|
| PH493518 | 170mm |
|----------|-------|



Cryers Elevator Right

| | |
|----------|-------|
| PH493519 | 170mm |
|----------|-------|



Mitchell Trimmer Double Ended

PH493520 160mm



Henderson Dissector (Periosteal Elevator)

PH493522 160mm



Intestinal Clamps

Denis Browne Intestinal Forceps Straight


PH513612 180mm



Doyen Intestinal Clamp Straight

PH513617 180mm

PH513621 230mm



Doyen Intestinal Clamp Curved

PH513618 180mm


PH513622 230mm



Doyen Intestinal Clamp Atraumatic Jaw Straight

PH513627 180mm

PH513631 230mm



Doyen Intestinal Clamp Atraumatic Jaw Curved

PH513628 180mm

PH513632 230mm




Kocher Intestinal Clamp Straight


PH513636 230mm
Blade Length 100mm

PH513637 250mm
Blade Length 125mm

PH513638 280mm
Blade Length 150mm



Lane Intestinal Clamp Straight



| | |
|--------------|-------|
| PH513647 | 180mm |
| Blade Length | |
| | 75mm |
| PH513648 | 230mm |
| Blade Length | |
| | 100mm |

Parker Kerr Intestinal Clamp Straight



| | |
|----------|-------|
| PH513671 | 230mm |
| PH513676 | 250mm |

Parker Kerr Intestinal Clamp Curved



| | |
|----------|-------|
| PH513680 | 230mm |
| PH513681 | 250mm |

Schoemaker Ogilvie Clamp with Hollow Blades



| | |
|----------|-------|
| PH513691 | 155mm |
| PH513692 | 170mm |
| PH513693 | 200mm |

Bone Instruments

Cairn Bone Rongeur

PH584147 150mm



Cleveland Bone Rongeur

PH584148 140mm



Daniel Miniature Bone Rongeur

PH584149 140mm



Glasgow Pattern Bone Rongeur

PH584150 150mm



Jansen Middleton Bone Rongeur Straight with Compound Action

PH584152 185mm



Jansen Middleton Bone Rongeur Curved with Compound Action

PH584153 185mm





Lempert Bone Rongeur Straight

PH584154 190mm



Luer Bone Rongeur Straight

PH584157 170mm



Luer Bone Rongeur Curved

PH584158 170mm



Leksell Bone Rongeur with Compound Action

PH584161 4mm Bite 230mm

PH584162 6mm Bite 230mm

PH584163 8mm Bite 230mm



Northfield Bone Rongeur Straight with 10mm Bite

PH584165 200mm



Olivecrona Bone Rongeur Angled with Compound Action

| | |
|----------|--------------------------|
| PH584170 | 5mm Bite 200mm |
| PH584171 | 6mm Bite 230mm Angled |



Pennybacker Bone Rongeur Straight with 4mm Bite

PH584175



Ruskin Bone Rongeur Curved 6mm Bite with Compound Action

| | |
|----------|-------|
| PH584178 | 180mm |
| PH584179 | 230mm |



Sauerburch Bone Rongeur with Compound Action

PH584180 300mm



Semb Bone Rongeur Curved with Double Action 20mm Wide Jaws with Compound Action

PH584185 240mm



Stille Luer Bone Rongeur Straight with Compound Action

PH584190 220mm





Stille Luer Bone Rongeur Curved
with Compound Action

PH584191 220mm



Stille Bone Rongeur Angled
with Compound Action

PH584194 230mm



Trotter Bone Rongeur Angled 10mm
Bite with Compound Action

PH584195 200mm



Ward Bone Rongeur Curved
with Compound Action

PH584197 180mm



Wilms Bone Rongeur Angled on Flat
with Narrow Jaws and a 4mm Bite
With A Simple Action

PH584198 200mm



Phoenix Bone Rongeur

PH584200 3mm Bite 230mm

PH584201 6mm Bite 230mm



Zielke Bone Rongeur 3mm Bite

PH584202 250mm



Freidman Bone Rongeur 3mm Bite

PH584250 130mm



Mini Bone Nibbler 1mm Bite

PH584253 140mm



Murphy Hip Skid Fully Insulated Double Ended

PH605502



Did you know we offer instrument repair services?


We offer flexible and comprehensive services including a rapid repair service. All our services are ISO 13485 accredited with warranties provided for all repairs.



Contact us today to find out more
+44 (0)3452 414818

Baby Hohman Elevator

PH252022 160mm



Hohman Bone Lever

PH252100 8mm Wide 240mm

PH252101 18mm Wide 240mm

PH252102 43mm Wide 240mm

PH252103 70mm Wide 240mm



Trethowen Bone Lever

PH614305 215mm



Trethowen Angled Bone Lever

PH614306 215mm



Bristow Bone Elevator

PH614310 230mm



Lane Bone Elevator with Slight Curve and Solid Handle

PH614314 240mm



Lane Bone Elevator with Slight Curve and Ring Handle

PH614315 265mm




Lane Bone Elevator Serrated with Slight Curve and Ring Handle

PH614316 260mm



Mitchel Bone Elevator

PH614317 200mm



Cobb Elevator 275mm

PH614321 10mm

PH614322 13mm

PH614323 19mm



Pennybacker Periosteal Elevator

| | |
|----------|----------------|
| PH624321 | 4mm Wide 170mm |
| PH624322 | 6mm Wide 170mm |



Lempert Periosteal Elevator

| | |
|----------|-------|
| PH624323 | 170mm |
|----------|-------|



Doyen Periosteal Elevator Curved to Left

| | |
|----------|-------------|
| PH624325 | Adult 165mm |
| PH624327 | Child 165mm |



Doyen Periosteal Elevator Curved to Right

| | |
|----------|-------------|
| PH624326 | Adult 165mm |
| PH624328 | Child 165mm |



Ward Periosteal Elevator

| | |
|----------|-------|
| PH624329 | 170mm |
|----------|-------|



Farabeuf Straight Bone Rugine with Chisel Edge

PH624330 150mm



Jones Rugine Pistol Shaped Handle Fan Shaped End and Convex Sharp Edge

PH624332 185mm



Miller Bone Rasp Large

PH654366 180mm



We also offer
decontamination and
sterilisation services.

With over 20 years' experience and more than 100 million instruments processed each year, we can offer a fully outsourced service or contingency for your sterile services department.



Contact us today to find out more
+44 (0)1332 387100

Martini Double Ended Bone Curette

PH624341 150mm



Stanley Boyd Bone Double Ended Curette Oval Scoops

PH624342 Small 150mm
PH624343 Medium 200mm
PH624344 Large 280mm



Volkman Spoon Double Ended Bone Curette

PH624346 Small 215mm
PH624347 Medium 215mm
PH624348 Large 215mm



Williger Bone Double Ended Curette

PH624350 Small 125mm
PH624351 Medium 125mm
PH624352 Large 125mm



Kocher Bone Hook

PH624354 250mm



Lucae Bone Mallet

PH624370 190mm



Heath Bone Mallet with Long Handle

PH624383 260mm



Heath Bone Mallet

PH624384 18 Oz 150mm

PH624385 28 Oz 150mm

PH624386 36 Oz 150mm



Did you know we offer instrument repair services?

We offer flexible and comprehensive services including a rapid repair service. All our services are ISO 13485 accredited with warranties provided for all repairs.

Call us today to find out more
(0)3452 414818



Light Weight Tufnol Mallet

PH624398 210mm



Heavy Tufnol Mallet

PH624399 270mm



STERIS IMS 03452 414818

Fickling Mallet

PH625001 200mm



Cottle Mallet 30mm Ø Head

PH625005 190mm





Mini Stille Bone Chisel 150mm

| | |
|----------|-----------|
| PH624435 | 3mm Wide |
| PH624436 | 5mm Wide |
| PH624437 | 7mm Wide |
| PH624438 | 9mm Wide |
| PH624439 | 12mm Wide |



Stille Bone Chisel 200mm

| | |
|----------|-----------|
| PH624440 | 10mm Wide |
| PH624441 | 15mm Wide |
| PH624442 | 20mm Wide |
| PH624443 | 25mm Wide |
| PH624444 | 30mm Wide |



Capener Tufnol Handled Chisel

| | |
|----------|-----------|
| PH626006 | 10mm Wide |
| PH626007 | 12mm Wide |
| PH626008 | 15mm Wide |
| PH626009 | 20mm Wide |
| PH626010 | 25mm Wide |
| PH626011 | 30mm Wide |



Jenkins Bone Gouge 150mm

| | |
|----------|-----------|
| PH624569 | 3mm Wide |
| PH624570 | 4mm Wide |
| PH624571 | 5mm Wide |
| PH624572 | 6mm Wide |
| PH624573 | 8mm Wide |
| PH624574 | 10mm Wide |
| PH624575 | 12mm Wide |
| PH624576 | 14mm Wide |



Stille Bone Gouge 200mm

| | |
|----------|-----------|
| PH624600 | 10mm Wide |
| PH624601 | 13mm Wide |
| PH624602 | 15mm Wide |
| PH624603 | 20mm Wide |
| PH624604 | 25mm Wide |
| PH624605 | 30mm Wide |



Coupland Dental Gouge 170mm

| | |
|----------|--------------|
| PH624635 | Small (3mm) |
| PH624636 | Medium (4mm) |
| PH624637 | Large (5mm) |



Tilley Gouge V Shaped

| | |
|----------|-------|
| PH624638 | 125mm |
|----------|-------|



Tufnol Handled Straight Gouge 270mm

| | |
|----------|-----------|
| PH628000 | 10mm Wide |
| PH628001 | 12mm Wide |
| PH628002 | 15mm Wide |
| PH628003 | 20mm Wide |
| PH628004 | 25mm Wide |
| PH628005 | 30mm Wide |



Capener Tufnol Handled Cranked Gouge 235mm

| | |
|----------|-----------|
| PH628010 | 6mm Wide |
| PH628011 | 8mm Wide |
| PH628012 | 11mm Wide |
| PH628013 | 16mm Wide |



Duray French Osteotome

| | |
|----------|-----------------|
| PH624496 | 140mm 5mm Wide |
| PH624497 | 140mm 7mm Wide |
| PH624498 | 140mm 9mm Wide |
| PH624499 | 140mm 11mm Wide |



Smith Peterson Bone Osteotomes

| | |
|----------|-----------------|
| PH624509 | 200mm 6mm Wide |
| PH624510 | 200mm 10mm Wide |
| PH624511 | 200mm 13mm Wide |
| PH624512 | 200mm 20mm Wide |
| PH624513 | 200mm 25mm Wide |
| PH624514 | 200mm 30mm Wide |




Stille Bone Osteotomes

| | |
|----------|-----------------|
| PH624515 | 200mm 2mm Wide |
| PH624516 | 200mm 5mm Wide |
| PH624517 | 200mm 6mm Wide |
| PH624518 | 200mm 7mm Wide |
| PH624519 | 200mm 8mm Wide |
| PH624520 | 200mm 9mm Wide |
| PH624521 | 200mm 10mm Wide |
| PH624522 | 200mm 11mm Wide |
| PH624523 | 200mm 13mm Wide |
| PH624524 | 200mm 15mm Wide |
| PH624525 | 200mm 20mm Wide |
| PH624526 | 200mm 25mm Wide |
| PH624527 | 200mm 30mm Wide |




Stille Curved Bone Osteotomes

| | |
|----------|-----------|
| PH624528 | 19mm Wide |
| PH624529 | 15mm Wide |
| PH624530 | 10mm Wide |




Lambotte Osteotome 150mm



| | | | |
|----------|----------------|----------|-----------------|
| PH624539 | 150mm 1mm Wide | PH624544 | 150mm 8mm Wide |
| PH624540 | 150mm 2mm Wide | PH624545 | 150mm 10mm Wide |
| PH624541 | 150mm 3mm Wide | PH624546 | 150mm 12mm Wide |
| PH624542 | 150mm 4mm Wide | PH664560 | 150mm 7mm Wide |
| PH624543 | 150mm 6mm Wide | | |



Lambotte Osteotome 200mm



| | |
|----------|-----------------|
| PH624547 | 200mm 15mm Wide |
| PH624548 | 200mm 20mm Wide |
| PH624549 | 200mm 25mm Wide |
| PH624550 | 200mm 30mm Wide |
| PH624552 | 200mm 10mm Wide |



Gillies (French) Osteotome



| | | | |
|------------|----------------|-------------|-----------------|
| PH624553/2 | 140mm 2mm Wide | PH624553/7 | 140mm 7mm Wide |
| PH624553/3 | 140mm 3mm Wide | PH624553/8 | 140mm 8mm Wide |
| PH624553/4 | 140mm 4mm Wide | PH624553/9 | 140mm 9mm Wide |
| PH624553/5 | 140mm 5mm Wide | PH624553/10 | 140mm 10mm Wide |
| PH624553/6 | 140mm 6mm Wide | PH624553/11 | 140mm 11mm Wide |




Jayes Osteotome Mini Lambotte



| | |
|----------|----------------|
| PH624554 | 125mm 1mm Wide |
| PH624555 | 125mm 2mm Wide |
| PH624556 | 125mm 4mm Wide |



Ward Osteotome




| | |
|----------|----------------|
| PH624557 | 180mm 3mm Wide |
| PH624559 | 180mm 5mm Wide |




Tufnol Handled Straight Osteotome

| | |
|----------|-----------------|
| PH627000 | 265mm 10mm Wide |
| PH627001 | 265mm 12mm Wide |
| PH627002 | 265mm 15mm Wide |
| PH627003 | 265mm 20mm Wide |
| PH627004 | 265mm 25mm Wide |
| PH627005 | 265mm 30mm Wide |




Tufnol Handled Curved Osteotome

| | |
|----------|-----------------|
| PH627006 | 265mm 10mm Wide |
| PH627007 | 265mm 12mm Wide |
| PH627008 | 265mm 15mm Wide |
| PH627009 | 265mm 20mm Wide |
| PH627010 | 265mm 25mm Wide |
| PH627011 | 265mm 30mm Wide |




Capener Tufnol Handled Straight Osteotome

| | |
|----------|-----------------|
| PH627015 | 250mm 10mm Wide |
| PH627016 | 250mm 12mm Wide |
| PH627017 | 250mm 15mm Wide |
| PH627018 | 250mm 20mm Wide |
| PH627019 | 250mm 25mm Wide |
| PH627020 | 250mm 30mm Wide |



Capener Tufnol Handled Curved Osteotome

| | |
|----------|-----------------|
| PH627021 | 250mm 10mm Wide |
| PH627022 | 250mm 12mm Wide |
| PH627023 | 250mm 15mm Wide |
| PH627024 | 250mm 20mm Wide |
| PH627025 | 250mm 25mm Wide |
| PH627026 | 250mm 30mm Wide |



Vehmehren Rib Shears

PH644687 270mm




Craford's Rib Approximator

PH644691 260mm



Sellor Rib Approximator and Spreader Combined

PH644692 200mm




Sellor Rib Spreader with 3 Pairs of Detachable Blades and Locking Key (Sternal Spreader)

PH644693 200mm



Burford Finochietto Rib Spreader with 2 Pairs of Detachable Blades

PH644694 290mm



Finochietto Rib Spreader

| | |
|----------|--|
| PH644696 | Baby size with 18mm x 21mm Blades and 75mm Spread |
| PH644698 | Small Size with 30mm x 45mm Blades and 165mm Spread |
| PH644699 | Medium Size with 45mm x 65mm Blades and 200mm Spread |
| PH644700 | Large Size with 70mm x 65mm Blades and 250mm Spread |



Watersons Rib Spreader Infant Size

| | |
|----------|---|
| PH644710 | 40mm Wide x 20mm Deep Blades and 105mm Spread |
|----------|---|



Plaster Cast

Bergmann Plaster Shears/Bandage Scissors with Angular Blades

PH654725 230mm



Bohler Plaster Shears

PH654727 230mm



Bruns Plaster Shears with Toothed Lower Blade

PH654729 240mm



Guys Heavy Pattern Plaster Shears

PH654733 300mm





Stille Lorenz Plaster Shears

| | |
|----------|-------|
| PH654737 | 230mm |
| PH654738 | 265mm |
| PH654739 | 300mm |
| PH654740 | 355mm |
| PH654741 | 380mm |

Henning Plaster Cast Spreaders

| | |
|----------|-------|
| PH654754 | 280mm |
|----------|-------|



We also offer
decontamination and
sterilisation services.

With over 20 years' experience and more than 100 million instruments processed each year, we can offer a fully outsourced service or contingency for your sterile services department.



Contact us today to find out more
+44 (0)1332 387100



Rongeurs

Caspar Rongeur Straight 150mm Shaft

| | |
|----------|----------|
| PH654767 | 2mm Bite |
|----------|----------|

| | |
|----------|----------|
| PH654768 | 3mm Bite |
|----------|----------|



Cushing Rongeur 130mm Effective Length Straight Jaw

| | |
|----------|----------|
| PH654770 | 3mm Bite |
|----------|----------|

| | |
|----------|----------|
| PH654771 | 4mm Bite |
|----------|----------|

| | |
|----------|----------|
| PH654772 | 5mm Bite |
|----------|----------|

| | |
|----------|----------|
| PH654773 | 7mm Bite |
|----------|----------|



Cushing Rongeur 130mm Effective Length Angled Upward Jaw

| | |
|----------|----------|
| PH654774 | 3mm Bite |
|----------|----------|

| | |
|----------|----------|
| PH654775 | 4mm Bite |
|----------|----------|

| | |
|----------|----------|
| PH654776 | 5mm Bite |
|----------|----------|

| | |
|----------|----------|
| PH654777 | 7mm Bite |
|----------|----------|



Did you know we offer instrument repair services?

We offer flexible and comprehensive services including a rapid repair service. All our services are ISO 13485 accredited with warranties provided for all repairs.




Contact us today to find out more
+44 (0)3452 414818



Beck Rongeur Straight
130mm Shaft Length

PH654783 2mm Bite




Clowards Rongeur 150mm
Effective Length Straight Jaw

PH654833 3mm Bite

PH654834 5mm Bite

PH654835 6mm Bite




Love Grunwald Pituitary Rongeur
200mm Straight Cut

PH654905 2mm Bite

PH654906 3mm Bite

PH654907 4mm Bite

PH654908 5mm Bite



Love Grunwald Pituitary Rongeur
200mm Upwards Cut

PH654910 2mm Bite

PH654911 3mm Bite

PH654912 4mm Bite

PH654913 5mm Bite



Pennybacker Rongeur Straight
130mm Shaft

PH654918 3mm Bite

Semi-Detachable Kerrison Rongeurs

The Semi-Detachable Kerrison rongeurs were designed with both the patient and surgeon in mind. The durable take-apart design enables the instrument to be completely cleaned, reducing the risk of bioburden buildup without sacrificing functionality.

Disassembly Instructions

- 1 Gently squeeze the handle
- 2 Press the catch downwards
- 3 Lift the top slide rail upwards

Please refer to Instructions for Use for more information.



Semi Detachable Kerrison Rongeur Upcut 180mm

| | | | | | |
|------------|----------|-----|------------|----------|-----|
| 45 Degrees | PH654950 | 1mm | 90 Degrees | PH654955 | 1mm |
| 45 Degrees | PH654951 | 2mm | 90 Degrees | PH654956 | 2mm |
| 45 Degrees | PH654952 | 3mm | 90 Degrees | PH654957 | 3mm |
| 45 Degrees | PH654953 | 4mm | 90 Degrees | PH654958 | 4mm |
| 45 Degrees | PH654954 | 5mm | 90 Degrees | PH654959 | 5mm |





The instrument is handcrafted from forge to finish, including a consistent cutting edge. Ejector bar and grooveless rails help prevent bone and soft tissue from sticking to the cutting edge. The Semi-Detachable Kerrison rongeurs can be disassembled for cleaning, allowing visualisation of the entire instrument.

Semi Detachable Kerrison Rongeur Upcut 200mm

| | | | | | |
|------------|----------|-----|------------|----------|-----|
| 45 Degrees | PH654960 | 1mm | 90 Degrees | PH654965 | 1mm |
| 45 Degrees | PH654961 | 2mm | 90 Degrees | PH654966 | 2mm |
| 45 Degrees | PH654962 | 3mm | 90 Degrees | PH654967 | 3mm |
| 45 Degrees | PH654963 | 4mm | 90 Degrees | PH654968 | 4mm |
| 45 Degrees | PH654964 | 5mm | 90 Degrees | PH654969 | 5mm |



Semi Detachable Kerrison Rongeur Upcut 210mm

| | | |
|------------|----------|-----|
| 45 Degrees | PH654970 | 1mm |
| 45 Degrees | PH654971 | 2mm |
| 45 Degrees | PH654972 | 3mm |
| 45 Degrees | PH654973 | 4mm |
| 45 Degrees | PH654974 | 5mm |

| | | |
|------------|----------|-----|
| 90 Degrees | PH654975 | 1mm |
| 90 Degrees | PH654976 | 2mm |
| 90 Degrees | PH654977 | 3mm |
| 90 Degrees | PH654978 | 4mm |
| 90 Degrees | PH654979 | 5mm |



Semi Detachable Kerrison Rongeur Upcut 180mm 90 Degrees (Black)

| | | |
|------------|----------|-----|
| 45 Degrees | PH654980 | 1mm |
| 45 Degrees | PH654981 | 2mm |
| 45 Degrees | PH654982 | 3mm |
| 45 Degrees | PH654983 | 4mm |
| 45 Degrees | PH654984 | 5mm |

| | | |
|------------|----------|-----|
| 90 Degrees | PH654985 | 1mm |
| 90 Degrees | PH654986 | 2mm |
| 90 Degrees | PH654987 | 3mm |
| 90 Degrees | PH654988 | 4mm |
| 90 Degrees | PH654989 | 5mm |



Semi Detach Kerrison Upcut 200mm (Black)

| | | | | | |
|------------|----------|-----|------------|----------|-----|
| 45 Degrees | PH654990 | 1mm | 90 Degrees | PH655000 | 1mm |
| 45 Degrees | PH654991 | 2mm | 90 Degrees | PH655001 | 2mm |
| 45 Degrees | PH654992 | 3mm | 90 Degrees | PH655002 | 3mm |
| 45 Degrees | PH654993 | 4mm | 90 Degrees | PH655003 | 4mm |
| 45 Degrees | PH654994 | 5mm | 90 Degrees | PH655004 | 5mm |



Semi Detachable Kerrison Rongeur 210mm (Black)

| | | | | | |
|------------|----------|-----|------------|----------|-----|
| 45 Degrees | PH654995 | 1mm | 90 Degrees | PH655005 | 1mm |
| 45 Degrees | PH654996 | 2mm | 90 Degrees | PH655006 | 2mm |
| 45 Degrees | PH654997 | 3mm | 90 Degrees | PH655007 | 3mm |
| 45 Degrees | PH654998 | 4mm | 90 Degrees | PH655008 | 4mm |
| 45 Degrees | PH654999 | 5mm | 90 Degrees | PH655009 | 5mm |

Diathermy

Monopolar (Diathermy) Lead 4mm
Machine End Standard 4.8mm Cutting
Socket 3 Metre Cable Length

PH757500



Monopolar (Diathermy) Lead 8mm
Bovie Machine End Standard 4.8mm
Cutting Socket 3 Metre Cable Length

PH757501



Laparoscopic Monopolar Lead with 4mm
Machine End and 4mm Laparoscopic
Cutting Socket to fit all Laparoscopic
Instruments 3 Metre Cable Length

PH757502




Laparoscopic Monopolar Lead with
8mm Bovie Machine End and 4mm
Laparoscopic Cutting Socket to fit
all Laparoscopic Instruments
3 Metre Cable Length

PH757503



Bipolar Lead Standard European
Instrument Block Fitting 2 x 4mm Plugs
Machine End 3 Metre Cable

PH757510




Diathermy Quiver Standard 60mm
Diameter 200mm Length

PH757560



Diathermy Quiver Telescopic 60mm
Diameter extending to 400mm Length

PH757561



Miscellaneous Instruments

| | | |
|-------|-------|----|
| XXXX | | XX |
| XXXXX | | XX |
| XXXX | | XX |

Ring Cutter

| | |
|----------|-------------------------|
| PH242100 | Complete with Saw 160mm |
| PH242101 | Forceps Only 160mm |
| PH242103 | Saw Blade 30mm diameter |



Inge Lamina Spreader with teeth

| | |
|----------|-------|
| PH594282 | 170mm |
| PH594284 | 270mm |



Inge Lamina Spreader without teeth

| | |
|----------|-------|
| PH594283 | 170mm |
| PH594285 | 270mm |



Rogers Hysterectomy Clamp Straight

PH726173



Rogers Hysterectomy Clamp Curved

PH726174



Anderson Obstetric Forceps

PH726660 380mm



Kielland Obstetric Forceps for Occipito-Posterior Position

PH726662 410mm




Wrigley Obstetric Forceps Short Model for Caesarean Section

PH726670 280mm



Barnes Neville Obstetric Forceps with Axis Traction Rod

PH726672 370mm



Barnes Neville Obstetric Forceps without Axis Traction Rod

PH726673 370mm



Sieves

- PH760528 150mm (6 inch) diameter
- PH760529 200mm (8 inch) diameter



Fomon Bone Rasp Tungsten Carbide Cut 205mm

| | |
|----------|------------|
| PH654358 | No. 1 & 2 |
| PH654359 | No. 3 & 4 |
| PH654360 | No. 5 & 6 |
| PH654361 | No. 7 & 8 |
| PH654362 | No. 9 & 10 |



Fomon Bone Rasp Fine

| | |
|----------|-------|
| PH654363 | 215mm |
|----------|-------|



Goligher Rectal Speculum with Third Blade

| | |
|----------|------|
| PH675310 | 76mm |
|----------|------|



Eisenhammer Rectal Speculum

| |
|----------|
| PH675341 |
|----------|



St Marks Rectal Dilator Hollow Stainless Steel



| | |
|----------|--------|
| PH675360 | Large |
| PH675362 | Medium |
| PH675364 | Small |

Mathieu Urethral Grasping Forceps with Crocodile Action and Serrated Jaws



| | |
|----------|-------------|
| PH695506 | 200mm Shaft |
| PH695507 | 280mm Shaft |

Thomson Walker Penile Clamp with Cross Action



| | |
|----------|-------|
| PH695510 | 115mm |
|----------|-------|

Clutton Urethral Bougies 267mm


| | | | |
|----------|------------|----------|------------|
| PH695529 | 6 - 10 FG | PH695535 | 18 - 22 FG |
| PH695530 | 8 - 12 FG | PH695536 | 20 - 24 FG |
| PH695531 | 10 - 14 FG | PH695537 | 22 - 26 FG |
| PH695532 | 12 - 16 FG | PH695538 | 24 - 28 FG |
| PH695533 | 14 - 18 FG | PH695539 | 26 - 30 FG |
| PH695534 | 16 - 20 FG | PH695540 | 28 - 32 FG |



Lister Urethral Bougies 280mm

PH695543 Full Set of 12 Sizes

| | | | | | |
|----------|------------|----------|-----------|----------|------------|
| PH695544 | 0/1 E.G. | PH695550 | 4/7 E.G. | PH695556 | 10/13 E.G. |
| PH695545 | 0/2 E.G. | PH695551 | 5/8 E.G. | PH695557 | 11/14 E.G. |
| PH695546 | 0.5/3 E.G. | PH695552 | 6/9 E.G. | PH695558 | 12/15 E.G. |
| PH695547 | 1/4 E.G. | PH695553 | 7/10 E.G. | PH695559 | 13/16 E.G. |
| PH695548 | 2/5 E.G. | PH695554 | 8/11 E.G. | PH695560 | 14/17 E.G. |
| PH695549 | 3/6 E.G. | PH695555 | 9/12 E.G. | PH695561 | 15/18 E.G. |



Whyndham Powell Urethral Bougies

PH695575 Full Set of 21 Sizes

| | | | | | | | |
|----------|-------|----------|-------|----------|-------|----------|-------|
| PH695573 | 8 FG | PH695581 | 17 FG | PH695589 | 25 FG | PH695596 | 32 FG |
| PH695574 | 10 FG | PH695582 | 18 FG | PH695590 | 26 FG | PH695597 | 36 FG |
| PH695576 | 12 FG | PH695583 | 19 FG | PH695591 | 27 FG | PH695598 | 38 FG |
| PH695577 | 13 FG | PH695584 | 20 FG | PH695592 | 28 FG | PH695599 | 40 FG |
| PH695578 | 14 FG | PH695586 | 21 FG | PH695593 | 29 FG | | |
| PH695579 | 15 FG | PH695587 | 22 FG | PH695594 | 30 FG | | |
| PH695580 | 16 FG | PH695588 | 24 FG | PH695595 | 31 FG | | |



Foley Catheter Introducer Curved

| | | |
|----------|------------------|-------|
| PH705592 | 2mm diameter tip | 430mm |
| PH705593 | 4mm diameter tip | 430mm |



Guyon Catheter Introducer Straight

| | | |
|----------|------------------|-------|
| PH705594 | 2mm diameter tip | 430mm |
| PH705595 | 4mm diameter tip | 430mm |



Guyon Catheter Introducer Curved

| | | |
|----------|------------------|-------|
| PH705596 | 2mm diameter tip | 430mm |
| PH705597 | 4mm diameter tip | 430mm |



We also offer decontamination and sterilisation services.

With over 20 years' experience and more than 100 million instruments processed each year, we can offer a fully outsourced service or contingency for your sterile services department.



Contact us today to find out more
+44 (0)3452 414818

Umbilical Dilator

PH715725 95mm



Desjardin Gall Stone Forceps

PH715746 Small Jaws 230mm
PH715747 Large Jaws 230mm



Hegar Uterine Dilator 185mm Single Ended

| | | | | | | | | |
|----------|-------|---|----------|--------|---|----------|--------|---|
| PH726503 | 1.5mm | Ø | PH726534 | 10mm | Ø | PH726551 | 16mm | Ø |
| PH726504 | 2.5mm | Ø | PH726535 | 7.5mm | Ø | PH726552 | 17mm | Ø |
| PH726505 | 3.5mm | Ø | PH726536 | 8.5mm | Ø | PH726553 | 18mm | Ø |
| PH726506 | 4.5mm | Ø | PH726537 | 9.5mm | Ø | PH726557 | 16.5mm | Ø |
| PH726507 | 5.5mm | Ø | PH726538 | 10.5mm | Ø | PH726558 | 17.5mm | Ø |
| PH726508 | 6.5mm | Ø | PH726540 | 11mm | Ø | PH726560 | 19mm | Ø |
| PH726520 | 1mm | Ø | PH726541 | 12mm | Ø | PH726561 | 20mm | Ø |
| PH726521 | 2mm | Ø | PH726542 | 13mm | Ø | PH726562 | 21mm | Ø |
| PH726522 | 3mm | Ø | PH726543 | 14mm | Ø | PH726563 | 22mm | Ø |
| PH726523 | 4mm | Ø | PH726545 | 11.5mm | Ø | PH726570 | 23mm | Ø |
| PH726524 | 5mm | Ø | PH726546 | 12.5mm | Ø | PH726571 | 24mm | Ø |
| PH726530 | 6mm | Ø | PH726547 | 13.5mm | Ø | PH726572 | 25mm | Ø |
| PH726531 | 7mm | Ø | PH726548 | 14.5mm | Ø | PH726573 | 26mm | Ø |
| PH726532 | 8mm | Ø | PH726549 | 15.5mm | Ø | PH726574 | 18.5mm | Ø |
| PH726533 | 9mm | Ø | PH726550 | 15mm | Ø | PH726575 | 19.5mm | Ø |



Hegar Uterine Dilator 220mm Double Ended

PH726580 Hegar Uterine Dilators
Set of 8 Double Ended Dilators Consisting of:
3mm and 4mm - 5mm and 6mm
7mm and 8mm - 9mm and 10mm
11mm and 12mm - 13mm and 14mm
15mm and 16mm - 17mm and 18mm

PH726590 Hegar Uterine Dilators
Set Of 12 Double Ended Dilators Consisting of:
1mm and 2mm - 3mm and 4mm
5mm and 6mm - 7mm and 8mm
9mm and 10mm - 11mm and 12mm
13mm and 14mm - 15mm and 16mm
17mm and 18mm - 19mm and 20mm
21mm and 22mm - 23mm and 24mm

| | | |
|----------|---------------------|-------|
| PH726600 | Sizes 1mm and 2mm | 220mm |
| PH726601 | Sizes 3mm and 4mm | 220mm |
| PH726602 | Sizes 5mm and 6mm | 220mm |
| PH726603 | Sizes 7mm and 8mm | 220mm |
| PH726605 | Sizes 9mm and 10mm | 220mm |
| PH726606 | Sizes 11mm and 12mm | 220mm |
| PH726607 | Sizes 13mm and 14mm | 220mm |
| PH726610 | Sizes 15mm and 16mm | 220mm |
| PH726611 | Sizes 17mm and 18mm | 220mm |
| PH726612 | Sizes 19mm and 20mm | 220mm |
| PH726615 | Sizes 21mm and 22mm | 220mm |
| PH726616 | Sizes 23mm and 24mm | 220mm |
| PH726617 | Sizes 25mm and 26mm | 220mm |



Hawkin Ambler Uterine Dilator Conical 170mm



| | | | | |
|----------|---------|---|------|---|
| PH726625 | Size 3 | = | 6mm | Ø |
| PH726626 | Size 4 | = | 7mm | Ø |
| PH726627 | Size 5 | = | 8mm | Ø |
| PH726628 | Size 6 | = | 9mm | Ø |
| PH726629 | Size 7 | = | 10mm | Ø |
| PH726630 | Size 8 | = | 11mm | Ø |
| PH726631 | Size 9 | = | 12mm | Ø |
| PH726632 | Size 10 | = | 13mm | Ø |
| PH726633 | Size 11 | = | 14mm | Ø |
| PH726634 | Size 12 | = | 15mm | Ø |
| PH726635 | Size 13 | = | 16mm | Ø |
| PH726636 | Size 14 | = | 17mm | Ø |
| PH726637 | Size 15 | = | 18mm | Ø |
| PH726638 | Size 16 | = | 19mm | Ø |
| PH726639 | Size 17 | = | 20mm | Ø |
| PH726640 | Size 18 | = | 21mm | Ø |

Parametrium Clamps Straight with Box Joint

| | |
|----------|-------|
| PH726075 | 210mm |
| PH726076 | 240mm |
| PH726077 | 300mm |
| PH726078 | 350mm |



Parametrium Clamps Slightly Curved with Box Joint

| | |
|----------|-------|
| PH726079 | 210mm |
| PH726080 | 240mm |
| PH726081 | 300mm |
| PH726082 | 350mm |



Parametrium Clamps Strongly Curved with Box Joint

| | |
|----------|-------|
| PH726083 | 210mm |
| PH726084 | 240mm |
| PH726085 | 300mm |
| PH726086 | 350mm |



Parametrium Clamps Angulated with Box Joint

| | |
|----------|-------|
| PH726087 | 210mm |
| PH726088 | 240mm |
| PH726089 | 300mm |
| PH726090 | 350mm |



McClintock Ovum Forceps Curved with Screw Joint Curved

| | |
|----------|-------|
| PH726118 | 240mm |
|----------|-------|



Gwilliam Hysterectomy Clamp Straight Longitudinal Serrated Jaws 1 into 2 Teeth Box Joint

| | |
|----------|-------|
| PH726152 | 200mm |
|----------|-------|



Gwilliam Hysterectomy Clamp CVD Longitudinal Serrated Jaws 1 into 2 Teeth Box Joint



PH726153 200mm

McCullagh Hysterectomy Clamp McCullogh 250mm Straight



PH726162 250mm

McCullagh Hysterectomy Clamp McCullogh Curved



PH726163 250mm

Doyen Myoma Screw with Small Spiral for use per Vagina with D Ring Handle



PH726180 175mm

Doyen Myoma Screw with Large Spiral for use per Abdomen with D Ring Handle



PH726182 150mm

Sims Uterine Curette Double Ended

| | |
|----------|--------|
| PH726222 | Small |
| PH726223 | Medium |
| PH726224 | Large |



Sims Uterine Curette Double Ended Sharp/Blunt 255mm

| | |
|----------|--------|
| PH726226 | Small |
| PH726227 | Medium |
| PH726228 | Large |



Sims Uterine Curette Single Ended Sharp 255mm

| | | | |
|----------|-----|----------|------|
| PH726248 | 3mm | PH726252 | 9mm |
| PH726249 | 5mm | PH726253 | 12mm |
| PH726250 | 6mm | PH726254 | 15mm |
| PH726251 | 7mm | PH726255 | 19mm |



Sims Uterine Sound Malleable with CM Graduated Shaft and Silver Plated Finish Sharp/Sharp 255mm

| | |
|----------|-------|
| PH726335 | 320mm |
|----------|-------|



Galabin Uterine Sound Malleable with CM Graduated Shaft and Silver Plated Finish

| | |
|----------|-------|
| PH726337 | 320mm |
|----------|-------|



Cusco Vaginal Speculum Stainless Steel with Holding Screw

| | | | |
|----------|-----------------------------------|----------|---|
| PH726370 | Virgin 86mm x 16mm at Distal End | PH726373 | Medium Long 115mm x 30mm at Distal End |
| PH726371 | Small 95mm x 25mm at Distal End | PH726374 | Large 110mm x 35mm at Distal End |
| PH726372 | Medium 102mm x 30mm at Distal End | PH726375 | Large Extra Long 115mm x 35mm at Distal End |



Cusco Vaginal Speculum with Holding Screw and Smoke Tube with Bead Finish

| | | | |
|----------|-----------------------------------|----------|--|
| PH726376 | Virgin 86mm x 16mm at Distal End | PH726379 | Medium Long 115mm x 30mm at Distal End |
| PH726377 | Small 95mm x 25mm at Distal End | PH726380 | Large 110mm x 35mm at Distal End |
| PH726378 | Medium 102mm x 30mm at Distal End | PH726381 | Extra Long 115mm x 35mm at Distal End |



Cusco Vaginal Speculum with Holding Screw and Smoke Tube with Black Finish

| | |
|----------|---|
| PH726382 | Virgin 86mm x 16mm at Distal End |
| PH726383 | Small 95mm x 25mm at Distal End |
| PH726384 | Medium 102mm x 30mm at Distal End |
| PH726385 | Medium Long 115mm x 30mm at Distal End |
| PH726386 | Large 110mm x 35mm at Distal End |
| PH726387 | Large Extra Long 115mm x 35mm at Distal End |



Cusco Vaginal Speculum with Holding Screw and Smoke Tube with Diathermy Insulated Finish

| | |
|----------|---|
| PH726388 | Virgin 86mm x 16mm at Distal End |
| PH726389 | Small 95mm x 25mm at Distal End |
| PH726390 | Medium 102mm x 30mm at Distal End |
| PH726391 | Medium Long 115mm x 30mm at Distal End |
| PH726392 | Large 110mm x 35mm at Distal End |
| PH726393 | Large Extra Long 115mm x 35mm at Distal End |



Graves (Brewer) Vaginal Speculum

| | |
|----------|--------------------|
| PH726395 | Small 75mm x 20mm |
| PH726396 | Medium 95mm x 35mm |
| PH726397 | Large 115mm x 32mm |



Kogan Endocervical Speculum without Ratchet (KOOGAN)

PH726400 240mm



Koogan Endocervical Speculum with Calibrated Adjustable Ratchet and Stop (KOOGAN)

PH726401 240mm



Sims Vaginal Speculum Double Ended

| | |
|----------|--------------------------------------|
| PH726420 | Small 60mm x 25mm/70mm x 29mm 130mm |
| PH726421 | Medium 70mm x 29mm/75mm x 33mm 155mm |
| PH726422 | Large 75mm x 33mm/85mm x 38mm 185mm |



Sims Vaginal Speculum Single Ended

| | |
|----------|--------------------------|
| PH726423 | Small 76mm X 21mm 180mm |
| PH726424 | Medium 86mm X 29mm 180mm |
| PH726425 | Large 115mm X 43mm 180mm |



Auvard Vaginal Speculum Tapered Blade with Fixed Weight and Chrome Plate Finish

| | | |
|----------|---------------------------|------------------|
| PH726439 | Tapered Blade 80mm x 45mm | Tapering to 40mm |
|----------|---------------------------|------------------|



Bailey Vaginal Speculum with Fixed Weight and Chrome Plate Finish

| | |
|----------|------------|
| PH726442 | 38mm Blade |
|----------|------------|



Winterton Vaginal Speculum

| |
|----------|
| PH726445 |
|----------|





Winterton Vaginal Speculum with
Smoke Tube and Bead Finish

PH726446



Winterton Vaginal Speculum with
Smoke Tube and Black Finish

PH726447



Winterton Vaginal Speculum Diathermy
Insulated with Smoke Tube

PH726448



Doyen Vaginal Speculum

| | |
|----------|--------------------------|
| PH726450 | Small 55mm x 45mm 240mm |
| PH726451 | Medium 85mm x 45mm 240mm |
| PH726452 | Large 115mm x 45mm 240mm |



Jessop Hospital Pattern Vaginal Speculum Double Ended

| | |
|----------|-------------------------------|
| PH726455 | 100mm x 28mm and 120mm x 30mm |
|----------|-------------------------------|



St Georges Pattern Vaginal Speculum Double Ended

| | |
|----------|-------------------------------|
| PH726465 | 100mm x 25mm and 120mm x 30mm |
|----------|-------------------------------|



Breisky Vaginal Speculum with Cranked Head Blade

| | | | | | |
|----------|---------------------------|-------|----------|---------------------------|-------|
| PH726480 | 100mm Length X 20mm Width | 300mm | PH726486 | 130mm Length X 25mm Width | 320mm |
| PH726481 | 100mm Length X 25mm Width | 300mm | PH726487 | 130mm Length X 30mm Width | 320mm |
| PH726482 | 100mm Length X 30mm Width | 300mm | PH726488 | 130mm Length X 35mm Width | 320mm |
| PH726483 | 100mm Length X 35mm Width | 300mm | PH726489 | 130mm Length X 40mm Width | 320mm |
| PH726484 | 100mm Length X 40mm Width | 300mm | PH726490 | 160mm Length X 40mm Width | 350mm |
| PH726485 | 130mm Length X 20mm Width | 320mm | PH726491 | 180mm Length X 40mm Width | 350mm |



Sims Diathermy Vaginal Speculum Double Ended

| | |
|----------|--------------------|
| PH757550 | Small 60mm X 30mm |
| PH757551 | Medium 70mm X 35mm |
| PH757552 | Large 80mm X 40mm |



Tischler Biopsy Punch Forceps Crocodile Action 210mm Effective Length with Single Ring Handle Grip Sharp/Sharp 255mm

PH727050 3 x 7mm 210mm



Baby Tischler Biopsy Punch Forceps Crocodile Action 210mm Effective Length with Single Ring Handle Grip

PH727051 2.3 x 5mm 210mm



Mini Townsend Biopsy Punch Forceps Crocodile Action 210mm Effective Length with Single Ring Handle Grip

PH727055 2.3 x 5mm 210mm



Uni Tischler Biopsy Punch Forceps Crocodile Action 230mm Effective Length with Rotating Punch Complete Straight Cut

PH727065 3 x 7mm 230mm



Uni Tischler Biopsy Punch Forceps Crocodile Action 230mm Effective Length with Rotating Punch Complete Upward Angle Cut

PH727066 3 x 7mm 230mm



Uni Eppendorf Biopsy Punch Forceps Crocodile Action 230mm Effective Length with Rotating Punch Complete Straight Cut

PH727079 3 x 6mm 230mm



Schumacher Biopsy Punch For Gynaecology Crocodile Action 230mm Effective Length Complete STR Cut

PH727084 4.5 x 7mm 230mm



Kevorkian Biopsy Punch for Gynaecology Crocodile Action 260mm Effective Length Complete Straight Cut

PH727086 3 x 10mm 230mm



Eppendorf Biopsy Punch for Gynaecology Crocodile Action 200mm Effective Length Complete Straight Cut

PH727089 4 x 8mm 200mm





Gigli Metal Handles (Supplied in a pair)

PH737700



Gigli Saw

| | |
|----------|-------|
| PH737701 | 300mm |
| PH737704 | 400mm |
| PH737705 | 500mm |
| PH737706 | 600mm |



Sterilisation Baskets & Trays


Urethroscope Baskets

PH760001 640mm x 150mm x 77mm



Endoscope Baskets Standard

PH760006 460mm x 80mm x 52mm



Arthroscope Baskets

For Two Arthroscopes

PH760004 290mm X 80mm X 52mm

For One Arthroscope

PH760005 290mm X 80mm X 52mm


PH760006 460mm X 80mm X 52mm



Double Scope Baskets

PH760003 290mm x 160mm x 52mm

PH760009 460mm x 160mm x 52mm




Endoscope Baskets (for two Endoscopes) Standard

PH760007 460mm x 80mm x 52mm




Endoscope Baskets Long

PH760008 670mm x 80mm x 52mm



Basket for Ophthalmic Instruments

PH760270 440mm x 250mm x 40mm



Universal Scope Basket

PH760002 310mm x 150mm x 77mm

PH760010 580mm x 160mm x 52mm


PH760011 290mm x 160mm x 52mm

PH760012 460mm x 160mm x 52mm




Flat Base Baskets with Perforated Sides

| | | | |
|----------|-----------------------|----------|-----------------------|
| PH760020 | 230mm x 130mm x 50mm | PH760057 | 480mm x 350mm x 60mm |
| PH760021 | 230mm x 130mm x 60mm | PH760058 | 480mm x 480mm x 50mm |
| PH760022 | 230mm x 130mm x 70mm | PH760059 | 480mm x 480mm x 60mm |
| PH760023 | 240mm x 130mm x 50mm | PH760060 | 480mm x 480mm x 70mm |
| PH760024 | 240mm x 130mm x 60mm | PH760061 | 480mm x 480mm x 100mm |
| PH760025 | 280mm x 200mm x 60mm | PH760062 | 480mm x 480mm x 120mm |
| PH760030 | 240mm x 250mm x 50mm | PH760065 | 500mm x 500mm x 50mm |
| PH760031 | 240mm x 250mm x 60mm | PH760066 | 520mm x 250mm x 50mm |
| PH760032 | 240mm x 250mm x 70mm | PH760067 | 520mm x 250mm x 60mm |
| PH760033 | 240mm x 250mm x 100mm | PH760068 | 520mm x 250mm x 75mm |
| PH760036 | 300mm x 300mm x 50mm | PH760069 | 520mm x 250mm x 95mm |
| PH760037 | 300mm x 300mm x 100mm | PH760070 | 520mm x 300mm x 60mm |
| PH760038 | 320mm x 320mm x 60mm | PH760075 | 540mm x 250mm x 95mm |
| PH760039 | 380mm x 240mm x 50mm | PH760076 | 550mm x 250mm x 60mm |
| PH760040 | 350mm x 350mm x 60mm | PH760077 | 550mm x 280mm x 50mm |
| PH760043 | 450mm x 230mm x 50mm | PH760078 | 550mm x 325mm x 50mm |
| PH760044 | 460mm x 320mm x 60mm | PH760079 | 550mm x 430mm x 50mm |
| PH760045 | 480mm x 120mm x 60mm | PH760085 | 600mm x 250mm x 150mm |
| PH760046 | 480mm x 150mm x 50mm | PH760086 | 600mm x 300mm x 60mm |
| PH760047 | 480mm x 250mm x 50mm | PH760087 | 600mm x 350mm x 60mm |
| PH760048 | 480mm x 250mm x 60mm | PH760088 | 670mm x 300mm x 70mm |
| PH760049 | 480mm x 250mm x 70mm | PH760089 | 700mm x 250mm x 60mm |
| PH760050 | 480mm x 250mm x 100mm | PH760090 | 700mm x 300mm x 90mm |
| PH760051 | 480mm x 250mm x 120mm | PH760091 | 900mm x 450mm x 50mm |
| PH760056 | 470mm x 400mm x 60mm | | |




Lids for Flat Base Baskets

| | |
|----------|---------------|
| PH760093 | 520mm X 250mm |
| PH760094 | 540mm X 250mm |
| PH760095 | 240mm x 250mm |
| PH760096 | 480mm x 250mm |
| PH760097 | 480mm x 480mm |
| PH760098 | 300mm x 300mm |
| PH760099 | 670mm x 330mm |




All Wire Mesh Boxes




| | |
|----------|---------------------|
| PH760140 | 240mm x 250mm x 50 |
| PH760141 | 240mm x 250mm x 60 |
| PH760142 | 240mm x 250mm x 70 |
| PH760143 | 480mm x 250mm x 50 |
| PH760144 | 480mm x 250mm x 60 |
| PH760145 | 480mm x 250mm x 70 |
| PH760146 | 480mm x 250mm x 100 |
| PH760147 | 480mm x 250mm x 120 |
| PH760148 | 480mm x 480mm x 50 |
| PH760149 | 480mm x 480mm x 60 |
| PH760150 | 480mm x 480mm x 70 |
| PH760151 | 300mm x 300mm x 60 |
| PH760152 | 470mm x 400mm x 60 |

Narrow Box Perforated Sides and Wire Mesh Lid



| | |
|----------|----------------------|
| PH760194 | 300mm x 80mm x 35mm |
| PH760195 | 360mm x 50mm x 25mm |
| PH760196 | 550mm x 50mm x 25mm |
| PH760197 | 720mm x 100mm x 25mm |
| PH760198 | 700mm x 50mm x 25mm |


All Fine Mesh Boxes



| | |
|----------|----------------------|
| PH760200 | 80mm x 40mm x 30mm |
| PH760201 | 80mm x 80mm x 35mm |
| PH760202 | 80mm x 80mm x 55mm |
| PH760203 | 100mm x 100mm x 55mm |
| PH760204 | 105mm x 70mm x 25mm |
| PH760206 | 150mm x 50mm x 35mm |
| PH760207 | 150mm x 70mm x 35mm |
| PH760208 | 150mm x 150mm x 35mm |
| PH760210 | 190mm x 85mm x 100mm |
| PH760211 | 190mm x 70mm x 70mm |
| PH760215 | 200mm x 100mm x 35mm |
| PH760216 | 220mm x 140mm x 50mm |
| PH760217 | 235mm x 80mm x 30mm |
| PH760218 | 235mm x 155mm x 40mm |
| PH760222 | 240mm x 90mm x 40mm |
| PH760224 | 240mm x 120mm x 35mm |
| PH760227 | 275mm x 175mm x 35mm |
| PH760228 | 280mm x 200mm x 60mm |


Fine Mesh Basket with Take Off Lid, Hinges and Locking Clip

| | |
|----------|----------------------|
| PH760255 | 120mm x 155mm x 40mm |
| PH760256 | 120mm x 155mm x 50mm |
| PH760259 | 235mm x 155mm x 40mm |
| PH760260 | 235mm x 155mm x 50mm |
| PH760264 | 380mm x 265mm x 45mm |
| PH760265 | 380mm x 265mm x 55mm |



Trinket Boxes Perforated

| | |
|----------|---------------------|
| PH760175 | 70mm x 30mm x 15mm |
| PH760176 | 90mm x 60mm x 25mm |
| PH760177 | 105mm x 70mm x 25mm |
| PH760178 | 120mm x 80mm x 25mm |



Trinket Boxes Perforated Sides and Base with Wire Mesh Lid

| | |
|----------|----------------------|
| PH760185 | 90mm x 60mm x 25mm |
| PH760186 | 105mm x 70mm x 25mm |
| PH760187 | 120mm x 80mm x 25mm |
| PH760189 | 150mm x 100mm x 25mm |
| PH760190 | 180mm x 130mm x 25mm |



Lid for All Wire Mesh Basket

| | |
|----------|---------------|
| PH760170 | 240mm X 250mm |
| PH760171 | 480mm X 250mm |
| PH760172 | 300mm X 300mm |
| PH760173 | 480mm X 250mm |
| PH760174 | 480mm X 480mm |

We also offer decontamination and sterilisation services.

With over 20 years' experience and more than 100 million instruments processed each year, we can offer a fully outsourced service or contingency for your sterile services department.



Contact us today to find out more
+44 (0)3452 414818

Receiver / Kidney Dish (without Lid)

| | | |
|-----------------|-------------|---------|
| Without Lid | PHS800049 | 8 inch |
| Without Lid | PHS800053 | 10 inch |
| Non Rolled Lips | PHS800049NR | 8 inch |
| Non Rolled Lips | PHS800053NR | 10 inch |



Gallipot Non Rolled Lips

| | |
|-------------|-------------|
| PHS800055NR | 80mm x 60mm |
|-------------|-------------|



Gallipot

| | |
|-----------|-------------|
| PHS800056 | 60mm x 50mm |
| PHS800056 | 4oz |



Product Codes Index

| Code | Description | Page |
|----------|--|------|
| PH100901 | Presbyterian Tubing Clamp Forceps With Serrated Jaws And Box Joint 165mm Straight | 76 |
| PH100902 | Presbyterian Tubing Clamp Forceps With Serrated Jaws And Box Joint 195mm Straight | 76 |
| PH100903 | St Bartholomews Forceps For Tubing Clamp With Square Jaws And Box Joint 150mm Straight | 76 |
| PH100904 | Square Jaw Tubing Clamp Straight 130mm | 76 |
| PH100905 | Square Jaw Tubing Clamp Angled 115mm | 76 |
| PH100906 | Square Jaw Tubing Clamp Straight 180mm | 76 |
| PH100907 | Square Jaw Tubing Clamp Angled 180mm | 76 |
| PH100908 | Square Jaw Tubing Clamp Straight 200mm | 76 |
| PH100915 | X-Ray Bag Clip 1inch | 77 |
| PH100917 | Zimmer Clip 2.5inch(63.5mm) Wide With Serrated Jaws. Suitable For Holding Towelling | 77 |
| PH100925 | Safety Pin For Holding Ring Handle Forceps 115mm | 77 |
| PH100926 | Safety Pin For Holding Ring Handle Forceps 130mm | 77 |
| PH100927 | Safety Pin For Holding Ring Handle Forceps 140mm | 77 |
| PH100928 | Safety Pin For Holding Ring Handle Forceps 150mm | 77 |
| PH110940 | Dissecting Forceps With Bayonet Shanks And Serrated Jaws 180mm | 52 |
| PH110941 | Dissecting Forceps With Bayonet Shanks And Serrated Jaws 200mm | 52 |
| PH110950 | Dissecting Forceps With Turned Over End And Serrated Jaws 130mm | 52 |
| PH110951 | Dissecting Forceps With Turned Over End And Serrated Jaws 150mm | 52 |
| PH110952 | Dissecting Forceps With Turned Over End And Serrated Jaws 180mm | 52 |
| PH110953 | Dissecting Forceps With Turned Over End And Serrated Jaws 200mm | 52 |
| PH110954 | Dissecting Forceps With Turned Over End And Serrated Jaws 230mm | 52 |
| PH110955 | Dissecting Forceps With Turned Over End And Serrated Jaws 250mm | 52 |
| PH110956 | Dissecting Forceps With Turned Over End And Serrated Jaws 300mm | 52 |
| PH110960 | Dissecting Forceps With Block End And Serrated Jaws With A Medium Point 130mm | 52 |
| PH110961 | Dissecting Forceps With Block End And Serrated Jaws With A Medium Point 150mm | 52 |
| PH110962 | Dissecting Forceps With Block End And Serrated Jaws With A Medium Point 180mm | 52 |
| PH110963 | Dissecting Forceps With Block End And Serrated Jaws With A Medium Point 200mm | 52 |
| PH110964 | Dissecting Forceps With Block End And Serrated Jaws With A Medium Point 230mm | 52 |
| PH110965 | Dissecting Forceps With Block End And Serrated Jaws With A Medium Point 250mm | 52 |
| PH110966 | Dissecting Forceps With Block End And Serrated Jaws With A Medium Point 300mm | 52 |
| PH110967 | Continental Pattern Dissecting Forceps With Serrated Jaws And Narrow End 130mm | 52 |
| PH110968 | Continental Pattern Dissecting Forceps With Serrated Jaws And Narrow End 140mm | 52 |

Product Codes Index

| Code | Description | Page |
|------------|--|------|
| PH110969 | Continental Pattern Dissecting Forceps With Serrated Jaws And Narrow End 160mm | 52 |
| PH110970 | Continental Pattern Dissecting Forceps With Serrated Jaws And Narrow End 180mm | 52 |
| PH110971 | Continental Pattern Dissecting Forceps With Serrated Jaws And Narrow End 200mm | 52 |
| PH110972 | Continental Pattern Dissecting Forceps With Serrated Jaws And Narrow End 250mm | 52 |
| PH110982 | Continental Pattern Dissecting Forceps With 1 Into 2 Teeth And Broad End 130mm | 52 |
| PH110983 | Continental Pattern Dissecting Forceps With 1 Into 2 Teeth And Broad End 140mm | 52 |
| PH110984 | Continental Pattern Dissecting Forceps With 1 Into 2 Teeth And Broad End 160mm | 52 |
| PH110985 | Continental Pattern Dissecting Forceps With 1 Into 2 Teeth And Broad End 180mm | 52 |
| PH110986 | Continental Pattern Dissecting Forceps With 1 Into 2 Teeth And Broad End 200mm | 52 |
| PH110987 | Continental Pattern Dissecting Forceps With 1 Into 2 Teeth And Broad End 250mm | 52 |
| PH110988 | Continental Pattern Dissecting Forceps With 1 Into 2 Teeth And Broad End 300mm | 52 |
| PH111016 | Adson Micro Dissecting Forceps With Serrated Jaws 130mm Straight | 52 |
| PH111017 | Adson Micro Dissecting Forceps With 1 Into 2 Teeth 130mm Straight | 52 |
| PH111019 | Adson Browne Dissecting Forceps With Serrated Jaws 130mm 130mm Straight | 53 |
| PH111020 | Adson Dissecting Forceps With Serrated Jaws 130mm Straight | 53 |
| PH111021 | Adson Dissecting Forceps With 1 Into 2 Teeth 130mm Straight | 53 |
| PH111021DA | Double Action Adson Dissecting Fcps 1-2 Teeth 13cm | 53 |
| PH111022 | Adson Dissecting Forceps With Serrated Jaws 180mm Straight | 53 |
| PH111023 | Adson Dissecting Forceps With 1 Into 2 Teeth 180mm Straight | 53 |
| PH111032 | Bickford Dissecting Forceps With Serrated Jaws 230mm | 53 |
| PH111033 | Bickford Dissecting Forceps With 1 Into 2 Teeth 230mm Straight | 54 |
| PH111035 | Bonney Dissecting Forceps With Serrated Jaws 180mm Straight | 54 |
| PH111036 | Bonney Dissecting Forceps With 1 Into 2 Teeth 180mm Straight | 54 |
| PH111045 | Cairn Dissecting Forceps With 2 Into 3 Teeth 165mm | 54 |
| PH111046 | Cairn Dissecting Forceps With 2 Into 3 Teeth 205mm | 54 |
| PH111051 | Charnley Suture Dissecting Forceps With Shutter And Slot 145mm | 54 |
| PH111052 | Charnley Suture Dissecting Forceps With Shutter And Slot 180mm 180mm Straight | 54 |
| PH111057 | College Pattern Forceps For Dental Tweezer 150mm Curved | 54 |
| PH111059 | Debakey Forceps For Tissue Dissecting With Atraumatic Jaws 150mm Straight | 55 |
| PH111060 | Debakey Forceps For Tissue Dissecting With Atraumatic Jaws 180mm Straight | 55 |
| PH111061 | Debakey Tissue Dissecting Fcps Atraumatic Jaw 20cm | 55 |
| PH111062 | Debakey Forceps For Tissue Dissecting With Atraumatic Jaws 230mm Straight | 55 |
| PH111064 | Dott Dissecting Forceps With Serrated Jaws 180mm Straight | 55 |

Product Codes Index

| Code | Description | Page |
|----------|--|------|
| PH111175 | Emmat Dissecting Forceps With Serrated Jaws 200mm Straight | 55 |
| PH111176 | Emmat Dissecting Forceps With Serrated Jaws 240mm Straight | 55 |
| PH111179 | Gillies Dissecting Forceps With Serrated Jaws 150mm Straight | 56 |
| PH111180 | Gillies Dissecting Forceps With 1 Into 2 Teeth 150mm Straight | 56 |
| PH111185 | Iris Dissecting Forceps With Serrated Jaws 100mm Str | 56 |
| PH111186 | Iris Dissecting Forceps With 1 Into 2 Teeth 100mm Straight | 56 |
| PH111187 | Iris Dissecting Forceps With Serrated Jaws 100mm Curved | 56 |
| PH111189 | Iris Dissecting Forceps With Serrated Jaws 115mm Straight | 56 |
| PH111190 | Iris Dissecting Forceps With 1 Into 2 Teeth 115mm Straight | 56 |
| PH111191 | Iris Dissecting Forceps With Serrated Jaws 115mm Curved | 56 |
| PH111195 | Jackson Burrows Dissecting Forceps Serrated With 1 Into 2 Teeth 140mm Straight | 56 |
| PH111202 | Jean Dissecting Forceps With 1 Into 2 Teeth 140mm Straight | 57 |
| PH111203 | Jean Dissecting Forceps With 1 Into 2 Teeth 180mm Straight | 57 |
| PH111211 | Jefferson Dissecting Forceps With 1 Into 2 Teeth 180mm Straight | 57 |
| PH111220 | Kurihara Endo Cervical Dissecting Forceps With Satin Finish 250mm Straight | 57 |
| PH111240 | Lane Dissecting Forceps With 1 Into 2 Teeth 130mm (Lanes) 130mm Straight | 57 |
| PH111241 | Lane Dissecting Forceps With 1 Into 2 Teeth (Lanes) 140mm Straight | 57 |
| PH111242 | Lane Dissecting Forceps With 1 Into 2 Teeth (Lanes) 150mm Straight | 57 |
| PH111244 | Lane Dissecting Forceps With 1 Into 2 Teeth (Lanes) 200mm Straight | 57 |
| PH111245 | Lane Dissecting Forceps With 2 Into 3 Teeth 140mm (Lanes) 140mm Straight | 57 |
| PH111246 | Lane Dissecting Forceps With 2 Into 3 Teeth 150mm (Lanes) 150mm Straight | 57 |
| PH111247 | Lane Dissecting Forceps With 2 Into 3 Teeth 180mm (Lanes) 180mm Straight | 57 |
| PH111248 | Lane Dissecting Forceps With 2 Into 3 Teeth 200mm (Lanes) 200mm Straight | 57 |
| PH111249 | Lane Dissecting Forceps With 3 Into 4 Teeth 150mm (Lanes) 150mm Straight | 57 |
| PH111250 | Lane Dissecting Forceps With 3 Into 4 Teeth (Lanes) 180mm Straight | 57 |
| PH111251 | Lane Dissecting Forceps With 3 Into 4 Teeth (Lanes) 200mm Straight | 57 |
| PH111264 | Mcindoe Dissecting Forceps With Serrated Jaws 150mm Straight | 58 |
| PH111265 | Mcindoe Dissecting Forceps With 1 Into 2 Teeth 150mm Straight | 58 |
| PH111268 | Moorfield Conjunctival Dissecting Forceps 115mm Straight | 58 |
| PH111280 | Officer Dissecting Forceps With 3 Into 4 Teeth 200mm Straight | 58 |
| PH111281 | Officer Dissecting Forceps With 3 Into 4 Teeth 250mm Straight | 58 |
| PH111290 | Ramsay Dissecting Forceps With Serrated Jaws (Ramsey) 180mm Straight | 58 |
| PH111291 | Ramsay Dissecting Forceps With 1 Into 2 Teeth (Ramsey) 180mm Straight | 58 |

Product Codes Index

| Code | Description | Page |
|----------|---|------|
| PH111292 | Ramsay Dissecting Forceps With 2 Into 3 Teeth (Ramsey) 180mm Straight | 59 |
| PH111295 | Russian Pattern Dissecting Forceps Oval Toothed End 150mm Straight | 59 |
| PH111296 | Russian Pattern Dissecting Forceps Oval Toothed End 200mm Straight | 59 |
| PH111297 | Russian Pattern Dissecting Forceps Oval Toothed End 250mm Straight | 59 |
| PH111298 | Russian Dissecting 180mm Forceps | 59 |
| PH111301 | St Marks Hospital Dissecting Forceps With 1 Into 2 Teeth 290mm Straight | 59 |
| PH111310 | Treves Dissecting Forceps With 1 Into 2 Teeth 130mm 130mm Straight | 59 |
| PH111311 | Treves Dissecting Forceps With 1 Into 2 Teeth 150mm Straight | 59 |
| PH111312 | Treves Dissecting Forceps With 1 Into 2 Teeth 180mm Straight | 59 |
| PH111313 | Treves Dissecting Forceps With 1 Into 2 Teeth 200mm Straight | 59 |
| PH111315 | Waughs Dissecting Forceps With Serrated Jaws 150mm (Tonsil) Straight | 59 |
| PH111316 | Waughs Dissecting Forceps With Serrated Jaws 180mm Straight | 59 |
| PH111317 | Waughs Dissecting Forceps With Serrated Jaws 200mm Straight | 59 |
| PH111318 | Waughs Dissecting Forceps With 1 Into 2 Teeth 150mm Straight | 59 |
| PH111319 | Waughs Dissecting Forceps With 1 Into 2 Teeth 180mm 180mm Straight | 59 |
| PH111320 | Waughs Dissecting Forceps With 1 Into 2 Teeth 200mm Straight | 59 |
| PH111330 | Adson Dissecting Forceps With Tungsten Carbide Jaws 120mm Straight | 60 |
| PH111331 | Adson Dissecting Forceps Tungsten Carbide Jaws 1:2 Teeth 120mm Straight | 60 |
| PH111350 | Debakey Dissecting Forceps With Tungsten Carbide Jaws 150mm Straight | 60 |
| PH111351 | Debakey Dissecting Forceps With Tungsten Carbide Jaws 200mm Straight | 60 |
| PH111352 | Debakey Dissecting Forceps With Tungsten Carbide Jaws 230mm Straight | 60 |
| PH121353 | Whitfield Epilation Forceps With Angled (Oblique) End 90mm Angled | |
| PH131405 | Martin Splinter Forceps With Box Joint 75mm Straight | 60 |
| PH131406 | Martin Splinter Forceps With Box Joint 115mm Straight | 60 |
| PH131407 | Martin Splinter Forceps With Box Joint 130mm Straight | 60 |
| PH141455 | Lister Sinus Forceps With Box Joint 130mm Straight | 60 |
| PH141456 | Lister Sinus Forceps With Box Joint 150mm Straight | 60 |
| PH141457 | Lister Sinus Forceps With Box Joint 180mm Straight | 60 |
| PH141458 | Lister Sinus Forceps With Box Joint 200mm Straight | 60 |
| PH161482 | Ramsey Dressing And Sponge Holding Forceps Box Joint 180mm Straight | 61 |
| PH161483 | Ramsey Dressing And Sponge Holding Forceps Box Joint 240mm Straight | 61 |
| PH161866 | Denis Browne Suture Forceps For Collar/Suture Crushing 145mm Straight | 61 |
| PH161867 | Guys Galabin / Suture Forceps 140mm Angled (Gallabin) | 61 |

Product Codes Index

| Code | Description | Page |
|----------|---|------|
| PH171550 | Backhaus Towel Clip With Box Joint 90mm Curved | 78 |
| PH171551 | Backhaus Towel Clip With Box Joint 115mm Curved | 78 |
| PH171552 | Backhaus Towel Clip With Box Joint 135mm Curved | 78 |
| PH171553 | Mayo (Backhaus) Towel Clip With Box Joint 135mm Curved | 78 |
| PH171558 | Towel Clip Non-Tearing With Ball And Socket. Box Joint 90mm Curved | 78 |
| PH171559 | Towel Clip Non-Tearing With Ball And Socket. Box Joint 100mm Curved | 78 |
| PH171560 | Towel Clip Non-Tearing With Ball And Socket. Box Joint 115mm Curved | 78 |
| PH171561 | Towel Clip Non-Tearing With Ball And Socket. Box Joint 140mm Curved | 78 |
| PH171562 | Towel Clip Non-Tearing With Serrated Jaws And Box Joint 100mm Curved | 78 |
| PH171563 | Towel Clip Non-Tearing With Serrated Jaws And Box Joint 115mm Curved | 78 |
| PH171564 | Towel Clip Non-Tearing With Serrated Jaws And Box Joint 140mm Curved | 78 |
| PH171565 | Towel Clip English Pattern With Cross Action 90mm | 78 |
| PH171567 | Schaedel Towel Clip With Cross Action 90mm | 78 |
| PH171574 | Lorna (Edna) Towel Clip With Box Joint 140mm Curved | 79 |
| PH171585 | Robin Anchoring Towel Clip With A 5mm Diameter Clip And Screw Joint 130mm Curved | 79 |
| PH171587 | Robin Anchoring Towel Clip With A 13mm Diameter Clip And Screw Joint 130mm Curved | 79 |
| PH171589 | Robin Anchoring Towel Clip With A 25mm Diameter Clip And Screw Joint 130mm Curved | 79 |
| PH181650 | Baumgartner Needle Holder With Box Joint 140mm Straight 140mm Straight | 82 |
| PH181684 | Crile Wood Needle Holder With Box Joint (Crilewood) 150mm Straight | 82 |
| PH181691 | Derf Needle Holder With Box Joint 120mm Straight | 82 |
| PH181706 | Foster Needle Holder With Scissor And Screw Joint 135mm Curved | 82 |
| PH181707 | Fry Needle Holder With Box Joint 115mm Straight | 82 |
| PH181710 | Gillies Needle Holder With Right Hand Scissor Screw Joint 160mm Curved | 82 |
| PH181712 | Halsey Needle Holder With Plain Jaws And Box Joint 130mm Straight | 83 |
| PH181713 | Halsey Needle Holder With Serrated Jaws And Box Joint 130mm Straight | 83 |
| PH181716 | Kilner Needle Holder With Box Joint With Curved Shanks 135mm Straight | 83 |
| PH181722 | Lawrence Needle Holder With Plain Jaws And Box Joint 150mm Straight | 83 |
| PH181724 | Lawson Tait Needle Holder With Box Joint 130mm Straight | 83 |
| PH181770 | Mayo Hegar Needle Holder With Box Joint 130mm Straight | 84 |
| PH181771 | Mayo Hegar Needle Holder With Box Joint 140mm Straight | 84 |
| PH181772 | Mayo Hegar Needle Holder With Box Joint 150mm Straight | 84 |
| PH181773 | Mayo Hegar Needle Holder With Box Joint 165mm Straight | 84 |
| PH181774 | Mayo Hegar Needle Holder With Box Joint 180mm Straight | 84 |

Product Codes Index

| Code | Description | Page |
|----------|--|------|
| PH181775 | Mayo Hegar Needle Holder With Box Joint 190mm Straight | 84 |
| PH181776 | Mayo Hegar Needle Holder With Box Joint 200mm Straight | 84 |
| PH181777 | Mayo Hegar Needle Holder With Box Joint 240mm Straight | 84 |
| PH181778 | Mayo Hegar Needle Holder With Box Joint 300mm Straight | 84 |
| PH181785 | Naunton-Morgan Needle Holder With Multi Action 180mm Straight | 84 |
| PH181786 | Naunton-Morgan Needle Holder With Multi Action 250mm Straight | 84 |
| PH181787 | Nievert Needle Holder With Offset Ring For Thumb And Box Joint 130mm Straight | 84 |
| PH181790 | Olsen - Hegar Needle Holder With Scissor And Screw Joint 135mm Straight | 84 |
| PH181791 | Olsen - Hegar Needle Holder With Scissors And Screw Joint 160mm Straight | 84 |
| PH181821 | Thomson Walker Needle Holder Curved Shanks And Box Joint 200mm Straight | 84 |
| PH181828 | Wright Needle Holder With Box Joint 130mm Straight | 85 |
| PH191550 | Baumgartner Needle Holder With Tungsten Carbide Jaws And Box Joint 135mm Straight | 85 |
| PH191555 | Bozemann Needle Holder With Tungsten Carbide Jaws And Box Joint With Curved Shanks 200mm | 85 |
| PH191560 | Crile Wood Needle Holder With Tungsten Carbide Jaws And Box Joint 140mm | 85 |
| PH191561 | Crile Wood Needle Holder With Tungsten Carbide Jaws And Box Joint 150mm | 85 |
| PH191562 | Crile Wood Needle Holder With Tungsten Carbide Jaws And Box Joint 180mm | 85 |
| PH191563 | Crile Wood Needle Holder With Tungsten Carbide Jaws And Box Joint 200mm | 85 |
| PH191564 | Crile Wood Needle Holder With Tungsten Carbide Jaws And Box Joint 230mm | 85 |
| PH191565 | Crile Wood Needle Holder With Tungsten Carbide Jaws And Box Joint 250mm | 85 |
| PH191566 | Crile Wood Needle Holder With Tungsten Carbide Jaws And Box Joint 300mm | 85 |
| PH191570 | Debakey Needle Holder With Tungsten Carbide Jaws And Box Joint 150mm Straight | 85 |
| PH191571 | Debakey Needle Holder With Tungsten Carbide Jaws And Box Joint 180mm Straight | 85 |
| PH191572 | Debakey Needle Holder With Tungsten Carbide Jaws And Box Joint 200mm Straight | 85 |
| PH191573 | Debakey Needle Holder With Tungsten Carbide Jaws And Box Joint 230mm Straight | 85 |
| PH191574 | Debakey Needle Holder With Tungsten Carbide Jaws And Box Joint 250mm Straight | 85 |
| PH191575 | Debakey Needle Holder With Tungsten Carbide Jaws And Box Joint 300mm Straight | 85 |
| PH191585 | Derf Needle Holder Tungsten Carbide Jaws And Box Joint 120mm Straight | 85 |
| PH191590 | Kilner Needle Holder Tungsten Carbide Jaws And Box Joint With Curved Shanks 135mm | 86 |
| PH191595 | Foster Gillies Needle Holder With Tungsten Carbide Jaws And Box Joint 135mm Angle | 86 |
| PH191598 | Fry Needle Holder With Tungsten Carbide Jaws And Box Joint 130mm Straight | 86 |
| PH191600 | Gillies Needle Holder Tungsten Carbide Jaws And Box Joint With Right Hand Scissors | 86 |
| PH191605 | Halsey Needle Holder Tungsten Carbide Jaws And Box Joint 130mm Straight | 86 |

Product Codes Index

| Code | Description | Page |
|----------|---|------|
| PH191610 | Heaney Needle Holder Tungsten Carbide Jaws And Box Joint 200mm Curved | 86 |
| PH191611 | Heaney Needle Holder Tungsten Carbide Jaws And Box Joint 250mm Curved | 86 |
| PH191613 | Lawrence Needle Holder Tungsten Carbide Jaws With Box Joint 150mm Straight | 87 |
| PH191615 | Mason Needle Holder Tungsten Carbide Jaws And Box Joint 265mm (Masson) Straight | 87 |
| PH191620 | Mayo Hegar Needle Holder With Tungsten Carbide Jaws And Box Joint 130mm Straight | 87 |
| PH191621 | Mayo Hegar Needle Holder With Tungsten Carbide Jaws And Box Joint 140mm Straight | 87 |
| PH191622 | Mayo Hegar Needle Holder With Tungsten Carbide Jaws And Box Joint 150mm Straight | 87 |
| PH191623 | Mayo Hegar Needle Holder With Tungsten Carbide Jaws And Box Joint 180mm Straight | 87 |
| PH191624 | Mayo Hegar Needle Holder With Tungsten Carbide Jaws And Box Joint 200mm Straight | 87 |
| PH191625 | Mayo Hegar Needle Holder With Tungsten Carbide Jaws And Box Joint 230mm Straight | 87 |
| PH191626 | Mayo Hegar Needle Holder With Tungsten Carbide Jaws And Box Joint 250mm Straight | 87 |
| PH191627 | Mayo Hegar Needle Holder With Tungsten Carbide Jaws And Box Joint 300mm Straight | 87 |
| PH191640 | Naunton Morgan Needle Holder With Tungsten Carbide Jaws And Box Joint 250mm Straight | 87 |
| PH191645 | Nievert Needle Holder With Tungsten Carbide Jaws And Offset Ring And Box Joint | 87 |
| PH191650 | Olsen Hegar Needle Holder With Tungsten Carbide Jaws And Scissors With Screw Joint 140mm Straight | 87 |
| PH191651 | Olsen Hegar Needle Holder With Tungsten Carbide Jaws And Scissors With Screw Joint 160mm Straight | 87 |
| PH191652 | Olsen Hegar Needle Holder With Tungsten Carbide Jaws And Scissors With Screw Joint 180mm Straight | 87 |
| PH191660 | Ryder Needle Holder With French Eye Tungsten Carbide Jaws And Box Joint 130mm Straight | 88 |
| PH191661 | Ryder Needle Holder With French Eye Tungsten Carbide Jaws And Box Joint 150mm Straight | 88 |
| PH191662 | Ryder Needle Holder With French Eye Tungsten Carbide Jaws And Box Joint 180mm Straight | 88 |
| PH191663 | Ryder Needle Holder With French Eye Tungsten Carbide Jaws And Box Joint 200mm Straight | 88 |
| PH191664 | Ryder Needle Holder With French Eye Tungsten Carbide Jaws And Box Joint 230mm Straight | 88 |
| PH191675 | Ryder Needle Holder With Micro French Eye Tungsten Carbide Jaws And Box Joint 130mm Straight | 88 |
| PH191676 | Ryder Needle Holder With Micro French Eye Tungsten Carbide Jaws And Box Joint 150mm Straight | 88 |
| PH191677 | Ryder Needle Holder With Micro French Eye Tungsten Carbide Jaws And Box Joint 180mm Straight | 88 |
| PH191678 | Ryder Needle Holder With Micro French Eye Tungsten Carbide Jaws And Box Joint 200mm Straight | 88 |
| PH191679 | Ryder Needle Holder With Micro French Eye Tungsten Carbide Jaws And Box Joint 230mm Straight | 88 |

Product Codes Index

| Code | Description | Page |
|----------|---|------|
| PH191705 | Sternum Needle Holder / Wire Twister With Heavy Tungsten Carbide Jaws And Box Joint | 88 |
| PH191706 | Sternum Needle Holder / Wire Twister With Heavy Tungsten Carbide Jaws And Box Joint | 88 |
| PH191707 | Sternum Needleholder Tc Jaws 200mm | 88 |
| PH191715 | Stratte Needle Holder With Tungsten Carbide Jaws And Curved Shanks 230mm | 88 |
| PH191750 | Webster Needle Holder With Tungsten Carbide Jaws And Box Joint 130mm Straight | 88 |
| PH201885 | Aneurysm Needle With Flat Handle And Curved Laterally 140mm | 92 |
| PH201886 | Aneurysm Needle With Oval Handle And Curved Laterally 140mm | 92 |
| PH201910 | Negus Knot Tier Or Pusher 200mm | 92 |
| PH201915 | Symes Aneurysm Needle Flat Handle And Curved Laterally 160mm | 92 |
| PH212000 | Probe With Eye 130mm Straight | 92 |
| PH212001 | Probe With Eye 150mm Straight | 92 |
| PH212002 | Probe With Eye 180mm Straight | 92 |
| PH212003 | Probe With Eye 200mm Straight | 92 |
| PH212004 | Probe With Eye 250mm Straight | 92 |
| PH212005 | Probe With Eye 300mm Straight | 92 |
| PH212010 | Probe Stainless Steel With Eye 130mm Straight | 92 |
| PH212011 | Probe Stainless Steel With Eye 150mm Straight | 92 |
| PH212012 | Probe Stainless Steel With Eye 180mm Straight | 92 |
| PH212013 | Probe Stainless Steel With Eye 200mm Straight | 92 |
| PH212014 | Probe Stainless Steel With Eye 250mm Straight | 92 |
| PH212015 | Probe Stainless Steel With Eye 300mm Straight | 92 |
| PH212020 | Double Ended Probe Stainless Steel 130mm Straight | 93 |
| PH212021 | Double Ended Probe Stainless Steel 150mm Straight | 93 |
| PH212022 | Double Ended Probe Stainless Steel 180mm Straight | 93 |
| PH212023 | Double Ended Probe Stainless Steel 200mm Straight | 93 |
| PH212024 | Double Ended Probe Stainless Steel 250mm Straight | 93 |
| PH212025 | St Bartholomews Wax Hook 180mm Angled | 93 |
| PH212029 | Jobson Horne Probe With Ring 140mm Straight | 93 |
| PH212030 | Jobson Horne Probe With Ring 180mm Straight | 93 |
| PH212031 | Leibrichs Probe Size 000 - 0000 | 93 |
| PH212032 | Leibrichs Probe Size 00 - 0 | 93 |
| PH212034 | Leibrichs Probe Size 2 - 3 | 93 |

Product Codes Index

| Code | Description | Page |
|----------|---|------|
| PH212035 | Leibrichs Probe Size 1 - 3 | 93 |
| PH212036 | Leibrichs Probe Size 4 - 5 | 93 |
| PH212037 | Leibrichs Probe Size 5 - 6 | 93 |
| PH212038 | Leibrichs Probe Size 3 - 4 | 93 |
| PH222050 | Brodie Fistula Director Grooved 130mm Straight | 94 |
| PH222051 | Brodie Fistula Director Grooved 140mm Straight | 94 |
| PH222052 | Brodie Fistula Director Grooved 160mm Straight | 94 |
| PH222053 | Brodie Fistula Director Grooved With Probe End (Ref Approx 3.5mm Probe Diameter) 130mm Straight | 94 |
| PH222054 | Brodie Fistula Director Grooved With Probe End (Ref Approx 3.5mm Probe Diameter) 140mm Straight | 94 |
| PH222055 | Brodie Fistula Director Grooved With Probe End (Ref Approx 3.5mm Probe Diameter) 160mm Straight | 94 |
| PH222059 | Lockhart Mummery Set Of 4 Fistula Directors | 94 |
| PH222060 | Lockhart Mummery Straight Fistula Director (No.1, Straight) 155mm Straight | 94 |
| PH222061 | Lockhart Mummery Half Curved Fistula Director (No.2, 45 Degree Curve) 150mm | 94 |
| PH222062 | Lockhart Mummery Curved Fistula Director (No.3, 90 Degree Curve) 140mm Curved | 94 |
| PH222063 | Lockhart Mummery Full Curved Fistula Director (No.4, Retrograde) 160mm Curved | 94 |
| PH232081 | Bowman Lachrymal Dilator Probe Double Ended Sizes 0.6mm And 0.7mm Diameter 100mm | 93 |
| PH232082 | Bowman Lachrymal Dilator Probe Double Ended Sizes 0.7mm And 0.8mm Diameter 100mm | 93 |
| PH232083 | Bowman Lachrymal Dilator Probe Double Ended Sizes 0.8mm And 0.9mm Diameter 100mm | 93 |
| PH232084 | Bowman Lachrymal Dilator Probe Double Ended Sizes 0.9mm And 1.1mm Diameter 100mm | 93 |
| PH232085 | Bowman Lachrymal Dilator Probe Double Ended Sizes 1.1mm And 1.3mm Diameter 100mm | 93 |
| PH232086 | Bowman Lachrymal Dilator Probe Double Ended Sizes 1.3mm And 1.4mm Diameter 100mm | 93 |
| PH232087 | Bowman Lachrymal Dilator Probe Double Ended Sizes 1.4mm And 1.5mm Diameter 100mm | 93 |
| PH232088 | Bowman Lachrymal Dilator Probe Double Ended Sizes 1.5mm And 1.6mm Diameter 100mm | 93 |
| PH232089 | Bowman Lachrymal Dilator Probe Double Ended Sizes 1.6mm And 1.8mm Diameter 100mm | 93 |
| PH232090 | Bowman Lachrymal Dilator Probe Double Ended Sizes 1.8mm And 1.9mm Diameter 100mm | 93 |
| PH232091 | Bowman Lachrymal Dilator Probe Double Ended Sizes 0.9mm And 1.0mm Diameter 100mm | 93 |
| PH232092 | Bowman Lachrymal Dilator Probe Double Ended Sizes 1.1mm And 1.2mm Diameter 100mm | 93 |

Product Codes Index

| Code | Description | Page |
|----------|--|------|
| PH232093 | Bowman Lachrymal Dilator Probe Double Ended Sizes 1.7mm And 1.8mm Diameter 100mm | 93 |
| PH242100 | Ring Cutter Complete With Saw 160mm | 182 |
| PH242101 | Ring Cutter Forceps Only 160mm | 182 |
| PH242103 | Ring Cutter Saw Blade. Standard Make (30mm Diameter Reference) | 182 |
| PH252007 | Charnley Pin Retractor | 98 |
| PH252008 | Small Teale Gorgets | 98 |
| PH252009 | Large Teale Gorget | 98 |
| PH252018 | Awty Channel Retractors 12mm Blade Width | 98 |
| PH252019 | Awty Channel Retractors 17mm Blade Width | 98 |
| PH252022 | Baby Hohman Bone Levers 6mm | 98 |
| PH252025 | Ronald Edwards Rib Spreaders But With A Rack Measuring 450mm And Blades Measuring 120mm Deep And 75mm Wide | 99 |
| PH252100 | Hohmann Muller Spike 8mm | 99 |
| PH252101 | Hohmann Muller Spike 18mm | 99 |
| PH252102 | Hohman Bone Lever 43mm Wide 240mm | 99 |
| PH252103 | Hohman Bone Lever 70mm Wide 240mm | 99 |
| PH252125 | Retractor Sharp With Single Hook 160mm Straight | 99 |
| PH252126 | Retractor Blunt With Single Hook 160mm Straight | 99 |
| PH252127 | Retractor Sharp With Double Hook 160mm Straight | 99 |
| PH252128 | Retractor Blunt With Double Hook 160mm Straight | 99 |
| PH252129 | Retractor Sharp With Hook And 1 Prong 160mm Straight | 99 |
| PH252130 | Retractor Sharp With Hook And 2 Prongs 160mm Straight | 100 |
| PH252131 | Retractor Sharp With Hook And 3 Prongs 160mm Straight | 100 |
| PH252132 | Retractor Sharp With Hook And 4 Prongs 160mm Straight | 100 |
| PH252133 | Retractor Blunt With Hook And 1 Prong 160mm Straight | 99 |
| PH252134 | Retractor Blunt With Hook And 2 Prongs 160mm Straight | 100 |
| PH252135 | Retractor Blunt With Hook And 3 Prongs 160mm Straight | 100 |
| PH252136 | Retractor Blunt With Hook And 4 Prongs 160mm Straight | 100 |
| PH252138 | Austin Kay Right Angled Retractor With Lip 150mm Angled | 100 |
| PH252141 | Allison Retractor Infant Size 280mm | 100 |
| PH252142 | Allison Retractor Child Size 280mm | 100 |
| PH252143 | Allison Retractor Adult Size 300mm | 100 |
| PH252144 | Aufricht Nasal Retractor 50mm Blade 160mm Angled | 101 |

Product Codes Index

| Code | Description | Page |
|----------|---|------|
| PH252145 | Cairn (Rake) Scalp Retractor With 3 Prongs 16mm Sharp 210mm Straight | 101 |
| PH252146 | Cairn (Rake) Scalp Retractor With 4 Prongs 19mm Sharp 210mm Straight | 101 |
| PH252160 | Czerny Retractor Double Ended 155mm | 101 |
| PH252170 | Deaver Retractor With Hollow Handle 18mm Blade x 180mm | 101 |
| PH252171 | Deaver Retractor Plain Handle With 25mm Blade 230mm | 101 |
| PH252172 | Deaver Retractor With Hollow Handle 25mm Blade x 250mm | 101 |
| PH252173 | Deaver Retractor Plain Handle With 38mm Blade 300mm | 101 |
| PH252174 | Deaver Retractor Plain Handle With 50mm Blade 300mm | 101 |
| PH252175 | Deaver Retractor Plain Handle With 75mm Blade 300mm | 101 |
| PH252176 | Deaver Retractor (Without Handle) Plain Strip 25mm Blade x 230mm | 101 |
| PH252177 | Deaver Retractor Without Handle Plain Strip With 25mm Blade 250mm | 101 |
| PH252178 | Deaver Retractor Without Handle Plain Strip With 38mm Blade 300mm | 101 |
| PH252179 | Deaver Retractor Without Handle Plain Strip With 50mm Blade 300mm | 101 |
| PH252180 | Deaver Retractor (Without Handle) Plain Strip 75mm Blade x 300mm | 101 |
| PH252185 | Doyen Abdominal Retractor Small With 90mm Width Blade 215mm | 102 |
| PH252186 | Doyen Abdominal Retractor Medium 100mm Width Blade 215mm | 102 |
| PH252187 | Doyen Abdominal Retractor Large With 115mm Width Blade 215mm | 102 |
| PH252188 | Doyen Retractor Small Ring Handles 95mm Width Blade 215mm | 102 |
| PH252189 | Doyen Retractor Medium Ring Handles 105mm Width Blade 215mm | 102 |
| PH252190 | Doyen Retractor Large Ring Handles 115mm Width Blade 215mm | 102 |
| PH252193 | Desmarres Retractor Small 9mm x 9mm (Eye Lid) 160mm | 102 |
| PH252194 | Desmarres Retractor Medium 9mm x 12mm (Eye Lid) 160mm | 102 |
| PH252195 | Desmarres Retractor Large 11mm x 14mm (Eye Lid) 160mm | 102 |
| PH252240 | Bowdler Henry Rake Retractor For Lower Wisdom Teeth 190mm | 102 |
| PH252251 | O'connell Double Ended Retractor For Nerve Root (Oconnel O Connel) 215mm | 102 |
| PH252256 | Ward Double Ended Cheek Retractor Used In The Removal Of Wisdom Teeth 150mm | 102 |
| PH252259 | Kelly Retractor With Small Size 38mm End Of Blade 190mm | 103 |
| PH252260 | Kelly Retractor With Large Size 57mm End Of Blade 190mm | 103 |
| PH252261 | Kilner (Cat Paw) Skin Retractor Double Ended Std Size 150mm | 103 |
| PH252263 | Kilner (Cheek) Small Retractor Double Ended Starting 25mm Tapering Down To 20mm | 103 |
| PH252264 | Kilner (Cheek) Large Retractor Double Ended Starting 35mm Tapering Down To 30mm | 103 |
| PH252265 | Kilner (Alae) Retractor Double Pronged 8mm Wide 85mm | 103 |
| PH252266 | Kilner (Alae) Retractor Double Pronged 10mm Wide 85mm | 103 |

Product Codes Index

| Code | Description | Page |
|-----------|---|------|
| PH252267 | Kilner (Alae) Retractor Double Pronged 13mm Wide 85mm | 103 |
| PH252269 | Kocher Retractor 40 x 10mm 215mm | 103 |
| PH252270 | Kocher Retractor 40 x 15mm 215mm | 103 |
| PH252271 | Kocher Retractor 40 x 20mm 215mm | 103 |
| PH252278 | Laster Cheek Retractor For Third Upper Molar 100mm | 103 |
| PH252285 | Langenbeck Small Retractor Single Blade Size 23mm x 7mm 215mm | 104 |
| PH252286 | Langenbeck Medium Retractor Single Blade Size 44mm x 13mm 215mm | 104 |
| PH252287 | Langenbeck Large Retractor Single Blade Size 65mm x 25mm 215mm | 104 |
| PH252288 | Langenbeck Extra Large Retractor Single Blade Size 90mm x 20mm 215mm | 104 |
| PH252300 | Morris Retractor Single Ended 38mm Wide x 50mm Deep 230mm | 104 |
| PH252301 | Morris Retractor Single Ended 50mm Wide x 52mm Deep 230mm | 104 |
| PH252302 | Morris Retractor Single Ended 65mm Wide x 50mm Deep 230mm | 104 |
| PH252303 | Morris Retractor Single Ended 50mm Wide x 50mm Deep 230mm | 104 |
| PH252304 | Morris Retractor Single Ended 50mm Wide x 60mm Deep 230mm | 104 |
| PH252305 | Morris Retractor Single Ended 50mm Wide x 65mm Deep 230mm | 104 |
| PH252306 | Morris Retractor Single Ended 50mm Wide x 70mm Deep 230mm | 104 |
| PH252307 | Morris Retractor Single Ended 65mm Wide x 70mm Deep 230mm | 104 |
| PH252308 | Morris Retractor Double Ended 32mm Wide x 38mm Wide And 57mm Deep 230mm | 104 |
| PH252309 | Morris Retractor Double Ended 44mm Wide x 50mm Wide And 57mm Deep 230mm | 104 |
| PH252310 | Morris Retractor Double Ended 57mm Wide x 64mm Wide And 57mm Deep 230mm | 104 |
| PH252310A | Morris Retractor Double Ended 30mm Wide x 40mm Wide And 70mm Deep 230mm | 104 |
| PH252310B | Morris Retractor Double Ended 45mm Wide x 50mm Wide And 70mm Deep 230mm | 104 |
| PH252310C | Morris Retractor Double Ended 60mm Wide x 65mm Wide And 70mm Deep 230mm | 104 |
| PH252333 | Park Anal And Sphincterotomy Combination Retractor Adult Size With Two Sets Of I | 105 |
| PH252338 | Ragnel Davis Double Ended Retractor 8mm x 4mm + 15mm x 5.5mm 150mm | 105 |
| PH252348 | St Marks / Lloyd Davis Retractor With Blade Length 130mm x 60mm Width With Lip | 105 |
| PH252349 | St Marks / Lloyd Davis Retractor With Blade Length 170mm Into Taper 64mm | 105 |
| PH252350 | St Marks / Lloyd Davis Retractor With Blade Length 130mm x 60mm Width Without Lip | 105 |
| PH252351 | St Marks / Lloyd Davis Retractor With Blade Length 170mm Into Taper 64mm | 105 |
| PH252352 | Volkman Retractor With 1 Sharp Prong 215mm | 106 |
| PH252353 | Volkman Retractor With 2 Sharp Prongs 215mm | 106 |
| PH252354 | Volkman Retractor With 3 Sharp Prongs 215mm | 106 |
| PH252355 | Volkman Retractor With 4 Sharp Prongs 215mm | 106 |

Product Codes Index

| Code | Description | Page |
|----------|--|------|
| PH252356 | Volkman Retractor With 5 Sharp Prongs 215mm | 106 |
| PH252357 | Volkman Retractor With 6 Sharp Prongs 215mm | 106 |
| PH252360 | Volkman Retractor With 1 Blunt Prong 215mm | 106 |
| PH252361 | Volkman Retractor With 2 Blunt Prongs 215mm | 106 |
| PH252362 | Volkman Retractor With 3 Blunt Prongs 215mm | 106 |
| PH252363 | Volkman Retractor With 4 Blunt Prongs 215mm | 106 |
| PH252364 | Volkman Retractor With 5 Blunt Prongs 215mm | 106 |
| PH252365 | Volkman Retractor With 6 Blunt Prongs 215mm | 106 |
| PH252368 | Baby Senn Miller Retractor Sharp Double Ended | 107 |
| PH252370 | Alms Skin Retractor With Sharp Prongs 70mm | 107 |
| PH252371 | Alms Skin Retractor With Sharp Prongs 100mm | 107 |
| PH252372 | Alms Skin Retractor With Blunt Prongs 70mm | 107 |
| PH252373 | Alms Skin Retractor With Blunt Prongs 100mm | 107 |
| PH252374 | Charnley Initial Incision Retractor Complete With 2 Blades | 108 |
| PH252375 | Balfour Abdominal Retractor With Centre Blade And Fenestrated End Blades 63mm x | 108 |
| PH252376 | Balfour Baby Abdominal Retractor With Centre Blade And Fenestrated End Blades 27mm Deep | 108 |
| PH252383 | Goligher Retractor Blades (2 Centre Blades And 2 Pairs Lateral Blades) | 109 |
| PH252404 | Joll Thyroid Retractor With 2 Self Retaining 2 Pincher Blades And Fluted Handle | 109 |
| PH252407 | Cone Laminectomy Retractor 3 Into 4 Teeth With Hinged Arms And Ratchet Sharp Teeth | 110 |
| PH252408 | Cone Laminectomy Retractor 3 Into 4 Teeth With Hinged Arms And Ratchet Blunt Teeth | 110 |
| PH252409 | Cone Laminectomy Retractor 3 Into 4 Teeth With Hinged Arms And Ratchet Sharp Teeth | 110 |
| PH252410 | Cone Laminectomy Retractor 3 Into 4 Teeth With Hinged Arms And Ratchet Blunt Teeth | 110 |
| PH252411 | Cone Laminectomy Retractor 3 Into 4 Teeth With Hinged Arms And Ratchet Sharp Teeth | 110 |
| PH252412 | Cone Laminectomy Retractor 3 Into 4 Teeth With Hinged Arms And Ratchet Blunt Teeth | 110 |
| PH252413 | Cone Laminectomy Retractor 3 Into 4 Teeth With Hinged Arms And Ratchet Sharp Teeth | 110 |
| PH252414 | Cone Laminectomy Retractor 3 Into 4 Teeth With Hinged Arms And Ratchet Blunt Teeth | 110 |
| PH252420 | Jackson Burrows Self Retaining Blunt Hook Jaws For Superficial Part Of The Knee 175mm | 110 |
| PH252424 | Norfolk And Norwich Self Retaining Retractor With 4 Into 5 Teeth And Cam Rack. B | 110 |
| PH252425 | Norfolk And Norwich Self Retaining Retractor Long Teeth With 4 Into 5 Teeth And Cam Rack. Blades 78mm Deep 220mm | 110 |
| PH252426 | Kocher Retractor Blade 45mm Wide x 85mm Deep 250mm | 113 |
| PH252427 | Kocher Retractor Blade 54mm Wide x 85mm Deep 250mm | 113 |
| PH252428 | Kocher Retractor Blade 62mm Wide x 85mm Deep 250mm | 113 |

Product Codes Index

| Code | Description | Page |
|----------|---|------|
| PH252437 | Gelpi Vulva Self Retaining Retractor Large Size With Ratchet 200mm | 110 |
| PH252438 | Gelpi Vulva Self Retaining Retractor Small Size With Ratchet 170mm | 110 |
| PH252440 | Irwin Retractor With 4 Into 4 Teeth And Cam Rack 180mm | 110 |
| PH252443 | Mollison Retractor Adult Size Self Retaining Curved With 4 Into 4 Teeth (Mollins) | 111 |
| PH252444 | Mollison Retractor Child Size Self Retaining Curved With 2 Into 2 Teeth 135mm | 111 |
| PH252446 | Travers Retractor With 4 Into 5 Teeth Blunt With Ratchet 210mm | 111 |
| PH252450 | West Retractor With 2 Into 3 Teeth Blunt 115mm | 111 |
| PH252451 | West Retractor With 2 Into 3 Teeth Sharp 115mm | 111 |
| PH252452 | West Retractor With 3 Into 4 Teeth Blunt 140mm | 111 |
| PH252453 | West Retractor With 3 Into 4 Teeth Sharp 140mm | 111 |
| PH252455 | Weitlander Retractor With 2 Into 3 Teeth Sharp (Weislander) 115mm | 111 |
| PH252456 | Weitlander Retractor With 3 Into 4 Teeth Sharp 140mm | 111 |
| PH252457 | Weitlander Retractor With 3 Into 4 Teeth Sharp 160mm | 111 |
| PH252458 | Weitlander Retractor With 2 Into 3 Teeth Blunt 115mm | 111 |
| PH252459 | Weitlander Retractor With 3 Into 4 Teeth Blunt 140mm | 111 |
| PH252460 | Weitlander Retractor With 3 Into 4 Teeth Blunt 160mm | 111 |
| PH252465 | Doyen Retractor Small 48mm x 60mm Rectangular Blade 250mm | 112 |
| PH252466 | Doyen Retractor Medium 48mm x 90mm Rectangular Blade 250mm | 112 |
| PH252467 | Doyen Retractor Large 48mm x 125mm Rectangular Blade 250mm | 112 |
| PH252470 | Landon Vaginal Retractor With 90mm x 25mm Blade And Hollow Handle (Langdon) 200mm | 112 |
| PH252499 | Tudor Edwards Scapula Retractor With 45mm Wide Blade x 50mm Deep 160mm | 112 |
| PH252508 | Large Purse String Applicator | 71 |
| PH252509 | Small Purse String Applicator | 71 |
| PH262500 | Gillies Small Skin Hook With Fine Point 155mm | 95 |
| PH262501 | Gillies Large Skin Hook With Fine Point 155mm | 95 |
| PH262502 | Kilner Detachable Skin Hook And Handle 155mm | 95 |
| PH262505 | Mcindoe Double Prong Skin Hooks With 9mm Wide End 150mm | 95 |
| PH262506 | Mcindoe Double Prong Skin Hooks With 10mm Wide End 150mm | 95 |
| PH262507 | Mcindoe Double Prong Skin Hooks With 11mm Wide End 150mm | 95 |
| PH262510 | Frazier Small Skin Hook 180mm | 95 |
| PH262549 | Gillies Medium Skin Hook With Fine Point 155mm | 95 |
| PH282535 | Lack Tongue Depressor Infant Size 12mm Wide | 116 |

Product Codes Index

| Code | Description | Page |
|-------------|--|------|
| PH282536 | Lack Tongue Depressor Child Size 19mm Wide | 116 |
| PH282537 | Lack Tongue Depressor Adult Size 22mm Wide | 116 |
| PH282538 | Mayo Tongue Depressor Double Ended Tapering Child Size 140mm | 116 |
| PH282539 | Mayo Tongue Depressor Double Ended Tapering Adult Size 200mm | 116 |
| PH292546 | Davis Boyle Mouth Gag Frame Only Child 135mm x 70mm 135mm | 116 |
| PH292547 | Davis Boyle Mouth Gag Frame Only Adult 150mm x 90mm 150mm | 116 |
| PH292552 | Davis Boyle Mouth Gag Tongue Plate Only Fig. 1 65mm | 116 |
| PH292553 | Davis Boyle Mouth Gag Tongue Plate Only Fig. 2 75mm | 116 |
| PH292554 | Davis Boyle Mouth Gag Tongue Plate Only Fig. 3 90mm | 116 |
| PH292555 | Davis Boyle Mouth Gag Tongue Plate Only Fig. 4 100mm | 116 |
| PH292556 | Davis Boyle Mouth Gag Tongue Plate Only Fig. 5 115mm | 116 |
| PH292562 | Daughty Mouth Gag Tongue Plate Only 64mm | 116 |
| PH292563 | Daughty Mouth Gag Tongue Plate Only 76mm | 116 |
| PH292564 | Daughty Mouth Gag Tongue Plate Only 89mm | 116 |
| PH292565 | Daughty Mouth Gag Tongue Plate Only 102mm | 116 |
| PH292566 | Daughty Mouth Gag Tongue Plate Only 114mm | 116 |
| PH292567AD | Draffin Bipod Fitted With 4 Ring For Varying Heights Adult 480mm | 117 |
| PH292567CH | Draffin Bipod Fitted With 4 Ring For Varying Heights Child 375mm | 117 |
| PH292568 | Maguaran Plate For Holding Draffin Bipod 407mm | 117 |
| PH302595 | Bowlby Tracheal Dilating Forceps Adult Size 140mm Curved | 118 |
| PH302596 | Bowlby Tracheal Dilating Forceps Child Size 115mm Curved | 118 |
| PH302601 | Magil Intra Tracheal Dilating Forceps Adult Size 240mm | 119 |
| PH302602 | Magil Intra Tracheal Dilating Forceps Child Size 200mm | 119 |
| PH322800-10 | Laryngeal Handle With 10mm Mirror Diameter s | 117 |
| PH322800-12 | Laryngeal Handle With 12mm Mirror Diameter s | 117 |
| PH322800-14 | Laryngeal Handle With 14mm Mirror Diameter s | 117 |
| PH322800-16 | Laryngeal Handle With 16mm Mirror Diameter s | 117 |
| PH322800-18 | Laryngeal Handle With 18mm Mirror Diameter s | 117 |
| PH322800-20 | Laryngeal Handle With 20mm Mirror Diameter s | 117 |
| PH322800-22 | Laryngeal Handle With 22mm Mirror Diameter s | 117 |
| PH322800-24 | Laryngeal Handle With 24mm Mirror Diameter s | 117 |
| PH322800-26 | Laryngeal Handle With 26mm Mirror Diameter s | 117 |
| PH322800-28 | Laryngeal Handle With 28mm Mirror Diameter s | 117 |

Product Codes Index

| Code | Description | Page |
|---------------|--|------|
| PH322800-30 | Laryngeal Handle With 30mm Mirror Diameter s | 117 |
| PH322800-8 | Laryngeal Handle With 8mm Mirror Diameter s | 117 |
| PH322801 | Laryngeal Handle Only 120mm | 117 |
| PH322810-10 | Dental Mirror With 10mm Mirror Diameter | 117 |
| PH322810-12 | Dental Mirror With 12mm Mirror Diameter | 117 |
| PH322810-14 | Dental Mirror With 14mm Mirror Diameter | 117 |
| PH322810-16 | Dental Mirror With 16mm Mirror Diameter | 117 |
| PH322810-18 | Dental Mirror With 18mm Mirror Diameter | 117 |
| PH322810-20 | Dental Mirror With 20mm Mirror Diameter | 117 |
| PH322810-20MO | Dental Mirror Head Only 20mm Mirror Diameter | 117 |
| PH322810-22 | Dental Mirror With 22mm Mirror Diameter | 117 |
| PH322810-24 | Dental Mirror With 24mm Mirror Diameter | 117 |
| PH322810-26 | Dental Mirror With 26mm Mirror Diameter | 117 |
| PH322810-28 | Dental Mirror With 28mm Mirror Diameter | 117 |
| PH322810-30 | Dental Mirror With 30mm Mirror Diameter | 117 |
| PH322810-8 | Dental Mirror With 8mm Mirror Diameter | 117 |
| PH342850 | Beale Aural Forceps Crocodile Action With Fine Serrated Jaw And Black Finish 70mm Shoulder Length | 62 |
| PH342851 | Cawthorne Aural Forceps Crocodile Action With Fine Smooth Jaw And Black Finish 70mm Shoulder Length | 62 |
| PH342852 | Hartmann Aural Forceps Crocodile Action With A Extra Fine Serrated Jaw And A Bla | 62 |
| PH342860 | Ormerod Aural Forceps Crocodile Action With Extra Fine Elongated Cup Jaw And Black Finish 70mm Shoulder Length (Ormerod) | 62 |
| PH342861 | Ormerod Aural Forceps Crocodile Action With A Extra Fine Cup Jaw And A Black Fin | 62 |
| PH342862 | Shea Aural Forceps Crocodile Action With Oval Cup Jaw Straight And Black Finish 70mm Shoulder Length | 62 |
| PH342863 | Shea Aural Forceps Crocodile Action With Oval Cup Jaw Curved Upwards And Black Finish 70mm Shoulder Length | 63 |
| PH342864 | Shea Aural Forceps Crocodile Action With Oval Cup Jaw Curved To The Right And Black Finish 70mm Shoulder Length | 63 |
| PH342879 | Tilley Aural Dressing Forceps With Serrated Jaws 150mm | 63 |
| PH342880 | Wilde Aural Spring Pattern Dressing Forceps With Angular Serrated Jaws 130mm | 63 |
| PH342883 | Hartmann Aural Forceps Crocodile Action With A Serrated Jaws And A 85mm Shoulder | 63 |
| PH342884 | Hartmann Aural Forceps Crocodile Action With 1 x 2 Teeth And 85mm Shoulder Length | 64 |
| PH342887 | Henckel Aural Forceps Crocodile Action With Fine Fenestrated Blades And 75mm Shoulder Length | 64 |
| PH362920 | House Aural Curette Double Ended Straight 1.0mm And 1.2mm Spoon 180mm Straight | 122 |

Product Codes Index

| Code | Description | Page |
|----------|--|------|
| PH362921 | House Aural Curette Double Ended Straight 1.5mm And 1.8mm Spoon 180mm Straight | 122 |
| PH362922 | House Aural Curette Double Ended Angled 1.0mm And 1.2mm Spoon 180mm Angled | 122 |
| PH362923 | House Aural Curette Double Ended Angled 1.5mm And 1.8mm Spoon 180mm Angled | 122 |
| PH362930 | Macewen Double Ended Curette Round Spoon And Seeker (Mcewan) 185mm | 122 |
| PH362931 | Macewen Double Ended Curette Oval Spoon And Seeker (Mcewan) 185mm | 122 |
| PH362940 | Cawthorne Hook Single End No.1 Slight Curve 165mm Curved | 123 |
| PH362941 | Cawthorne Hook Single End No.2 Medium Curve 165mm Curved | 123 |
| PH362942 | Cawthorne Hook Single End No.3 Full Curve 165mm Curved | 123 |
| PH362943 | Cawthorne Hook Single End No.4 Double Curve 165mm Curved | 123 |
| PH362944 | Cawthorne Hook Single End No.5 90 Deg Angle 0.5mm Point 165mm Angled | 123 |
| PH362945 | Cawthorne Hook Single End No.6 90 Deg Angle 1.0mm Point 165mm Angled | 123 |
| PH362946 | Cawthorne Hook Single Ended No.7 Retrograde 165mm Curved | 124 |
| PH362947 | Cawthorne Hook Single End No.8 45 Deg Angle 0.5mm Point 165mm Angled | 124 |
| PH362948 | Cawthorne Hook No.9 Straight 0.5mm Point 165mm | 124 |
| PH362960 | Volkmann Aural Small Curette Double Ended Scoops 215mm Straight | 124 |
| PH362961 | Volkmann Aural Medium Curette Double Ended Scoops 215mm Straight | 124 |
| PH362962 | Volkmann Aural Large Curette Double Ended Scoops 215mm Straight | 124 |
| PH362963 | Volkmann Aural Curette Single End Size 00 Hollow Handle 180mm Straight | 124 |
| PH362964 | Volkmann Aural Curette Single End Size 0 Hollow Handle 180mm Straight | 124 |
| PH362965 | Volkmann Aural Curette Single End Size 1 Hollow Handle 180mm Straight | 124 |
| PH362966 | Volkmann Aural Curette Single End Size 2 Hollow Handle 180mm Straight | 124 |
| PH362967 | Volkmann Aural Curette Single End Size 3 Hollow Handle 180mm Straight | 124 |
| PH362968 | Volkmann Aural Curette Single End Size 4 Hollow Handle 180mm Straight | 124 |
| PH362969 | Volkmann Aural Curette Single End Size 5 Hollow Handle 180mm Straight | 124 |
| PH362970 | Volkmann Aural Curette Single End Size 6 Hollow Handle 180mm Straight | 124 |
| PH362981 | Musgrove Aural Double Ended Wool Carrier And Cerumen Pick 170mm | 125 |
| PH362982 | St Barts Blunt Hook Wool Carrier 170mm Angled | 125 |
| PH362983 | Formby Cerumen Scoop And Hook Double Ended 180mm | 125 |
| PH362986 | Gruber Aural Speculum 4mm Diameter Round 38mm | 126 |
| PH362987 | Gruber Aural Speculum 5mm Diameter Round 38mm | 126 |
| PH362988 | Gruber Aural Speculum 6mm Diameter Round 38mm | 126 |
| PH362989 | Gruber Aural Speculum 7mm Diameter Round 38mm | 126 |
| PH362991 | Gruber Aural Speculum 4mm Diameter Oval | 126 |

Product Codes Index

| Code | Description | Page |
|----------|--|------|
| PH362992 | Gruber Aural Speculum 5mm Diameter Oval | 126 |
| PH362993 | Gruber Aural Speculum 6mm Diameter Oval | 126 |
| PH362994 | Gruber Aural Speculum 7mm Diameter Oval | 126 |
| PH362995 | Shea Speculum Oblique End Set Of Five | 126 |
| PH362996 | Shea Speculum Oblique End Size No.1 = 3.7mm | 126 |
| PH362997 | Shea Speculum Oblique End Size No.2 = 4.0mm | 126 |
| PH362998 | Shea Speculum Oblique End Size No.3 = 5.0mm | 126 |
| PH362999 | Shea Speculum Oblique End Size No.4 = 6.0mm | 126 |
| PH363000 | Shea Speculum Oblique End Size No.5 = 7.0mm | 126 |
| PH363021 | Rosen Aural Speculum 4mm Diameter Round 38mm | 126 |
| PH363022 | Rosen Aural Speculum 5mm Diameter Round 38mm | 126 |
| PH363023 | Rosen Aural Speculum 6mm Diameter Round 38mm | 126 |
| PH363024 | Rosen Aural Speculum 7mm Diameter Round 38mm | 126 |
| PH363026 | Rosen Aural Speculum 4mm Diameter Oval 38mm | 127 |
| PH363027 | Rosen Aural Speculum 5mm Diameter Oval 38mm | 127 |
| PH363028 | Rosen Aural Speculum 6mm Diameter Oval 38mm | 127 |
| PH363029 | Rosen Aural Speculum 7mm Diameter Oval 38mm | 127 |
| PH363031 | Toynbee Aural Speculum 4mm Diameter Round 41mm | 127 |
| PH363032 | Toynbee Aural Speculum 5mm Diameter Round 41mm | 127 |
| PH363033 | Toynbee Aural Speculum 6mm Diameter Round 41mm | 127 |
| PH363034 | Toynbee Aural Speculum 7mm Diameter Round 41mm | 127 |
| PH373051 | Agnew Myringotome Ear Knife Angled 190mm Angled | 127 |
| PH383083 | Dental Cartridge Syringe 2.2cc Straight | 64 |
| PH413176 | Tilley Nasal Dressing Forceps With Serrated Jaws And Angled Shanks 165mm | 64 |
| PH413179 | Weil Straight Blakley Forcep | 65 |
| PH413180 | Weil Upwards Blakley Forcep | 65 |
| PH413181 | Weil Blakley Forcep 90deg Up | 65 |
| PH423211 | Irwin Moore Nasal Turbinate Forceps With Screw Joint And Angled Shanks 190mm | 65 |
| PH423215 | Tilley Henckel Nasal Turbinate Forceps 4.5mm Ethmoid Punches | 65 |
| PH423216 | Tilley Henckel Nasal Turbinate Forceps 5.5mm Ethmoid Punches | 65 |
| PH423217 | Tilley Henckel Nasal Turbinate Forceps 6.5mm Ethmoid Punches | 65 |
| PH423218 | Tilley Henckel Nasal Turbinate Forceps 8mm Ethmoid Punches | 65 |
| PH423220 | Jansen Middleton Nasal Turbinate Septum Forceps Double Action 190mm | 66 |

Product Codes Index

| Code | Description | Page |
|----------|---|------|
| PH423221 | Luc Nasal Turbinate Forceps Oval Small With Screw Joint And Angled Shanks (Lucs) | 66 |
| PH423222 | Luc Nasal Turbinate Forceps Oval Medium With Screw Joint And Angled Shanks 180mm | 66 |
| PH423223 | Luc Nasal Turbinate Forceps Oval Large With Screw Joint And Angled Shanks 180mm | 66 |
| PH433255 | Walsham Nasal Septum Forceps For Straightening Box Joint 220mm | 66 |
| PH433256 | Walsham Nasal Septum Forceps For Redressing Left With Box Joint 220mm | 66 |
| PH433257 | Walsham Nasal Septum Forceps For Redressing Right With Box Joint 220mm | 66 |
| PH443272 | St Clair Thomson Adenoid Curette Size 0 With Cage And Metal Handle 8mm 220mm | 130 |
| PH443273 | St Clair Thomson Adenoid Curette Size 1 With Cage And Metal Handle 10mm 220mm | 130 |
| PH443274 | St Clair Thomson Adenoid Curette Size 1 1/2 With Cage And Metal Handle 12mm 220mm | 130 |
| PH443275 | St Clair Thomson Adenoid Curette Size 2 With Cage And Metal Handle 14mm 220mm | 130 |
| PH443276 | St Clair Thomson Adenoid Curette Size 3 With Cage And Metal Handle 16mm 220mm | 130 |
| PH443277 | St Clair Thomson Adenoid Curette Size 4 With Cage And Metal Handle 18mm 220mm | 130 |
| PH443279 | Beckmann Adenoid Curette Size 0 With Hollow Handle 8mm 210mm | 130 |
| PH443280 | Beckmann Adenoid Curette Size 1 With Hollow Handle 10mm 210mm | 130 |
| PH443281 | Beckmann Adenoid Curette Size 1 1/2 With Hollow Handle 12mm 210mm | 130 |
| PH443282 | Beckmann Adenoid Curette Size 2 With Hollow Handle 14mm 210mm | 130 |
| PH443283 | Beckmann Adenoid Curette Size 3 With Hollow Handle 16mm 210mm | 130 |
| PH443284 | Beckmann Adenoid Curette Size 4 With Hollow Handle 18mm 210mm | 130 |
| PH453310 | Cottle Nasal Speculum No.1 Blade 20mm. Screw Joint With Side Locking Screw | 130 |
| PH453311 | Cottle Nasal Speculum No.2 Blade 32mm. Screw Joint With Side Locking Screw | 130 |
| PH453312 | Cottle Nasal Speculum No.3 Blade 41mm. Screw Joint With Side Locking Screw | 130 |
| PH453313 | Cottle Nasal Speculum No.4 Blade 60mm. Screw Joint With Side Locking Screw | 130 |
| PH453314 | Killian Nasal Speculum Blade 38mm Length. Screw Joint And Side Locking Screw | 130 |
| PH453315 | Killian Nasal Speculum Blade 64mm Length. Screw Joint And Side Locking Screw | 130 |
| PH453316 | Killian Nasal Speculum Blade 76mm Length. Screw Joint And Side Locking Screw | 130 |
| PH453317 | Killian Nasal Speculum Blade 89mm Length. Screw Joint And Side Locking Screw | 130 |
| PH453320 | St Clair Thomson Nasal Speculum Blade 38mm Length. Small | 131 |
| PH453321 | St Clair Thomson Nasal Speculum Blade 51mm Length. Medium | 131 |
| PH453322 | St Clair Thomson Nasal Speculum Blade 64mm Length. Large | 131 |
| PH453323 | St Clair Thomson Nasal Speculum Blade 89mm Length. Extra Large | 131 |
| PH453325 | Thudichum Nasal Speculum Size 1 Blade (Smallest) 70mm | 131 |
| PH453326 | Thudichum Nasal Speculum Size 2 Blade 70mm | 131 |
| PH453327 | Thudichum Nasal Speculum Size 3 Blade 70mm | 131 |

Product Codes Index

| Code | Description | Page |
|----------|--|------|
| PH453328 | Thudichum Nasal Speculum Size 4 Blade (Medium) 70mm | 131 |
| PH453329 | Thudichum Nasal Speculum Size 5 Blade 70mm | 131 |
| PH453330 | Thudichum Nasal Speculum Size 6 Blade 70mm | 131 |
| PH453331 | Thudichum Nasal Speculum Size 7 Blade (Biggest) 70mm | 131 |
| PH453335 | Tumarkin Aural Speculum Tapered Slot Set Of 8 | 127 |
| PH453336 | Tumarkin Aural Speculum With Tapered Slot And A Black Finish. Individual Size 0 4.5mm 35mm | 127 |
| PH453337 | Tumarkin Tapered Spec Slotted Size 1 = 6.5mm Black | 127 |
| PH453338 | Tumarkin Tapered Spec Slotted Size 2 = 7.0mm Black | 127 |
| PH453339 | Tumarkin Tapered Spec Slotted Size 3 = 7.5mm Black | 127 |
| PH453340 | Tumarkin Tapered Spec Slotted Size 4 = 8.0mm Black | 127 |
| PH453341 | Tumarkin Nasal Speculum With Tapered Slot And A Black Finish. | 127 |
| PH453342 | Tumarkin Nasal Speculum With Tapered Slot And A Black Finish. | 127 |
| PH453343 | Tumarkin Nasal Speculum With Tapered Slot And A Black Finish. | 127 |
| PH463377 | Eve Nasal Tonsil Snare Sliding Type 280mm Straight | 131 |
| PH463389 | Ballenger Nasal Cartilage Knife 3mm Wide Blade With Straight Body And A Revolving Blade 200mm Straight | 132 |
| PH463390 | Ballenger Nasal Cartilage Knife 5mm Wide Blade With Straight Body And A Revolving Blade 200mm Straight | 132 |
| PH463391 | Ballenger Nasal Cartilage Knife 8mm Wide Blade With Straight Body And A Revolving Blade 200mm Straight | 132 |
| PH473401 | Denis Browne Tonsil Seizing Forceps With Box Joint Tonsil Holding Forceps And Angled Shanks 215mm | 66 |
| PH483450 | Silver Nasal Chisel Straight With Side Probe Guard 180mm | 133 |
| PH483451 | Silver Nasal Chisel Curved To Right With Side Probe Guard 180mm | 133 |
| PH483452 | Silver Nasal Chisel Curved To Left With Side Probe Guard 180mm | 133 |
| PH483458 | Joseph Nasal Saw Vertical Shaft And Bayonet Crank Right 185mm | 133 |
| PH483460 | Joseph Nasal Saw Vertical Shaft And Bayonet Crank Left 185mm | 133 |
| PH483470 | Collingwood Stewart Hernia Forceps Ring Bladed Screw Joint 130mm | 38 |
| PH493482 | Freer Double Ended Separator And Elevator 200mm | 133 |
| PH493486 | Hurd Double Ended Pillar Retractor And Tonsil Dissector 200mm | 133 |
| PH493495 | Gwynne Evans Double Ended 6mm Wide Dissector 190mm | 136 |
| PH493497 | Hill Double Ended Elevator Curved To Side With Rugine End 5mm Wide (Hills) 220mm | 136 |
| PH493498 | Howarth Double Ended Elevator Curved On Flat With Rugine End 220mm Curved | 136 |
| PH493505 | Macdonald Double Ended Dissector (Mcdonald) 190mm | 136 |
| PH493508 | Syme Double Ended Dissector One Blunt End And One Semi Sharp End 180mm | 136 |

Product Codes Index

| Code | Description | Page |
|----------|---|------|
| PH493510 | Watson Cheyne Double Ended Elevator Fine With Probe 130mm | 137 |
| PH493511 | Watson Cheyne Double Ended Elevator Fine With Probe 170mm 170mm | 137 |
| PH493512 | Watson Cheyne Double Ended Elevator Fine With Probe 230mm | 137 |
| PH493515 | Warwick James Single End Elevator Straight Light Model 170mm Straight | 137 |
| PH493516 | Warwick James Single End Elevator Curved To Right Light Model 170mm Curved | 137 |
| PH493517 | Warwick James Single End Elevator Curved To Left Light Model 170mm Curved | 137 |
| PH493518 | Cryers Elevator Left 170mm | 137 |
| PH493519 | Cryers Elevator Right 170mm | 137 |
| PH493520 | Mitchell Trimmer Double Ended 160mm | 138 |
| PH493521 | Henderson Dissector (Periosteal Elevator) 210mm | 138 |
| PH503550 | Allis Tissue Forceps With 3 Into 4 Teeth And Box Joint 150mm | 38 |
| PH503551 | Allis Tissue Forceps With 4 Into 5 Teeth And Box Joint 150mm | 38 |
| PH503552 | Allis Tissue Forceps With 5 Into 6 Teeth And Box Joint 150mm | 38 |
| PH503553 | Allis Tissue Forceps With 5 Into 6 Teeth And Box Joint 180mm | 38 |
| PH503554 | Allis Tissue Forceps With 5 Into 6 Teeth And Box Joint 200mm | 38 |
| PH503555 | Allis Tissue Forceps With 5 Into 6 Teeth And Box Joint 250mm | 38 |
| PH503556 | Allis Tissue Forceps With 5 Into 6 Teeth And Box Joint 280mm | 38 |
| PH503558 | Babcock Tissue Forceps With Box Joint 160mm | 38 |
| PH503559 | Babcock Tissue Forceps With Box Joint 180mm | 38 |
| PH503560 | Babcock Tissue Forceps With Box Joint 200mm | 38 |
| PH503561 | Babcock Tissue Forceps With Box Joint 230mm | 38 |
| PH503567 | Collin Duval Tissue Forceps With Box Joint 180mm | 38 |
| PH503568 | Collin Duval Tissue Forceps With Box Joint 200mm | 38 |
| PH503575 | Judd Allis Tissue Forceps With 3 Into 4 Teeth And Box Joint 150mm | 39 |
| PH503576 | Judd Allis Tissue Forceps With 3 Into 4 Teeth And Box Joint 190mm | 39 |
| PH503577 | Lahey Tissue Forceps With 3 Into 3 Teeth And Box Joint Thyroid Vulsellum 150mm | 39 |
| PH503578 | Lahey Tissue Forceps With 3 Into 3 Teeth And Box Joint Thyroid Vulsellum 200mm | 39 |
| PH503580 | Lane Tissue Forceps With 1 Into 2 Teeth And Box Joint 130mm | 39 |
| PH503581 | Lane Tissue Forceps With 1 Into 2 Teeth And Box Joint 150mm | 39 |
| PH503582 | Lane Tissue Forceps With 1 Into 2 Teeth And Box Joint 195mm | 39 |
| PH503583 | Lane Tissue Forceps With 1 Into 2 Teeth And Box Joint 200mm | 39 |
| PH503590 | Littlewood Tissue Forceps With 2 Into 3 Teeth And Box Joint (Littlewoods) 190mm | 39 |
| PH503600 | Poirier Tissue Forceps With 2 Into 3 Teeth And Screw Joint 130mm | 39 |

Product Codes Index

| Code | Description | Page |
|----------|---|------|
| PH503605 | Rutherford Morison Tissue Forceps 4 Into 5 Teeth Box Joint 160mm | 39 |
| PH503607 | Stiles Tissue Forceps With 6 Into 7 Teeth And Box Joint 200mm | 40 |
| PH503611 | Craft Vasectomy Forceps With Sharp Hooked Jaws 150mm | 40 |
| PH503630 | Ring Type Vasectomy Clamps 3.0mm 140mm | 40 |
| PH503631 | Ring Type Vasectomy Clamps 3.5mm 140mm | 40 |
| PH503632 | Ring Type Vasectomy Clamps 4.0mm 140mm | 40 |
| PH503635 | Scalpelless Vasectomy Forceps Curved Sharp 140mm | 40 |
| PH513612 | Denis Browne Intestinal Forceps Spring Pattern With Serrated Grips And Fenestrated Jaws 180mm Straight | 140 |
| PH513617 | Doyen Intestinal Clamp With Longitudinal Serrated Blades And Box Joint 180mm Straight | 140 |
| PH513618 | Doyen Intestinal Clamp With Longitudinal Serrated Blades And Box Joint 180mm Cur | 140 |
| PH513621 | Doyen Intestinal Clamp With 100mm Longitudinal Serrated Blades And Box Joint 230 | 140 |
| PH513622 | Doyen Intestinal Clamp With 100mm Longitudinal Serrated Blades And Box Joint 230 | 140 |
| PH513627 | Doyen Intestinal Clamp With Atraumatic Jaw Blades And Box Joint 180mm Straight | 140 |
| PH513628 | Doyen Intestinal Clamp With Atraumatic Jaw Blades And Box Joint 180mm Curved | 140 |
| PH513631 | Doyen Intestinal Clamp With Atraumatic Jaw Blades And Box Joint 230mm Straight | 140 |
| PH513632 | Doyen Intestinal Clamp Atraumatic Jaw Bj 23cm Cvd | 140 |
| PH513636 | Kocher Intestinal Clamp With 100mm Longitudinal Serrated Light Spring Blades And Box Joint 230mm Straight | 140 |
| PH513637 | Kocher Intestinal Clamp With 125mm Longitudinal Serrated Light Spring Blades And Box Joint 250mm Straight | 140 |
| PH513638 | Kocher Intestinal Clamp With 150mm Longitudinal Serrated Light Spring Blades And Box Joint 280mm Straight | 140 |
| PH513647 | Lane Intestinal Clamp With 75mm Longitudinal Serrated Blades And Box Joint 180mm Straight | 141 |
| PH513648 | Lane Intestinal Clamp With 100mm Longitudinal Serrated Blades And Box Joint 230mm Straight | 141 |
| PH513671 | Parker Kerr Intestinal Clamp With Screw Joint 230mm Straight | 141 |
| PH513676 | Parker Kerr Intestinal Clamp With Screw Joint 250mm Straight | 141 |
| PH513680 | Parker Kerr Intestinal Clamp With Screw Joint 230mm Curved | 141 |
| PH513681 | Parker Kerr Intestinal Clamp With Screw Joint 250mm Curved | 141 |
| PH513691 | Schoemaker Ogilvie Clamp For Hemisphincterectomy With Hollow Blades And Box Joint 155mm Straight | 141 |
| PH513692 | Schoemaker Ogilvie Clamp For Hemisphincterectomy With Hollow Blades And Box Joint 170mm Straight | 141 |
| PH513693 | Schoemaker Ogilvie Clamp For Hemisphincterectomy With Hollow Blades And Box Join | 141 |
| PH523745 | Lahey Cholecystectomy Curved Forceps With Longitudinal Serrations And A Box Join | 67 |
| PH523746 | Lahey Cholecystectomy Curved Forceps With Longitudinal Serrations And A Box Join | 67 |

Product Codes Index

| Code | Description | Page |
|----------|---|------|
| PH523759 | Mixer Cholecystectomy Forceps With A Box Joint 160mm Curved | 67 |
| PH523760 | Mixer Cholecystectomy Forceps With A Box Joint 190mm Curved | 67 |
| PH523763 | Moynihan Cholecystectomy Large Curved Forceps With Box Joint 230mm Curved | 67 |
| PH523765 | O'shaughnessy Cholecystectomy Forceps With A Box Joint 200mm Curved | 67 |
| PH533801 | Duval Lung Forceps With 13mm Wide Jaws And A Box Joint 180mm | 67 |
| PH533802 | Duval Lung Forceps With 13mm Wide Jaws And A Box Joint 200mm | 67 |
| PH533803 | Duval Lung Forceps With 19mm Wide Jaws Screw Joint 200mm | 67 |
| PH533804 | Duval Lung Forceps With 19mm Wide Jaws And A Box Joint 230mm | 67 |
| PH533805 | Duval Lung Forceps With 25mm Wide Jaws Screw Joint 200mm | 67 |
| PH533806 | Duval Lung Forceps With 25mm Wide Jaws Screw Joint 230mm | 67 |
| PH533807 | Duval Lung Forceps With 14mm Wide Jaws Screw Joint 230mm | 67 |
| PH574070 | Liston Bone Cutting Forceps With A Box Joint 130mm Straight | 68 |
| PH574071 | Liston Bone Cutting Forceps With A Box Joint 150mm Straight | 68 |
| PH574072 | Liston Bone Cutting Forceps With A Box Joint 180mm Straight | 68 |
| PH574073 | Liston Bone Cutting Forceps With A Box Joint 200mm Straight | 68 |
| PH574074 | Liston Bone Cutting Forceps With A Box Joint 230mm Straight | 68 |
| PH574112 | Mcindoe Bone Cutting Forceps With Compound Action And Angled On Flat 170mm | 68 |
| PH574113 | Mcindoe Bone Cutting Forceps With Compound Action And Angled On Flat 190mm | 68 |
| PH574115 | Ruskin Liston Bone Cutting Forceps With Compound Action 150mm Straight | 68 |
| PH574116 | Ruskin Liston Bone Cutting Forceps With Compound Action 180mm Straight | 68 |
| PH574118 | Ruskin Liston Bone Cutting Forceps With Compound Action 200mm Straight | 68 |
| PH574119 | Ruskin Liston Bone Cutting Forceps With Compound Action 150mm Curved | 68 |
| PH574120 | Ruskin Liston Bone Cutting Forceps With Compound Action 180mm Curved | 68 |
| PH574121 | Ruskin Liston Bone Cutting Forceps Curved Compound Action 200mm Curved | 68 |
| PH574126 | Stamm Bone Cutting Forceps 150mm | 68 |
| PH574136 | Horsley Bone Cutting Forceps With Compound Action 270mm Curved | 68 |
| PH584147 | Cairn Bone Rongeur With Fine Jaws And Angle On Flat 170mm | 144 |
| PH584148 | Cleveland Bone Rongeur Curved With A Box Joint 140mm Curved | 144 |
| PH584149 | Daniel Miniature Bone Rongeur With 1.5mm Bite And Curved Handles With A Screw Joint | 144 |
| PH584150 | Glasgow Pattern Bone Rongeur With 4mm Bite And Slightly Curved On Flat | 144 |
| PH584152 | Jansen Middleton Bone Rongeur With Compound Action 185mm Straight | 144 |
| PH584153 | Jansen Middleton Bone Rongeur With Compound Action 185mm Curved | 144 |
| PH584154 | Lempert Bone Rongeur With Simple Action And A Box Joint 190mm Straight | 145 |

Product Codes Index

| Code | Description | Page |
|----------|--|------|
| PH584157 | Luer Bone Rongeur With Simple Action And A Box Joint 170mm Straight | 145 |
| PH584158 | Luer Bone Rongeur With Simple Action And A Box Joint 170mm Curved | 145 |
| PH584161 | Leksell Bone Rongeur With 4mm Bite And Compound Action 230mm Curved | 145 |
| PH584162 | Leksell Bone Rongeur With 6mm Bite And Compound Action 230mm Curved | 145 |
| PH584163 | Leksell Bone Rongeur With 8mm Bite And Compound Action 230mm Curved | 145 |
| PH584165 | Northfield Bone Rongeur With 10mm Bite And Simple Action With A Box Joint 200mm | 145 |
| PH584170 | Olivecrona Bone Rongeur Angled To Side With 5mm Bite And A Compound Action 200mm | 146 |
| PH584171 | Olivecrona Bone Rongeur Angled To Side With 6mm Bite And A Compound Action 230mm Angled | 146 |
| PH584175 | Pennybacker Bone Rongeur Straight Fine Jaws With 4mm Bite And A Simple Action Wi | 146 |
| PH584178 | Ruskin Bone Rongeur Curved On Flat With 6mm Bite And A Compound Action 180mm Curved | 146 |
| PH584179 | Ruskin Bone Rongeur Curved On Flat With 6mm Bite And A Compound Action 230mm Curved | 146 |
| PH584180 | Sauerburch Bone Rongeur With Round Cutting Jaws And A Compound Action 300mm Straight | 146 |
| PH584185 | Semb Bone Rongeur Slightly Curved With Double Action 20mm Wide Jaws And A Compound Action 240mm Curved | 146 |
| PH584190 | Stille Luer Bone Rongeur Straight With A Compound Action (Or Luer Stille) 220mm | 146 |
| PH584191 | Stille Luer Bone Rongeur Curved With A Compound Action (Or Luer Stille) 220mm Curved | 147 |
| PH584194 | Stille Bone Rongeur Angled To One Side With A Compound Double Action 230mm Angle | 147 |
| PH584195 | Trotter Bone Rongeur Angled To One Side With A 10mm Bite And Simple Action With | 147 |
| PH584197 | Ward Bone Rongeur Small Slightly Curved On Flat With A Compound Action 180mm Curved | 147 |
| PH584198 | Wilms Bone Rongeur Angled On Flat With Narrow Jaws And A 4mm Bite With A Simple Action 200mm Angled | 147 |
| PH584200 | Phoenix Bone Nibbler 3mm Bite With A Compound Action 230mm 230mm Straight | 147 |
| PH584201 | Phoenix Bone Nibbler 6mm Bite With A Compound Action 230mm Straight | 147 |
| PH584202 | Zielke Bone Nibbler 250mm x 3mm Bite 250mm Curved | 148 |
| PH584250 | Freidman Bone Nibbler Single Action Small (Friedman) 125mm Straight | 148 |
| PH584253 | Mini Bone Nibbler 1mm Jaw | 148 |
| PH614305 | Trethowen Bone Lever 215mm | 149 |
| PH614306 | Trethowen Angled Bone Lever Angled | 149 |
| PH614310 | Bristow Bone Elevator Sharp Blade With A Hollow Handle 230mm Curved | 149 |
| PH614314 | Lane Bone Elevator Slight Curve With A Solid Handle 240mm Curved | 149 |
| PH614315 | Lane Bone Elevator With A Slight Curve And A Ring Handle 265mm Curved | 150 |
| PH614316 | Lane Bone Elevator With A Serrated Slight Curved End And A Finger Ring Handle 260mm Curved | 150 |

Product Codes Index

| Code | Description | Page |
|----------|---|------|
| PH614317 | Mitchel Bone Elevator Curved Blade With A Hollow Handle 200mm Curved | 150 |
| PH614321 | Cobb Elevators 10mm 275mm | 150 |
| PH614322 | Cobb Elevators 13mm 275mm | 150 |
| PH614323 | Cobb Elevators 19mm 275mm | 150 |
| PH624321 | Pennybacker Elevator 4mm Blade Width Straight 170mm | 151 |
| PH624322 | Pennybacker Elevator 6mm Blade Width 170mm | 151 |
| PH624323 | Lempert Periosteal Elevator 170mm | 151 |
| PH624325 | Doyen Periosteal Elevator Curved To Left Adult 165mm Curved | 151 |
| PH624326 | Doyen Periosteal Elevator Curved To Right Adult 165mm Curved | 151 |
| PH624327 | Doyen Periosteal Elevator Curved To Left Child 165mm Curved | 151 |
| PH624328 | Doyen Periosteal Elevator Curved To Right Child 165mm Curved | 151 |
| PH624329 | Ph654329 Ward Periosteal Elevator (Plastic And Oral Surgery) 170mm | 151 |
| PH624330 | Farabeuf Bone Rugine With Chisel Edge With A Flat Handle (Elevator) 150mm Straigt | 152 |
| PH624332 | Jones Rugine Pistol Shaped Handle Fan Shaped End And Convex Sharp Edge 185mm | 152 |
| PH624341 | Martini Double Ended Bone Curette 150mm | 153 |
| PH624342 | Stanley Boyd Bone Curette Small Double End Oval Scoops 150mm | 153 |
| PH624343 | Stanley Boyd Bone Curette Medium Double End Oval Scoops 200mm | 153 |
| PH624344 | Stanley Boyd Bone Curette Large Double End Oval Scoops 280mm | 153 |
| PH624346 | Volkman Spoon Bone Curette Small Double Ended 215mm | 153 |
| PH624347 | Volkman Spoon Bone Curette Medium Double Ended 215mm | 153 |
| PH624348 | Volkman Spoon Bone Curette Large Double Ended 215mm | 153 |
| PH624350 | Williger Bone Curette Small Double Ended 125mm | 153 |
| PH624351 | Williger Bone Curette Medium Double Ended 125mm | 153 |
| PH624352 | Williger Bone Curette Large Double Ended 125mm | 153 |
| PH624354 | Kocher Bone Hook Sharp Pointed With A Fenestrated Handle 250mm | 153 |
| PH624370 | Lucae Bone Mallet 8 Ounces Light Pattern Toffee Mallet 190mm | 154 |
| PH624383 | Heath Bone Mallet 35mm Diameter Head & Long Handle | 154 |
| PH624384 | Heath Bone Mallet 18 Oz | 154 |
| PH624385 | Heath Bone Mallet 28 Oz | 154 |
| PH624386 | Heath Bone Mallet 36 Oz | 154 |
| PH624398 | Light Weight Tufnol Mallet | 154 |
| PH624399 | Heavy Tufnol Mallet 270mm | 154 |
| PH624435 | Mini Stille Bone Chisel 3mm Wide 150mm | 156 |

Product Codes Index

| Code | Description | Page |
|----------|---|------|
| PH624436 | Mini Stille Bone Chisel 5mm Wide 150mm | 156 |
| PH624437 | Mini Stille Bone Chisel 7mm Wide 150mm | 156 |
| PH624438 | Mini Stille Bone Chisel 9mm Wide 150mm | 156 |
| PH624439 | Mini Stille Bone Chisel 12mm Wide 150mm | 156 |
| PH624440 | Stille Bone Chisel 10mm Wide x 200mm | 156 |
| PH624441 | Stille Bone Chisel 15mm Wide 200mm | 156 |
| PH624442 | Stille Bone Chisel 20mm Wide 200mm | 156 |
| PH624443 | Stille Bone Chisel 25mm Wide 200mm | 156 |
| PH624444 | Stille Bone Chisel 30mm Wide 200mm | 156 |
| PH624496 | Duray French Osteotome 5mm Wide 140mm | 158 |
| PH624497 | Duray French Osteotome 7mm Wide 140mm | 158 |
| PH624498 | Duray French Osteotome 9mm Wide 140mm | 158 |
| PH624499 | Duray French Osteotome 11mm Wide 140mm | 158 |
| PH624509 | Smith Peterson Bone Osteotomes 6mm Wide 200mm | 158 |
| PH624510 | Smith Peterson Bone Osteotomes 10mm Wide 200mm | 158 |
| PH624511 | Smith Peterson Bone Osteotomes 13mm Wide 200mm | 158 |
| PH624512 | Smith Peterson Bone Osteotomes 20mm Wide 200mm | 158 |
| PH624513 | Smith Peterson Bone Osteotomes 25mm Wide 200mm | 158 |
| PH624514 | Smith Peterson Bone Osteotomes 30mm Wide 200mm | 158 |
| PH624515 | Stille Bone Osteotomes 2mm Wide 200mm | 158 |
| PH624516 | Stille Bone Osteotomes 5mm Wide 200mm | 158 |
| PH624517 | Stille Bone Osteotomes 6mm Wide 200mm | 158 |
| PH624518 | Stille Bone Osteotomes 7mm Wide 200mm | 158 |
| PH624519 | Stille Bone Osteotomes 8mm Wide 200mm | 158 |
| PH624520 | Stille Bone Osteotomes 9mm Wide 200mm | 158 |
| PH624521 | Stille Bone Osteotomes 10mm Wide 200mm | 158 |
| PH624522 | Stille Bone Osteotomes 11mm Wide 200mm | 158 |
| PH624523 | Stille Bone Osteotomes 13mm Wide 200mm | 158 |
| PH624524 | Stille Bone Osteotomes 15mm Wide 200mm | 158 |
| PH624525 | Stille Bone Osteotomes 20mm Wide 200mm | 158 |
| PH624526 | Stille Bone Osteotomes 25mm Wide 200mm | 158 |
| PH624527 | Stille Bone Osteotomes 30mm Wide 200mm | 158 |
| PH624528 | Stille Bone Osteotomes 19mm Wide (3/4inch) Curved | 158 |

Product Codes Index

| Code | Description | Page |
|-------------|---|------|
| PH624529 | Stille Bone Osteotomes 15mm Wide Curved | 158 |
| PH624530 | Stille Bone Osteotomes 10mm Wide Curved 200mm | 158 |
| PH624539 | Lambotte Osteotome 1mm Wide 150mm | 159 |
| PH624540 | Lambotte Osteotome 2mm Wide 150mm | 159 |
| PH624541 | Lambotte Osteotome 3mm Wide 150mm | 159 |
| PH624542 | Lambotte Osteotome 4mm Wide 150mm | 159 |
| PH624543 | Lambotte Osteotome 6mm Wide 150mm | 159 |
| PH624544 | Lambotte Osteotome 8mm Wide 150mm | 159 |
| PH624545 | Lambotte Osteotome 10mm Wide 150mm | 159 |
| PH624546 | Lambotte Osteotome 12mm Wide 150mm | 159 |
| PH624547 | Lambotte Osteotome 15mm Wide 200mm | 159 |
| PH624548 | Lambotte Osteotome 20mm Wide 200mm | 159 |
| PH624549 | Lambotte Osteotome 25mm Wide 200mm | 159 |
| PH624550 | Lambotte Osteotome 30mm Wide 200mm | 159 |
| PH624552 | Lambotte Osteotome 10mm Wide x 200mm | 159 |
| PH624553/10 | Gillies (French) Osteotome 10mm Wide 140mm | 159 |
| PH624553/11 | Gillies (French) Osteotome 11mm Wide 140mm | 159 |
| PH624553/2 | Gillies (French) Osteotome 2mm Wide 140mm | 159 |
| PH624553/3 | Gillies (French) Osteotome 3mm Wide 140mm | 159 |
| PH624553/4 | Gillies (French) Osteotome 4mm Wide 140mm | 159 |
| PH624553/5 | Gillies (French) Osteotome 5mm Wide 140mm | 159 |
| PH624553/6 | Gillies (French) Osteotome 6mm Wide 140mm | 159 |
| PH624553/7 | Gillies (French) Osteotome 7mm Wide 140mm | 159 |
| PH624553/8 | Gillies (French) Osteotome 8mm Wide 140mm | 159 |
| PH624553/9 | Gillies (French) Osteotome 9mm Wide 140mm | 159 |
| PH624554 | Jayes Osteotome 1mm Mini Lambotte 125mm | 159 |
| PH624555 | Jayes Osteotome 2mm Mini Lambotte 125mm | 159 |
| PH624556 | Jayes Osteotome 4mm Mini Lambotte 125mm | 159 |
| PH624557 | Ward Osteotome 3mm Wide 180mm | 159 |
| PH624559 | Ward Osteotome 5mm Wide 180mm | 159 |
| PH624569 | Jenkins Bone Gouge 3mm Wide 150mm | 157 |
| PH624570 | Jenkins Bone Gouge 4mm Wide 150mm | 157 |
| PH624571 | Jenkins Bone Gouge 5mm Wide 150mm | 157 |

Product Codes Index

| Code | Description | Page |
|----------|--|------|
| PH624572 | Jenkins Bone Gouge 6mm Wide 150mm | 157 |
| PH624573 | Jenkins Bone Gouge 8mm Wide 150mm | 157 |
| PH624574 | Jenkins Bone Gouge 10mm Wide 150mm | 157 |
| PH624575 | Jenkins Bone Gouge 12mm Wide 150mm | 157 |
| PH624576 | Jenkins Bone Gouge 14mm Wide 150mm | 157 |
| PH624600 | Stille Bone Gouge 10mm Wide 200mm | 157 |
| PH624601 | Stille Bone Gouge 13mm Wide 200mm | 157 |
| PH624602 | Stille Bone Gouge 15mm Wide 200mm | 157 |
| PH624603 | Stille Bone Gouge 20mm Wide 200mm | 157 |
| PH624604 | Stille Bone Gouge 25mm Wide 200mm | 157 |
| PH624605 | Stille Bone Gouge 30mm Wide 200mm | 157 |
| PH624635 | Coupland Dental Gouge Small (3mm) 170mm | 157 |
| PH624636 | Coupland Dental Gouge Medium (4mm) 170mm | 157 |
| PH624637 | Coupland Dental Gouge Large (5mm) 170mm | 157 |
| PH624638 | Tilley Gouge V Shaped 125mm | 157 |
| PH625001 | Fickling Mallet 200mm | 155 |
| PH625005 | Cottle Mallet 30mm Diameter Head 190mm | 155 |
| PH626006 | Capener Tufnol Handled Chisel 10mm Wide (Capner) 250mm | 156 |
| PH626007 | Capener Tufnol Handled Chisel 12mm Wide x 250mm | 156 |
| PH626008 | Capener Tufnol Handled Chisel 15mm Wide 250mm | 156 |
| PH626009 | Capener Tufnol Handled Chisel 20mm Wide 250mm | 156 |
| PH626010 | Capener Tufnol Handled Chisel 25mm Wide x 250mm | 156 |
| PH626011 | Capener Tufnol Handled Chisel 30mm Wide x 250mm | 156 |
| PH627000 | Tufnol Handled Straight Osteotome 10mm Wide 265mm Straight | 160 |
| PH627001 | Tufnol Handled Straight Osteotome 12mm Wide 265mm Straight | 160 |
| PH627002 | Tufnol Handled Straight Osteotome 15mm Wide 265mm Straight | 160 |
| PH627003 | Tufnol Handled Straight Osteotome 20mm Wide 265mm Straight | 160 |
| PH627004 | Tufnol Handled Straight Osteotome 25mm Wide 265mm Straight | 160 |
| PH627005 | Tufnol Handled Straight Osteotome 30mm Wide 265mm Straight | 160 |
| PH627006 | Tufnol Handled Curved Osteotome 10mm Wide 265mm Curved | 160 |
| PH627007 | Tufnol Handled Curved Osteotome 12mm Wide 265mm Curved | 160 |
| PH627008 | Tufnol Handled Curved Osteotome 15mm Wide 265mm Curved | 160 |
| PH627009 | Tufnol Handled Curved Osteotome 20mm Wide 265mm Curved | 160 |

Product Codes Index

| Code | Description | Page |
|----------|--|------|
| PH627010 | Tufnol Handled Curved Osteotome 25mm Wide 265mm Curved | 160 |
| PH627011 | Tufnol Handled Curved Osteotome 30mm Wide 265mm Curved | 160 |
| PH627015 | Capener Tufnol Handled Straight Osteotome 10mm Wide 250mm Straight | 160 |
| PH627016 | Capener Tufnol Handled Straight Osteotome 12mm Wide 250mm Straight | 160 |
| PH627017 | Capener Tufnol Handled Straight Osteotome 15mm Wide 250mm Straight | 160 |
| PH627018 | Capener Tufnol Handled Straight Osteotome 20mm Wide 250mm Straight | 160 |
| PH627019 | Capener Tufnol Handled Straight Osteotome 25mm Wide 250mm Straight | 160 |
| PH627020 | Capener Tufnol Handled Straight Osteotome 30mm Wide 250mm Straight | 160 |
| PH627021 | Capener Tufnol Handled Curved Osteotome 10mm Wide 250mm Curved | 160 |
| PH627022 | Capener Tufnol Handled Curved Osteotome 12mm Wide 250mm Curved | 160 |
| PH627023 | Capener Tufnol Handled Curved Osteotome 15mm Wide 250mm Curved | 160 |
| PH627024 | Capener Tufnol Handled Curved Osteotome 20mm Wide 250mm Curved | 160 |
| PH627025 | Capener Tufnol Handled Curved Osteotome 25mm Wide 250mm Curved | 160 |
| PH627026 | Capener Tufnol Handled Curved Osteotome 30mm Wide 250mm Curved | 160 |
| PH628000 | Tufnol Handled Straight Gouge 10mm Wide 270mm | 157 |
| PH628001 | Tufnol Handled Straight Gouge 12mm Wide 270mm | 157 |
| PH628002 | Tufnol Handled Straight Gouge 15mm Wide 270mm | 157 |
| PH628003 | Tufnol Handled Straight Gouge 20mm Wide 270mm | 157 |
| PH628004 | Tufnol Handled Straight Gouge 25mm Wide 270mm | 157 |
| PH628005 | Tufnol Handled Straight Gouge 30mm Wide 270mm | 157 |
| PH628010 | Capener Tufnol Handled Cranked Gouge 6mm Wide (Capner) 235mm | 157 |
| PH628011 | Capener Tufnol Handled Cranked Gouge 8mm Wide 235mm | 157 |
| PH628012 | Capener Tufnol Handled Cranked Gouge 11mm Wide 235mm | 157 |
| PH628013 | Capener Tufnol Handled Cranked Gouge 16mm Wide 235mm | 157 |
| PH644687 | Vehmhren Rib Shears With Sliding Blade Action | 161 |
| PH644691 | Craford's Rib Approximator 260mm | 161 |
| PH644692 | Sellor Rib Approximator And Spreader Combined 200mm | 161 |
| PH644693 | Sellor Rib Spreader With 3 Pairs Of Detachable Blades And Locking Key (Sternal Spreader) 200mm | 162 |
| PH644694 | Burford Finochietto Rib Spreader With 2 Pairs Of Detachable Blades 290mm | 162 |
| PH644696 | Finochietto Rib Spreader Baby Size With 18mm x 21mm Blades And A 75mm Spread | 163 |
| PH644698 | Finochietto Rib Spreader Small Size With 30mm x 45mm Blades And A 165mm Spread | 163 |
| PH644699 | Finochietto Rib Spreader Medium Size With 45mm x 65mm Blades And A 200mm Spread | 163 |

Product Codes Index

| Code | Description | Page |
|----------|---|------|
| PH644700 | Finochietto Rib Spreader Large Size With 70mm x 65mm Blades And A 250mm Spread | 163 |
| PH644710 | Watersons Rib Spreader Infant Size With 40mm Wide x 20mm Deep Blades And 105mm Spread | 163 |
| PH654358 | Fomon Bone Rasp No. 1 & 2 Tc Cut 205mm | 186 |
| PH654359 | Fomon Bone Rasp No. 3 & 4 Tc Cut 205mm | 186 |
| PH654360 | Fomon Bone Rasp No. 5 & 6 Tc Cut 205mm | 186 |
| PH654361 | Fomon Bone Rasp No. 7 & 8 Tc Cut 205mm | 186 |
| PH654362 | Fomon Bone Rasp No. 9 & 10 Tc Cut 205mm | 186 |
| PH654363 | Fomon Rasp Fine 215mm | 186 |
| PH654366 | Miller Bone Rasp 180mm Large | 152 |
| PH654725 | Bergmann Plaster Shears/ Bandage Scissors With Angular Blades 230mm | 166 |
| PH654727 | Bohler Plaster Shears With A Screw Joint 230mm | 166 |
| PH654729 | Bruns Plaster Shears With A Toothed Lower Blade 240mm | 166 |
| PH654733 | Guys Heavy Pattern Plaster Shears 300mm | 166 |
| PH654737 | Stille Lorenz Plaster Shears 230mm | 167 |
| PH654738 | Stille Lorenz Plaster Shears 265mm | 167 |
| PH654739 | Stille Lorenz Plaster Shears 300mm | 167 |
| PH654740 | Stille Lorenz Plaster Shears 355mm | 167 |
| PH654741 | Stille Lorenz Plaster Shears 380mm | 167 |
| PH654754 | Henning Plaster Cast Spreaders With A Screw Joint 280mm | 167 |
| PH654767 | Caspar Rongeur Straight 2mm Straight | 170 |
| PH654768 | Caspar Rongeur Straight 150mm Length Shaft 3mm Straight | 170 |
| PH654770 | Cushing Rongeur Straight Jaw 3mm Bite 130mm Effective Length Straight | 170 |
| PH654771 | Cushing Rongeur Straight Jaw 4mm Bite 130mm Effective Length Straight | 170 |
| PH654772 | Cushing Rongeur Straight Jaw 5mm Bite 130mm Effective Length Straight | 170 |
| PH654773 | Cushing Rongeur Straight Jaw 7mm Bite 130mm Effective Length Straight | 170 |
| PH654774 | Cushing Rongeur Angled Upward Jaw 3mm Bite 130mm Effective Length Angled | 170 |
| PH654775 | Cushing Rongeur Angled Upward Jaw 4mm Bite 130mm Effective Length Angled | 170 |
| PH654776 | Cushing Rongeur Angled Upward Jaw 5mm Bite 130mm Effective Length Angled | 170 |
| PH654777 | Cushing Rongeur Angled Upward Jaw 7mm Bite 130mm Effective Length Angled | 170 |
| PH654783 | Beck Rongeur Straight 130mm Shaft Length With A 2mm Bite For Intervertebral Disc Straight | 171 |
| PH654833 | Clowards Rongeur 3mm Bite Straight Jaw 150mm Effective Length | 171 |
| PH654834 | Clowards Rongeur 5mm Bite Straight Jaw 150mm Effective Length | 171 |

Product Codes Index

| Code | Description | Page |
|----------|---|------|
| PH654835 | Clowards Rongeur 6mm Bite Straight Jaw 150mm Effective Length | 171 |
| PH654905 | Love Grunwald Pituitary Rongeur 200mm Shaft 2mm Straight Cut | 171 |
| PH654906 | Love Grunwald Pituitary Rongeur 200mm Shaft 3mm Straight Cut | 171 |
| PH654907 | Love Grunwald Pituitary Rongeur 200mm Shaft 4mm Straight Cut | 171 |
| PH654908 | Love Grunwald Pituitary Rongeur 200mm Shaft 5mm Straight Cut | 171 |
| PH654910 | Love Grunwald Pituitary Rongeur 200mm Shaft 2mm Upwards Cut | 171 |
| PH654911 | Love Grunwald Pituitary Rongeur 200mm Shaft 3mm Upwards Cut | 171 |
| PH654912 | Love Grunwald Pituitary Rongeur 200mm Shaft 4mm Upwards Cut | 171 |
| PH654913 | Love Grunwald Pituitary Rongeur 200mm Shaft 5mm Upwards Cut | 171 |
| PH654918 | Pennybacker Ronguer Straight 130mm Shaft 3mm Bite | 171 |
| PH654950 | Semi Detach Kerrison Rong Upcut 180mm 45 Deg x 1mm | 172 |
| PH654951 | Semi Detach Kerrison Rong Upcut 180mm 45 Deg x 2mm | 172 |
| PH654952 | Semi Detach Kerrison Rong Upcut 180mm 45 Deg x 3mm | 172 |
| PH654953 | Semi Detach Kerrison Rong Upcut 180mm 45 Deg x 4mm | 172 |
| PH654954 | Semi Detach Kerrison Rong Upcut 180mm 45 Deg x 5mm | 172 |
| PH654955 | Semi Detach Kerrison Rong Upcut 180mm 90 Deg x 1mm | 172 |
| PH654956 | Semi Detach Kerrison Rong Upcut 180mm 90 Deg x 2mm | 172 |
| PH654957 | Semi Detach Kerrison Rong Upcut 180mm 90 Deg x 3mm | 172 |
| PH654958 | Semi Detach Kerrison Rong Upcut 180mm 90 Deg x 4mm | 172 |
| PH654959 | Semi Detach Kerrison Rong Upcut 180mm 90 Deg x 5mm | 172 |
| PH654960 | Semi Detach Kerrison Rong Upcut 200mm 45 Deg x 1mm | 173 |
| PH654961 | Semi Detach Kerrison Rong Upcut 200mm 45 Deg x 2mm | 173 |
| PH654962 | Semi Detach Kerrison Rong Upcut 200mm 45 Deg x 3mm | 173 |
| PH654963 | Semi Detach Kerrison Rong Upcut 200mm 45 Deg x 4mm | 173 |
| PH654964 | Semi Detach Kerrison Rong Upcut 200mm 45 Deg x 5mm | 173 |
| PH654965 | Semi Detach Kerrison Rong Upcut 200mm 90 Deg x 1mm | 173 |
| PH654966 | Semi Detach Kerrison Rong Upcut 200mm 90 Deg x 2mm | 173 |
| PH654967 | Semi Detach Kerrison Rong Upcut 200mm 90 Deg x 3mm | 173 |
| PH654968 | Semi Detach Kerrison Rong Upcut 200mm 90 Deg x 4mm | 173 |
| PH654969 | Semi Detach Kerrison Rong Upcut 200mm 90 Deg x 5mm | 173 |
| PH654970 | Semi Detach Kerrison Rong Upcut 210mm 45 Deg x 1mm | 174 |
| PH654971 | Semi Detach Kerrison Rong Upcut 210mm 45 Deg x 2mm | 174 |

Product Codes Index

| Code | Description | Page |
|----------|--|------|
| PH654972 | Semi Detach Kerrison Rong Upcut 210mm 45 Deg x 3mm | 174 |
| PH654973 | Semi Detach Kerrison Rong Upcut 210mm 45 Deg x 4mm | 174 |
| PH654974 | Semi Detach Kerrison Rong Upcut 210mm 45 Deg x 5mm | 174 |
| PH654975 | Semi Detach Kerrison Rong Upcut 210mm 90 Deg x 1mm | 174 |
| PH654976 | Semi Detach Kerrison Rong Upcut 210mm 90 Deg x 2mm | 174 |
| PH654977 | Semi Detach Kerrison Rong Upcut 210mm 90 Deg x 3mm | 174 |
| PH654978 | Semi Detach Kerrison Rong Upcut 210mm 90 Deg x 4mm | 174 |
| PH654979 | Semi Detach Kerrison Rong Upcut 210mm 90 Deg x 5mm | 174 |
| PH654980 | Semi Detach Kerrison Upcut 18cm 45 Deg x 1mm/Black | 174 |
| PH654981 | Semi Detach Kerrison Upcut 18cm 45 Deg x 2mm/Black | 174 |
| PH654982 | Semi Detach Kerrison Upcut 18cm 45 Deg x 3mm/Black | 174 |
| PH654983 | Semi Detach Kerrison Upcut 18cm 45 Deg x 4mm/Black | 174 |
| PH654984 | Semi Detach Kerrison Upcut 18cm 45 Deg x 5mm/Black | 174 |
| PH654985 | Semi Detach Kerrison Upcut 18cm 90 Deg x 1mm/Black | 174 |
| PH654986 | Semi Detach Kerrison Upcut 18cm 90 Deg x 2mm/Black | 174 |
| PH654987 | Semi Detach Kerrison Upcut 18cm 90 Deg x 3mm/Black | 174 |
| PH654988 | Semi Detach Kerrison Upcut 18cm 90 Deg x 4mm/Black | 174 |
| PH654989 | Semi Detach Kerrison Upcut 18cm 90 Deg x 5mm/Black | 174 |
| PH654990 | Semi Detach Kerrison Upcut 20cm 45 Deg x 1mm/Black | 175 |
| PH654991 | Semi Detach Kerrison Upcut 20cm 45 Deg x 2mm/Black | 175 |
| PH654992 | Semi Detach Kerrison Upcut 20cm 45 Deg x 3mm/Black | 175 |
| PH654993 | Semi Detach Kerrison Upcut 20cm 45 Deg x 4mm/Black | 175 |
| PH654994 | Semi Detach Kerrison Upcut 20cm 45 Deg x 5mm/Black | 175 |
| PH654995 | Semi Detach Kerrison Upcut 21cm 45 Deg x 1mm/Black | 175 |
| PH654996 | Semi Detach Kerrison Upcut 21cm 45 Deg x 2mm/Black | 175 |
| PH654997 | Semi Detach Kerrison Upcut 21cm 45 Deg x 3mm/Black | 175 |
| PH654998 | Semi Detach Kerrison Upcut 21cm 45 Deg x 4mm/Black | 175 |
| PH654999 | Semi Detach Kerrison Upcut 21cm 45 Deg x 5mm/Black | 175 |
| PH655000 | Semi Detach Kerrison Upcut 20cm 90 Deg x 1mm/Black | 175 |
| PH655001 | Semi Detach Kerrison Upcut 20cm 90 Deg x 2mm/Black | 175 |
| PH655002 | Semi Detach Kerrison Upcut 20cm 90 Deg x 3mm/Black | 175 |
| PH655003 | Semi Detach Kerrison Upcut 20cm 90 Deg x 4mm/Black | 175 |
| PH655004 | Semi Detach Kerrison Upcut 20cm 90 Deg x 5mm/Black | 175 |

Product Codes Index

| Code | Description | Page |
|----------|--|------|
| PH655005 | Semi Detach Kerrison Upcut 21cm 90 Deg x 1mm/Black | 175 |
| PH655006 | Semi Detach Kerrison Upcut 21cm 90 Deg x 2mm/Black | 175 |
| PH655007 | Semi Detach Kerrison Upcut 21cm 90 Deg x 3mm/Black | 175 |
| PH655008 | Semi Detach Kerrison Upcut 21cm 90 Deg x 4mm/Black | 175 |
| PH655009 | Semi Detach Kerrison Upcut 21cm 90 Deg x 5mm/Black | 175 |
| PH664560 | Lambotte Osteotome 7mm Wide | 159 |
| PH665257 | Castroviejo Needle Holder Tungsten Carbide Smooth Straight Jaws With A Spring Gold Plated Handle And Centre Lock With A Screw Joint 140mm | 89 |
| PH665258 | Castroviejo Needle Holder Tungsten Carbide Smooth Curved Jaws With A Spring Gold Plated Handle And Centre Lock With A Screw Joint 140mm | 89 |
| PH665259 | Castroviejo Needle Holder Tungsten Carbide Cross Serrated Straight Jaws With A Spring Gold Plated Handle And Centre Lock With A Screw Joint 140mm Straight | 89 |
| PH665260 | Castroviejo Needle Holder Tungsten Carbide Cross Serrated Curved Jaws With A Spring Gold Plated Handle And Centre Lock With A Screw Joint 140mm | 89 |
| PH675000 | Split Proctoscope | N/A |
| PH675310 | Goligher Rectal Speculum With Third Blade 76mm | 186 |
| PH675341 | Eisenhammer Rectal Speculum | 186 |
| PH675360 | St Marks Rectal Dilator Large Hollow Stainless Steel 75mm | 187 |
| PH675362 | St Marks Rectal Dilator Medium Hollow Stainless Steel 75mm | 187 |
| PH675364 | St Marks Rectal Dilator Small Hollow Stainless Steel 75mm | 187 |
| PH695506 | Mathieu Urethral Grasping Forceps With Crocodile Action And Serrated Jaws With 200mm Shaft | 187 |
| PH695507 | Mathieu Urethral Grasping Forceps With Crocodile Action And Serrated Jaws With 280mm Shaft | 187 |
| PH695510 | Thomson Walker Penile Clamp With Cross Action 115mm | 187 |
| PH695529 | Clutton Urethral Bougies (Sounds) Individual Size 6 - 10 Fg 267mm | 188 |
| PH695530 | Clutton Urethral Bougies (Sounds) Individual Size 8 - 12 Fg 267mm | 188 |
| PH695531 | Clutton Urethral Bougies (Sounds) Individual Size 10 - 14 Fg 267mm | 188 |
| PH695532 | Clutton Urethral Bougies (Sounds) Individual Size 12 - 16 Fg 267mm | 188 |
| PH695533 | Clutton Urethral Bougies (Sounds) Individual Size 14 - 18 Fg 267mm | 188 |
| PH695534 | Clutton Urethral Bougies (Sounds) Individual Size 16 - 20 Fg 267mm | 188 |
| PH695535 | Clutton Urethral Bougies (Sounds) Individual Size 18 - 22 Fg 267mm | 188 |
| PH695536 | Clutton Urethral Bougies (Sounds) Individual Size 20 - 24 Fg 267mm | 188 |
| PH695537 | Clutton Urethral Bougies (Sounds) Individual Size 22 - 26 Fg 267mm | 188 |
| PH695538 | Clutton Urethral Bougies (Sounds) Individual Size 24 - 28 Fg 267mm | 188 |
| PH695539 | Clutton Urethral Bougies (Sounds) Individual Size 26 - 30 Fg 267mm | 188 |

Product Codes Index

| Code | Description | Page |
|----------|--|------|
| PH695540 | Clutton Urethral Bougies (Sounds) Individual Size 28 - 32 Fg 267mm | 188 |
| PH695543 | Lister Urethral Bougies (Sounds) Set Of 12 Sizes 0 / 1 Eg To 9 / 12 Eg | 188 |
| PH695544 | Lister Urethral Bougies (Sounds) Individual Size 0 / 1 E.g. 280mm | 188 |
| PH695545 | Lister Urethral Bougies (Sounds) Individual Size 0 / 2 E.g. 280mm | 188 |
| PH695546 | Lister Urethral Bougies (Sounds) Individual Size 0.5 / 3 E.g. 280mm | 188 |
| PH695547 | Lister Urethral Bougies (Sounds) Individual Size 1 / 4 E.g. 280mm | 188 |
| PH695548 | Lister Urethral Bougies (Sounds) Individual Size 2 / 5 E.g. 280mm | 188 |
| PH695549 | Lister Urethral Bougies (Sounds) Individual Size 3 / 6 E.g. 280mm | 188 |
| PH695550 | Lister Urethral Bougies (Sounds) Individual Size 4 / 7 E.g. 280mm | 188 |
| PH695551 | Lister Urethral Bougies (Sounds) Individual Size 5 / 8 E.g. 280mm | 188 |
| PH695552 | Lister Urethral Bougies (Sounds) Individual Size 6 / 9 E.g. 280mm | 188 |
| PH695553 | Lister Urethral Bougies (Sounds) Individual Size 7 / 10 E.g. 280mm | 188 |
| PH695554 | Lister Urethral Bougies (Sounds) Individual Size 8 / 11 E.g. 280mm | 188 |
| PH695555 | Lister Urethral Bougies (Sounds) Individual Size 9 / 12 E.g. 280mm | 188 |
| PH695556 | Lister Urethral Bougies (Sounds) Individual Size 10 / 13 E.g. 280mm | 188 |
| PH695557 | Lister Urethral Bougies (Sounds) Individual Size 11 / 14 E.g. 280mm | 188 |
| PH695558 | Lister Urethral Bougies (Sounds) Individual Size 12 / 15 E.g. 280mm | 188 |
| PH695559 | Lister Urethral Bougies (Sounds) Individual Size 13 / 16 E.g. 280mm | 188 |
| PH695560 | Lister Urethral Bougies (Sounds) Individual Size 14 / 17 E.g. 280mm | 188 |
| PH695561 | Lister Urethral Bougies (Sounds) Individual Size 15 / 18 E.g. 280mm | 188 |
| PH695573 | Whyndham Powell Urethral Bougies (Sound) Size 8fg | 188 |
| PH695574 | Whyndham Powell Urethral Bougies (Sound) Size 10fg | 188 |
| PH695575 | Whyndham Powell Urethral Bougies (Sounds) | 188 |
| PH695576 | Whyndham Powell Urethral Bougies (Sounds) Individual Size 12 Fg | 188 |
| PH695577 | Whyndham Powell Urethral Bougies (Sounds) Individual Size 13 Fg | 188 |
| PH695578 | Whyndham Powell Urethral Bougies (Sounds) Individual Size 14 Fg | 188 |
| PH695579 | Whyndham Powell Urethral Bougies (Sounds) Individual Size 15 Fg | 188 |
| PH695580 | Whyndham Powell Urethral Bougies (Sounds) Individual Size 16 Fg | 188 |
| PH695581 | Whyndham Powell Urethral Bougies (Sounds) Individual Size 17 Fg | 188 |
| PH695582 | Whyndham Powell Urethral Bougies (Sounds) Individual Size 18 Fg | 188 |
| PH695583 | Whyndham Powell Urethral Bougies (Sounds) Individual Size 19 Fg | 188 |
| PH695584 | Whyndham Powell Urethral Bougies (Sounds) Individual Size 20 Fg | 188 |
| PH695585 | Whyndham Powell Urethral Bougies (Sounds) Individual Size 21 Fg | 188 |

Product Codes Index

| Code | Description | Page |
|----------|---|------|
| PH695586 | Whyndham Powell Urethral Bougies (Sounds) Individual Size 22 Fg | 188 |
| PH695587 | Whyndham Powell Urethral Bougies (Sounds) Individual Size 23 Fg | 188 |
| PH695588 | Whyndham Powell Urethral Bougies (Sounds) Individual Size 24 Fg | 188 |
| PH695589 | Whyndham Powell Urethral Bougies (Sounds) Individual Size 25 Fg | 188 |
| PH695590 | Whyndham Powell Urethral Bougies (Sounds) Individual Size 26 Fg | 188 |
| PH695591 | Whyndham Powell Urethral Bougies (Sounds) Individual Size 27 Fg | 188 |
| PH695592 | Whyndham Powell Urethral Bougies (Sounds) Individual Size 28 Fg | 188 |
| PH695593 | Whyndham Powell Urethral Bougies (Sounds) Individual Size 29 Fg | 188 |
| PH695594 | Whyndham Powell Urethral Bougies (Sounds) Individual Size 30 Fg | 188 |
| PH695595 | Whyndham Powell Urethral Bougies (Sounds) Individual Size 31 Fg | 188 |
| PH695596 | Whyndham Powell Urethral Bougies (Sounds) Individual Size 32 Fg | 188 |
| PH695597 | Whyndham Powell Urethral Bougies (Sounds) Individual Size 36 Fg | 188 |
| PH695598 | Whyndham Powell Urethral Bougies (Sounds) Individual Size 38 Fg | 188 |
| PH695599 | Whyndham Powell Urethral Bougies (Sounds) Individual Size 40 Fg | 188 |
| PH70200 | Dressing And Stitch Scissors Sharp / Sharp 100mm Straight | 12 |
| PH70201 | Dressing And Stitch Scissors Sharp / Sharp 115mm Straight | 12 |
| PH70202 | Dressing And Stitch Scissors Sharp / Sharp 130mm Straight | 12 |
| PH70203 | Dressing And Stitch Scissors Sharp / Sharp 140mm Straight | 12 |
| PH70204 | Dressing And Stitch Scissors Sharp / Sharp 150mm Straight | 12 |
| PH70205 | Dressing And Stitch Scissors Sharp / Sharp 180mm Straight | 12 |
| PH70208 | Dressing And Stitch Scissors Sharp / Sharp 100mm Curved | 12 |
| PH70209 | Dressing And Stitch Scissors Sharp / Sharp 115mm Curved | 12 |
| PH70210 | Dressing And Stitch Scissors Sharp / Sharp 130mm Curved | 12 |
| PH70211 | Dressing And Stitch Scissors Sharp / Sharp 140mm Curved | 12 |
| PH70212 | Dressing And Stitch Scissors Sharp / Sharp 150mm Curved | 12 |
| PH70213 | Dressing And Stitch Scissors Sharp / Sharp 180mm Curved | 12 |
| PH70216 | Dressing Scissors Sharp / Blunt 100mm Straight | 12 |
| PH70217 | Dressing Scissors Sharp / Blunt 115mm Straight | 12 |
| PH70218 | Dressing Scissors Sharp / Blunt 130mm Straight | 12 |
| PH70219 | Dressing Scissors Sharp / Blunt 140mm Straight | 12 |
| PH70220 | Dressing Scissors Sharp / Blunt 150mm Straight | 12 |
| PH70221 | Dressing Scissors Sharp / Blunt 180mm Straight | 12 |
| PH70224 | Dressing Scissors Sharp / Blunt 100mm Curved | 12 |

Product Codes Index

| Code | Description | Page |
|---------|---|------|
| PH70225 | Dressing Scissors Sharp / Blunt 115mm Curved | 12 |
| PH70226 | Dressing Scissors Sharp / Blunt 130mm Curved | 12 |
| PH70227 | Dressing Scissors Sharp / Blunt 140mm Curved | 12 |
| PH70228 | Dressing Scissors Sharp / Blunt 150mm Curved | 12 |
| PH70229 | Dressing Scissors Sharp / Blunt 180mm Curved | 12 |
| PH70232 | Dressing Scissors Blunt / Blunt 100mm Straight | 12 |
| PH70233 | Dressing Scissors Blunt / Blunt 115mm Straight | 12 |
| PH70234 | Dressing Scissors Blunt / Blunt 130mm Straight | 12 |
| PH70235 | Dressing Scissors Blunt / Blunt 140mm Straight | 12 |
| PH70236 | Dressing Scissors Blunt / Blunt 150mm Straight | 12 |
| PH70237 | Dressing Scissors Blunt / Blunt 180mm Straight | 12 |
| PH70240 | Dressing Scissors Blunt / Blunt 100mm Curved | 12 |
| PH70241 | Dressing Scissors Blunt / Blunt 115mm Curved | 12 |
| PH70242 | Dressing Scissors Blunt / Blunt 130mm Curved | 12 |
| PH70243 | Dressing Scissors Blunt / Blunt 140mm Curved | 12 |
| PH70244 | Dressing Scissors Blunt / Blunt 150mm Curved | 12 |
| PH70245 | Dressing Scissors Blunt / Blunt 180mm Curved | 12 |
| PH70260 | Dissecting Scissors With Open Shanks 115mm Straight | 13 |
| PH70261 | Dissecting Scissors With Open Shanks 130mm Straight | 13 |
| PH70262 | Dissecting Scissors With Open Shanks 150mm Straight | 13 |
| PH70275 | Aufricht Scissors 140mm Straight | 13 |
| PH70276 | Aufricht Scissors 150mm Curved | 13 |
| PH70281 | Iris Scissors Fine Pointed 100mm Straight | 13 |
| PH70282 | Iris Scissors Fine Pointed 115mm Straight | 13 |
| PH70283 | Iris Scissors Fine Pointed 130mm Straight | 13 |
| PH70284 | Iris Scissors Fine Pointed 100mm Curved | 13 |
| PH70285 | Iris Scissors Fine Pointed 115mm Curved | 13 |
| PH70286 | Iris Scissors Fine Pointed 130mm Curved | 13 |
| PH70289 | Lister Bandage Scissors 115mm Angled | 13 |
| PH70290 | Lister Bandage Scissors 140mm Angled | 13 |
| PH70291 | Lister Bandage Scissors 185mm Angled | 13 |
| PH70292 | Lister Bandage Scissors 200mm Angled | 13 |
| PH70293 | Lister Bandage Scissors 230mm Angled | 13 |

Product Codes Index

| Code | Description | Page |
|---------|---|------|
| PH70300 | Mayo Scissors 140mm Straight | 14 |
| PH70301 | Mayo Scissors 165mm Straight | 14 |
| PH70302 | Mayo Scissors 190mm Straight | 14 |
| PH70303 | Mayo Scissors 215mm Straight | 14 |
| PH70304 | Mayo Scissors 230mm Straight | 14 |
| PH70305 | Mayo Scissors 140mm Curved | 14 |
| PH70306 | Mayo Scissors 165mm Curved | 14 |
| PH70307 | Mayo Scissors 190mm Curved | 14 |
| PH70308 | Mayo Scissors 215mm Curved | 14 |
| PH70309 | Mayo Scissors 230mm Curved | 14 |
| PH70350 | Metzenbaum Scissors 140mm Straight | 14 |
| PH70351 | Metzenbaum Scissors 165mm Straight | 14 |
| PH70352 | Metzenbaum Scissors 190mm Straight | 14 |
| PH70353 | Metzenbaum Scissors 215mm Straight | 14 |
| PH70354 | Metzenbaum Scissors 240mm Straight | 14 |
| PH70360 | Metzenbaum Scissors 140mm Curved | 14 |
| PH70361 | Metzenbaum Scissors 165mm Curved | 14 |
| PH70362 | Metzenbaum Scissors 190mm Curved | 14 |
| PH70363 | Metzenbaum Scissors 215mm Curved | 14 |
| PH70364 | Metzenbaum Scissors 240mm Curved | 14 |
| PH70370 | Mcindoe Scissors 115mm Curved | 14 |
| PH70371 | Mcindoe Scissors 190mm Curved | 14 |
| PH70380 | Kilner Scissors 115mm Straight | 14 |
| PH70381 | Kilner Scissors 115mm Curved | 15 |
| PH70382 | Strabismus Scissors 115mm Straight | 15 |
| PH70383 | Strabismus Scissors 115mm Curved | 15 |
| PH70402 | Universal Scissors For Ligature And Wire Cutting 115mm Angled | 15 |
| PH70420 | Nelson Scissors 230mm Straight (Thoracic) | 15 |
| PH70421 | Nelson Scissors 280mm Straight (Thoracic) | 15 |
| PH70422 | Nelson Scissors 230mm Curved (Thoracic) | 15 |
| PH70430 | Wilson Tonsil Scissors 200mm Straight | 16 |
| PH70431 | Wilson Tonsil Scissors 200mm Curved | 16 |
| PH70454 | Carless Light Weight Scissors 165mm Straight | 16 |

Product Codes Index

| Code | Description | Page |
|----------|--|------|
| PH70455 | Carless Light Weight Scissors 215mm Straight | 16 |
| PH70457 | Carless Light Weight Scissors 165mm Curved | 16 |
| PH70458 | Carless Light Weight Scissors 215mm Curved | 16 |
| PH70465 | Pott Smith Scissors Angled To Side At 25 Deg 190mm | 16 |
| PH70466 | Pott Smith Scissors Angled To Side At 45 Deg 190mm Angled | 16 |
| PH70467 | Pott Smith Scissors Angled On Flat 190mm Angled | 17 |
| PH70470 | Lawson Tait Scissors (Barnes) (Episiotomy) 130mm Angled | 17 |
| PH70474 | Braun Stadler Scissors 150mm (Episiotomy) Angled | 17 |
| PH70475 | Braun Stadler Scissors 190mm (Episiotomy) Angled | 17 |
| PH70476 | Braun Stadler Scissors 200mm (Episiotomy) Angled | 17 |
| PH70478 | Doyen Scissors (Uterine) 170mm Straight | 17 |
| PH70479 | Doyen Scissors (Uterine) 170mm Curved | 17 |
| PH70480 | Doyen Scissors (Uterine) 200mm Straight | 17 |
| PH70481 | Doyen Scissors (Uterine) 200mm Curved | 17 |
| PH70483 | Fergusson Scissors (Uterine) 210mm Straight | 17 |
| PH70484 | Fergusson Scissors (Uterine) Angled To Flat 210mm Angled | 18 |
| PH70488 | Sims Scissors (Uterine) 200mm Straight | 18 |
| PH70489 | Sims Scissors (Uterine) 230mm Straight | 18 |
| PH70490 | Sims Scissors (Uterine) 200mm Curved | 18 |
| PH70491 | Sims Scissors (Uterine) 230mm Curved | 18 |
| PH70492 | Sims Scissors (Uterine) Angled To One Side 200mm | 18 |
| PH70493 | Sims Scissors (Uterine) Angled To One Side 230mm | 18 |
| PH70500 | American Pattern Umbilical Cord Cutter 115mm Straight | 18 |
| PH70520 | Universal Tough Cut Scissors With Plastic Handles 140mm Angled (Tubing Scissors) | 19 |
| PH70521 | Universal Tough Cut Scissors With Plastic Handles 190mm Angled (Tubing Scissors) | 19 |
| PH70522 | Capsular Scissors Angled On Flat 250mm Angled | 19 |
| PH70530 | American Pattern Fistula Scissors 180mm Curved | 19 |
| PH70538 | Mckissock Dural Scissor Straight 150mm | 19 |
| PH70539 | Mckissock Dural Scissors Curved 150mm | 19 |
| PH705592 | Foley Catheter Introducer Curved With A 2mm Diameter Tip 430mm | 19 |
| PH705593 | Foley Catheter Introducer Curved With A 4mm Diameter Tip 430mm | 19 |
| PH705594 | Guyon Catheter Introducer Straight With A 2mm Diameter Tip 430mm | 19 |
| PH705595 | Guyon Catheter Introducer Straight With A 4mm Diameter Tip 430mm | 19 |

Product Codes Index

| Code | Description | Page |
|----------|--|------|
| PH705596 | Guyon Catheter Introducer Curved With A 2mm Diameter Tip 430mm | 19 |
| PH705597 | Guyon Catheter Introducer Curved With A 4mm Diameter Tip 430mm | 19 |
| PH70570 | Lahey Scissors Curved 140mm | 19 |
| PH70600 | Dressing / Stitch Scissors With Tungsten Carbide Jaws Sharp / Sharp 150mm Straight | 20 |
| PH70601 | Dressing / Stitch Scissors With Tungsten Carbide Jaws Sharp / Sharp 180mm Straight | 20 |
| PH70605 | Dressing / Stitch Scissors With Tungsten Carbide Jaws Sharp / Sharp 150mm Curved | 20 |
| PH70606 | Dressing / Stitch Scissors With Tungsten Carbide Jaws Sharp / Sharp 180mm Curved | 20 |
| PH70610 | Dressing / Stitch Scissors With Tungsten Carbide Jaws Sharp / Blunt 150mm Straight | 20 |
| PH70612 | Dressing / Stitch Scissors With Tungsten Carbide Jaws Sharp / Blunt 180mm Straight | 20 |
| PH70615 | Dressing / Stitch Scissors With Tungsten Carbide Jaws Sharp / Blunt 150mm Curved | 20 |
| PH70616 | Dressing / Stitch Scissors With Tungsten Carbide Jaws Sharp / Blunt 180mm Curved | 20 |
| PH70620 | Dressing / Stitch Scissors With Tungsten Carbide Jaws Blunt / Blunt 150mm Straight | 20 |
| PH70621 | Dressing / Stitch Scissors With Tungsten Carbide Jaws Blunt / Blunt 180mm Straight | 20 |
| PH70625 | Iris Scissors With Tungsten Carbide Jaws Fine Pointed 115mm Straight | 20 |
| PH70626 | Iris Scissors With Tungsten Carbide Jaws Fine Pointed 115mm Curved | 21 |
| PH70630 | Kilner Scissors With Tungsten Carbide Jaws 115mm Straight | 21 |
| PH70631 | Kilner Scissors With Tungsten Carbide Jaws 115mm Curved | 21 |
| PH70635 | Strabismus Scissors With Tungsten Carbide Jaws 115mm Straight | 21 |
| PH70636 | Strabismus Scissors With Tungsten Carbide Jaws 115mm Curved | 21 |
| PH70640 | Mayo Scissors With Tungsten Carbide Jaws 140mm Straight | 21 |
| PH70641 | Mayo Scissors With Tungsten Carbide Jaws 150mm Straight | 21 |
| PH70642 | Mayo Scissors With Tungsten Carbide Jaws 180mm Straight | 21 |
| PH70643 | Mayo Scissors With Tungsten Carbide Jaws 200mm Straight | 21 |
| PH70644 | Mayo Scissors With Tungsten Carbide Jaws 230mm Straight | 21 |
| PH70650 | Mayo Scissors With Tungsten Carbide Jaws 140mm Curved | 22 |
| PH70651 | Mayo Scissors With Tungsten Carbide Jaws 150mm Curved | 22 |
| PH70652 | Mayo Scissors With Tungsten Carbide Jaws 180mm Curved | 22 |
| PH70653 | Mayo Scissors With Tungsten Carbide Jaws 200mm Curved | 22 |
| PH70654 | Mayo Scissors With Tungsten Carbide Jaws 230mm Curved | 22 |
| PH70655 | Lister Bandage Scissors With Tungsten Carbide Jaws 115mm Angled | 22 |
| PH70656 | Lister Bandage Scissors With Tungsten Carbide Jaws 140mm Angled | 22 |
| PH70657 | Lister Bandage Scissors With Tungsten Carbide Jaws 185mm Angled | 22 |
| PH70658 | Lister Bandage Scissors With Tungsten Carbide Jaws 200mm Angled | 22 |

Product Codes Index

| Code | Description | Page |
|---------|---|------|
| PH70659 | Lister Bandage Scissors With Tungsten Carbide Jaws 230mm Angled | 22 |
| PH70660 | Metzenbaum Scissors With Tungsten Carbide Jaws 140mm Straight | 22 |
| PH70661 | Metzenbaum Scissors With Tungsten Carbide Jaws 180mm Straight | 22 |
| PH70662 | Metzenbaum Scissors With Tungsten Carbide Jaws 200mm Straight | 22 |
| PH70663 | Metzenbaum Scissors With Tungsten Carbide Jaws 230mm Straight | 22 |
| PH70670 | Metzenbaum Scissors With Tungsten Carbide Jaws 140mm Curved | 22 |
| PH70671 | Metzenbaum Scissors With Tungsten Carbide Jaws 180mm Curved | 22 |
| PH70672 | Metzenbaum Scissors With Tungsten Carbide Jaws 200mm Curved | 22 |
| PH70673 | Metzenbaum Scissors With Tungsten Carbide Jaws 230mm Curved | 22 |
| PH70679 | Dental Plier Cutter | 22 |
| PH70680 | Universal Scissors With Tungsten Carbide Jaws For Ligature And Wire Cutting 115m | 22 |
| PH70685 | Small Wire Cutters Tungsten Carbide Jaws 125mm Straight | 28 |
| PH70686 | Top Cutters Tc Jaws 150mm | 28 |
| PH70690 | Compound Action Wire Cutter With Tungsten Carbide Jaws 7inch | 29 |
| PH70691 | Wire Cutter Tc 200mm | 29 |
| PH70692 | Compound Action Top Wire Cutter Tungsten Carbide Jaws 200mm Straight | 29 |
| PH70693 | Compound Action Side Wire Cutter Tungsten Carbide Jaws 230mm Straight | 30 |
| PH70694 | Compound Action Large Wire Cutter To Cut 5.5mm Diameter Max With Tungsten Carbide Jaw | 30 |
| PH70695 | Orthopaedic Plier Cutter With Serrated Jaw Hammer End And Side Cutter 180mm | 30 |
| PH70696 | Power Cutter Tc Jaws | 31 |
| PH70698 | Small Pliers | 31 |
| PH70700 | Large Pliers | 31 |
| PH70701 | Iris Supercut Scissors Fine Pointed 115mm Straight | 23 |
| PH70702 | Iris Supercut Scissors Fine Pointed 115mm Curved | 23 |
| PH70705 | Kilner Supercut Scissors 115mm Straight | 23 |
| PH70706 | Kilner Supercut Scissors 115mm Curved | 23 |
| PH70710 | Stevens Supercut Tenotomy Scissors 115mm Straight | 23 |
| PH70711 | Stevens Supercut Tenotomy Scissors 115mm Curved | 23 |
| PH70713 | Stevens Supercut Tenotomy Scissors 135mm Curved | 23 |
| PH70715 | Mayo Supercut Scissors With Black Handles 140mm Straight | 24 |
| PH70716 | Mayo Supercut Scissors With Black Handles 150mm Straight | 24 |
| PH70717 | Mayo Supercut Scissors With Black Handles 180mm Straight | 24 |
| PH70718 | Mayo Supercut Scissors With Black Handles 200mm Straight | 24 |

Product Codes Index

| Code | Description | Page |
|----------|--|------|
| PH70720 | Mayo Supercut Scissors With Black Handles 140mm Curved | 24 |
| PH70721 | Mayo Supercut Scissors With Black Handles 150mm Curved | 24 |
| PH70722 | Mayo Supercut Scissors With Black Handles 180mm Curved | 24 |
| PH70723 | Mayo Supercut Scissors With Black Handles 200mm Curved | 24 |
| PH70730 | Metzenbaum Supercut Scissors With Black Handles 150mm Straight | 24 |
| PH70731 | Metzenbaum Supercut Scissors With Black Handles 180mm Straight | 24 |
| PH70732 | Metzenbaum Supercut Scissors With Black Handles 200mm Straight | 24 |
| PH70735 | Metzenbaum Supercut Scissors With Black Handles 150mm Curved | 24 |
| PH70736 | Metzenbaum Supercut Scissors With Black Handles 180mm Curved | 24 |
| PH70737 | Metzenbaum Supercut Scissors With Black Handles 200mm Curved | 24 |
| PH70738 | Jameson Reynolds Supercut Scissors With Black Handles 160mm Curved | 24 |
| PH70739 | Mcindoe Supercut Scissors With Black Handles 180mm Curved | 24 |
| PH70745 | Strabismus Supercut Scissors Straight | 25 |
| PH70746 | Strabismus Supercut Scissors Curved | 25 |
| PH70755 | Gorney Freeman Facelift Scissors Tc Supercut | 25 |
| PH70756 | Gorney Freeman Facelift Scissors Tc Supercut | 25 |
| PH70800 | Mcindoe Scissors Tc Jaws 140mm Straight | 25 |
| PH70801 | Mcindoe Scissors Tc Jaws 180mm Straight | 25 |
| PH70802 | Mcindoe Scissors Tc Jaws 200mm Straight | 25 |
| PH70803 | Mcindoe Scissors Tc Jaws 230mm Straight | 25 |
| PH70804 | Mcindoe Scissors Tc Jaws 140mm Curved | 25 |
| PH70806 | Mcindoe Scissors Tc Jaws 200mm Curved | 25 |
| PH70807 | Mcindoe Scissors Tc Jaws 230mm Curved | 25 |
| PH70808 | Mcindoe Scissors Tc Jaws 165mm Curved | 25 |
| PH70809 | Mcindoe Scissors Tc Jaws 165mm Straight | 25 |
| PH70810 | Mcindoe Scissors Tc Jaws 280mm Curved | 25 |
| PH70811 | Mcindoe Scissors Tc Jaws 190mm Curved | 25 |
| PH715725 | Umbilical Dilator 95mm | 190 |
| PH715746 | Desjardin Gall Stone Forceps Small Jaws With A Screw Joint 230mm | 190 |
| PH715747 | Desjardin Gall Stone Forceps Large Jaws With A Screw Joint 230mm | 190 |
| PH725903 | Teale Uterine Vulsellum Forceps Curved With 3 Into 4 Teeth And A Box Joint 180mm | 72 |
| PH725904 | Teale Uterine Vulsellum Forceps Curved With 3 Into 4 Teeth And A Box Joint 230mm | 72 |

Product Codes Index

| Code | Description | Page |
|----------|--|------|
| PH725906 | Schroeder Uterine Vulsellum Forceps With 2 Into 2 Teeth And A Box Joint 240mm Straight | 72 |
| PH725911 | Luer Uterine Vulsellum Forceps With 1 Into 1 Teeth And A Box Joint 240mm Straight | 72 |
| PH725925 | Pratt Vulsellum Forceps With Teeth And A Screw Joint 280mm Straight | 72 |
| PH725932 | Lewer Vulsellum Forceps With 2 Into 2 Teeth And A Box Joint 240mm Curved | 72 |
| PH725934 | Rotunda Hosp Pattern Vulsellum Fcps 1 Into 1 Teeth | 72 |
| PH725964 | Schroeder Tenaculum Forceps With 1 Into 1 Teeth And A Box Joint 250mm Straight | 73 |
| PH725984 | Green Armytage Uterine Dressing Forceps With Wide Serrated Jaws And A Screw Join | 73 |
| PH725986 | Green Armytage Uterine Dressing Forceps Angled Ends With Wide Serrated Jaws | 73 |
| PH726000 | Bonney Uterine Polypus Forceps With Fenestrated Serrated Jaws And A Box Joint | 73 |
| PH726075 | Parametrium Clamps Straight With Box Joint 210mm | 192 |
| PH726076 | Parametrium Clamps Straight With Box Joint 240mm | 192 |
| PH726077 | Parametrium Clamps Straight With Box Joint 300mm | 192 |
| PH726078 | Parametrium Clamps Straight With Box Joint 350mm | 192 |
| PH726079 | Parametrium Clamps Slightly Curved With Box Joint 210mm | 192 |
| PH726080 | Parametrium Clamps Slightly Curved With Box Joint 240mm | 192 |
| PH726081 | Parametrium Clamps Slightly Curved With Box Joint 300mm | 192 |
| PH726082 | Parametrium Clamps Slightly Curved With Box Joint 350mm | 192 |
| PH726083 | Parametrium Clamps Strongly Curved With Box Joint 210mm | 192 |
| PH726084 | Parametrium Clamps Strongly Curved With Box Joint 240mm | 192 |
| PH726085 | Parametrium Clamps Strongly Curved With Box Joint 300mm | 192 |
| PH726086 | Parametrium Clamps Strongly Curved With Box Joint 350mm | 192 |
| PH726087 | Parametrium Clamps Angulated With Box Joint 210mm | 192 |
| PH726088 | Parametrium Clamps Angulated With Box Joint 240mm | 192 |
| PH726089 | Parametrium Clamps Angulated With Box Joint 300mm | 192 |
| PH726090 | Parametrium Clamps Angulated With Box Joint 350mm | 192 |
| PH726118 | Mcclintock Ovum Forceps Curved With A Screw Joint 240mm Curved | 192 |
| PH726152 | Gwilliam Hysterectomy Clamp Straight Longitudinal Serrated Jaws With 1 Into 2 Teeth | 192 |
| PH726153 | Gwilliam Hysterectomy Clamp Curved Longitudinal Serrated Jaws With 1 Into 2 Teeth | 193 |
| PH726162 | Mccullagh Hysterectomy Clamp Mccullogh 250mm Straight | 193 |
| PH726163 | Mccullagh Hysterectomy Clamp Mccullogh 250mm Curved | 193 |
| PH726173 | Rogers Hysterectomy Clamp Straight | 193 |
| PH726174 | Rogers Hysterectomy Clamp Curved | 193 |

Product Codes Index

| Code | Description | Page |
|----------|--|------|
| PH726180 | Doyen Myoma Screw With Small Spiral For Use Per Vagina With D Ring Handle 175mm | 193 |
| PH726182 | Doyen Myoma Screw With Large Spiral For Use Per Abdomen With D Ring Handle 150mm | 193 |
| PH726222 | Sims Uterine Curette Small Double Ended Sharp / Sharp 255mm | 194 |
| PH726223 | Sims Uterine Curette Medium Double Ended Sharp / Sharp 255mm | 194 |
| PH726224 | Sims Uterine Curette Large Double Ended Sharp / Sharp 255mm | 194 |
| PH726226 | Sims Uterine Curette Small Double Ended Sharp / Blunt 255mm | 194 |
| PH726227 | Sims Uterine Curette Medium Double Ended Sharp / Blunt 255mm | 194 |
| PH726228 | Sims Uterine Curette Large Double Ended Sharp / Blunt 255mm | 194 |
| PH726248 | Sims Uterine Curette 3mm Single Ended Sharp 255mm | 194 |
| PH726249 | Sims Uterine Curette 5mm Single Ended Sharp 255mm | 194 |
| PH726250 | Sims Uterine Curette 6mm Single Ended Sharp 255mm | 194 |
| PH726251 | Sims Uterine Curette 7mm Single Ended Sharp 255mm | 194 |
| PH726252 | Sims Uterine Curette 9mm Single Ended Sharp 255mm | 194 |
| PH726253 | Sims Uterine Curette 12mm Single Ended Sharp 255mm | 194 |
| PH726254 | Sims Uterine Curette 15mm Single Ended Sharp 255mm | 194 |
| PH726255 | Sims Uterine Curette 19mm Single Ended Sharp 255mm | 194 |
| PH726335 | Sims Uterine Sound Malleable Cm Graduated | 194 |
| PH726337 | Galabin Uterine Sound Malleable Cm Graduated | 194 |
| PH726370 | Cusco Virgin 86mm x 16mm | 195 |
| PH726371 | Cusco Small 95mm x 25mm | 195 |
| PH726372 | Cusco Medium 102mm x 30mm | 195 |
| PH726373 | Cusco Medium Long 115mm x 30mm | 195 |
| PH726374 | Cusco Large 110mm x 35mm | 195 |
| PH726375 | Cusco Large Extra Long 115 x 35mm | 195 |
| PH726376 | Cusco Virgin 86 x 16mm Smoke Tube Bead | 195 |
| PH726377 | Cusco Small 95 x 25mm Smoke Tube Bead | 195 |
| PH726378 | Cusco Medium 102 x 30mm Smoke Tube Bead | 195 |
| PH726379 | Cusco Medium Long 115 x 30mm Smoke Tube Bead | 195 |
| PH726380 | Cusco Large 110 x 35mm Smoke Tube Bead | 195 |
| PH726381 | Cusco Large Ex Long 115 x 35mm Smoke Tube Bead | 195 |
| PH726382 | Cusco Virgin 86 x 16mm Smoke Tube Black | 196 |
| PH726383 | Cusco Small 95 x 25mm Smoke Tube Black | 196 |

Product Codes Index

| Code | Description | Page |
|----------|---|------|
| PH726384 | Cusco Medium 102 x 30mm Smoke Tube Black | 196 |
| PH726385 | Cusco Medium Long 115 x 30mm Smoke Tube Black | 196 |
| PH726386 | Cusco Large 110 x 35mm Smoke Tube Black | 196 |
| PH726387 | Cusco Large Ex Long 115 x 35mm Smoke Tube Black | 196 |
| PH726388 | Cusco Virgin 86 x 16mm Smoke Tube Insulated | 196 |
| PH726389 | Cusco Small 95 x 25mm Smoke Tube Insulated | 196 |
| PH726390 | Cusco Medium 102 x 30mm Smoke Tube Insulated | 196 |
| PH726391 | Cusco Medium Long 115 x 30mm Smoke Tube Insulated | 196 |
| PH726392 | Cusco Large 110 x 35mm Smoke Tube Insulated | 196 |
| PH726393 | Cusco Lge Ex Long 115 x 35mm Smoke Tube Insulated | 196 |
| PH726395 | Graves (Brewer) Vaginal Speculum Small 75mm x 20mm | 197 |
| PH726396 | Graves (Brewer) Vaginal Speculum Medium 95mm x 35mm | 197 |
| PH726397 | Graves (Brewer) Vaginal Speculum Large 115mm x 32mm | 197 |
| PH726400 | Kogan Endocervical Speculum Without Ratchet (Koogan) 240mm | 197 |
| PH726401 | Koogan Endocervical Speculum & Adjustable Ratchet | 197 |
| PH726420 | Sims Vaginal Speculum Small Double Ended 60mm x 30mm / 70mm x 30mm 130mm | 198 |
| PH726421 | Sims Vaginal Speculum Medium Double Ended 70mm x 30mm / 75mm x 30mm 155mm | 198 |
| PH726422 | Sims Vaginal Speculum Large Double Ended 75mm x 30mm / 85mm x 30mm 185mm | 198 |
| PH726423 | Sims Vaginal Speculum Small Single Ended 76mm x 21mm 180mm | 198 |
| PH726424 | Sims Vaginal Speculum Medium Single End 86mm x 29mm 180mm | 198 |
| PH726425 | Sims Vaginal Speculum Large Single End 115mm x 43mm 180mm | 198 |
| PH726439 | Auvard Speculum 80 x 45 Tapered 40mm Fixed Weight | 198 |
| PH726442 | Bailey Vaginal Speculum 38mm Blade & Fixed Weight | 198 |
| PH726445 | Winterton Vaginal Speculum | 198 |
| PH726446 | Winterton Vaginal Speculum With Smoke Tube And Bead Finish | 199 |
| PH726447 | Winterton Vaginal Speculum With Smoke Tube And Black Finish | 199 |
| PH726448 | Winterton Vaginal Speculum Diameterthermy Insulated With Smoke Tube | 199 |
| PH726450 | Doyen Vaginal Speculum Small 55mm x 45mm 240mm | 200 |
| PH726451 | Doyen Vaginal Speculum Medium 85mm x 45mm 240mm | 200 |
| PH726452 | Doyen Vaginal Speculum Large 115mm x 45mm 240mm | 200 |
| PH726455 | Jessop Hospital Pattern Vaginal Speculum Double Ended 100mm x 28mm And 120mm x 30mm | 200 |
| PH726465 | St Georges Pattern Vaginal Speculum Double Ended 100mm x 25mm And 120mm x 30mm | 200 |

Product Codes Index

| Code | Description | Page |
|----------|--|------|
| PH726480 | Breisky Vaginal Speculum With Cranked Head Blade 100mm Length x 20mm Width 300mm | 200 |
| PH726481 | Breisky Vaginal Speculum With Cranked Head Blade 100mm Length x 25mm Width 300mm | 200 |
| PH726482 | Breisky Vaginal Speculum With Cranked Head Blade 100mm Length x 30mm Width 300mm | 200 |
| PH726483 | Breisky Vaginal Speculum With Cranked Head Blade 100mm Length x 35mm Width 300mm | 200 |
| PH726484 | Breisky Vaginal Speculum With Cranked Head Blade 100mm Length x 40mm Width 300mm | 200 |
| PH726485 | Breisky Vaginal Speculum With Cranked Head Blade 130mm Length x 20mm Width 320mm | 200 |
| PH726486 | Breisky Vaginal Speculum With Cranked Head Blade 130mm Length x 25mm Width 320mm | 200 |
| PH726487 | Breisky Vaginal Speculum With Cranked Head Blade 130mm Length x 30mm Width 320mm | 200 |
| PH726488 | Breisky Vaginal Speculum With Cranked Head Blade 130mm Length x 35mm Width 320mm | 200 |
| PH726489 | Breisky Vaginal Speculum With Cranked Head Blade 130mm Length x 40mm Width 320mm | 200 |
| PH726490 | Breisky Vaginal Speculum With Cranked Head Blade 160mm Length x 40mm Width 350mm | 200 |
| PH726491 | Breisky Vaginal Speculum With Cranked Head Blade 180mm Length x 40mm Width 350mm | 200 |
| PH726503 | Hegar Uterine Dilator 1.5mm Diameter Single Ended 185mm | 190 |
| PH726504 | Hegar Uterine Dilator S/E 185mm 2.5mm Diameter | 190 |
| PH726505 | Hegar Uterine Dilator 3.5mm Diameter Single Ended 185mm | 190 |
| PH726506 | Hegar Uterine Dilator 4.5mm Diameter Single Ended 185mm | 190 |
| PH726507 | Hegar Uterine Dilator 5.5mm Diameter Single Ended 185mm | 190 |
| PH726508 | Hegar Uterine Dilator 6.5mm Diameter Single Ended 185mm | 190 |
| PH726520 | Hegar Uterine Dilator 1mm Diameter Single Ended 185mm | 190 |
| PH726521 | Hegar Uterine Dilator 2mm Diameter Single Ended 185mm | 190 |
| PH726522 | Hegar Uterine Dilator 3mm Diameter Single Ended 185mm | 190 |
| PH726523 | Hegar Uterine Dilator 4mm Diameter Single Ended 185mm | 190 |
| PH726524 | Hegar Uterine Dilator 5mm Diameter Single Ended 185mm | 190 |
| PH726530 | Hegar Uterine Dilator 6mm Diameter Single Ended 185mm | 190 |
| PH726531 | Hegar Uterine Dilator 7mm Diameter Single Ended 185mm | 190 |
| PH726532 | Hegar Uterine Dilator 8mm Diameter Single Ended 185mm | 190 |
| PH726533 | Hegar Uterine Dilator 9mm Diameter Single Ended 185mm | 190 |
| PH726534 | Hegar Uterine Dilator 10mm Diameter Single Ended 185mm | 190 |

Product Codes Index

| Code | Description | Page |
|----------|---|------|
| PH726535 | Hegar Uterine Dilator 7.5mm Diameter Single Ended 185mm | 190 |
| PH726536 | Hegar Uterine Dilator 8.5mm Diameter Single Ended 185mm | 190 |
| PH726537 | Hegar Uterine Dilator 9.5mm Diameter Single Ended 185mm | 190 |
| PH726538 | Hegar Uterine Dilator 10.5mm Diameter Single Ended 185mm | 190 |
| PH726540 | Hegar Uterine Dilator 11mm Diameter Single Ended 185mm | 190 |
| PH726541 | Hegar Uterine Dilator 12mm Diameter Single Ended 185mm | 190 |
| PH726542 | Hegar Uterine Dilator 13mm Diameter Single Ended 185mm | 190 |
| PH726543 | Hegar Uterine Dilator 14mm Diameter Single Ended 185mm | 190 |
| PH726545 | Hegar Uterine Dilator 11.5mm Diameter Single Ended 185mm | 190 |
| PH726546 | Hegar Uterine Dilator 12.5mm Diameter Single Ended 185mm | 190 |
| PH726547 | Hegar Uterine Dilator 13.5mm Diameter Single Ended 185mm | 190 |
| PH726548 | Hegar Uterine Dilator 14.5mm Diameter Single Ended 185mm | 190 |
| PH726549 | Hegar Uterine Dilator S/E 185mm 15.5mm Diameter | 190 |
| PH726550 | Hegar Uterine Dilator 15mm Diameter Single Ended 185mm | 190 |
| PH726551 | Hegar Uterine Dilator 16mm Diameter Single Ended 185mm | 190 |
| PH726552 | Hegar Uterine Dilator 17mm Diameter Single Ended 185mm | 190 |
| PH726553 | Hegar Uterine Dilator 18mm Diameter Single Ended 185mm | 190 |
| PH726557 | Hegar Uterine Dilator S/E 185mm 16.5mm Diameter | 190 |
| PH726558 | Hegar Uterine Dilator S/E 185mm 17.5mm Diameter | 190 |
| PH726560 | Hegar Uterine Dilator 19mm Diameter Single Ended 185mm | 190 |
| PH726561 | Hegar Uterine Dilator 20mm Diameter Single Ended 185mm | 190 |
| PH726562 | Hegar Uterine Dilator 21mm Diameter Single Ended 185mm | 190 |
| PH726563 | Hegar Uterine Dilator 22mm Diameter Single Ended 185mm | 190 |
| PH726570 | Hegar Uterine Dilator 23mm Diameter Single Ended 185mm | 190 |
| PH726571 | Hegar Uterine Dilator 24mm Diameter Single Ended 185mm | 190 |
| PH726572 | Hegar Uterine Dilator 25mm Diameter Single Ended 185mm | 190 |
| PH726573 | Hegar Uterine Dilator 26mm Diameter Single Ended 185mm | 190 |
| PH726574 | Hegar Uterine Dilator S/E 185mm 18.5mm Diameter | 190 |
| PH726575 | Hegar Uterine Dilator S/E 185mm 19.5mm Diameter | 190 |
| PH726580 | Hegar Uterine Dilators. Set Of 8 Double Ended Dilators | 191 |
| PH726590 | Hegar Uterine Dilators, Set Of 12 Double Ended Dilators | 191 |
| PH726600 | Hegar Uterine Dilator Double Ended Sizes 1mm And 2mm Diameter 220mm | 191 |
| PH726601 | Hegar Uterine Dilator Double Ended Sizes 3mm And 4mm Diameter 220mm | 191 |

Product Codes Index

| Code | Description | Page |
|----------|--|------|
| PH726602 | Hegar Uterine Dilator Double Ended Sizes 5mm And 6mm Diameter 220mm | 191 |
| PH726603 | Hegar Uterine Dilator Double Ended Sizes 7mm And 8mm Diameter 220mm | 191 |
| PH726605 | Hegar Uterine Dilator Double Ended Sizes 9mm And 10mm Diameter 220mm | 191 |
| PH726606 | Hegar Uterine Dilator Double Ended Sizes 11mm And 12mm Diameter 220mm | 191 |
| PH726607 | Hegar Uterine Dilator Double Ended Sizes 13mm And 14mm Diameter 220mm | 191 |
| PH726610 | Hegar Uterine Dilator Double Ended Sizes 15mm And 16mm Diameter 220mm | 191 |
| PH726611 | Hegar Uterine Dilator Double Ended Sizes 17mm And 18mm Diameter 220mm | 191 |
| PH726612 | Hegar Uterine Dilator Double Ended Sizes 19mm And 20mm Diameter 220mm | 191 |
| PH726615 | Hegar Uterine Dilator Double Ended Sizes 21mm And 22mm Diameter 220mm | 191 |
| PH726616 | Hegar Uterine Dilator Double Ended Sizes 23mm And 24mm Diameter 220mm | 191 |
| PH726617 | Hegar Uterine Dilator Double Ended Sizes 25mm And 26mm Diameter 220mm | 191 |
| PH726625 | Hawkin Ambler Uterine Dilator Conical Size 3 6mm Diameter Single Ended 170mm | 191 |
| PH726626 | Hawkin Ambler Uterine Dilator Conical Size 4 7mm Diameter Single Ended 170mm | 191 |
| PH726627 | Hawkin Ambler Uterine Dilator Conical Size 5 8mm Diameter Single Ended 170mm | 191 |
| PH726628 | Hawkin Ambler Uterine Dilator Conical Size 6 9mm Diameter Single Ended 170mm | 191 |
| PH726629 | Hawkin Ambler Uterine Dilator Conical Size 7 10mm Diameter Single Ended 170mm | 191 |
| PH726630 | Hawkin Ambler Uterine Dilator Conical Size 8 11mm Diameter Single Ended 170mm | 191 |
| PH726631 | Hawkin Ambler Uterine Dilator Conical Size 9 12mm Diameter Single Ended 170mm | 191 |
| PH726632 | Hawkin Ambler Uterine Dilator Conical Size 10 13mm Diameter Single Ended 170mm | 191 |
| PH726633 | Hawkin Ambler Uterine Dilator Conical Size 11 14mm Diameter Single Ended 170mm | 191 |
| PH726634 | Hawkin Ambler Uterine Dilator Conical Size 12 15mm Diameter Single Ended 170mm | 191 |
| PH726635 | Hawkin Ambler Uterine Dilator Conical Size 13 16mm Diameter Single Ended 170mm | 191 |
| PH726636 | Hawkin Ambler Uterine Dilator Conical Size 14 17mm Diameter Single Ended 170mm | 191 |
| PH726637 | Hawkin Ambler Uterine Dilator Conical Size 15 18mm Diameter Single Ended 170mm | 191 |
| PH726638 | Hawkin Ambler Uterine Dilator Conical Size 16 19mm Diameter Single Ended 170mm | 191 |
| PH726639 | Hawkin Ambler Uterine Dilator Conical Size 17 20mm Diameter Single Ended 170mm | 191 |
| PH726640 | Hawkin Ambler Uterine Dilator Conical Size 18 21mm Diameter Single Ended 170mm | 191 |
| PH726660 | Anderson Obstetric Forceps 380mm | 183 |
| PH726662 | Kielland Obstetric Forceps For Occipito-Posterior Position 410mm | 183 |
| PH726670 | Wrigley Obstetric Forceps Short Model For Caesarian Section 280mm | 184 |
| PH726672 | Barnes Neville Obstetric Forceps With Axis Traction Rod (Neville Barnes) 370mm | 184 |
| PH726673 | Barnes Neville Obstetric Forceps Without Axis Traction Rod (Neville Barnes) 370 | 185 |
| PH727050 | Tischler Biopsy Punch Forceps Crocodile Action 210mm Effective Length With Singl | 201 |

Product Codes Index

| Code | Description | Page |
|----------|---|------|
| PH727051 | Baby Tischler Biopsy Punch Forceps Crocodile Action 210mm Effective Length With | 201 |
| PH727055 | Mini Townsend Biopsy Punch Forceps Crocodile Action 210mm Effective Length With | 201 |
| PH727065 | Uni Tischler Biopsy Punch Forceps Crocodile Action 230mm Effective Length With R | 202 |
| PH727066 | Uni Tischler Biopsy Punch Forceps Crocodile Action 230mm Effective Length With Rotating Punch Complete Upward Angle Cut 3 x 7mm 230mm | 202 |
| PH727079 | Uni Eppendorf Biopsy Punch Forceps Crocodile Action 230mm Effective Length With Rotating Punch Complete Straight Cut 3 x 6mm 230mm | 202 |
| PH727084 | Schumacher Biopsy Punch For Gynaecology Crocodile Action 230mm Effective Length Complete Straight Cut 4.5 x 7mm 230mm | 203 |
| PH727086 | Kevorkian Biopsy Punch For Gynaecology Crocodile Action 260mm Effective Length Complete Straight Cut 3 x 10mm 260mm | 203 |
| PH727089 | Eppendorf Biopsy Punch For Gynaecology Crocodile Action 200mm Effective Length Complete Straight Cut 4 x 8mm 200mm | 203 |
| PH737700 | Gigli Metal Handle (Supplied In A Pair) | 204 |
| PH737701 | Gigli Saw 30cm Length | 204 |
| PH737704 | Gigli Saw 40cm Length 40cm | 204 |
| PH737705 | Gigli Saw 50cm Length 50cm | 204 |
| PH737706 | Gigli Saw 60cm Length | 204 |
| PH757500 | Monopolar (Diameterthermy) Lead With 4mm Machine End And Std 4.8mm Cutting Socket | 178 |
| PH757501 | Monopolar (Diameterthermy) Lead With 8mm Bovie Machine End And Std 4.8mm Cutting Sock | 178 |
| PH757502 | Laparoscopic Monopolar Lead With 4mm Machine End And 4mm Laparoscopic Cutting Sock | 178 |
| PH757503 | Laparoscopic Monopolar Lead With 8mm Bovie Machine End And 4mm Laparoscopic Cutting Sock | 178 |
| PH757510 | Bipolar Lead With Std European Instrument Block Fitting And 2 x 4mm Plugs On Mac | 178 |
| PH757550 | Sims Diameterthermy Vaginal Speculum Double Ended Small 60mm x 30mm | 200 |
| PH757551 | Sims Diameterthermy Vaginal Speculum Double Ended Medium 70mm x 35mm | 200 |
| PH757552 | Sims Diameterthermy Vaginal Speculum Double Ended Large 80mm x 40mm | 200 |
| PH757560 | Diameterthermy Quiver Standard 60mm Diameter 230mm | 179 |
| PH757561 | Diameterthermy Quiver Telescopic 60mm Diameter Extending 400mm | 179 |
| PH760001 | Endoscope / Arthroscope Basket With Small Lens Box 640mm x 150mm x 77mm | 206 |
| PH760002 | Universal Scope Basket 310mm x 150mm x 77mm | 207 |
| PH760003 | Double Scope Basket 290mm x 160mm x 52mm | 206 |
| PH760004 | Endoscope / Arthroscope Basket Double Scope Inserts 290mm x 80mm x 52mm | 206 |
| PH760005 | Endoscope / Arthroscope Basket Single Scope Inserts 290mm x 80mm x 52mm | 206 |
| PH760006 | Endoscope / Arthroscope Basket Single Scope Inserts 460mm x 80mm x 52mm | 206 |
| PH760007 | Endoscope / Arthroscope Basket Double Scope Inserts 460mm x 80mm x 52mm | 206 |

Product Codes Index

| Code | Description | Page |
|----------|---|------|
| PH760008 | Endoscope / Arthroscope Basket 670mm x 80mm x 52mm | 207 |
| PH760009 | Double Scope Basket 460mm x 160mm x 52mm | 206 |
| PH760010 | Universal Scope Basket 580mm x 160mm x 52mm | 207 |
| PH760011 | Universal Scope Basket 290mm x 160mm x 52mm | 207 |
| PH760012 | Universal Scope Basket 460mm x 160mm x 52mm | 207 |
| PH760020 | Wire Mesh Base Basket With Perforated Sides 230mm x 130mm x 50mm | 208 |
| PH760021 | Wire Mesh Base Basket With Perforated Sides 230mm x 130mm x 60mm | 208 |
| PH760022 | Wire Mesh Base Basket With Perforated Sides 230mm x 130mm x 70mm | 208 |
| PH760023 | Wire Mesh Base Basket With Perforated Sides 240mm x 130mm x 50mm | 208 |
| PH760024 | Wire Mesh Base Basket With Perforated Sides 240mm x 130mm x 60mm | 208 |
| PH760025 | Wire Mesh Base Basket With Perforated Sides 280mm x 200mm x 60mm | 208 |
| PH760030 | Wire Mesh Base Basket With Perforated Sides 240mm x 250mm x 50mm Half Din | 208 |
| PH760031 | Wire Mesh Base Basket With Perforated Sides 240mm x 250mm x 60mm | 208 |
| PH760032 | Wire Mesh Base Basket With Perforated Sides 240mm x 250mm x 70mm | 208 |
| PH760033 | Wire Mesh Base Basket With Perforated Sides 240mm x 250mm x 100mm | 208 |
| PH760036 | Wire Mesh Base Basket With Perforated Sides 300mm x 300mm x 50mm | 208 |
| PH760037 | Wire Mesh Base Basket With Perforated Sides 300mm x 300mm x 100mm | 208 |
| PH760038 | Wire Mesh Base Basket With Perforated Sides 320mm x 320mm x 60mm | 208 |
| PH760039 | Wire Mesh Base Basket With Perforated Sides 380mm x 240mm x 50mm | 208 |
| PH760040 | Wire Mesh Base Basket With Perforated Sides 350mm x 350mm x 60mm | 208 |
| PH760043 | Wire Mesh Base Basket With Perforated Sides 450mm x 230mm x 50mm | 208 |
| PH760044 | Wire Mesh Base Basket With Perforated Sides 460mm x 320mm x 60mm | 208 |
| PH760045 | Wire Mesh Base Basket With Perforated Sides 480mm x 120mm x 60mm | 208 |
| PH760046 | Wire Mesh Base Basket With Perforated Sides 480mm x 150mm x 50mm | 208 |
| PH760047 | Wire Mesh Base Basket With Perforated Sides 480mm x 250mm x 50mm One Din | 208 |
| PH760048 | Wire Mesh Base Basket With Perforated Sides 480mm x 250mm x 60mm | 208 |
| PH760049 | Wire Mesh Base Basket With Perforated Sides 480mm x 250mm x 70mm | 208 |
| PH760050 | Wire Mesh Base Basket With Perforated Sides 480mm x 250mm x 100mm | 208 |
| PH760051 | Wire Mesh Base Basket With Perforated Sides 480mm x 250mm x 120mm | 208 |
| PH760056 | Wire Mesh Base Basket With Perforated Sides 470mm x 400mm x 60mm | 208 |
| PH760057 | Wire Mesh Base Basket With Perforated Sides 480mm x 350mm x 60mm | 208 |
| PH760058 | Wire Mesh Base Basket With Perforated Sides 480mm x 480mm x 50mm Two Din | 208 |

Product Codes Index

| Code | Description | Page |
|----------|---|------|
| PH760059 | Wire Mesh Base Basket With Perforated Sides 480mm x 480mm x 60mm | 208 |
| PH760060 | Wire Mesh Base Basket With Perforated Sides 480mm x 480mm x 70mm | 208 |
| PH760061 | Wire Mesh Base Basket With Perforated Sides 480mm x 480mm x 100mm | 208 |
| PH760062 | Wire Mesh Base Basket With Perforated Sides 480mm x 480mm x 120mm | 208 |
| PH760065 | Wire Mesh Base Basket With Perforated Sides 500mm x 500mm x 50mm | 208 |
| PH760066 | Wire Mesh Base Basket With Perforated Sides 520mm x 250mm x 50mm | 208 |
| PH760067 | Wire Mesh Base Basket With Perforated Sides 520mm x 250mm x 60mm | 208 |
| PH760068 | Wire Mesh Base Basket With Perforated Sides 520mm x 250mm x 75mm | 208 |
| PH760069 | Wire Mesh Base Basket With Perforated Sides 520mm x 250mm x 95mm | 208 |
| PH760070 | Wire Mesh Base Basket With Perforated Sides 520mm x 300mm x 60mm | 208 |
| PH760075 | Wire Mesh Base Basket With Perforated Sides 540mm x 250mm x 95mm | 208 |
| PH760076 | Wire Mesh Base Basket With Perforated Sides 550mm x 250mm x 60mm | 208 |
| PH760077 | Wire Mesh Base Basket With Perforated Sides 550mm x 280mm x 50mm | 208 |
| PH760078 | Wire Mesh Base Basket With Perforated Sides 550mm x 325mm x 50mm | 208 |
| PH760079 | Wire Mesh Base Basket With Perforated Sides 550mm x 430mm x 50mm | 208 |
| PH760085 | Wire Mesh Base Basket With Perforated Sides 600mm x 250mm x 150mm | 208 |
| PH760086 | Wire Mesh Base Basket With Perforated Sides 600mm x 300mm x 60mm | 208 |
| PH760087 | Wire Mesh Base Basket With Perforated Sides 600mm x 350mm x 60mm | 208 |
| PH760088 | Wire Mesh Base Basket With Perforated Sides 670mm x 300mm x 70mm | 208 |
| PH760089 | Wire Mesh Base Basket With Perforated Sides 700mm x 250mm x 60mm | 208 |
| PH760090 | Wire Mesh Base Basket With Perforated Sides 700mm x 300mm x 90mm | 208 |
| PH760091 | Wire Mesh Base Basket With Perforated Sides 900mm x 450mm x 50mm | 208 |
| PH760093 | Lid For Flat Base Basket 520 x 250mm | 208 |
| PH760094 | Lid For Flat Base Basket 540mm x 250mm | 208 |
| PH760095 | Lid For Flat Base Basket 240mm x 250mm | 208 |
| PH760096 | Lid For Flat Base Basket 480mm x 250mm | 208 |
| PH760097 | Lid For Flat Base Basket 480mm x 480mm | 208 |
| PH760098 | Lid For Flat Base Basket 300mm x 300mm | 208 |
| PH760099 | Lid For Flat Base Basket 670mm x 300mm | 208 |
| PH760140 | All Wire Mesh Basket 240mm x 250mm x 50mm | 209 |
| PH760141 | All Wire Mesh Basket 240mm x 250mm x 60mm | 209 |
| PH760142 | All Wire Mesh Basket 240mm x 250mm x 70mm | 209 |

Product Codes Index

| Code | Description | Page |
|----------|---|------|
| PH760143 | All Wire Mesh Basket 480mm x 250mm x 50mm | 209 |
| PH760144 | All Wire Mesh Basket 480mm x 250mm x 60mm | 209 |
| PH760145 | All Wire Mesh Basket 480mm x 250mm x 70mm | 209 |
| PH760146 | All Wire Mesh Basket 480mm x 250mm x 100mm | 209 |
| PH760147 | All Wire Mesh Basket 480mm x 250mm x 120mm | 209 |
| PH760148 | All Wire Mesh Basket 480mm x 480mm x 50mm | 209 |
| PH760149 | All Wire Mesh Basket 480mm x 480mm x 60mm | 209 |
| PH760150 | All Wire Mesh Basket 480mm x 480mm x 70mm | 209 |
| PH760151 | All Wire Mesh Basket 300mm x 300mm x 60mm | 209 |
| PH760152 | All Wire Mesh Basket 470mm x 400mm x 60mm | 209 |
| PH760170 | Lid For All Wire Mesh Basket 240mm x 250mm | 210 |
| PH760171 | Lid For All Wire Mesh Basket 480mm x 250mm | 210 |
| PH760172 | Lid For All Wire Mesh Basket 300mm x 300mm | 210 |
| PH760173 | Lid For All Wire Mesh Basket 480mm x 250mm | 210 |
| PH760174 | Lid For All Wire Mesh Basket 480mm x 480mm | 210 |
| PH760175 | Trinket Box Perforated 70mm x 30mm x 15mm | 210 |
| PH760176 | Trinket Box Perforated 90mm x 60mm x 25mm | 210 |
| PH760177 | Trinket Box Perforated 105mm x 70mm x 25mm | 210 |
| PH760178 | Trinket Box Perforated 120mm x 80mm x 25mm | 210 |
| PH760185 | Trinket Box Perforated Sides And Wire Mesh Lid 90mm x 60mm x 25mm | 210 |
| PH760186 | Trinket Box Perforated Sides And Wire Mesh Lid 105mm x 70mm x 25mm | 210 |
| PH760187 | Trinket Box Perforated Sides And Wire Mesh Lid 120mm x 80mm x 25mm | 210 |
| PH760189 | Trinket Box Perforated Sides And Wire Mesh Lid 150mm x 100mm x 25mm | 210 |
| PH760190 | Trinket Box Perforated Sides And Wire Mesh Lid 180mm x 130mm x 25mm | 210 |
| PH760194 | Narrow Box Perforated Sides And Wire Mesh Lid 300mm x 80mm x 35mm | 209 |
| PH760195 | Narrow Box Perforated Sides And Wire Mesh Lid 360mm x 50mm x 25mm | 209 |
| PH760196 | Narrow Box Perforated Sides And Wire Mesh Lid 550mm x 50mm x 25mm | 209 |
| PH760197 | Narrow Box Perf Side Wire Mesh Lid 720 X100 x 25mm | 209 |
| PH760198 | Narrow Box Perf Sides & Mesh Lid 700 x 50 x 25mm | 209 |
| PH760200 | All Fine Mesh Box 80mm x 40mm x 30mm | 209 |
| PH760201 | All Fine Mesh Box 80mm x 80mm x 35mm | 209 |
| PH760202 | All Fine Mesh Box 80mm x 80mm x 55mm | 209 |

Product Codes Index

| Code | Description | Page |
|----------|--|------|
| PH760203 | All Fine Mesh Box 100mm x 100mm x 55mm | 209 |
| PH760204 | All Fine Mesh Box 105mm x 70mm x 25mm | 209 |
| PH760206 | All Fine Mesh Box 150mm x 50mm x 35mm | 209 |
| PH760207 | All Fine Mesh Box 150mm x 70mm x 35mm | 209 |
| PH760208 | All Fine Mesh Box 150mm x 150mm x 35mm | 209 |
| PH760210 | All Fine Mesh Box 190mm x 85mm x 100mm | 209 |
| PH760211 | All Fine Mesh Box 190mm x 70mm x 70mm | 209 |
| PH760215 | All Fine Mesh Box 200mm x 100mm x 35mm | 209 |
| PH760216 | All Fine Mesh Box 220mm x 140mm x 50mm | 209 |
| PH760217 | All Fine Mesh Box 235mm x 80mm x 30mm | 209 |
| PH760218 | All Fine Mesh Box 235mm x 155mm x 40mm | 209 |
| PH760222 | All Fine Mesh Box 240mm x 90mm x 40mm | 209 |
| PH760224 | All Fine Mesh Box 240mm x 120mm x 35mm | 209 |
| PH760227 | All Fine Mesh Box 275mm x 175mm x 35mm | 209 |
| PH760228 | All Fine Mesh Box 280mm x 200mm x 60mm | 209 |
| PH760255 | Fine Mesh Basket With Take Off Lockable Lid 120mm x 155mm x 40mm | 210 |
| PH760256 | Fine Mesh Basket With Take Off Lockable Lid 120mm x 155mm x 50mm | 210 |
| PH760259 | Fine Mesh Basket With Take Off Lockable Lid 235mm x 155mm x 40mm | 210 |
| PH760260 | Fine Mesh Basket With Take Off Lockable Lid 235mm x 155mm x 50mm | 210 |
| PH760264 | Fine Mesh Basket With Take Off Lockable Lid 380mm x 265mm x 45mm | 210 |
| PH760265 | Fine Mesh Basket With Take Off Lockable Lid 380mm x 265mm x 55mm | 210 |
| PH760270 | Basket For Ophthalmic Instruments 440mm x 250mm x 40mm | 207 |
| PH760528 | Sieves 150mm (6inch) Diameter | 185 |
| PH760529 | Sieves 200mm (8inch) Diameter | 185 |
| PH80526 | Nail Nipper Straight Double Spring With Lock 140mm Straight | 34 |
| PH80527 | Nail Nipper Concave Double Spring With Lock 140mm Curved | 34 |
| PH80528 | Nail Nipper Concave Double Spring With Lock 160mm Curved | 34 |
| PH80534 | Straight Nail Nipper Compound Action 150mm Straight | 34 |
| PH80536 | Dovo German Pattern Nail Nipper Straight 140mm Straight | 34 |
| PH80537 | Ingrown Toe Nail Nipper Sharp Pointed Straight 130mm Straight | 34 |
| PH80541 | Thwaites Nail Nipper Pistol Grip Single Spring 140mm | 34 |
| PH80546 | Double Ended Spatula And Packer With Solid Handle 155mm Straight | 35 |

Product Codes Index

| Code | Description | Page |
|---------|--|------|
| PH80549 | Fickling Nail Elevator 115mm Straight | 35 |
| PH80553 | Blacks File With A Double Ended Course Cut 140mm | 35 |
| PH90555 | Adson Artery Forceps With Box Joint 180mm Straight | 41 |
| PH90556 | Adson Artery Forceps With Box Joint 180mm Curved | 41 |
| PH90558 | Birkett Artery Forceps With Box Joint 180mm Straight | 41 |
| PH90559 | Birkett Artery Forceps With Box Joint 180mm Curved | 41 |
| PH90561 | Bland Sutton Artery Forceps With Screw Joint 190mm Curved | 41 |
| PH90570 | Cairn Artery Forceps With Box Joint 145mm Straight | 41 |
| PH90571 | Cairn Artery Forceps With Box Joint 145mm Curved | 42 |
| PH90572 | Cairn Curved To One Side Artery Forceps With Box Joint 145mm Curved | 42 |
| PH90582 | Crile Artery Forceps With Box Joint 140mm Straight | 42 |
| PH90583 | Crile Artery Forceps With Box Joint 160mm Straight | 42 |
| PH90584 | Crile Artery Forceps With Box Joint 180mm Straight | 42 |
| PH90585 | Crile Artery Forceps With Box Joint 140mm Curved | 42 |
| PH90586 | Crile Artery Forceps With Box Joint 160mm Curved | 42 |
| PH90587 | Crile Artery Forceps With Box Joint 180mm Curved | 42 |
| PH90588 | Crile Rankin Artery Forceps With Box Joint 140mm Straight | 42 |
| PH90589 | Crile Rankin Artery Forceps With Box Joint 160mm Straight | 42 |
| PH90590 | Crile Rankin Artery Forceps With Box Joint 140mm Curved | 42 |
| PH90591 | Crile Rankin Artery Forceps With Box Joint 160mm Curved | 42 |
| PH90592 | Cushing Artery Forceps With Box Joint 145mm Straight | 43 |
| PH90593 | Cushing Artery Forceps With Box Joint 145mm Curved | 43 |
| PH90594 | Cushing Curved To One Side Artery Forceps With Box Joint 145mm Curved | 43 |
| PH90596 | Denis Browne Tonsil Forceps Angled To One Side 200mm (Dennis Brown) | 43 |
| PH90605 | Dunhill Artery Forceps With Box Joint 130mm Straight | 43 |
| PH90607 | Dunhill Artery Forceps With Box Joint 130mm Curved | 43 |
| PH90608 | Dunhill Artery Forceps With Box Joint 180mm Curved | 43 |
| PH90612 | Fickling Artery Forceps Angled And Plain With Box Joint 180mm Angled | 44 |
| PH90613 | Fickling Artery Forceps Angled 1 Into 2 Teeth And Box Joint 180mm Angled | 44 |
| PH90618 | Grey Turner Artery Forceps With Box Joint 145mm Straight | 44 |
| PH90619 | Grey Turner Artery Forceps With Box Joint 140mm Curved | 44 |
| PH90620 | Grey Turner Artery Forceps With Box Joint 180mm Straight | 44 |
| PH90621 | Grey Turner Artery Forceps With Box Joint 180mm Curved | 44 |

Product Codes Index

| Code | Description | Page |
|---------|--|------|
| PH90622 | Hartmann Artery Forceps Extra Fine With Box Joint 90mm Straight | 44 |
| PH90623 | Hartmann Artery Forceps Extra Fine With Box Joint 90mm Curved | 44 |
| PH90624 | Hartmann Artery Forceps Extra Fine With Box Joint 100mm Straight | 45 |
| PH90626 | Halsted Mosquito Artery Forceps With Box Joint 130mm Straight | 45 |
| PH90627 | Halstead Mosquito Artery Forceps With Box Joint 200mm Straight | |
| PH90628 | Halsted Mosquito Artery Forceps With Box Joint 130mm Curved | 45 |
| PH90629 | Halsted Mosquito Artery Forceps 1 Into 2 Teeth With Box Joint 130mm Straight | 45 |
| PH90630 | Halsted Mosquito Artery Forceps 1 Into 2 Teeth With Box Joint 130mm Curved | 46 |
| PH90631 | Micro Halsted Mosquito Artery Forceps With Box Joint 130mm Straight | 46 |
| PH90632 | Micro Halsted Mosquito Artery Forceps With Box Joint 130mm Curved | 46 |
| PH90636 | Heiss Half Curve Artery Forceps With Box Joint 200mm | 46 |
| PH90637 | Heiss Full Curve Artery Forceps With Box Joint 200mm | 46 |
| PH90639 | Halsted Mosquito Artery Forceps Curved With Atraumatic Jaws 130mm | 46 |
| PH90645 | Jolls Artery Forceps With Box Joint 145mm Straight | 47 |
| PH90646 | Jolls Artery Forceps With Box Joint 145mm Curved | 47 |
| PH90650 | Kelly Artery Forceps With Box Joint 140mm Straight | 47 |
| PH90651 | Kelly Artery Forceps With Box Joint 160mm Straight | 47 |
| PH90652 | Kelly Artery Forceps With Box Joint 140mm Curved | 47 |
| PH90653 | Kelly Artery Forceps With Box Joint 160mm Curved | 47 |
| PH90655 | Kelly Rankin Artery Forceps With Box Joint 160mm Curved | 47 |
| PH90657 | Kelly Frazier Artery Forceps With Box Joint 180mm Curved | 47 |
| PH90682 | Kocher Artery Forceps With Box Joint 1 Into 2 Teeth 140mm Straight | 48 |
| PH90683 | Kocher Artery Forceps With Box Joint 1 Into 2 Teeth 150mm Straight | 48 |
| PH90684 | Kocher Artery Forceps With Box Joint 1 Into 2 Teeth 180mm Straight | 48 |
| PH90685 | Kocher Artery Forceps With Box Joint 1 Into 2 Teeth 200mm Straight | 48 |
| PH90686 | Kocher Artery Forceps With Box Joint 1 Into 2 Teeth 230mm Straight | 48 |
| PH90691 | Kocher Artery Forceps With Box Joint 1 Into 2 Teeth 130mm Curved | 48 |
| PH90692 | Kocher Artery Forceps With Box Joint 1 Into 2 Teeth 140mm Curved | 48 |
| PH90693 | Kocher Artery Forceps With Box Joint 1 Into 2 Teeth 150mm Curved | 48 |
| PH90694 | Kocher Artery Forceps With Box Joint 1 Into 2 Teeth 180mm Curved | 48 |
| PH90695 | Kocher Artery Forceps With Box Joint 1 Into 2 Teeth 200mm Curved | 48 |
| PH90696 | Kocher Artery Forceps With Box Joint 1 Into 2 Teeth 230mm Curved | 48 |
| PH90708 | Baby Mixer Artery Forceps With Box Joint 140mm Curved | 48 |

Product Codes Index

| Code | Description | Page |
|---------|--|------|
| PH90709 | Baby Mixer Artery Forceps With Box Joint 180mm Curved | 48 |
| PH90712 | Mixer Artery Forceps With Box Joint 200mm Curved | 48 |
| PH90713 | Mixer Artery Forceps With Box Joint 230mm Curved | 48 |
| PH90714 | Mixer Artery Forceps With Box Joint 280mm Curved | 48 |
| PH90716 | Mixer O'shaughnessy Artery Forceps With Box Joint 150mm Curved | 48 |
| PH90717 | Mixer O'shaughnessy Artery Forceps With Box Joint (O Shaughnessy) 190mm Curved | 48 |
| PH90722 | Moynihan Artery Forceps With Box Joint 145mm Curved | 48 |
| PH90723 | Moynihan Artery Forceps With Box Joint 180mm Curved | 48 |
| PH90730 | Negus Short Curve Artery Forceps With Screw Joint 190mm | 49 |
| PH90731 | Negus Long Curve Artery Forceps With Screw Joint 190mm Curved | 49 |
| PH90760 | Rochester Pean Artery Forceps With Box Joint 130mm Straight | 49 |
| PH90761 | Rochester Pean Artery Forceps With Box Joint 140mm Straight | 49 |
| PH90762 | Rochester Pean Artery Forceps With Box Joint 160mm Straight | 49 |
| PH90763 | Rochester Pean Artery Forceps With Box Joint 180mm Straight | 49 |
| PH90764 | Rochester Pean Artery Forceps With Box Joint 190mm Straight | 49 |
| PH90765 | Rochester Pean Artery Forceps With Box Joint 200mm Straight | 49 |
| PH90766 | Rochester Pean Artery Forceps With Box Joint 230mm Straight | 49 |
| PH90770 | Rochester Pean Artery Forceps With Box Joint 130mm Curved | 49 |
| PH90771 | Rochester Pean Artery Forceps With Box Joint 140mm Curved | 49 |
| PH90772 | Rochester Pean Artery Forceps With Box Joint 160mm Curved | 49 |
| PH90773 | Rochester Pean Artery Forceps With Box Joint 180mm Curved | 49 |
| PH90774 | Rochester Pean Artery Forceps With Box Joint 190mm Curved | 49 |
| PH90775 | Rochester Pean Artery Forceps With Box Joint 200mm Curved | 49 |
| PH90776 | Rochester Pean Artery Forceps With Box Joint 230mm Curved | 49 |
| PH90804 | Robert Artery Forceps With Box Joint 230mm (Roberts) 230mm Curved | 49 |
| PH90805 | Robert Artery Forceps With Screw Joint 230mm Curved | 49 |
| PH90807 | Sawtell Artery Forceps With Box Joint 180mm Curved | 50 |
| PH90809 | Wilson Tonsil Artery Forceps With Angled Shanks Screw Joint 200mm Curved | 50 |
| PH90815 | Robert Artery Forceps With Screw Joint (300mm Curved | 49 |
| PH90836 | Spencer Wells Artery Forceps With Box Joint 130mm Straight | 50 |
| PH90837 | Spencer Wells Artery Forceps With Box Joint 150mm Straight | 50 |
| PH90838 | Spencer Wells Artery Forceps With Box Joint 180mm Straight | 50 |

Product Codes Index

| Code | Description | Page |
|-------------|--|------|
| PH90839 | Spencer Wells Artery Forceps With Box Joint 200mm Straight | 50 |
| PH90840 | Spencer Wells Artery Forceps With Box Joint 230mm Straight | 50 |
| PH90841 | Spencer Wells Artery Forceps With Box Joint 300mm Straight | 50 |
| PH90847 | Spencer Wells Artery Forceps With Box Joint 150mm Curved | 50 |
| PH90848 | Spencer Wells Artery Forceps With Box Joint 180mm Curved | 50 |
| PH90849 | Spencer Wells Artery Forceps With Box Joint 200mm Curved | 50 |
| PH90850 | Spencer Wells Artery Forceps With Box Joint 230mm Curved | 50 |
| PH90851 | Spencer Wells Artery Forceps With Box Joint 300mm Curved | 50 |
| PH90860 | Schmidt Artery Forceps With Box Joint 190mm Curved | 50 |
| PH90861 | Sarrot Artery Forceps With Box Joint 240mm Straight | 51 |
| PH90862 | Sarrot Artery Forceps With Box Joint 240mm Curved | 51 |
| PH90870 | Belcher Gigli Saw Passing Forceps 300mm | 51 |
| PH90880 | Lloyd Davies Artery Forceps 250mm Straight | 51 |
| PH90890 | Kilner Artery Forceps Straight 145mm | 51 |
| PH90891 | Kilner Artery Forceps Curved 145mm 145mm | 51 |
| PHS800049 | 8" Receiver / Kidney Dish 203.2mm | 211 |
| PHS800049NR | Receivers SS 8" (Without Lid) Non Rolled Lips | 211 |
| PHS800053 | Receivers SS 10" (Without Lid) | 211 |
| PHS800053NR | Receivers SS 10" (Without Lid) Non Rolled Lips | 211 |
| PHS800055NR | Gallipots SS (80mm x 60mm) Non Rolled Lips | 211 |
| PHS800056 | Gallipots SS (60mm x 50mm) | 211 |

Product Name A-Z

| Name | Page |
|---|------|
| 8" Receiver / Kidney Dish 203.2mm | 211 |
| Adson Artery Forceps With Box Joint 180mm Curved | 41 |
| Adson Artery Forceps With Box Joint 180mm Straight | 41 |
| Adson Browne Dissecting Forceps With Serrated Jaws 130mm 130mm Straight | 53 |
| Adson Dissecting Forceps Tungsten Carbide Jaws 1:2 Teeth 120mm Straight | 60 |
| Adson Dissecting Forceps With 1 Into 2 Teeth 130mm Straight | 53 |
| Adson Dissecting Forceps With 1 Into 2 Teeth 180mm Straight | 53 |
| Adson Dissecting Forceps With Serrated Jaws 130mm Straight | 53 |
| Adson Dissecting Forceps With Serrated Jaws 180mm Straight | 53 |
| Adson Dissecting Forceps With Tungsten Carbide Jaws 120mm Straight | 60 |
| Adson Micro Dissecting Forceps With 1 Into 2 Teeth 130mm Straight | 52 |
| Adson Micro Dissecting Forceps With Serrated Jaws 130mm Straight | 52 |
| Agnew Myringotome Ear Knife Angled 190mm Angled | 127 |
| All Fine Mesh Box 100mm x 100mm x 55mm | 209 |
| All Fine Mesh Box 105mm x 70mm x 25mm | 209 |
| All Fine Mesh Box 150mm x 150mm x 35mm | 209 |
| All Fine Mesh Box 150mm x 50mm x 35mm | 209 |
| All Fine Mesh Box 150mm x 70mm x 35mm | 209 |
| All Fine Mesh Box 190mm x 70mm x 70mm | 209 |
| All Fine Mesh Box 190mm x 85mm x 100mm | 209 |
| All Fine Mesh Box 200mm x 100mm x 35mm | 209 |
| All Fine Mesh Box 220mm x 140mm x 50mm | 209 |
| All Fine Mesh Box 235mm x 155mm x 40mm | 209 |
| All Fine Mesh Box 235mm x 80mm x 30mm | 209 |
| All Fine Mesh Box 240mm x 120mm x 35mm | 209 |
| All Fine Mesh Box 240mm x 90mm x 40mm | 209 |

| Name | Page |
|--|------|
| All Fine Mesh Box 275mm x 175mm x 35mm | 209 |
| All Fine Mesh Box 280mm x 200mm x 60mm | 209 |
| All Fine Mesh Box 80mm x 40mm x 30mm | 209 |
| All Fine Mesh Box 80mm x 80mm x 35mm | 209 |
| All Fine Mesh Box 80mm x 80mm x 55mm | 209 |
| All Wire Mesh Basket 240mm x 250mm x 50mm | 209 |
| All Wire Mesh Basket 240mm x 250mm x 60mm | 209 |
| All Wire Mesh Basket 240mm x 250mm x 70mm | 209 |
| All Wire Mesh Basket 300mm x 300mm x 60mm | 209 |
| All Wire Mesh Basket 470mm x 400mm x 60mm | 209 |
| All Wire Mesh Basket 480mm x 250mm x 100mm | 209 |
| All Wire Mesh Basket 480mm x 250mm x 120mm | 209 |
| All Wire Mesh Basket 480mm x 250mm x 50mm | 209 |
| All Wire Mesh Basket 480mm x 250mm x 60mm | 209 |
| All Wire Mesh Basket 480mm x 250mm x 70mm | 209 |
| All Wire Mesh Basket 480mm x 480mm x 50mm | 209 |
| All Wire Mesh Basket 480mm x 480mm x 60mm | 209 |
| All Wire Mesh Basket 480mm x 480mm x 70mm | 209 |
| Allis Tissue Forceps With 3 Into 4 Teeth And Box Joint 150mm | 38 |
| Allis Tissue Forceps With 4 Into 5 Teeth And Box Joint 150mm | 38 |
| Allis Tissue Forceps With 5 Into 6 Teeth And Box Joint 150mm | 38 |
| Allis Tissue Forceps With 5 Into 6 Teeth And Box Joint 180mm | 38 |
| Allis Tissue Forceps With 5 Into 6 Teeth And Box Joint 200mm | 38 |
| Allis Tissue Forceps With 5 Into 6 Teeth And Box Joint 250mm | 38 |
| Allis Tissue Forceps With 5 Into 6 Teeth And Box Joint 280mm | 38 |

Product Name A-Z

| Name | Page |
|---|------|
| Allison Retractor Adult Size 300mm | 100 |
| Allison Retractor Child Size 280mm | 100 |
| Allison Retractor Infant Size 280mm | 100 |
| Alms Skin Retractor With Blunt Prongs 100mm | 107 |
| Alms Skin Retractor With Blunt Prongs 70mm | 107 |
| Alms Skin Retractor With Sharp Prongs 100mm | 107 |
| Alms Skin Retractor With Sharp Prongs 70mm | 107 |
| American Pattern Fistula Scissors 180mm Curved | 19 |
| American Pattern Umbilical Cord Cutter 115mm Straight | 18 |
| Anderson Obstetric Forceps 380mm | 183 |
| Aneurysm Needle With Flat Handle And Curved Laterally 140mm | 92 |
| Aneurysm Needle With Oval Handle And Curved Laterally 140mm | 92 |
| Aufricht Nasal Retractor 50mm Blade 160mm Angled | 101 |
| Aufricht Scissors 140mm Straight | 13 |
| Aufricht Scissors 150mm Curved | 13 |
| Austin Kay Right Angled Retractor With Lip 150mm Angled | 100 |
| Auvarde Speculum 80 x 45 Tapered 40mm Fixed Weight | 198 |
| Awty Channel Retractors 12mm Blade Width | 98 |
| Awty Channel Retractors 17mm Blade Width | 98 |
| Babcock Tissue Forceps With Box Joint 160mm | 38 |
| Babcock Tissue Forceps With Box Joint 180mm | 38 |
| Babcock Tissue Forceps With Box Joint 200mm | 38 |
| Babcock Tissue Forceps With Box Joint 230mm | 38 |
| Baby Hohman Bone Levers 6mm | 98 |
| Baby Mixter Artery Forceps With Box Joint 140mm Curved | 48 |
| Baby Mixter Artery Forceps With Box Joint 180mm Curved | 48 |

| Name | Page |
|--|------|
| Baby Senn Miller Retractor Sharp Double Ended | 107 |
| Baby Tischler Biopsy Punch Forceps Crocodile Action 210mm Effective Length With | 201 |
| Backhaus Towel Clip With Box Joint 115mm Curved | 78 |
| Backhaus Towel Clip With Box Joint 135mm Curved | 78 |
| Backhaus Towel Clip With Box Joint 90mm Curved | 78 |
| Bailey Vaginal Speculum 38mm Blade & Fixed Weight | 198 |
| Balfour Abdominal Retractor With Centre Blade And Fenestrated End Blades 63mm X | 108 |
| Balfour Baby Abdominal Retractor With Centre Blade And Fenestrated End Blades 27mm Deep | 108 |
| Ballenger Nasal Cartilage Knife 3mm Wide Blade With Straight Body And A Revolving Blade 200mm Straight | 132 |
| Ballenger Nasal Cartilage Knife 5mm Wide Blade With Straight Body And A Revolving Blade 200mm Straight | 132 |
| Ballenger Nasal Cartilage Knife 8mm Wide Blade With Straight Body And A Revolving Blade 200mm Straight | 132 |
| Barnes Neville Obstetric Forceps With Axis Traction Rod (Neville Barnes) 370mm | 184 |
| Barnes Neville Obstetric Forceps Without Axis Traction Rod (Neville Barnes) 370 | 185 |
| Basket For Ophthalmic Instruments 440mm x 250mm x 40mm | 207 |
| Baumgartner Needle Holder With Box Joint 140mm Straight 140mm Straight | 82 |
| Baumgartner Needle Holder With Tungsten Carbide Jaws And Box Joint 135mm Straight | 85 |
| Beale Aural Forceps Crocodile Action With Fine Serrated Jaw And Black Finish 70mm Shoulder Length | 62 |
| Beck Rongeur Straight 130mm Shaft Length With A 2mm Bite For Intervertebral Disc Straight | 171 |
| Beckmann Adenoid Curette Size 0 With Hollow Handle 8mm 210mm | 130 |
| Beckmann Adenoid Curette Size 1 1/2 With Hollow Handle 12mm 210mm | 130 |
| Beckmann Adenoid Curette Size 1 With Hollow Handle 10mm 210mm | 130 |

Product Name A-Z

| Name | Page |
|--|------|
| Beckmann Adenoid Curette Size 2 With Hollow Handle 14mm 210mm | 130 |
| Beckmann Adenoid Curette Size 3 With Hollow Handle 16mm 210mm | 130 |
| Beckmann Adenoid Curette Size 4 With Hollow Handle 18mm 210mm | 130 |
| Belcher Gigli Saw Passing Forceps 300mm | 51 |
| Bergmann Plaster Shears/ Bandage Scissors With Angular Blades 230mm | 166 |
| Bickford Dissecting Forceps With 1 Into 2 Teeth 230mm Straight | 54 |
| Bickford Dissecting Forceps With Serrated Jaws 230mm | 53 |
| Bipolar Lead With Std European Instrument Block Fitting And 2 x 4mm Plugs On Mac | 178 |
| Birkett Artery Forceps With Box Joint 180mm Curved | 41 |
| Birkett Artery Forceps With Box Joint 180mm Straight | 41 |
| Blacks File With A Double Ended Course Cut 140mm | 35 |
| Bland Sutton Artery Forceps With Screw Joint 190mm Curved | 41 |
| Bohler Plaster Shears With A Screw Joint 230mm | 166 |
| Bone Holding Forcep 180mm | N/A |
| Bone Holding Forcep 200mm | N/A |
| Bone Holding Forcep 230mm | N/A |
| Bonney Dissecting Forceps With 1 Into 2 Teeth 180mm Straight | 54 |
| Bonney Dissecting Forceps With Serrated Jaws 180mm Straight | 54 |
| Bonney Uterine Polypus Forceps With Fenestrated Serrated Jaws And A Box Joint | 73 |
| Bowdler Henry Rake Retractor For Lower Wisdom Teeth 190mm | 102 |
| Bowlby Tracheal Dilating Forceps Adult Size 140mm Curved | |
| Bowlby Tracheal Dilating Forceps Child Size 115mm Curved | |
| Bowman Lachrymal Dilator Probe Double Ended Sizes 0.6mm And 0.7mm Dia 100mm | 93 |
| Bowman Lachrymal Dilator Probe Double Ended Sizes 0.7mm And 0.8mm Dia 100mm | 93 |
| Bowman Lachrymal Dilator Probe Double Ended Sizes 0.8mm And 0.9mm Dia 100mm | 93 |

| Name | Page |
|--|------|
| Bowman Lachrymal Dilator Probe Double Ended Sizes 0.9mm And 1.0mm Dia 100mm | 93 |
| Bowman Lachrymal Dilator Probe Double Ended Sizes 0.9mm And 1.1mm Dia 100mm | 93 |
| Bowman Lachrymal Dilator Probe Double Ended Sizes 1.1mm And 1.2mm Dia 100mm | 93 |
| Bowman Lachrymal Dilator Probe Double Ended Sizes 1.1mm And 1.3mm Dia 100mm | 93 |
| Bowman Lachrymal Dilator Probe Double Ended Sizes 1.3mm And 1.4mm Dia 100mm | 93 |
| Bowman Lachrymal Dilator Probe Double Ended Sizes 1.4mm And 1.5mm Dia 100mm | 93 |
| Bowman Lachrymal Dilator Probe Double Ended Sizes 1.5mm And 1.6mm Dia 100mm | 93 |
| Bowman Lachrymal Dilator Probe Double Ended Sizes 1.6mm And 1.8mm Dia 100mm | 93 |
| Bowman Lachrymal Dilator Probe Double Ended Sizes 1.7mm And 1.8mm Dia 100mm | 93 |
| Bowman Lachrymal Dilator Probe Double Ended Sizes 1.8mm And 1.9mm Dia 100mm | 93 |
| Bozemann Needle Holder With Tungsten Carbide Jaws And Box Joint With Curved Shanks 200mm | 85 |
| Braun Stadler Scissors 150mm (Episiotomy) Angled | 17 |
| Braun Stadler Scissors 190mm (Episiotomy) Angled | 17 |
| Braun Stadler Scissors 200mm (Episiotomy) Angled | 17 |
| Breisky Vaginal Speculum With Cranked Head Blade 100mm Length x 20mm Width 300mm | 200 |
| Breisky Vaginal Speculum With Cranked Head Blade 100mm Length x 25mm Width 300mm | 200 |
| Breisky Vaginal Speculum With Cranked Head Blade 100mm Length x 30mm Width 300mm | 200 |
| Breisky Vaginal Speculum With Cranked Head Blade 100mm Length x 35mm Width 300mm | 200 |
| Breisky Vaginal Speculum With Cranked Head Blade 100mm Length x 40mm Width 300mm | 200 |
| Breisky Vaginal Speculum With Cranked Head Blade 130mm Length x 20mm Width 320mm | 200 |
| Breisky Vaginal Speculum With Cranked Head Blade 130mm Length x 25mm Width 320mm | 200 |

Product Name A-Z

| Name | Page |
|--|------|
| Breisky Vaginal Speculum With Cranked Head Blade 130mm Length x 30mm Width 320mm | 200 |
| Breisky Vaginal Speculum With Cranked Head Blade 130mm Length x 35mm Width 320mm | 200 |
| Breisky Vaginal Speculum With Cranked Head Blade 130mm Length x 40mm Width 320mm | 200 |
| Breisky Vaginal Speculum With Cranked Head Blade 160mm Length x 40mm Width 350mm | 200 |
| Breisky Vaginal Speculum With Cranked Head Blade 180mm Length x 40mm Width 350mm | 200 |
| Bristow Bone Elevator Sharp Blade With A Hollow Handle 230mm Curved | 149 |
| Brodie Fistula Director Grooved 130mm Straight | 94 |
| Brodie Fistula Director Grooved 140mm Straight | 94 |
| Brodie Fistula Director Grooved 160mm Straight | 94 |
| Brodie Fistula Director Grooved With Probe End (Ref Approx 3.5mm Probe Dia) 130mm Straight | 94 |
| Brodie Fistula Director Grooved With Probe End (Ref Approx 3.5mm Probe Dia) 140mm Straight | 94 |
| Brodie Fistula Director Grooved With Probe End (Ref Approx 3.5mm Probe Dia) 160mm Straight | 94 |
| Bruns Plaster Shears With A Toothed Lower Blade 240mm | 166 |
| Burford Finochietto Rib Spreader With 2 Pairs Of Detachable Blades 290mm | 162 |
| Cairn (Rake) Scalp Retractor With 3 Prongs 16mm Sharp 210mm Straight | 101 |
| Cairn (Rake) Scalp Retractor With 4 Prongs 19mm Sharp 210mm Straight | 101 |
| Cairn Artery Forceps With Box Joint 145mm Curved | 42 |
| Cairn Artery Forceps With Box Joint 145mm Straight | 41 |
| Cairn Bone Rongeur With Fine Jaws And Angle On Flat 170mm | 144 |
| Cairn Curved To One Side Artery Forceps With Box Joint 145mm Curved | 42 |
| Cairn Dissecting Forceps With 2 Into 3 Teeth 165mm | 54 |

| Name | Page |
|--|------|
| Cairn Dissecting Forceps With 2 Into 3 Teeth 205mm | 54 |
| Capener Tufnol Handled Chisel 10mm Wide (Capner) 250mm | 156 |
| Capener Tufnol Handled Chisel 12mm Wide x 250mm | 156 |
| Capener Tufnol Handled Chisel 15mm Wide 250mm | 156 |
| Capener Tufnol Handled Chisel 20mm Wide 250mm | 156 |
| Capener Tufnol Handled Chisel 25mm Wide x 250mm | 156 |
| Capener Tufnol Handled Chisel 30mm Wide x 250mm | 156 |
| Capener Tufnol Handled Cranked Gouge 11mm Wide 235mm | 157 |
| Capener Tufnol Handled Cranked Gouge 16mm Wide 235mm | 157 |
| Capener Tufnol Handled Cranked Gouge 6mm Wide (Capner) 235mm | 157 |
| Capener Tufnol Handled Cranked Gouge 8mm Wide 235mm | 157 |
| Capener Tufnol Handled Curved Osteotome 10mm Wide 250mm Curved | 160 |
| Capener Tufnol Handled Curved Osteotome 12mm Wide 250mm Curved | 160 |
| Capener Tufnol Handled Curved Osteotome 15mm Wide 250mm Curved | 160 |
| Capener Tufnol Handled Curved Osteotome 20mm Wide 250mm Curved | 160 |
| Capener Tufnol Handled Curved Osteotome 25mm Wide 250mm Curved | 160 |
| Capener Tufnol Handled Curved Osteotome 30mm Wide 250mm Curved | 160 |
| Capener Tufnol Handled Straight Osteotome 10mm Wide 250mm Straight | 160 |
| Capener Tufnol Handled Straight Osteotome 12mm Wide 250mm Straight | 160 |
| Capener Tufnol Handled Straight Osteotome 15mm Wide 250mm Straight | 160 |
| Capener Tufnol Handled Straight Osteotome 20mm Wide 250mm Straight | 160 |
| Capener Tufnol Handled Straight Osteotome 25mm Wide 250mm Straight | 160 |
| Capener Tufnol Handled Straight Osteotome 30mm Wide 250mm Straight | 160 |
| Capsular Scissors Angled On Flat 250mm Angled | 19 |

Product Name A-Z

| Name | Page |
|--|------|
| Carless Light Weight Scissors 165mm Straight | 16 |
| Carless Light Weight Scissors 165mm Curved | 16 |
| Carless Light Weight Scissors 215mm Curved | 16 |
| Carless Light Weight Scissors 215mm Straight | 16 |
| Caspar Rongeur Straight 150mm Length Shaft 3mm Straight | 170 |
| Caspar Rongeur Straight 2mm Straight | 170 |
| Castroviejo Needle Holder Tungsten Carbide Cross Serrated Curved Jaws With A Spring Gold Plated Handle And Centre Lock With A Screw Joint 140mm | 89 |
| Castroviejo Needle Holder Tungsten Carbide Cross Serrated Straight Jaws With A Spring Gold Plated Handle And Centre Lock With A Screw Joint 140mm Straight | 89 |
| Castroviejo Needle Holder Tungsten Carbide Smooth Curved Jaws With A Spring Gold Plated Handle And Centre Lock With A Screw Joint 140mm | 89 |
| Castroviejo Needle Holder Tungsten Carbide Smooth Straight Jaws With A Spring Gold Plated Handle And Centre Lock With A Screw Joint 140mm | 89 |
| Cawthorne Aural Forceps Crocodile Action With Fine Smooth Jaw And Black Finish 70mm Shoulder Length | 62 |
| Cawthorne Hook No.9 Straight 0.5mm Point 165mm | 124 |
| Cawthorne Hook Single End No.1 Slight Curve 165mm Curved | 123 |
| Cawthorne Hook Single End No.2 Medium Curve 165mm Curved | 123 |
| Cawthorne Hook Single End No.3 Full Curve 165mm Curved | 123 |
| Cawthorne Hook Single End No.4 Double Curve 165mm Curved | 123 |
| Cawthorne Hook Single End No.5 90 Deg Angle 0.5mm Point 165mm Angled | 123 |
| Cawthorne Hook Single End No.6 90 Deg Angle 1.0mm Point 165mm Angled | 123 |
| Cawthorne Hook Single End No.8 45 Deg Angle 0.5mm Point 165mm Angled | 124 |
| Cawthorne Hook Single Ended No.7 Retrograde 165mm Curved | 124 |
| Charnley Initial Incision Retractor Complete With 2 Blades | 108 |

| Name | Page |
|---|------|
| Charnley Pin Retractor | 98 |
| Charnley Suture Dissecting Forceps With Shutter And Slot 145mm | 54 |
| Charnley Suture Dissecting Forceps With Shutter And Slot 180mm 180mm Straight | 54 |
| Cleveland Bone Rongeur Curved With A Box Joint 140mm Curved | 144 |
| Clowards Rongeur 3mm Bite Straight Jaw 150mm Effective Length | 171 |
| Clowards Rongeur 5mm Bite Straight Jaw 150mm Effective Length | 171 |
| Clowards Rongeur 6mm Bite Straight Jaw 150mm Effective Length | 171 |
| Clutton Urethral Bougies (Sounds) Individual Size 10 - 14 Fg 267mm | 188 |
| Clutton Urethral Bougies (Sounds) Individual Size 12 - 16 Fg 267mm | 188 |
| Clutton Urethral Bougies (Sounds) Individual Size 14 - 18 Fg 267mm | 188 |
| Clutton Urethral Bougies (Sounds) Individual Size 16 - 20 Fg 267mm | 188 |
| Clutton Urethral Bougies (Sounds) Individual Size 18 - 22 Fg 267mm | 188 |
| Clutton Urethral Bougies (Sounds) Individual Size 20 - 24 Fg 267mm | 188 |
| Clutton Urethral Bougies (Sounds) Individual Size 22 - 26 Fg 267mm | 188 |
| Clutton Urethral Bougies (Sounds) Individual Size 24 - 28 Fg 267mm | 188 |
| Clutton Urethral Bougies (Sounds) Individual Size 26 - 30 Fg 267mm | 188 |
| Clutton Urethral Bougies (Sounds) Individual Size 28 - 32 Fg 267mm | 188 |
| Clutton Urethral Bougies (Sounds) Individual Size 6 - 10 Fg 267mm | 188 |
| Clutton Urethral Bougies (Sounds) Individual Size 8 - 12 Fg 267mm | 188 |
| Cobb Elevators 10mm 275mm | 150 |
| Cobb Elevators 13mm 275mm | 150 |
| Cobb Elevators 19mm 275mm | 150 |
| College Pattern Forceps For Dental Tweezer 150mm Curved | 54 |
| Collin Duval Tissue Forceps With Box Joint 180mm | 38 |
| Collin Duval Tissue Forceps With Box Joint 200mm | 38 |

Product Name A-Z

| Name | Page |
|--|------|
| Collingwood Stewart Hernia Forceps Ring Bladed Screw Joint 130mm | 38 |
| Compound Action Large Wire Cutter To Cut 5.5mm Dia Max With Tungsten Carbide Jaw | 30 |
| Compound Action Side Wire Cutter Tungsten Carbide Jaws 230mm Straight | 30 |
| Compound Action Top Wire Cutter Tungsten Carbide Jaws 200mm Straight | 29 |
| Compound Action Wire Cutter With Tungsten Carbide Jaws 7inch | 29 |
| Cone Laminectomy Retractor 3 Into 4 Teeth With Hinged Arms And Ratchet Blunt Teeth | 110 |
| Cone Laminectomy Retractor 3 Into 4 Teeth With Hinged Arms And Ratchet Blunt Teeth | 110 |
| Cone Laminectomy Retractor 3 Into 4 Teeth With Hinged Arms And Ratchet Blunt Teeth | 110 |
| Cone Laminectomy Retractor 3 Into 4 Teeth With Hinged Arms And Ratchet Blunt Teeth | 110 |
| Cone Laminectomy Retractor 3 Into 4 Teeth With Hinged Arms And Ratchet Sharp Teeth | 110 |
| Cone Laminectomy Retractor 3 Into 4 Teeth With Hinged Arms And Ratchet Sharp Teeth | 110 |
| Cone Laminectomy Retractor 3 Into 4 Teeth With Hinged Arms And Ratchet Sharp Teeth | 110 |
| Cone Laminectomy Retractor 3 Into 4 Teeth With Hinged Arms And Ratchet Sharp Teeth | 110 |
| Cone Laminectomy Retractor 3 Into 4 Teeth With Hinged Arms And Ratchet Sharp Teeth | 110 |
| Continental Pattern Dissecting Forceps With 1 Into 2 Teeth And Broad End 130mm | 52 |
| Continental Pattern Dissecting Forceps With 1 Into 2 Teeth And Broad End 140mm | 52 |
| Continental Pattern Dissecting Forceps With 1 Into 2 Teeth And Broad End 160mm | 52 |
| Continental Pattern Dissecting Forceps With 1 Into 2 Teeth And Broad End 180mm | 52 |
| Continental Pattern Dissecting Forceps With 1 Into 2 Teeth And Broad End 200mm | 52 |
| Continental Pattern Dissecting Forceps With 1 Into 2 Teeth And Broad End 250mm | 52 |
| Continental Pattern Dissecting Forceps With 1 Into 2 Teeth And Broad End 300mm | 52 |
| Continental Pattern Dissecting Forceps With Serrated Jaws And Narrow End 130mm | 52 |
| Continental Pattern Dissecting Forceps With Serrated Jaws And Narrow End 140mm | 52 |
| Continental Pattern Dissecting Forceps With Serrated Jaws And Narrow End 160mm | 52 |
| Continental Pattern Dissecting Forceps With Serrated Jaws And Narrow End 180mm | 52 |

| Name | Page |
|--|------|
| Continental Pattern Dissecting Forceps With Serrated Jaws And Narrow End 200mm | 52 |
| Continental Pattern Dissecting Forceps With Serrated Jaws And Narrow End 250mm | 52 |
| Cottle Mallet 30mm Diameter Head 190mm | 155 |
| Cottle Nasal Speculum No.1 Blade 20mm. Screw Joint With Side Locking Screw | 130 |
| Cottle Nasal Speculum No.2 Blade 32mm. Screw Joint With Side Locking Screw | 130 |
| Cottle Nasal Speculum No.3 Blade 41mm. Screw Joint With Side Locking Screw | 130 |
| Cottle Nasal Speculum No.4 Blade 60mm. Screw Joint With Side Locking Screw | 130 |
| Coupland Dental Gouge Large (5mm) 170mm | 157 |
| Coupland Dental Gouge Medium (4mm) 170mm | 157 |
| Coupland Dental Gouge Small (3mm) 170mm | 157 |
| Craforde Rib Approximator 260mm | 161 |
| Craft Vasectomy Forceps With Sharp Hooked Jaws 150mm | 40 |
| Crile Artery Forceps With Box Joint 140mm Curved | 42 |
| Crile Artery Forceps With Box Joint 140mm Straight | 42 |
| Crile Artery Forceps With Box Joint 160mm Curved | 42 |
| Crile Artery Forceps With Box Joint 160mm Straight | 42 |
| Crile Artery Forceps With Box Joint 180mm Curved | 42 |
| Crile Artery Forceps With Box Joint 180mm Straight | 42 |
| Crile Rankin Artery Forceps With Box Joint 140mm Curved | 42 |
| Crile Rankin Artery Forceps With Box Joint 140mm Straight | 42 |
| Crile Rankin Artery Forceps With Box Joint 160mm Curved | 42 |
| Crile Rankin Artery Forceps With Box Joint 160mm Straight | 42 |
| Crile Wood Needle Holder With Box Joint (Crilewood) 150mm Straight | 82 |
| Crile Wood Needle Holder With Tungsten Carbide Jaws And Box Joint 140mm | 85 |

Product Name A-Z

| Name | Page |
|---|------|
| Crile Wood Needle Holder With Tungsten Carbide Jaws And Box Joint 150mm | 85 |
| Crile Wood Needle Holder With Tungsten Carbide Jaws And Box Joint 180mm | 85 |
| Crile Wood Needle Holder With Tungsten Carbide Jaws And Box Joint 200mm | 85 |
| Crile Wood Needle Holder With Tungsten Carbide Jaws And Box Joint 230mm | 85 |
| Crile Wood Needle Holder With Tungsten Carbide Jaws And Box Joint 250mm | 85 |
| Crile Wood Needle Holder With Tungsten Carbide Jaws And Box Joint 300mm | 85 |
| Cryers Elevator Left 170mm | 137 |
| Cryers Elevator Right 170mm | 137 |
| Cusco Large 110 x 35mm Smoke Tube Bead | 195 |
| Cusco Large 110 x 35mm Smoke Tube Insulated | 196 |
| Cusco Large 110mm x 35mm | 195 |
| Cusco Large Ex Long 115 x 35mm Smoke Tube Bead | 195 |
| Cusco Large Ex Long 115 x 35mm Smoke Tube Black | 196 |
| Cusco Large Extra Long 115 x 35mm | 195 |
| Cusco Lge Ex Long 115 x 35mm Smoke Tube Insulated | 196 |
| Cusco Medium 102 x 30mm Smoke Tube Bead | 195 |
| Cusco Medium 102 x 30mm Smoke Tube Black | 196 |
| Cusco Medium 102 x 30mm Smoke Tube Insulated | 196 |
| Cusco Medium 102mm x 30mm | 195 |
| Cusco Medium Long 115 x 30mm Smoke Tube Bead | 195 |
| Cusco Medium Long 115 x 30mm Smoke Tube Black | 196 |
| Cusco Medium Long 115 x 30mm Smoke Tube Insulated | 196 |
| Cusco Medium Long 115mm x 30mm | 195 |
| Cusco Small 95 x 25mm Smoke Tube Bead | 195 |
| Cusco Small 95 x 25mm Smoke Tube Black | 196 |
| Cusco Small 95 x 25mm Smoke Tube Insulated | 196 |

| Name | Page |
|---|------|
| Cusco Small 95mm x 25mm | 195 |
| Cusco Virgin 86 x 16mm Smoke Tube Black | 196 |
| Cusco Virgin 86 x 16mm Smoke Tube Insulated | 196 |
| Cusco Virgin 86 x 16mm Smoke Tube Bead | 195 |
| Cusco Virgin 86mm x 16mm | 195 |
| Cushing Artery Forceps With Box Joint 145mm Curved | 43 |
| Cushing Artery Forceps With Box Joint 145mm Straight | 43 |
| Cushing Curved To One Side Artery Forceps With Box Joint 145mm Curved | 43 |
| Cushing Rongeur Angled Upward Jaw 3mm Bite 130mm Effective Length Angled | 170 |
| Cushing Rongeur Angled Upward Jaw 4mm Bite 130mm Effective Length Angled | 170 |
| Cushing Rongeur Angled Upward Jaw 5mm Bite 130mm Effective Length Angled | 170 |
| Cushing Rongeur Angled Upward Jaw 7mm Bite 130mm Effective Length Angled | 170 |
| Cushing Rongeur Straight Jaw 3mm Bite 130mm Effective Length Straight | 170 |
| Cushing Rongeur Straight Jaw 4mm Bite 130mm Effective Length Straight | 170 |
| Cushing Rongeur Straight Jaw 5mm Bite 130mm Effective Length Straight | 170 |
| Cushing Rongeur Straight Jaw 7mm Bite 130mm Effective Length Straight | 170 |
| Czerny Retractor Double Ended 155mm | 101 |
| Daniel Miniature Bone Rongeur With 1.5mm Bite And Curved Handles With A Screw Joint | 144 |
| Daughty Mouth Gag Tongue Plate Only 102mm | 116 |
| Daughty Mouth Gag Tongue Plate Only 114mm | 116 |
| Daughty Mouth Gag Tongue Plate Only 64mm | 116 |
| Daughty Mouth Gag Tongue Plate Only 76mm | 116 |
| Daughty Mouth Gag Tongue Plate Only 89mm | 116 |
| Davis Boyle Mouth Gag Frame Only Adult 150mm x 90mm 150mm | 116 |
| Davis Boyle Mouth Gag Frame Only Child 135mm x 70mm 135mm | 116 |

Product Name A-Z

| Name | Page |
|---|------|
| Davis Boyle Mouth Gag Tongue Plate Only Fig. 1 65mm | 116 |
| Davis Boyle Mouth Gag Tongue Plate Only Fig. 2 75mm | 116 |
| Davis Boyle Mouth Gag Tongue Plate Only Fig. 3 90mm | 116 |
| Davis Boyle Mouth Gag Tongue Plate Only Fig. 4 100mm | 116 |
| Davis Boyle Mouth Gag Tongue Plate Only Fig. 5 115mm | 116 |
| Deaver Retractor (Without Handle) Plain Strip 25mm Blade x 230mm | 101 |
| Deaver Retractor (Without Handle) Plain Strip 75mm Blade x 300mm | 101 |
| Deaver Retractor Plain Handle With 25mm Blade 230mm | 101 |
| Deaver Retractor Plain Handle With 38mm Blade 300mm | 101 |
| Deaver Retractor Plain Handle With 50mm Blade 300mm | 101 |
| Deaver Retractor Plain Handle With 75mm Blade 300mm | 101 |
| Deaver Retractor With Hollow Handle 18mm Blade x 180mm | 101 |
| Deaver Retractor With Hollow Handle 25mm Blade x 250mm | 101 |
| Deaver Retractor Without Handle Plain Strip With 25mm Blade 250mm | 101 |
| Deaver Retractor Without Handle Plain Strip With 38mm Blade 300mm | 101 |
| Deaver Retractor Without Handle Plain Strip With 50mm Blade 300mm | 101 |
| Debakey Dissecting Forceps With Tungsten Carbide Jaws 150mm Straight | 60 |
| Debakey Dissecting Forceps With Tungsten Carbide Jaws 200mm Straight | 60 |
| Debakey Dissecting Forceps With Tungsten Carbide Jaws 230mm Straight | 60 |
| Debakey Forceps For Tissue Dissecting With Atraumatic Jaws 150mm Straight | 55 |
| Debakey Forceps For Tissue Dissecting With Atraumatic Jaws 180mm Straight | 55 |
| Debakey Forceps For Tissue Dissecting With Atraumatic Jaws 230mm Straight | 55 |
| Debakey Needle Holder With Tungsten Carbide Jaws And Box Joint 150mm Straight | 85 |

| Name | Page |
|--|------|
| Debakey Needle Holder With Tungsten Carbide Jaws And Box Joint 180mm Straight | 85 |
| Debakey Needle Holder With Tungsten Carbide Jaws And Box Joint 200mm Straight | 85 |
| Debakey Needle Holder With Tungsten Carbide Jaws And Box Joint 230mm Straight | 85 |
| Debakey Needle Holder With Tungsten Carbide Jaws And Box Joint 250mm Straight | 85 |
| Debakey Needle Holder With Tungsten Carbide Jaws And Box Joint 300mm Straight | 85 |
| Debakey Tissue Dissecting Fcps Atraumatic Jaw 20cm | 55 |
| Denis Browne Intestinal Forceps Spring Pattern With Serrated Grips And Fenestrated Jaws 180mm Straight | 140 |
| Denis Browne Suture Forceps For Collar/ Suture Crushing 145mm Straight | 61 |
| Denis Browne Tonsil Forceps Angled To One Side 200mm (Dennis Brown) | 43 |
| Denis Browne Tonsil Seizing Forceps With Box Joint Tonsil Holding Forceps And Angled Shanks 215mm | 66 |
| Dental Cartridge Syringe 2.2cc Straight | 64 |
| Dental Mirror Head Only 20mm Mirror Diameter | 117 |
| Dental Mirror With 10mm Mirror Diameter | 117 |
| Dental Mirror With 12mm Mirror Diameter | 117 |
| Dental Mirror With 14mm Mirror Diameter | 117 |
| Dental Mirror With 16mm Mirror Diameter | 117 |
| Dental Mirror With 18mm Mirror Diameter | 117 |
| Dental Mirror With 20mm Mirror Diameter | 117 |
| Dental Mirror With 22mm Mirror Diameter | 117 |
| Dental Mirror With 24mm Mirror Diameter | 117 |
| Dental Mirror With 26mm Mirror Diameter | 117 |
| Dental Mirror With 28mm Mirror Diameter | 117 |
| Dental Mirror With 30mm Mirror Diameter | 117 |
| Dental Mirror With 8mm Mirror Diameter | 117 |
| Dental Plier Cutter | 22 |

Product Name A-Z

| Name | Page |
|---|------|
| Derf Needle Holder Tungsten Carbide Jaws And Box Joint 120mm Straight | 85 |
| Derf Needle Holder With Box Joint 120mm Straight | 82 |
| Desjardin Gall Stone Forceps Large Jaws With A Screw Joint 230mm | 190 |
| Desjardin Gall Stone Forceps Small Jaws With A Screw Joint 230mm | 190 |
| Desmarres Retractor Large 11mm x 14mm (Eye Lid) 160mm | 102 |
| Desmarres Retractor Medium 9mm x 12mm (Eye Lid) 160mm | 102 |
| Desmarres Retractor Small 9mm x 9mm (Eye Lid) 160mm | 102 |
| Diathermy Quiver Standard 60mm Dia 230mm | 179 |
| Diathermy Quiver Telescopic 60mm Dia Extending 400mm | 179 |
| Dissecting Forceps With Bayonet Shanks And Serrated Jaws 180mm | 52 |
| Dissecting Forceps With Bayonet Shanks And Serrated Jaws 200mm | 52 |
| Dissecting Forceps With Block End And Serrated Jaws With A Medium Point 130mm | 52 |
| Dissecting Forceps With Block End And Serrated Jaws With A Medium Point 150mm | 52 |
| Dissecting Forceps With Block End And Serrated Jaws With A Medium Point 180mm | 52 |
| Dissecting Forceps With Block End And Serrated Jaws With A Medium Point 200mm | 52 |
| Dissecting Forceps With Block End And Serrated Jaws With A Medium Point 230mm | 52 |
| Dissecting Forceps With Block End And Serrated Jaws With A Medium Point 250mm | 52 |
| Dissecting Forceps With Block End And Serrated Jaws With A Medium Point 300mm | 52 |
| Dissecting Forceps With Turned Over End And Serrated Jaws 130mm | 52 |
| Dissecting Forceps With Turned Over End And Serrated Jaws 150mm | 52 |
| Dissecting Forceps With Turned Over End And Serrated Jaws 180mm | 52 |
| Dissecting Forceps With Turned Over End And Serrated Jaws 200mm | 52 |
| Dissecting Forceps With Turned Over End And Serrated Jaws 230mm | 52 |
| Dissecting Forceps With Turned Over End And Serrated Jaws 250mm | 52 |

| Name | Page |
|--|------|
| Dissecting Forceps With Turned Over End And Serrated Jaws 300mm | 52 |
| Dissecting Scissors With Open Shanks 115mm Straight | 13 |
| Dissecting Scissors With Open Shanks 130mm Straight | 13 |
| Dissecting Scissors With Open Shanks 150mm Straight | 13 |
| Dott Dissecting Forceps With Serrated Jaws 180mm Straight | 55 |
| Double Action Adson Dissecting Fcps 1-2 Teeth 13cm | 53 |
| Double Ended Probe Stainless Steel 130mm Straight | 93 |
| Double Ended Probe Stainless Steel 150mm Straight | 93 |
| Double Ended Probe Stainless Steel 180mm Straight | 93 |
| Double Ended Probe Stainless Steel 200mm Straight | 93 |
| Double Ended Probe Stainless Steel 250mm Straight | 93 |
| Double Ended Spatula And Packer With Solid Handle 155mm Straight | 35 |
| Double Scope Basket 290mm x 160mm x 52mm | 206 |
| Double Scope Basket 460mm x 160mm x 52mm | 206 |
| Dovo German Pattern Nail Nipper Straight 140mm Straight | 34 |
| Doyen Abdominal Retractor Large With 115mm Width Blade 215mm | 102 |
| Doyen Abdominal Retractor Medium 100mm Width Blade 215mm | 102 |
| Doyen Abdominal Retractor Small With 90mm Width Blade 215mm | 102 |
| Doyen Intestinal Clamp Atraumatic Jaw Bj 23cm Cvd | 140 |
| Doyen Intestinal Clamp With 100mm Longitudinal Serrated Blades And Box Joint 230 | 140 |
| Doyen Intestinal Clamp With 100mm Longitudinal Serrated Blades And Box Joint 230 | 140 |
| Doyen Intestinal Clamp With Atraumatic Jaw Blades And Box Joint 180mm Curved | 140 |
| Doyen Intestinal Clamp With Atraumatic Jaw Blades And Box Joint 180mm Straight | 140 |

Product Name A-Z

| Name | Page |
|---|------|
| Doyen Intestinal Clamp With Atraumatic Jaw Blades And Box Joint 230mm Straight | 140 |
| Doyen Intestinal Clamp With Longitudinal Serrated Blades And Box Joint 180mm Cur | 140 |
| Doyen Intestinal Clamp With Longitudinal Serrated Blades And Box Joint 180mm Straight | 140 |
| Doyen Myoma Screw With Large Spiral For Use Per Abdomen With D Ring Handle 150mm | 193 |
| Doyen Myoma Screw With Small Spiral For Use Per Vagina With D Ring Handle 175mm | 193 |
| Doyen Periosteal Elevator Curved To Left Child 165mm Curved | 151 |
| Doyen Periosteal Elevator Curved To Left Adult 165mm Curved | 151 |
| Doyen Periosteal Elevator Curved To Right Child 165mm Curved | 151 |
| Doyen Periosteal Elevator Curved To Right Adult 165mm Curved | 151 |
| Doyen Retractor Large 48mm x 125mm Rectangular Blade 250mm | 112 |
| Doyen Retractor Large Ring Handles 115mm Width Blade 215mm | 102 |
| Doyen Retractor Medium 48mm x 90mm Rectangular Blade 250mm | 112 |
| Doyen Retractor Medium Ring Handles 105mm Width Blade 215mm | 102 |
| Doyen Retractor Small 48mm x 60mm Rectangular Blade 250mm | 112 |
| Doyen Retractor Small Ring Handles 95mm Width Blade 215mm | 102 |
| Doyen Scissors (Uterine) 170mm Curved | 17 |
| Doyen Scissors (Uterine) 200mm Curved | 17 |
| Doyen Scissors (Uterine) 200mm Straight | 17 |
| Doyen Vaginal Speculum Large 115mm x 45mm 240mm | 200 |
| Doyen Vaginal Speculum Medium 85mm x 45mm 240mm | 200 |
| Doyen Vaginal Speculum Small 55mm x 45mm 240mm | 200 |
| Draffin Bipod Fitted With 4 Ring For Varying Heights Adult 480mm | 117 |
| Draffin Bipod Fitted With 4 Ring For Varying Heights Child 375mm | 117 |
| Dressing / Stitch Scissors With Tungsten Carbide Jaws Blunt / Blunt 150mm Straight | 20 |

| Name | Page |
|--|------|
| Dressing / Stitch Scissors With Tungsten Carbide Jaws Blunt / Blunt 180mm Straight | 20 |
| Dressing / Stitch Scissors With Tungsten Carbide Jaws Sharp / Blunt 150mm Curved | 20 |
| Dressing / Stitch Scissors With Tungsten Carbide Jaws Sharp / Blunt 150mm Straight | 20 |
| Dressing / Stitch Scissors With Tungsten Carbide Jaws Sharp / Blunt 180mm Curved | 20 |
| Dressing / Stitch Scissors With Tungsten Carbide Jaws Sharp / Blunt 180mm Straight | 20 |
| Dressing / Stitch Scissors With Tungsten Carbide Jaws Sharp / Sharp 150mm Curved | 20 |
| Dressing / Stitch Scissors With Tungsten Carbide Jaws Sharp / Sharp 150mm Straight | 20 |
| Dressing / Stitch Scissors With Tungsten Carbide Jaws Sharp / Sharp 180mm Curved | 20 |
| Dressing / Stitch Scissors With Tungsten Carbide Jaws Sharp / Sharp 180mm Straight | 20 |
| Dressing And Stitch Scissors Sharp / Sharp 100mm Curved | 12 |
| Dressing And Stitch Scissors Sharp / Sharp 100mm Straight | 12 |
| Dressing And Stitch Scissors Sharp / Sharp 115mm Curved | 12 |
| Dressing And Stitch Scissors Sharp / Sharp 115mm Straight | 12 |
| Dressing And Stitch Scissors Sharp / Sharp 130mm Curved | 12 |
| Dressing And Stitch Scissors Sharp / Sharp 130mm Straight | 12 |
| Dressing And Stitch Scissors Sharp / Sharp 140mm Curved | 12 |
| Dressing And Stitch Scissors Sharp / Sharp 140mm Straight | 12 |
| Dressing And Stitch Scissors Sharp / Sharp 150mm Curved | 12 |
| Dressing And Stitch Scissors Sharp / Sharp 150mm Straight | 12 |
| Dressing And Stitch Scissors Sharp / Sharp 180mm Curved | 12 |
| Dressing And Stitch Scissors Sharp / Sharp 180mm Straight | 12 |
| Dressing Scissors Blunt / Blunt 100mm Curved | 12 |
| Dressing Scissors Blunt / Blunt 100mm Straight | 12 |

Product Name A-Z

| Name | Page |
|--|------|
| Dressing Scissors Blunt / Blunt 115mm Curved | 12 |
| Dressing Scissors Blunt / Blunt 115mm Straight | 12 |
| Dressing Scissors Blunt / Blunt 130mm Curved | 12 |
| Dressing Scissors Blunt / Blunt 130mm Straight | 12 |
| Dressing Scissors Blunt / Blunt 140mm Curved | 12 |
| Dressing Scissors Blunt / Blunt 140mm Straight | 12 |
| Dressing Scissors Blunt / Blunt 150mm Curved | 12 |
| Dressing Scissors Blunt / Blunt 150mm Straight | 12 |
| Dressing Scissors Blunt / Blunt 180mm Curved | 12 |
| Dressing Scissors Blunt / Blunt 180mm Straight | 12 |
| Dressing Scissors Sharp / Blunt 100mm Curved | 12 |
| Dressing Scissors Sharp / Blunt 100mm Straight | 12 |
| Dressing Scissors Sharp / Blunt 115mm Curved | 12 |
| Dressing Scissors Sharp / Blunt 115mm Straight | 12 |
| Dressing Scissors Sharp / Blunt 130mm Curved | 12 |
| Dressing Scissors Sharp / Blunt 130mm Straight | 12 |
| Dressing Scissors Sharp / Blunt 140mm Curved | 12 |
| Dressing Scissors Sharp / Blunt 140mm Straight | 12 |
| Dressing Scissors Sharp / Blunt 150mm Curved | 12 |
| Dressing Scissors Sharp / Blunt 150mm Straight | 12 |
| Dressing Scissors Sharp / Blunt 180mm Curved | 12 |
| Dressing Scissors Sharp / Blunt 180mm Straight | 12 |
| Dunhill Artery Forceps With Box Joint 130mm Curved | 43 |
| Dunhill Artery Forceps With Box Joint 130mm Straight | 43 |

| Name | Page |
|--|------|
| Dunhill Artery Forceps With Box Joint 180mm Curved | 43 |
| Duray French Osteotome 11mm Wide 140mm | 158 |
| Duray French Osteotome 5mm Wide 140mm | 158 |
| Duray French Osteotome 7mm Wide 140mm | 158 |
| Duray French Osteotome 9mm Wide 140mm | 158 |
| Duval Lung Forceps With 13mm Wide Jaws And A Box Joint 180mm | 67 |
| Duval Lung Forceps With 13mm Wide Jaws And A Box Joint 200mm | 67 |
| Duval Lung Forceps With 14mm Wide Jaws Screw Joint 230mm | 67 |
| Duval Lung Forceps With 19mm Wide Jaws And A Box Joint 230mm | 67 |
| Duval Lung Forceps With 19mm Wide Jaws Screw Joint 200mm | 67 |
| Duval Lung Forceps With 25mm Wide Jaws Screw Joint 200mm | 67 |
| Duval Lung Forceps With 25mm Wide Jaws Screw Joint 230mm | 67 |
| Eisenhammer Rectal Speculum | 186 |
| Emmat Dissecting Forceps With Serrated Jaws 200mm Straight | 55 |
| Emmat Dissecting Forceps With Serrated Jaws 240mm Straight | 55 |
| Endoscope / Arthroscope Basket 670mm x 80mm x 52mm | 207 |
| Endoscope / Arthroscope Basket Double Scope Inserts 290mm x 80mm x 52mm | 206 |
| Endoscope / Arthroscope Basket Double Scope Inserts 460mm x 80mm x 52mm | 206 |
| Endoscope / Arthroscope Basket Single Scope Inserts 290mm x 80mm x 52mm | 206 |
| Endoscope / Arthroscope Basket Single Scope Inserts 460mm x 80mm x 52mm | 206 |
| Endoscope / Arthroscope Basket With Small Lens Box 640mm x 150mm x 77mm | 206 |
| Eppendorf Biopsy Punch For Gynaecology Crocodile Action 200mm Effective Length Complete Straight Cut 4 x 8mm 200mm | 203 |
| Eve Nasal Tonsil Snare Sliding Type 280mm Straight | 131 |

Product Name A-Z

| Name | Page |
|--|------|
| Farabeuf Bone Ruge With Chisel Edge With A Flat Handle (Elevator) 150mm Straight | 152 |
| Fergusson Bone Holding Forceps With Screw Joint 150mm | N/A |
| Fergusson Bone Holding Forceps With Screw Joint 200mm | N/A |
| Fergusson Scissors (Uterine) 210mm Straight | 17 |
| Fergusson Scissors (Uterine) Angled To Flat 210mm Angled | 18 |
| Fickling Artery Forceps Angled 1 Into 2 Teeth And Box Joint 180mm Angled | 44 |
| Fickling Artery Forceps Angled And Plain With Box Joint 180mm Angled | 44 |
| Fickling Mallet 200mm | 155 |
| Fickling Nail Elevator 115mm Straight | 35 |
| Fine Mesh Basket With Take Off Lockable Lid 120mm x 155mm x 40mm | 210 |
| Fine Mesh Basket With Take Off Lockable Lid 120mm x 155mm x 50mm | 210 |
| Fine Mesh Basket With Take Off Lockable Lid 235mm x 155mm x 40mm | 210 |
| Fine Mesh Basket With Take Off Lockable Lid 235mm x 155mm x 50mm | 210 |
| Fine Mesh Basket With Take Off Lockable Lid 380mm x 265mm x 45mm | 210 |
| Fine Mesh Basket With Take Off Lockable Lid 380mm x 265mm x 55mm | 210 |
| Finochietto Rib Spreader Baby Size With 18mm x 21mm Blades And A 75mm Spread | 163 |
| Finochietto Rib Spreader Large Size With 70mm x 65mm Blades And A 250mm Spread | 163 |
| Finochietto Rib Spreader Medium Size With 45mm x 65mm Blades And A 200mm Spread | 163 |
| Finochietto Rib Spreader Small Size With 30mm x 45mm Blades And A 165mm Spread | 163 |
| Foley Catheter Introducer Curved With A 2mm Diameter Tip 430mm | 19 |
| Foley Catheter Introducer Curved With A 4mm Diameter Tip 430mm | 19 |
| Fomon Bone Rasp No. 1 & 2 Tc Cut 205mm | 186 |
| Fomon Bone Rasp No. 3 & 4 Tc Cut 205mm | 186 |
| Fomon Bone Rasp No. 5 & 6 Tc Cut 205mm | 186 |

| Name | Page |
|--|------|
| Fomon Bone Rasp No. 7 & 8 Tc Cut 205mm | 186 |
| Fomon Bone Rasp No. 9 & 10 Tc Cut 205mm | 186 |
| Fomon Rasp Fine 215mm | 186 |
| Formby Cerumen Scoop And Hook Double Ended 180mm | 125 |
| Foster Gillies Needle Holder With Tungsten Carbide Jaws And Box Joint 135mm Angl | 86 |
| Foster Needle Holder With Scissor And Screw Joint 135mm Curved | 82 |
| Frazier Small Skin Hook 180mm | 95 |
| Freer Double Ended Separator And Elevator 200mm | 133 |
| Freidman Bone Nibbler Single Action Small (Friedman) 125mm Straight | 148 |
| Fry Needle Holder With Box Joint 115mm Straight | 82 |
| Fry Needle Holder With Tungsten Carbide Jaws And Box Joint 130mm Straight | 86 |
| Galabin Uterine Sound Malleable Cm Graduated | 194 |
| Gallipots Ss (60mm x 50mm) | 211 |
| Gallipots Ss (80mm x 60mm) Non Rolled Lips | 211 |
| Gelpi Vulva Self Retaining Retractor Large Size With Ratchet 200mm | 110 |
| Gelpi Vulva Self Retaining Retractor Small Size With Ratchet 170mm | 110 |
| Gigli Metal Handle (Supplied In A Pair) | 204 |
| Gigli Saw 30cm Length | 204 |
| Gigli Saw 40cm Length 40cm | 204 |
| Gigli Saw 50cm Length 50cm | 204 |
| Gigli Saw 60cm Length | 204 |
| Gillies (French) Osteotome 10mm Wide 140mm | 159 |
| Gillies (French) Osteotome 11mm Wide 140mm | 159 |
| Gillies (French) Osteotome 2mm Wide 140mm | 159 |
| Gillies (French) Osteotome 3mm Wide 140mm | 159 |
| Gillies (French) Osteotome 4mm Wide 140mm | 159 |

Product Name A-Z

| Name | Page |
|--|------|
| Gillies (French) Osteotome 5mm Wide 140mm | 159 |
| Gillies (French) Osteotome 6mm Wide 140mm | 159 |
| Gillies (French) Osteotome 7mm Wide 140mm | 159 |
| Gillies (French) Osteotome 8mm Wide 140mm | 159 |
| Gillies (French) Osteotome 9mm Wide 140mm | 159 |
| Gillies Dissecting Forceps With 1 Into 2 Teeth 150mm Straight | 56 |
| Gillies Dissecting Forceps With Serrated Jaws 150mm Straight | 56 |
| Gillies Large Skin Hook With Fine Point 155mm | 95 |
| Gillies Medium Skin Hook With Fine Point 155mm | 95 |
| Gillies Needle Holder Tungsten Carbide Jaws And Box Joint With Right Hand Scissors | 86 |
| Gillies Needle Holder With Right Hand Scissor Screw Joint 160mm Curved | 82 |
| Gillies Small Skin Hook With Fine Point 155mm | 95 |
| Glasgow Pattern Bone Rongeur With 4mm Bite And Slightly Curved On Flat | 144 |
| Goligher Rectal Speculum With Third Blade 76mm | 186 |
| Goligher Retractor Blades (2 Centre Blades And 2 Pairs Lateral Blades) | 109 |
| Gorney Freeman Facelift Scissors Tc Supercut | 25 |
| Gorney Freeman Facelift Scissors Tc Supercut | 25 |
| Graves (Brewer) Vaginal Speculum Large 115mm x 32mm | 197 |
| Graves (Brewer) Vaginal Speculum Medium 95mm x 35mm | 197 |
| Graves (Brewer) Vaginal Speculum Small 75mm x 20mm | 197 |
| Green Armytage Uterine Dressing Forceps Angled Ends With Wide Serrated Jaws | 73 |
| Green Armytage Uterine Dressing Forceps With Wide Serrated Jaws And A Screw Join | 73 |
| Grey Turner Artery Forceps With Box Joint 140mm Curved | 44 |
| Grey Turner Artery Forceps With Box Joint 145mm Straight | 44 |

| Name | Page |
|---|------|
| Grey Turner Artery Forceps With Box Joint 180mm Curved | 44 |
| Grey Turner Artery Forceps With Box Joint 180mm Straight | 44 |
| Gruber Aural Speculum 4mm Dia Oval | 126 |
| Gruber Aural Speculum 4mm Dia Round 38mm | 126 |
| Gruber Aural Speculum 5mm Dia Oval | 126 |
| Gruber Aural Speculum 5mm Dia Round 38mm | 126 |
| Gruber Aural Speculum 6mm Dia Oval | 126 |
| Gruber Aural Speculum 6mm Dia Round 38mm | 126 |
| Gruber Aural Speculum 7mm Dia Oval | 126 |
| Gruber Aural Speculum 7mm Dia Round 38mm | 126 |
| Guyon Catheter Introducer Curved With A 2mm Diameter Tip 430mm | 19 |
| Guyon Catheter Introducer Curved With A 4mm Diameter Tip 430mm | 19 |
| Guyon Catheter Introducer Straight With A 2mm Diameter Tip 430mm | 19 |
| Guyon Catheter Introducer Straight With A 4mm Diameter Tip 430mm | 19 |
| Guys Galabin / Suture Forceps 140mm Angled (Gallabin) | 61 |
| Guys Heavy Pattern Plaster Shears 300mm | 166 |
| Gwilliam Hysterectomy Clamp Curved Longitudinal Serrated Jaws With 1 Into 2 Teeth | 193 |
| Gwilliam Hysterectomy Clamp Straight Longitudinal Serrated Jaws With 1 Into 2 Teeth | 192 |
| Gwynne Evans Double Ended 6mm Wide Dissector 190mm | 136 |
| Halsey Needle Holder Tungsten Carbide Jaws And Box Joint 130mm Straight | 86 |
| Halsey Needle Holder With Plain Jaws And Box Joint 130mm Straight | 83 |
| Halsey Needle Holder With Serrated Jaws And Box Joint 130mm Straight | 83 |
| Halstead Mosquito Artery Forceps With Box Joint 200mm Straight | |
| Halsted Mosquito Artery Forceps 1 Into 2 Teeth With Box Joint 130mm Curved | 46 |

Product Name A-Z

| Name | Page |
|---|------|
| Halsted Mosquito Artery Forceps 1 Into 2 Teeth With Box Joint 130mm Straight | 45 |
| Halsted Mosquito Artery Forceps Curved With Atraumatic Jaws 130mm | 46 |
| Halsted Mosquito Artery Forceps With Box Joint 130mm Curved | 45 |
| Halsted Mosquito Artery Forceps With Box Joint 130mm Straight | 45 |
| Hartmann Artery Forceps Extra Fine With Box Joint 100mm Straight | 45 |
| Hartmann Artery Forceps Extra Fine With Box Joint 90mm Curved | 44 |
| Hartmann Artery Forceps Extra Fine With Box Joint 90mm Straight | 44 |
| Hartmann Aural Forceps Crocodile Action With 1 x 2 Teeth And 85mm Shoulder Length | 64 |
| Hartmann Aural Forceps Crocodile Action With A Extra Fine Serrated Jaw And A Bla | 62 |
| Hartmann Aural Forceps Crocodile Action With A Serrated Jaws And A 85mm Shoulder | 63 |
| Hawkin Ambler Uterine Dilator Conical Size 10 13mm Diameter Single Ended 170mm | 191 |
| Hawkin Ambler Uterine Dilator Conical Size 11 14mm Diameter Single Ended 170mm | 191 |
| Hawkin Ambler Uterine Dilator Conical Size 13 16mm Diameter Single Ended 170mm | 191 |
| Hawkin Ambler Uterine Dilator Conical Size 14 17mm Diameter Single Ended 170mm | 191 |
| Hawkin Ambler Uterine Dilator Conical Size 15 18mm Diameter Single Ended 170mm | 191 |
| Hawkin Ambler Uterine Dilator Conical Size 16 19mm Diameter Single Ended 170mm | 191 |
| Hawkin Ambler Uterine Dilator Conical Size 17 20mm Diameter Single Ended 170mm | 191 |
| Hawkin Ambler Uterine Dilator Conical Size 18 21mm Diameter Single Ended 170mm | 191 |
| Hawkin Ambler Uterine Dilator Conical Size 3 6mm Diameter Single Ended 170mm | 191 |
| Hawkin Ambler Uterine Dilator Conical Size 4 7mm Diameter Single Ended 170mm | 191 |
| Hawkin Ambler Uterine Dilator Conical Size 5 8mm Diameter Single Ended 170mm | 191 |
| Hawkin Ambler Uterine Dilator Conical Size 6 9mm Diameter Single Ended 170mm | 191 |
| Hawkin Ambler Uterine Dilator Conical Size 7 10mm Diameter Single Ended 170mm | 191 |

| Name | Page |
|---|------|
| Hawkin Ambler Uterine Dilator Conical Size 8 11mm Diameter Single Ended 170mm | 191 |
| Hawkin Ambler Uterine Dilator Conical Size 9 12mm Diameter Single Ended 170mm | 191 |
| Heaney Needle Holder Tungsten Carbide Jaws And Box Joint 200mm Curved | 86 |
| Heaney Needle Holder Tungsten Carbide Jaws And Box Joint 250mm Curved | 86 |
| Heath Bone Mallet 18 Oz | 154 |
| Heath Bone Mallet 28 Oz | 154 |
| Heath Bone Mallet 35mm Dia Head & Long Handle | 154 |
| Heath Bone Mallet 36 Oz | 154 |
| Heavy Tufnol Mallet 270mm | 154 |
| Hegar Uterine Dilator 1.5mm Diameter Single Ended 185mm | 190 |
| Hegar Uterine Dilator 10.5mm Diameter Single Ended 185mm | 190 |
| Hegar Uterine Dilator 10mm Diameter Single Ended 185mm | 190 |
| Hegar Uterine Dilator 11.5mm Diameter Single Ended 185mm | 190 |
| Hegar Uterine Dilator 11mm Diameter Single Ended 185mm | 190 |
| Hegar Uterine Dilator 12.5mm Diameter Single Ended 185mm | 190 |
| Hegar Uterine Dilator 12mm Diameter Single Ended 185mm | 190 |
| Hegar Uterine Dilator 13.5mm Diameter Single Ended 185mm | 190 |
| Hegar Uterine Dilator 13mm Diameter Single Ended 185mm | 190 |
| Hegar Uterine Dilator 14.5mm Diameter Single Ended 185mm | 190 |
| Hegar Uterine Dilator 14mm Diameter Single Ended 185mm | 190 |
| Hegar Uterine Dilator 15mm Diameter Single Ended 185mm | 190 |
| Hegar Uterine Dilator 16mm Diameter Single Ended 185mm | 190 |
| Hegar Uterine Dilator 17mm Diameter Single Ended 185mm | 190 |
| Hegar Uterine Dilator 18mm Diameter Single Ended 185mm | 190 |
| Hegar Uterine Dilator 19mm Diameter Single Ended 185mm | 190 |

Product Name A-Z

| Name | Page |
|---|------|
| Hegar Uterine Dilator 1mm Diameter Single Ended 185mm | 190 |
| Hegar Uterine Dilator 20mm Diameter Single Ended 185mm | 190 |
| Hegar Uterine Dilator 21mm Diameter Single Ended 185mm | 190 |
| Hegar Uterine Dilator 22mm Diameter Single Ended 185mm | 190 |
| Hegar Uterine Dilator 23mm Diameter Single Ended 185mm | 190 |
| Hegar Uterine Dilator 24mm Diameter Single Ended 185mm | 190 |
| Hegar Uterine Dilator 25mm Diameter Single Ended 185mm | 190 |
| Hegar Uterine Dilator 26mm Diameter Single Ended 185mm | 190 |
| Hegar Uterine Dilator 2mm Diameter Single Ended 185mm | 190 |
| Hegar Uterine Dilator 3.5mm Diameter Single Ended 185mm | 190 |
| Hegar Uterine Dilator 3mm Diameter Single Ended 185mm | 190 |
| Hegar Uterine Dilator 4.5mm Diameter Single Ended 185mm | 190 |
| Hegar Uterine Dilator 4mm Diameter Single Ended 185mm | 190 |
| Hegar Uterine Dilator 5.5mm Diameter Single Ended 185mm | 190 |
| Hegar Uterine Dilator 5mm Diameter Single Ended 185mm | 190 |
| Hegar Uterine Dilator 6.5mm Diameter Single Ended 185mm | 190 |
| Hegar Uterine Dilator 6mm Diameter Single Ended 185mm | 190 |
| Hegar Uterine Dilator 7.5mm Diameter Single Ended 185mm | 190 |
| Hegar Uterine Dilator 7mm Diameter Single Ended 185mm | 190 |
| Hegar Uterine Dilator 8.5mm Diameter Single Ended 185mm | 190 |
| Hegar Uterine Dilator 8mm Diameter Single Ended 185mm | 190 |
| Hegar Uterine Dilator 9.5mm Diameter Single Ended 185mm | 190 |
| Hegar Uterine Dilator 9mm Diameter Single Ended 185mm | 190 |
| Hegar Uterine Dilator Double Ended Sizes 11mm And 12mm Diameter 220mm | 191 |

| Name | Page |
|--|------|
| Hegar Uterine Dilator Double Ended Sizes 13mm And 14mm Diameter 220mm | 191 |
| Hegar Uterine Dilator Double Ended Sizes 15mm And 16mm Diameter 220mm | 191 |
| Hegar Uterine Dilator Double Ended Sizes 17mm And 18mm Diameter 220mm | 191 |
| Hegar Uterine Dilator Double Ended Sizes 19mm And 20mm Diameter 220mm | 191 |
| Hegar Uterine Dilator Double Ended Sizes 1mm And 2mm Diameter 220mm | 191 |
| Hegar Uterine Dilator Double Ended Sizes 21mm And 22mm Diameter 220mm | 191 |
| Hegar Uterine Dilator Double Ended Sizes 23mm And 24mm Diameter 220mm | 191 |
| Hegar Uterine Dilator Double Ended Sizes 25mm And 26mm Diameter 220mm | 191 |
| Hegar Uterine Dilator Double Ended Sizes 3mm And 4mm Diameter 220mm | 191 |
| Hegar Uterine Dilator Double Ended Sizes 5mm And 6mm Diameter 220mm | 191 |
| Hegar Uterine Dilator Double Ended Sizes 7mm And 8mm Diameter 220mm | 191 |
| Hegar Uterine Dilator Double Ended Sizes 9mm And 10mm Diameter 220mm | 191 |
| Hegar Uterine Dilator S/E 185mm 15.5mm Diameter | 190 |
| Hegar Uterine Dilator S/E 185mm 16.5mm Diameter | 190 |
| Hegar Uterine Dilator S/E 185mm 17.5mm Diameter | 190 |
| Hegar Uterine Dilator S/E 185mm 18.5mm Diameter | 190 |
| Hegar Uterine Dilator S/E 185mm 19.5mm Diameter | 190 |
| Hegar Uterine Dilator S/E 185mm 2.5mm Diameter | 190 |
| Hegar Uterine Dilators, Set Of 12 Double Ended Dilators | 191 |
| Hegar Uterine Dilators. Set Of 8 Double Ended Dilators | 191 |
| Heiss Full Curve Artery Forceps With Box Joint 200mm | 46 |
| Heiss Half Curve Artery Forceps With Box Joint 200mm | 46 |
| Henckel Aural Forceps Crocodile Action With Fine Fenestrated Blades And 75mm Shoulder Length | 64 |

Product Name A-Z

| Name | Page |
|--|------|
| Henderson Dissector (Periosteal Elevator) 210mm | 138 |
| Henning Plaster Cast Spreaders With A Screw Joint 280mm | 167 |
| Hill Double Ended Elevator Curved To Side With Rugine End 5mm Wide (Hills) 220mm | 136 |
| Hohman Bone Lever 43mm Wide 240mm | 99 |
| Hohman Bone Lever 70mm Wide 240mm | 99 |
| Hohmann Muller Spike 18mm | 99 |
| Hohmann Muller Spike 8mm | 99 |
| Horsley Bone Cutting Forceps With Compound Action 270mm Curved | 68 |
| House Aural Curette Double Ended Angled 1.0mm And 1.2mm Spoon 180mm Angled | 122 |
| House Aural Curette Double Ended Angled 1.5mm And 1.8mm Spoon 180mm Angled | 122 |
| House Aural Curette Double Ended Straight 1.0mm And 1.2mm Spoon 180mm Straight | 122 |
| House Aural Curette Double Ended Straight 1.5mm And 1.8mm Spoon 180mm Straight | 122 |
| Howarth Double Ended Elevator Curved On Flat With Rugine End 220mm Curved | 136 |
| Hurd Double Ended Pillar Retractor And Tonsil Dissector 200mm | 133 |
| Ingrown Toe Nail Nipper Sharp Pointed Straight 130mm Straight | 34 |
| Iris Dissecting Forceps With 1 Into 2 Teeth 100mm Straight | 56 |
| Iris Dissecting Forceps With 1 Into 2 Teeth 115mm Straight | 56 |
| Iris Dissecting Forceps With Serrated Jaws 100mm Curved | 56 |
| Iris Dissecting Forceps With Serrated Jaws 100mm Str | 56 |
| Iris Dissecting Forceps With Serrated Jaws 115mm Curved | 56 |
| Iris Dissecting Forceps With Serrated Jaws 115mm Straight | 56 |
| Iris Scissors Fine Pointed 100mm Curved | 13 |
| Iris Scissors Fine Pointed 100mm Straight | 13 |
| Iris Scissors Fine Pointed 115mm Curved | 13 |
| Iris Scissors Fine Pointed 115mm Straight | 13 |
| Iris Scissors Fine Pointed 130mm Curved | 13 |

| Name | Page |
|---|------|
| Iris Scissors Fine Pointed 130mm Straight | 13 |
| Iris Scissors With Tungsten Carbide Jaws Fine Pointed 115mm Curved | 21 |
| Iris Scissors With Tungsten Carbide Jaws Fine Pointed 115mm Straight | 20 |
| Iris Supercut Scissors Fine Pointed 115mm Curved | 23 |
| Iris Supercut Scissors Fine Pointed 115mm Straight | 23 |
| Irwin Moore Nasal Turbinate Forceps With Screw Joint And Angled Shanks 190mm | 65 |
| Irwin Retractor With 4 Into 4 Teeth And Cam Rack 180mm | 110 |
| Jackson Burrows Dissecting Forceps Serrated With 1 Into 2 Teeth 140mm Straight | 56 |
| Jackson Burrows Self Retaining Blunt Hook Jaws For Superficial Part Of The Knee 175mm | 110 |
| Jameson Reynolds Supercut Scissors With Black Handles 160mm Curved | 24 |
| Jansen Middleton Bone Rongeur With Compound Action 185mm Curved | 144 |
| Jansen Middleton Bone Rongeur With Compound Action 185mm Straight | 144 |
| Jansen Middleton Nasal Turbinate Septum Forceps Double Action 190mm | 66 |
| Jayes Osteotome 1mm Mini Lambotte 125mm | 159 |
| Jayes Osteotome 2mm Mini Lambotte 125mm | 159 |
| Jayes Osteotome 4mm Mini Lambotte 125mm | 159 |
| Jean Dissecting Forceps With 1 Into 2 Teeth 140mm Straight | 57 |
| Jean Dissecting Forceps With 1 Into 2 Teeth 180mm Straight | 57 |
| Jefferson Dissecting Forceps With 1 Into 2 Teeth 180mm Straight | 57 |
| Jenkins Bone Gouge 10mm Wide 150mm | 157 |
| Jenkins Bone Gouge 12mm Wide 150mm | 157 |
| Jenkins Bone Gouge 14mm Wide 150mm | 157 |
| Jenkins Bone Gouge 3mm Wide 150mm | 157 |
| Jenkins Bone Gouge 4mm Wide 150mm | 157 |
| Jenkins Bone Gouge 5mm Wide 150mm | 157 |

Product Name A-Z

| Name | Page |
|---|------|
| Jenkins Bone Gouge 6mm Wide 150mm | 157 |
| Jenkins Bone Gouge 8mm Wide 150mm | 157 |
| Jessop Hospital Pattern Vaginal Speculum Double Ended 100mm x 28mm And 120mm x 30mm | 200 |
| Jobson Horne Probe With Ring 140mm Straight | 93 |
| Jobson Horne Probe With Ring 180mm Straight | 93 |
| Joll Thyroid Retractor With 2 Self Retaining 2 Pincher Blades And Fluted Handle | 109 |
| Jolls Artery Forceps With Box Joint 145mm Curved | 47 |
| Jolls Artery Forceps With Box Joint 145mm Straight | 47 |
| Jones Rugine Pistol Shaped Handle Fan Shaped End And Convex Sharp Edge 185mm | 152 |
| Joseph Nasal Saw Vertical Shaft And Bayonet Crank Left 185mm | 133 |
| Joseph Nasal Saw Vertical Shaft And Bayonet Crank Right 185mm | 133 |
| Judd Allis Tissue Forceps With 3 Into 4 Teeth And Box Joint 150mm | 39 |
| Judd Allis Tissue Forceps With 3 Into 4 Teeth And Box Joint 190mm | 39 |
| Kelly Artery Forceps With Box Joint 140mm Curved | 47 |
| Kelly Artery Forceps With Box Joint 140mm Straight | 47 |
| Kelly Artery Forceps With Box Joint 160mm Curved | 47 |
| Kelly Artery Forceps With Box Joint 160mm Straight | 47 |
| Kelly Frazier Artery Forceps With Box Joint 180mm Curved | 47 |
| Kelly Rankin Artery Forceps With Box Joint 160mm Curved | 47 |
| Kelly Retractor With Large Size 57mm End Of Blade 190mm | 103 |
| Kelly Retractor With Small Size 38mm End Of Blade 190mm | 103 |
| Kevorkian Biopsy Punch For Gynaecology Crocodile Action 260mm Effective Length Complete Straight Cut 3 x 10mm 260mm | 203 |
| Kielland Obstetric Forceps For Occipito-Posterior Position 410mm | 183 |
| Killian Nasal Speculum Blade 38mm Length. Screw Joint And Side Locking Screw | 130 |

| Name | Page |
|---|------|
| Killian Nasal Speculum Blade 64mm Length. Screw Joint And Side Locking Screw | 130 |
| Killian Nasal Speculum Blade 76mm Length. Screw Joint And Side Locking Screw | 130 |
| Killian Nasal Speculum Blade 89mm Length. Screw Joint And Side Locking Screw | 130 |
| Kilner (Alae) Retractor Double Pronged 10mm Wide 85mm | 103 |
| Kilner (Alae) Retractor Double Pronged 13mm Wide 85mm | 103 |
| Kilner (Alae) Retractor Double Pronged 8mm Wide 85mm | 103 |
| Kilner (Cat Paw) Skin Retractor Double Ended Std Size 150mm | 103 |
| Kilner (Cheek) Large Retractor Double Ended Starting 35mm Tapering Down To 30mm | 103 |
| Kilner (Cheek) Small Retractor Double Ended Starting 25mm Tapering Down To 20mm | 103 |
| Kilner Artery Forceps Curved 145mm 145mm | 51 |
| Kilner Artery Forceps Straight 145mm | 51 |
| Kilner Detachable Skin Hook And Handle 155mm | 95 |
| Kilner Needle Holder Tungsten Carbide Jaws And Box Joint With Curved Shanks 135m | 86 |
| Kilner Needle Holder With Box Joint With Curved Shanks 135mm Straight | 83 |
| Kilner Scissors 115mm Curved | 15 |
| Kilner Scissors 115mm Straight | 14 |
| Kilner Scissors With Tungsten Carbide Jaws 115mm Curved | 21 |
| Kilner Scissors With Tungsten Carbide Jaws 115mm Straight | 21 |
| Kilner Supercut Scissors 115mm Curved | 23 |
| Kilner Supercut Scissors 115mm Straight | 23 |
| Kocher Artery Forceps With Box Joint 1 Into 2 Teeth 130mm Curved | 48 |
| Kocher Artery Forceps With Box Joint 1 Into 2 Teeth 140mm Curved | 48 |
| Kocher Artery Forceps With Box Joint 1 Into 2 Teeth 140mm Straight | 48 |
| Kocher Artery Forceps With Box Joint 1 Into 2 Teeth 150mm Curved | 48 |

Product Name A-Z

| Name | Page |
|---|------|
| Kocher Artery Forceps With Box Joint 1 Into 2 Teeth 150mm Straight | 48 |
| Kocher Artery Forceps With Box Joint 1 Into 2 Teeth 180mm Curved | 48 |
| Kocher Artery Forceps With Box Joint 1 Into 2 Teeth 180mm Straight | 48 |
| Kocher Artery Forceps With Box Joint 1 Into 2 Teeth 200mm Curved | 48 |
| Kocher Artery Forceps With Box Joint 1 Into 2 Teeth 200mm Straight | 48 |
| Kocher Artery Forceps With Box Joint 1 Into 2 Teeth 230mm Curved | 48 |
| Kocher Artery Forceps With Box Joint 1 Into 2 Teeth 230mm Straight | 48 |
| Kocher Bone Hook Sharp Pointed With A Fenestrated Handle 250mm | 153 |
| Kocher Intestinal Clamp With 100mm Longitudinal Serrated Light Spring Blades And Box Joint 230mm Straight | 140 |
| Kocher Intestinal Clamp With 125mm Longitudinal Serrated Light Spring Blades And Box Joint 250mm Straight | 140 |
| Kocher Intestinal Clamp With 150mm Longitudinal Serrated Light Spring Blades And Box Joint 280mm Straight | 140 |
| Kocher Retractor 40 x 10mm 215mm | 103 |
| Kocher Retractor 40 x 15mm 215mm | 103 |
| Kocher Retractor 40 x 20mm 215mm | 103 |
| Kocher Retractor Blade 45mm Wide x 85mm Deep 250mm | 113 |
| Kocher Retractor Blade 54mm Wide x 85mm Deep 250mm | 113 |
| Kocher Retractor Blade 62mm Wide x 85mm Deep 250mm | 113 |
| Kogan Endocervical Speculum Without Ratchet (Koogan) 240mm | 197 |
| Koogan Endocervical Speculum & Adjustable Ratchet | 197 |
| Kurihara Endo Cervical Dissecting Forceps With Satin Finish 250mm Straight | 57 |
| Lack Tongue Depressor Adult Size 22mm Wide | 116 |
| Lack Tongue Depressor Child Size 19mm Wide | 116 |
| Lack Tongue Depressor Infant Size 12mm Wide | 116 |
| Lahey Cholecystectomy Curved Forceps With Longitudinal Serrations And A Box Join | 67 |

| Name | Page |
|--|------|
| Lahey Cholecystectomy Curved Forceps With Longitudinal Serrations And A Box Join | 67 |
| Lahey Scissors Curved 140mm | 19 |
| Lahey Tissue Forceps With 3 Into 3 Teeth And Box Joint Thyroid Vulsellum 150mm | 39 |
| Lahey Tissue Forceps With 3 Into 3 Teeth And Box Joint Thyroid Vulsellum 200mm | 39 |
| Lambotte Osteotome 10mm Wide 150mm | 159 |
| Lambotte Osteotome 10mm Wide x 200mm | 159 |
| Lambotte Osteotome 12mm Wide 150mm | 159 |
| Lambotte Osteotome 15mm Wide 200mm | 159 |
| Lambotte Osteotome 1mm Wide 150mm | 159 |
| Lambotte Osteotome 20mm Wide 200mm | 159 |
| Lambotte Osteotome 25mm Wide 200mm | 159 |
| Lambotte Osteotome 2mm Wide 150mm | 159 |
| Lambotte Osteotome 30mm Wide 200mm | 159 |
| Lambotte Osteotome 3mm Wide 150mm | 159 |
| Lambotte Osteotome 4mm Wide 150mm | 159 |
| Lambotte Osteotome 6mm Wide 150mm | 159 |
| Lambotte Osteotome 7mm Wide | 159 |
| Lambotte Osteotome 8mm Wide 150mm | 159 |
| Landon Vaginal Retractor With 90mm x 25mm Blade And Hollow Handle (Langdon) 200mm | 112 |
| Lane Bone Elevator Slight Curve With A Solid Handle 240mm Curved | 149 |
| Lane Bone Elevator With A Serrated Slight Curved End And A Finger Ring Handle 260mm Curved | 150 |
| Lane Bone Elevator With A Slight Curve And A Ring Handle 265mm Curved | 150 |
| Lane Dissecting Forceps With 1 Into 2 Teeth (Lanes) 140mm Straight | 57 |
| Lane Dissecting Forceps With 1 Into 2 Teeth (Lanes) 150mm Straight | 57 |
| Lane Dissecting Forceps With 1 Into 2 Teeth (Lanes) 180mm Straight | 57 |
| Lane Dissecting Forceps With 1 Into 2 Teeth (Lanes) 200mm Straight | 57 |
| Lane Dissecting Forceps With 1 Into 2 Teeth 130mm (Lanes) 130mm Straight | 57 |

Product Name A-Z

| Name | Page |
|--|------|
| Lane Dissecting Forceps With 2 Into 3 Teeth 140mm (Lanes) 140mm Straight | 57 |
| Lane Dissecting Forceps With 2 Into 3 Teeth 150mm (Lanes) 150mm Straight | 57 |
| Lane Dissecting Forceps With 2 Into 3 Teeth 180mm (Lanes) 180mm Straight | 57 |
| Lane Dissecting Forceps With 2 Into 3 Teeth 200mm (Lanes) 200mm Straight | 57 |
| Lane Dissecting Forceps With 3 Into 4 Teeth (Lanes) 180mm Straight | 57 |
| Lane Dissecting Forceps With 3 Into 4 Teeth (Lanes) 200mm Straight | 57 |
| Lane Dissecting Forceps With 3 Into 4 Teeth 150mm (Lanes) 150mm Straight | 57 |
| Lane Intestinal Clamp With 100mm Longitudinal Serrated Blades And Box Joint 230mm Straight | 141 |
| Lane Intestinal Clamp With 75mm Longitudinal Serrated Blades And Box Joint 180mm Straight | 141 |
| Lane Tissue Forceps With 1 Into 2 Teeth And Box Joint 130mm | 39 |
| Lane Tissue Forceps With 1 Into 2 Teeth And Box Joint 150mm | 39 |
| Lane Tissue Forceps With 1 Into 2 Teeth And Box Joint 195mm | 39 |
| Lane Tissue Forceps With 1 Into 2 Teeth And Box Joint 200mm | 39 |
| Langenbeck Extra Large Retractor Single Blade Size 90mm x 20mm 215mm | 104 |
| Langenbeck Large Retractor Single Blade Size 65mm x 25mm 215mm | 104 |
| Langenbeck Medium Retractor Single Blade Size 44mm x 13mm 215mm | 104 |
| Langenbeck Small Retractor Single Blade Size 23mm x 7mm 215mm | 104 |
| Laparoscopic Monopolar Lead With 4mm Machine End And 4mm Laparoscopic Cutting Sock | 178 |
| Laparoscopic Monopolar Lead With 8mm Bovie Machine End And 4mm Laparoscopic Cutting Sock | 178 |
| Large Pliers | 31 |
| Large Purse String Applicator | 71 |
| Large Teale Gorget | 98 |
| Laryngeal Handle Only 120mm | 117 |

| Name | Page |
|--|------|
| Laryngeal Handle With 10mm Mirror Diameters | 117 |
| Laryngeal Handle With 12mm Mirror Diameters | 117 |
| Laryngeal Handle With 14mm Mirror Diameters | 117 |
| Laryngeal Handle With 16mm Mirror Diameters | 117 |
| Laryngeal Handle With 18mm Mirror Diameters | 117 |
| Laryngeal Handle With 20mm Mirror Diameters | 117 |
| Laryngeal Handle With 22mm Mirror Diameters | 117 |
| Laryngeal Handle With 24mm Mirror Diameters | 117 |
| Laryngeal Handle With 26mm Mirror Diameters | 117 |
| Laryngeal Handle With 28mm Mirror Diameters | 117 |
| Laryngeal Handle With 30mm Mirror Diameters | 117 |
| Laryngeal Handle With 8mm Mirror Diameters | 117 |
| Laster Cheek Retractor For Third Upper Molar 100mm | 103 |
| Lawrence Needle Holder Tungsten Carbide Jaws With Box Joint 150mm Straight | 87 |
| Lawrence Needle Holder With Plain Jaws And Box Joint 150mm Straight | 83 |
| Lawson Tait Needle Holder With Box Joint 130mm Straight | 83 |
| Lawson Tait Scissors (Barnes) (Episiotomy) 130mm Angled | 17 |
| Leibrichs Probe Size 1 - 2 | 93 |
| Leibrichs Probe Size 1 - 3 | 93 |
| Leibrichs Probe Size 2 - 3 | 93 |
| Leibrichs Probe Size 3 - 4 | 93 |
| Leibrichs Probe Size 4 - 5 | 93 |
| Leibrichs Probe Size 5 - 6 | 93 |
| Leibrichs Probe Size 00 - 0 | 93 |
| Leibrichs Probe Size 000 - 0000 | 93 |
| Leksell Bone Rongeur With 4mm Bite And Compound Action 230mm Curved | 145 |

Product Name A-Z

| Name | Page |
|--|------|
| Leksell Bone Rongeur With 6mm Bite And Compound Action 230mm Curved | 145 |
| Leksell Bone Rongeur With 8mm Bite And Compound Action 230mm Curved | 145 |
| Lempert Bone Rongeur With Simple Action And A Box Joint 190mm Straight | 145 |
| Lempert Periosteal Elevator 170mm | 151 |
| Lewer Vulsellum Forceps With 2 Into 2 Teeth And A Box Joint 240mm Curved | 72 |
| Lid For All Wire Mesh Basket 240mm x 250mm | 210 |
| Lid For All Wire Mesh Basket 300mm x 300mm | 210 |
| Lid For All Wire Mesh Basket 480mm x 250mm | 210 |
| Lid For All Wire Mesh Basket 480mm x 250mm | 210 |
| Lid For All Wire Mesh Basket 480mm x 250mm | 210 |
| Lid For All Wire Mesh Basket 480mm x 480mm | 210 |
| Lid For Flat Base Basket 240mm x 250mm | 208 |
| Lid For Flat Base Basket 300mm x 300mm | 208 |
| Lid For Flat Base Basket 480mm x 250mm | 208 |
| Lid For Flat Base Basket 480mm x 480mm | 208 |
| Lid For Flat Base Basket 520 x 250mm | 208 |
| Lid For Flat Base Basket 540mm x 250mm | 208 |
| Lid For Flat Base Basket 670mm x 300mm | 208 |
| Light Weight Tufnol Mallet | 154 |
| Lister Bandage Scissors 115mm Angled | 13 |
| Lister Bandage Scissors 140mm Angled | 13 |
| Lister Bandage Scissors 185mm Angled | 13 |
| Lister Bandage Scissors 200mm Angled | 13 |
| Lister Bandage Scissors 230mm Angled | 13 |
| Lister Bandage Scissors With Tungsten Carbide Jaws 115mm Angled | 22 |
| Lister Bandage Scissors With Tungsten Carbide Jaws 140mm Angled | 22 |
| Lister Bandage Scissors With Tungsten Carbide Jaws 185mm Angled | 22 |
| Lister Bandage Scissors With Tungsten Carbide Jaws 200mm Angled | 22 |
| Lister Bandage Scissors With Tungsten Carbide Jaws 230mm Angled | 22 |

| Name | Page |
|--|------|
| Lister Sinus Forceps With Box Joint 130mm Straight | 60 |
| Lister Sinus Forceps With Box Joint 150mm Straight | 60 |
| Lister Sinus Forceps With Box Joint 180mm Straight | 60 |
| Lister Sinus Forceps With Box Joint 200mm Straight | 60 |
| Lister Urethral Bougies (Sounds) Individual Size 0 / 1 E.g. 280mm | 188 |
| Lister Urethral Bougies (Sounds) Individual Size 0 / 2 E.g. 280mm | 188 |
| Lister Urethral Bougies (Sounds) Individual Size 0.5 / 3 E.g. 280mm | 188 |
| Lister Urethral Bougies (Sounds) Individual Size 1 / 4 E.g. 280mm | 188 |
| Lister Urethral Bougies (Sounds) Individual Size 10 / 13 E.g. 280mm | 188 |
| Lister Urethral Bougies (Sounds) Individual Size 11 / 14 E.g. 280mm | 188 |
| Lister Urethral Bougies (Sounds) Individual Size 12 / 15 E.g. 280mm | 188 |
| Lister Urethral Bougies (Sounds) Individual Size 13 / 16 E.g. 280mm | 188 |
| Lister Urethral Bougies (Sounds) Individual Size 14 / 17 E.g. 280mm | 188 |
| Lister Urethral Bougies (Sounds) Individual Size 15 / 18 E.g. 280mm | 188 |
| Lister Urethral Bougies (Sounds) Individual Size 2 / 5 E.g. 280mm | 188 |
| Lister Urethral Bougies (Sounds) Individual Size 3 / 6 E.g. 280mm | 188 |
| Lister Urethral Bougies (Sounds) Individual Size 4 / 7 E.g. 280mm | 188 |
| Lister Urethral Bougies (Sounds) Individual Size 5 / 8 E.g. 280mm | 188 |
| Lister Urethral Bougies (Sounds) Individual Size 6 / 9 E.g. 280mm | 188 |
| Lister Urethral Bougies (Sounds) Individual Size 7 / 10 E.g. 280mm | 188 |
| Lister Urethral Bougies (Sounds) Individual Size 8 / 11 E.g. 280mm | 188 |
| Lister Urethral Bougies (Sounds) Individual Size 9 / 12 E.g. 280mm | 188 |
| Lister Urethral Bougies (Sounds) Set Of 12 Sizes 0 / 1 Eg To 9 / 12 Eg | 188 |
| Liston Bone Cutting Forceps With A Box Joint 130mm Straight | 68 |

Product Name A-Z

| Name | Page | Name | Page |
|--|------|--|------|
| Liston Bone Cutting Forceps With A Box Joint 150mm Straight | 68 | Lucae Bone Mallet 8 Ounces Light Pattern Toffee Mallet 190mm | 154 |
| Liston Bone Cutting Forceps With A Box Joint 180mm Straight | 68 | Luer Bone Rongeur With Simple Action And A Box Joint 170mm Curved | 145 |
| Liston Bone Cutting Forceps With A Box Joint 200mm Straight | 68 | Luer Bone Rongeur With Simple Action And A Box Joint 170mm Straight | 145 |
| Liston Bone Cutting Forceps With A Box Joint 230mm Straight | 68 | Luer Uterine Vulsellum Forceps With 1 Into 1 Teeth And A Box Joint 240mm Straight | 72 |
| Littlewood Tissue Forceps With 2 Into 3 Teeth And Box Joint (Littlewoods) 190mm | 39 | Macdonald Double Ended Dissector (Mcdonald) 190mm | 136 |
| Lloyd Davies Artery Forceps 250mm Straight | 51 | Macewen Double Ended Curette Oval Spoon And Seeker (Mcewan) 185mm | 122 |
| Lockhart Mummery Curved Fistula Director (No.3, 90 Degree Curve) 140mm Curved | 94 | Macewen Double Ended Curette Round Spoon And Seeker (Mcewan) 185mm | 122 |
| Lockhart Mummery Full Curved Fistula Director (No.4, Retrograde) 160mm Curved | 94 | Magil Intra Tracheal Dilating Forceps Adult Size 240mm | 119 |
| Lockhart Mummery Half Curved Fistula Director (No.2, 45 Degree Curve) 150mm | 94 | Magil Intra Tracheal Dilating Forceps Child Size 200mm | 119 |
| Lockhart Mummery Set Of 4 Fistula Directors | 94 | Maguaran Plate For Holding Draffin Bipod 407mm | 117 |
| Lockhart Mummery Straight Fistula Director (No.1, Straight) 155mm Straight | 94 | Martin Splinter Forceps With Box Joint 115mm Straight | 60 |
| Lorna (Edna) Towel Clip With Box Joint 140mm Curved | 79 | Martin Splinter Forceps With Box Joint 130mm Straight | 60 |
| Love Grunwald Pituitary Rongeur 200mm Shaft 2mm Straight Cut | 171 | Martin Splinter Forceps With Box Joint 75mm Straight | 60 |
| Love Grunwald Pituitary Rongeur 200mm Shaft 2mm Upwards Cut | 171 | Martini Double Ended Bone Curette 150mm | 153 |
| Love Grunwald Pituitary Rongeur 200mm Shaft 3mm Straight Cut | 171 | Mason Needle Holder Tungsten Carbide Jaws And Box Joint 265mm (Masson) Straight | 87 |
| Love Grunwald Pituitary Rongeur 200mm Shaft 3mm Upwards Cut | 171 | Mathieu Urethral Grasping Forceps With Crocodile Action And Serrated Jaws With 200mm Shaft | 187 |
| Love Grunwald Pituitary Rongeur 200mm Shaft 4mm Straight Cut | 171 | Mathieu Urethral Grasping Forceps With Crocodile Action And Serrated Jaws With 280mm Shaft | 187 |
| Love Grunwald Pituitary Rongeur 200mm Shaft 4mm Upwards Cut | 171 | Mayo (Backhaus) Towel Clip With Box Joint 135mm Curved | 78 |
| Love Grunwald Pituitary Rongeur 200mm Shaft 5mm Straight Cut | 171 | Mayo Hegar Needle Holder With Box Joint 130mm Straight | 84 |
| Love Grunwald Pituitary Rongeur 200mm Shaft 5mm Upwards Cut | 171 | Mayo Hegar Needle Holder With Box Joint 140mm Straight | 84 |
| Luc Nasal Turbinate Forceps Oval Large With Screw Joint And Angled Shanks 180mm | 66 | Mayo Hegar Needle Holder With Box Joint 150mm Straight | 84 |
| Luc Nasal Turbinate Forceps Oval Medium With Screw Joint And Angled Shanks 180mm | 66 | Mayo Hegar Needle Holder With Box Joint 165mm Straight | 84 |
| Luc Nasal Turbinate Forceps Oval Small With Screw Joint And Angled Shanks (Lucs) | 66 | Mayo Hegar Needle Holder With Box Joint 180mm Straight | 84 |

Product Name A-Z

| Name | Page |
|--|------|
| Mayo Hegar Needle Holder With Box Joint 190mm Straight | 84 |
| Mayo Hegar Needle Holder With Box Joint 200mm Straight | 84 |
| Mayo Hegar Needle Holder With Box Joint 240mm Straight | 84 |
| Mayo Hegar Needle Holder With Box Joint 300mm Straight | 84 |
| Mayo Hegar Needle Holder With Tungsten Carbide Jaws And Box Joint 130mm Straight | 87 |
| Mayo Hegar Needle Holder With Tungsten Carbide Jaws And Box Joint 140mm Straight | 87 |
| Mayo Hegar Needle Holder With Tungsten Carbide Jaws And Box Joint 150mm Straight | 87 |
| Mayo Hegar Needle Holder With Tungsten Carbide Jaws And Box Joint 180mm Straight | 87 |
| Mayo Hegar Needle Holder With Tungsten Carbide Jaws And Box Joint 200mm Straight | 87 |
| Mayo Hegar Needle Holder With Tungsten Carbide Jaws And Box Joint 230mm Straight | 87 |
| Mayo Hegar Needle Holder With Tungsten Carbide Jaws And Box Joint 250mm Straight | 87 |
| Mayo Hegar Needle Holder With Tungsten Carbide Jaws And Box Joint 300mm Straight | 87 |
| Mayo Scissors 140mm Curved | 14 |
| Mayo Scissors 140mm Straight | 14 |
| Mayo Scissors 165mm Curved | 14 |
| Mayo Scissors 165mm Straight | 14 |
| Mayo Scissors 190mm Curved | 14 |
| Mayo Scissors 190mm Straight | 14 |
| Mayo Scissors 215mm Curved | 14 |
| Mayo Scissors 215mm Straight | 14 |
| Mayo Scissors 230mm Curved | 14 |
| Mayo Scissors 230mm Straight | 14 |
| Mayo Scissors With Tungsten Carbide Jaws 140mm Curved | 22 |

| Name | Page |
|--|------|
| Mayo Scissors With Tungsten Carbide Jaws 140mm Straight | 21 |
| Mayo Scissors With Tungsten Carbide Jaws 150mm Curved | 22 |
| Mayo Scissors With Tungsten Carbide Jaws 150mm Straight | 21 |
| Mayo Scissors With Tungsten Carbide Jaws 180mm Curved | 22 |
| Mayo Scissors With Tungsten Carbide Jaws 180mm Straight | 21 |
| Mayo Scissors With Tungsten Carbide Jaws 200mm Curved | 22 |
| Mayo Scissors With Tungsten Carbide Jaws 200mm Straight | 21 |
| Mayo Scissors With Tungsten Carbide Jaws 230mm Curved | 22 |
| Mayo Scissors With Tungsten Carbide Jaws 230mm Straight | 21 |
| Mayo Supercut Scissors With Black Handles 140mm Curved | 24 |
| Mayo Supercut Scissors With Black Handles 140mm Straight | 24 |
| Mayo Supercut Scissors With Black Handles 150mm Curved | 24 |
| Mayo Supercut Scissors With Black Handles 150mm Straight | 24 |
| Mayo Supercut Scissors With Black Handles 180mm Curved | 24 |
| Mayo Supercut Scissors With Black Handles 180mm Straight | 24 |
| Mayo Supercut Scissors With Black Handles 200mm Curved | 24 |
| Mayo Supercut Scissors With Black Handles 200mm Straight | 24 |
| Mayo Tongue Depressor Double Ended Tapering Adult Size 200mm | 116 |
| Mayo Tongue Depressor Double Ended Tapering Child Size 140mm | 116 |
| Mcclintock Ovum Forceps Curved With A Screw Joint 240mm Curved | 192 |
| Mccullagh Hysterectomy Clamp Mccullogh 250mm Curved | 193 |
| Mccullagh Hysterectomy Clamp Mccullogh 250mm Straight | 193 |
| Mcindoe Bone Cutting Forceps With Compound Action And Angled On Flat 170mm | 68 |

Product Name A-Z

| Name | Page |
|--|------|
| Mcindoe Bone Cutting Forceps With Compound Action And Angled On Flat 190mm | 68 |
| Mcindoe Dissecting Forceps With 1 Into 2 Teeth 150mm Straight | 58 |
| Mcindoe Dissecting Forceps With Serrated Jaws 150mm Straight | 58 |
| Mcindoe Double Prong Skin Hooks With 10mm Wide End 150mm | 95 |
| Mcindoe Double Prong Skin Hooks With 11mm Wide End 150mm | 95 |
| Mcindoe Double Prong Skin Hooks With 9mm Wide End 150mm | 95 |
| Mcindoe Scissors 115mm Curved | 14 |
| Mcindoe Scissors 190mm Curved | 14 |
| Mcindoe Scissors Tc Jaws 140mm Curved | 25 |
| Mcindoe Scissors Tc Jaws 140mm Straight | 25 |
| Mcindoe Scissors Tc Jaws 165mm Curved | 25 |
| Mcindoe Scissors Tc Jaws 165mm Straight | 25 |
| Mcindoe Scissors Tc Jaws 180mm Curved | 25 |
| Mcindoe Scissors Tc Jaws 180mm Straight | 25 |
| Mcindoe Scissors Tc Jaws 190mm Curved | 25 |
| Mcindoe Scissors Tc Jaws 200mm Curved | 25 |
| Mcindoe Scissors Tc Jaws 200mm Straight | 25 |
| Mcindoe Scissors Tc Jaws 230mm Curved | 25 |
| Mcindoe Scissors Tc Jaws 230mm Straight | 25 |
| Mcindoe Scissors Tc Jaws 280mm Curved | 25 |
| Mcindoe Supercut Scissors With Black Handles 180mm Curved | 24 |
| Mckissock Dural Scissor Straight 150mm | 19 |
| Mckissock Dural Scissors Curved 150mm | 19 |
| Metzenbaum Scissors 140mm Curved | 14 |
| Metzenbaum Scissors 140mm Straight | 14 |
| Metzenbaum Scissors 165mm Curved | 14 |
| Metzenbaum Scissors 165mm Straight | 14 |
| Metzenbaum Scissors 190mm Curved | 14 |
| Metzenbaum Scissors 190mm Straight | 14 |

| Name | Page |
|---|------|
| Metzenbaum Scissors 215mm Curved | 14 |
| Metzenbaum Scissors 215mm Straight | 14 |
| Metzenbaum Scissors 240mm Curved | 14 |
| Metzenbaum Scissors 240mm Straight | 14 |
| Metzenbaum Scissors With Tungsten Carbide Jaws 140mm Curved | 22 |
| Metzenbaum Scissors With Tungsten Carbide Jaws 140mm Straight | 22 |
| Metzenbaum Scissors With Tungsten Carbide Jaws 180mm Curved | 22 |
| Metzenbaum Scissors With Tungsten Carbide Jaws 180mm Straight | 22 |
| Metzenbaum Scissors With Tungsten Carbide Jaws 200mm Curved | 22 |
| Metzenbaum Scissors With Tungsten Carbide Jaws 200mm Straight | 22 |
| Metzenbaum Scissors With Tungsten Carbide Jaws 230mm Curved | 22 |
| Metzenbaum Scissors With Tungsten Carbide Jaws 230mm Straight | 22 |
| Metzenbaum Supercut Scissors With Black Handles 150mm Curved | 24 |
| Metzenbaum Supercut Scissors With Black Handles 150mm Straight | 24 |
| Metzenbaum Supercut Scissors With Black Handles 180mm Curved | 24 |
| Metzenbaum Supercut Scissors With Black Handles 180mm Straight | 24 |
| Metzenbaum Supercut Scissors With Black Handles 200mm Curved | 24 |
| Metzenbaum Supercut Scissors With Black Handles 200mm Straight | 24 |
| Micro Halsted Mosquito Artery Forceps With Box Joint 130mm Curved | 46 |
| Micro Halsted Mosquito Artery Forceps With Box Joint 130mm Straight | 46 |
| Miller Bone Rasp 180mm Large | 152 |
| Mini Bone Nibbler 1mm Jaw | 148 |
| Mini Stille Bone Chisel 12mm Wide 150mm | 156 |
| Mini Stille Bone Chisel 3mm Wide 150mm | 156 |
| Mini Stille Bone Chisel 5mm Wide 150mm | 156 |
| Mini Stille Bone Chisel 7mm Wide 150mm | 156 |
| Mini Stille Bone Chisel 9mm Wide 150mm | 156 |

Product Name A-Z

| Name | Page |
|--|------|
| Mini Townsend Biopsy Punch Forceps Crocodile Action 210mm Effective Length With | 201 |
| Mitchel Bone Elevator Curved Blade With A Hollow Handle 200mm Curved | 150 |
| Mitchell Trimmer Double Ended 160mm | 138 |
| Mixer Artery Forceps With Box Joint 200mm Curved | 48 |
| Mixer Artery Forceps With Box Joint 230mm Curved | 48 |
| Mixer Artery Forceps With Box Joint 280mm Curved | 48 |
| Mixer Cholecystectomy Forceps With A Box Joint 160mm Curved | 67 |
| Mixer Cholecystectomy Forceps With A Box Joint 190mm Curved | 67 |
| Mixer O'shaughnessy Artery Forceps With Box Joint (O Shaughnessy) 190mm Curved | 48 |
| Mixer O'shaughnessy Artery Forceps With Box Joint 150mm Curved | 48 |
| Mollison Retractor Adult Size Self Retaining Curved With 4 Into 4 Teeth (Mollins) | 111 |
| Mollison Retractor Child Size Self Retaining Curved With 2 Into 2 Teeth 135mm | 111 |
| Monopolar (Diathermy) Lead With 4mm Machine End And Std 4.8mm Cutting Socket | 178 |
| Monopolar (Diathermy) Lead With 8mm Bovie Machine End And Std 4.8mm Cutting Sock | 178 |
| Moorfield Conjunctival Dissecting Forceps 115mm Straight | 58 |
| Morris Retractor Double Ended 30mm Wide x 40mm Wide And 70mm Deep 230mm | 104 |
| Morris Retractor Double Ended 32mm Wide x 38mm Wide And 57mm Deep 230mm | 104 |
| Morris Retractor Double Ended 44mm Wide x 50mm Wide And 57mm Deep 230mm | 104 |
| Morris Retractor Double Ended 45mm Wide x 50mm Wide And 70mm Deep 230mm | 104 |
| Morris Retractor Double Ended 57mm Wide x 64mm Wide And 57mm Deep 230mm | 104 |
| Morris Retractor Double Ended 60mm Wide x 65mm Wide And 70mm Deep 230mm | 104 |
| Morris Retractor Single Ended 38mm Wide x 50mm Deep 230mm | 104 |
| Morris Retractor Single Ended 50mm Wide x 50mm Deep 230mm | 104 |

| Name | Page |
|--|------|
| Morris Retractor Single Ended 50mm Wide x 52mm Deep 230mm | 104 |
| Morris Retractor Single Ended 50mm Wide x 60mm Deep 230mm | 104 |
| Morris Retractor Single Ended 50mm Wide x 65mm Deep 230mm | 104 |
| Morris Retractor Single Ended 50mm Wide x 70mm Deep 230mm | 104 |
| Morris Retractor Single Ended 65mm Wide x 50mm Deep 230mm | 104 |
| Morris Retractor Single Ended 65mm Wide x 70mm Deep 230mm | 104 |
| Moynihan Artery Forceps With Box Joint 145mm Curved | 48 |
| Moynihan Artery Forceps With Box Joint 180mm Curved | 48 |
| Moynihan Cholecystectomy Large Curved Forceps With Box Joint 230mm Curved | 67 |
| Murphy Hip Skid Fully Insulated Double Ended | 148 |
| Musgrove Aural Double Ended Wool Carrier And Cerumen Pick 170mm | 125 |
| Nail Nipper Concave Double Spring With Lock 140mm Curved | 34 |
| Nail Nipper Concave Double Spring With Lock 160mm Curved | 34 |
| Nail Nipper Straight Double Spring With Lock 140mm Straight | 34 |
| Narrow Box Perf Side Wire Mesh Lid 720 X100 x 25mm | 209 |
| Narrow Box Perf Sides & Mesh Lid 700 x 50 x 25mm | 209 |
| Narrow Box Perforated Sides And Wire Mesh Lid 300mm x 80mm x 35mm | 209 |
| Narrow Box Perforated Sides And Wire Mesh Lid 360mm x 50mm x 25mm | 209 |
| Narrow Box Perforated Sides And Wire Mesh Lid 550mm x 50mm x 25mm | 209 |
| Naunton Morgan Needle Holder With Tungsten Carbide Jaws And Box Joint 250mm Straight | 87 |
| Naunton-Morgan Needle Holder With Multi Action 180mm Straight | 84 |
| Naunton-Morgan Needle Holder With Multi Action 250mm Straight | 84 |
| Negus Knot Tier Or Pusher 200mm | 92 |
| Negus Long Curve Artery Forceps With Screw Joint 190mm Curved | 49 |

Product Name A-Z

| Name | Page |
|---|------|
| Negus Short Curve Artery Forceps With Screw Joint 190mm | 49 |
| Nelson Scissors 230mm Curved (Thoracic) | 15 |
| Nelson Scissors 230mm Straight (Thoracic) | 15 |
| Nelson Scissors 280mm Straight (Thoracic) | 15 |
| Nievert Needle Holder With Offset Ring For Thumb And Box Joint 130mm Straight | 84 |
| Nievert Needle Holder With Tungsten Carbide Jaws And Offset Ring And Box Joint | 87 |
| Norfolk And Norwich Self Retaining Retractor Long Teeth With 4 Into 5 Teeth And Cam Rack. Blades 78mm Deep 220mm | 110 |
| Norfolk And Norwich Self Retaining Retractor With 4 Into 5 Teeth And Cam Rack. B | 110 |
| Northfield Bone Rongeur With 10mm Bite And Simple Action With A Box Joint 200mm | 145 |
| O'connell Double Ended Retractor For Nerve Root (Oconnel O Connel) 215mm | 102 |
| Officer Dissecting Forceps With 3 Into 4 Teeth 200mm Straight | 58 |
| Officer Dissecting Forceps With 3 Into 4 Teeth 250mm Straight | 58 |
| Olivecrona Bone Rongeur Angled To Side With 5mm Bite And A Compound Action 200mm | 146 |
| Olivecrona Bone Rongeur Angled To Side With 6mm Bite And A Compound Action 230mm Angled | 146 |
| Olsen - Hegar Needle Holder With Scissor And Screw Joint 135mm Straight | 84 |
| Olsen - Hegar Needle Holder With Scissors And Screw Joint 160mm Straight | 84 |
| Olsen Hegar Needle Holder With Tungsten Carbide Jaws And Scissors With Screw Joint 140mm Straight | 87 |
| Olsen Hegar Needle Holder With Tungsten Carbide Jaws And Scissors With Screw Joint 160mm Straight | 87 |
| Olsen Hegar Needle Holder With Tungsten Carbide Jaws And Scissors With Screw Joint 180mm Straight | 87 |
| Ormerod Aural Forceps Crocodile Action With A Extra Fine Cup Jaw And A Black Fin | 62 |
| Ormerod Aural Forceps Crocodile Action With Extra Fine Elongated Cup Jaw And Black Finish 70mm Shoulder Length (Omerod) | 62 |

| Name | Page |
|--|------|
| Orthopaedic Plier Cutter With Serrated Jaw Hammer End And Side Cutter 180mm | 30 |
| O'shaughnessy Cholecystectomy Forceps With A Box Joint 200mm Curved | 67 |
| Parametrium Clamps Angulated With Box Joint 240mm | 192 |
| Parametrium Clamps Angulated With Box Joint 210mm | 192 |
| Parametrium Clamps Angulated With Box Joint 300mm | 192 |
| Parametrium Clamps Angulated With Box Joint 350mm | 192 |
| Parametrium Clamps Slightly Curved With Box Joint 210mm | 192 |
| Parametrium Clamps Slightly Curved With Box Joint 240mm | 192 |
| Parametrium Clamps Slightly Curved With Box Joint 300mm | 192 |
| Parametrium Clamps Slightly Curved With Box Joint 350mm | 192 |
| Parametrium Clamps Straight With Box Joint 210mm | 192 |
| Parametrium Clamps Straight With Box Joint 240mm | 192 |
| Parametrium Clamps Straight With Box Joint 300mm | 192 |
| Parametrium Clamps Straight With Box Joint 350mm | 192 |
| Parametrium Clamps Strongly Curved With Box Joint 240mm | 192 |
| Parametrium Clamps Strongly Curved With Box Joint 210mm | 192 |
| Parametrium Clamps Strongly Curved With Box Joint 300mm | 192 |
| Parametrium Clamps Strongly Curved With Box Joint 350mm | 192 |
| Park Anal And Sphincterotomy Combination Retractor Adult Size With Two Sets Of I | 105 |
| Parker Kerr Intestinal Clamp With Screw Joint 230mm Curved | 141 |
| Parker Kerr Intestinal Clamp With Screw Joint 230mm Straight | 141 |
| Parker Kerr Intestinal Clamp With Screw Joint 250mm Curved | 141 |
| Parker Kerr Intestinal Clamp With Screw Joint 250mm Straight | 141 |
| Pennybacker Bone Rongeur Straight Fine Jaws With 4mm Bite And A Simple Action Wi | 146 |

Product Name A-Z

| Name | Page |
|---|------|
| Pennybacker Elevator 4mm Blade Width Straight 170mm | 151 |
| Pennybacker Elevator 6mm Blade Width 170mm | 151 |
| Pennybacker Ronguer Straight 130mm Shaft 3mm Bite | 171 |
| Ph654329 Ward Periosteal Elevator (Plastic And Oral Surgery) 170mm | 151 |
| Phoenix Bone Nibbler 3mm Bite With A Compound Action 230mm 230mm Straight | 147 |
| Phoenix Bone Nibbler 6mm Bite With A Compound Action 230mm Straight | 147 |
| Poirier Tissue Forceps With 2 Into 3 Teeth And Screw Joint 130mm | 39 |
| Pott Smith Scissors Angled On Flat 190mm Angled | 17 |
| Pott Smith Scissors Angled To Side At 25 Deg 190mm | 16 |
| Pott Smith Scissors Angled To Side At 45 Deg 190mm Angled | 16 |
| Power Cutter Tc Jaws | 31 |
| Pratt Vulsellum Forceps With Teeth And A Screw Joint 280mm Straight | 72 |
| Presbyterian Tubing Clamp Forceps With Serrated Jaws And Box Joint 165mm Straight | 76 |
| Presbyterian Tubing Clamp Forceps With Serrated Jaws And Box Joint 195mm Straight | 76 |
| Probe Stainless Steel With Eye 130mm Straight | 92 |
| Probe Stainless Steel With Eye 150mm Straight | 92 |
| Probe Stainless Steel With Eye 180mm Straight | 92 |
| Probe Stainless Steel With Eye 200mm Straight | 92 |
| Probe Stainless Steel With Eye 250mm Straight | 92 |
| Probe Stainless Steel With Eye 300mm Straight | 92 |
| Probe With Eye 130mm Straight | 92 |
| Probe With Eye 150mm Straight | 92 |
| Probe With Eye 180mm Straight | 92 |
| Probe With Eye 200mm Straight | 92 |
| Probe With Eye 250mm Straight | 92 |

| Name | Page |
|---|------|
| Probe With Eye 300mm Straight | 92 |
| Ragnel Davis Double Ended Retractor 8mm x 4mm + 15mm x 5.5mm 150mm | 105 |
| Ramsey Dressing And Sponge Holding Forceps Box Joint 180mm Straight | 61 |
| Ramsey Dressing And Sponge Holding Forceps Box Joint 240mm Straight | 61 |
| Ramsay Dissecting Forceps With 1 Into 2 Teeth (Ramsey) 180mm Straight | 58 |
| Ramsay Dissecting Forceps With 2 Into 3 Teeth (Ramsey) 180mm Straight | 59 |
| Ramsay Dissecting Forceps With Serrated Jaws (Ramsey) 180mm Straight | 58 |
| Receivers Ss 10" (Without Lid) | 211 |
| Receivers Ss 10" (Without Lid) Non Rolled Lips | 211 |
| Receivers Ss 8" (Without Lid) Non Rolled Lips | 211 |
| Retractor Blunt With Double Hook 160mm Straight | 99 |
| Retractor Blunt With Hook And 1 Prong 160mm Straight | 99 |
| Retractor Blunt With Hook And 2 Prongs 160mm Straight | 100 |
| Retractor Blunt With Hook And 3 Prongs 160mm Straight | 100 |
| Retractor Blunt With Hook And 4 Prongs 160mm Straight | 100 |
| Retractor Blunt With Single Hook 160mm Straight | 99 |
| Retractor Sharp With Double Hook 160mm Straight | 99 |
| Retractor Sharp With Hook And 1 Prong 160mm Straight | 99 |
| Retractor Sharp With Hook And 2 Prongs 160mm Straight | 100 |
| Retractor Sharp With Hook And 3 Prongs 160mm Straight | 100 |
| Retractor Sharp With Hook And 4 Prongs 160mm Straight | 100 |
| Retractor Sharp With Single Hook 160mm Straight | 99 |
| Ring Cutter Complete With Saw 160mm | 182 |
| Ring Cutter Forceps Only 160mm | 182 |
| Ring Cutter Saw Blade. Standard Make (30mm Diameter Reference) | 182 |

Product Name A-Z

| Name | Page |
|---|------|
| Ring Type Vasectomy Clamps 3.0mm 140mm | 40 |
| Ring Type Vasectomy Clamps 4.0mm 140mm | 40 |
| Ring Type Vasectomy Clamps 3.5mm 140mm | 40 |
| Robert Artery Forceps With Box Joint 230mm (Roberts) 230mm Curved | 49 |
| Robert Artery Forceps With Screw Joint (300mm Curved | 49 |
| Robert Artery Forceps With Screw Joint 230mm Curved | 49 |
| Robin Anchoring Towel Clip With A 13mm Diameter Clip And Screw Joint 130mm Curved | 79 |
| Robin Anchoring Towel Clip With A 25mm Diameter Clip And Screw Joint 130mm Curved | 79 |
| Robin Anchoring Towel Clip With A 5mm Dia Diameter | 79 |
| Rochester Pean Artery Forceps With Box Joint 130mm Curved | 49 |
| Rochester Pean Artery Forceps With Box Joint 130mm Straight | 49 |
| Rochester Pean Artery Forceps With Box Joint 140mm Curved | 49 |
| Rochester Pean Artery Forceps With Box Joint 140mm Straight | 49 |
| Rochester Pean Artery Forceps With Box Joint 160mm Curved | 49 |
| Rochester Pean Artery Forceps With Box Joint 160mm Straight | 49 |
| Rochester Pean Artery Forceps With Box Joint 180mm Curved | 49 |
| Rochester Pean Artery Forceps With Box Joint 180mm Straight | 49 |
| Rochester Pean Artery Forceps With Box Joint 190mm Curved | 49 |
| Rochester Pean Artery Forceps With Box Joint 190mm Straight | 49 |
| Rochester Pean Artery Forceps With Box Joint 200mm Curved | 49 |
| Rochester Pean Artery Forceps With Box Joint 200mm Straight | 49 |
| Rochester Pean Artery Forceps With Box Joint 230mm Curved | 49 |
| Rochester Pean Artery Forceps With Box Joint 230mm Straight | 49 |

| Name | Page |
|--|------|
| Rogers Hysterectomy Clamp Curved | 193 |
| Rogers Hysterectomy Clamp Straight | 193 |
| Ronald Edwards Rib Spreaders But With A Rack Measuring 450mm And Blades Measuring 120mm Deep And 75mm Wide | 99 |
| Rosen Aural Speculum 4mm Diameter Oval 38mm | 127 |
| Rosen Aural Speculum 4mm Diameter Round 38mm | 126 |
| Rosen Aural Speculum 5mm Diameter Oval 38mm | 127 |
| Rosen Aural Speculum 5mm Diameter Round 38mm | 126 |
| Rosen Aural Speculum 6mm Diameter Oval 38mm | 127 |
| Rosen Aural Speculum 6mm Diameter Round 38mm | 126 |
| Rosen Aural Speculum 7mm Diameter Oval 38mm | 127 |
| Rosen Aural Speculum 7mm Dia Round 38mm | 126 |
| Rotunda Hosp Pattern Vulsellum Fcps 1 Into 1 Teeth | 72 |
| Ruskin Bone Rongeur Curved On Flat With 6mm Bite And A Compound Action 180mm Curved | 146 |
| Ruskin Bone Rongeur Curved On Flat With 6mm Bite And A Compound Action 230mm Cur | 146 |
| Ruskin Liston Bone Cutting Forceps Curved Compound Action 200mm Curved | 68 |
| Ruskin Liston Bone Cutting Forceps With Compound Action 150mm Curved | 68 |
| Ruskin Liston Bone Cutting Forceps With Compound Action 150mm Straight | 68 |
| Ruskin Liston Bone Cutting Forceps With Compound Action 180mm Curved | 68 |
| Ruskin Liston Bone Cutting Forceps With Compound Action 180mm Straight | 68 |
| Ruskin Liston Bone Cutting Forceps With Compound Action 200mm Straight | 68 |
| Russian Dissecting 180mm Forceps | 59 |
| Russian Pattern Dissecting Forceps Oval Toothed End 150mm Straight | 59 |
| Russian Pattern Dissecting Forceps Oval Toothed End 200mm Straight | 59 |
| Russian Pattern Dissecting Forceps Oval Toothed End 250mm Straight | 59 |

Product Name A-Z

| Name | Page |
|--|------|
| Rutherford Morison Tissue Forceps 4 Into 5 Teeth Box Joint 160mm | 39 |
| Ryder Needle Holder With French Eye Tungsten Carbide Jaws And Box Joint 130mm Straight | 88 |
| Ryder Needle Holder With French Eye Tungsten Carbide Jaws And Box Joint 150mm Straight | 88 |
| Ryder Needle Holder With French Eye Tungsten Carbide Jaws And Box Joint 180mm Straight | 88 |
| Ryder Needle Holder With French Eye Tungsten Carbide Jaws And Box Joint 200mm Straight | 88 |
| Ryder Needle Holder With French Eye Tungsten Carbide Jaws And Box Joint 230mm Straight | 88 |
| Ryder Needle Holder With Micro French Eye Tungsten Carbide Jaws And Box Joint 130mm Straight | 88 |
| Ryder Needle Holder With Micro French Eye Tungsten Carbide Jaws And Box Joint 150mm Straight | 88 |
| Ryder Needle Holder With Micro French Eye Tungsten Carbide Jaws And Box Joint 180mm Straight | 88 |
| Ryder Needle Holder With Micro French Eye Tungsten Carbide Jaws And Box Joint 200mm Straight | 88 |
| Ryder Needle Holder With Micro French Eye Tungsten Carbide Jaws And Box Joint 230mm Straight | 88 |
| Safety Pin For Holding Ring Handle Forceps 115mm | 77 |
| Safety Pin For Holding Ring Handle Forceps 130mm | 77 |
| Safety Pin For Holding Ring Handle Forceps 140mm | 77 |
| Safety Pin For Holding Ring Handle Forceps 150mm | 77 |
| Sarrot Artery Forceps With Box Joint 240mm Curved | 51 |
| Sarrot Artery Forceps With Box Joint 240mm Straight | 51 |
| Sauerburch Bone Rongeur With Round Cutting Jaws And A Compound Action 300mm Straight | 146 |
| Sawtell Artery Forceps With Box Joint 180mm Curved | 50 |
| Scalpelless Vasectomy Forceps Curved Sharp 140mm | 40 |

| Name | Page |
|---|------|
| Schaedel Towel Clip With Cross Action 90mm | 78 |
| Schmidt Artery Forceps With Box Joint 190mm Curved | 50 |
| Schoemaker Ogilvie Clamp For Hemisphincterectomy With Hollow Blades And Box Join | 141 |
| Schoemaker Ogilvie Clamp For Hemisphincterectomy With Hollow Blades And Box Joint 155mm Straight | 141 |
| Schoemaker Ogilvie Clamp For Hemisphincterectomy With Hollow Blades And Box Joint 170mm Straight | 141 |
| Schroeder Tenaculum Forceps With 1 Into 1 Teeth And A Box Joint 250mm Straight | 73 |
| Schroeder Uterine Vulsellum Forceps With 2 Into 2 Teeth And A Box Joint 240mm Straight | 72 |
| Schumacher Biopsy Punch For Gynaecology Crocodile Action 230mm Effective Length Complete Straight Cut 4.5 x 7mm 230mm | 203 |
| Sellor Rib Approximator And Spreader Combined 200mm | 161 |
| Sellor Rib Spreader With 3 Pairs Of Detachable Blades And Locking Key (Sternal Spreader) 200mm | 162 |
| Semb Bone Rongeur Slightly Curved With Double Action 20mm Wide Jaws And A Compound Action 240mm Curved | 146 |
| Semi Detach Kerrison Rong Upcut 180mm 45 Deg x 1mm | 172 |
| Semi Detach Kerrison Rong Upcut 180mm 45 Deg x 2mm | 172 |
| Semi Detach Kerrison Rong Upcut 180mm 45 Deg x 3mm | 172 |
| Semi Detach Kerrison Rong Upcut 180mm 45 Deg x 4mm | 172 |
| Semi Detach Kerrison Rong Upcut 180mm 45 Deg x 5mm | 172 |
| Semi Detach Kerrison Rong Upcut 180mm 90 Deg x 1mm | 172 |
| Semi Detach Kerrison Rong Upcut 180mm 90 Deg x 2mm | 172 |
| Semi Detach Kerrison Rong Upcut 180mm 90 Deg x 3mm | 172 |
| Semi Detach Kerrison Rong Upcut 180mm 90 Deg x 4mm | 172 |

Product Name A-Z

| Name | Page |
|--|------|
| Semi Detach Kerrison Rong Upcut 180mm 90 Deg x 5mm | 172 |
| Semi Detach Kerrison Rong Upcut 200mm 45 Deg x 1mm | 173 |
| Semi Detach Kerrison Rong Upcut 200mm 45 Deg x 2mm | 173 |
| Semi Detach Kerrison Rong Upcut 200mm 45 Deg x 3mm | 173 |
| Semi Detach Kerrison Rong Upcut 200mm 45 Deg x 4mm | 173 |
| Semi Detach Kerrison Rong Upcut 200mm 45 Deg x 5mm | 173 |
| Semi Detach Kerrison Rong Upcut 200mm 90 Deg x 1mm | 173 |
| Semi Detach Kerrison Rong Upcut 200mm 90 Deg x 2mm | 173 |
| Semi Detach Kerrison Rong Upcut 200mm 90 Deg x 3mm | 173 |
| Semi Detach Kerrison Rong Upcut 200mm 90 Deg x 4mm | 173 |
| Semi Detach Kerrison Rong Upcut 200mm 90 Deg x 5mm | 173 |
| Semi Detach Kerrison Rong Upcut 210mm 45 Deg x 1mm | 174 |
| Semi Detach Kerrison Rong Upcut 210mm 45 Deg x 2mm | 174 |
| Semi Detach Kerrison Rong Upcut 210mm 45 Deg x 3mm | 174 |
| Semi Detach Kerrison Rong Upcut 210mm 45 Deg x 4mm | 174 |
| Semi Detach Kerrison Rong Upcut 210mm 45 Deg x 5mm | 174 |
| Semi Detach Kerrison Rong Upcut 210mm 90 Deg x 1mm | 174 |
| Semi Detach Kerrison Rong Upcut 210mm 90 Deg x 2mm | 174 |
| Semi Detach Kerrison Rong Upcut 210mm 90 Deg x 3mm | 174 |
| Semi Detach Kerrison Rong Upcut 210mm 90 Deg x 4mm | 174 |
| Semi Detach Kerrison Rong Upcut 210mm 90 Deg x 5mm | 174 |
| Semi Detach Kerrison Upcut 18cm 45 Deg x 1mm/Black | 174 |
| Semi Detach Kerrison Upcut 18cm 45 Deg x 2mm/Black | 174 |
| Semi Detach Kerrison Upcut 18cm 45 Deg x 3mm/Black | 174 |

| Name | Page |
|--|------|
| Semi Detach Kerrison Upcut 18cm 45 Deg x 4mm/Black | 174 |
| Semi Detach Kerrison Upcut 18cm 45 Deg x 5mm/Black | 174 |
| Semi Detach Kerrison Upcut 18cm 90 Deg x 1mm/Black | 174 |
| Semi Detach Kerrison Upcut 18cm 90 Deg x 2mm/Black | 174 |
| Semi Detach Kerrison Upcut 18cm 90 Deg x 3mm/Black | 174 |
| Semi Detach Kerrison Upcut 18cm 90 Deg x 4mm/Black | 174 |
| Semi Detach Kerrison Upcut 18cm 90 Deg x 5mm/Black | 174 |
| Semi Detach Kerrison Upcut 20cm 45 Deg x 1mm/Black | 175 |
| Semi Detach Kerrison Upcut 20cm 45 Deg x 2mm/Black | 175 |
| Semi Detach Kerrison Upcut 20cm 45 Deg x 3mm/Black | 175 |
| Semi Detach Kerrison Upcut 20cm 45 Deg x 4mm/Black | 175 |
| Semi Detach Kerrison Upcut 20cm 45 Deg x 5mm/Black | 175 |
| Semi Detach Kerrison Upcut 20cm 90 Deg x 1mm/Black | 175 |
| Semi Detach Kerrison Upcut 20cm 90 Deg x 2mm/Black | 175 |
| Semi Detach Kerrison Upcut 20cm 90 Deg x 3mm/Black | 175 |
| Semi Detach Kerrison Upcut 20cm 90 Deg x 4mm/Black | 175 |
| Semi Detach Kerrison Upcut 20cm 90 Deg x 5mm/Black | 175 |
| Semi Detach Kerrison Upcut 21cm 45 Deg x 1mm/Black | 175 |
| Semi Detach Kerrison Upcut 21cm 45 Deg x 2mm/Black | 175 |
| Semi Detach Kerrison Upcut 21cm 45 Deg x 3mm/Black | 175 |
| Semi Detach Kerrison Upcut 21cm 45 Deg x 4mm/Black | 175 |
| Semi Detach Kerrison Upcut 21cm 45 Deg x 5mm/Black | 175 |
| Semi Detach Kerrison Upcut 21cm 90 Deg x 1mm/Black | 175 |
| Semi Detach Kerrison Upcut 21cm 90 Deg x 2mm/Black | 175 |

Product Name A-Z

| Name | Page |
|---|------|
| Semi Detach Kerrison Upcut 21cm 90 Deg x 3mm/Black | 175 |
| Semi Detach Kerrison Upcut 21cm 90 Deg x 4mm/Black | 175 |
| Semi Detach Kerrison Upcut 21cm 90 Deg x 5mm/Black | 175 |
| Shea Aural Forceps Crocodile Action With Oval Cup Jaw Curved To The Right And Black Finish 70mm Shoulder Length | 63 |
| Shea Aural Forceps Crocodile Action With Oval Cup Jaw Curved Upwards And Black Finish 70mm Shoulder Length | 63 |
| Shea Aural Forceps Crocodile Action With Oval Cup Jaw Straight And Black Finish 70mm Shoulder Length | 62 |
| Shea Speculum Oblique End Set Of Five | 126 |
| Shea Speculum Oblique End Size No.1 = 3.7mm | 126 |
| Shea Speculum Oblique End Size No.2 = 4.0mm | 126 |
| Shea Speculum Oblique End Size No.3 = 5.0mm | 126 |
| Shea Speculum Oblique End Size No.4 = 6.0mm | 126 |
| Shea Speculum Oblique End Size No.5 = 7.0mm | 126 |
| Sieves 150mm (6inch) Diameter | 185 |
| Sieves 200mm (8inch) Diameter | 185 |
| Silver Nasal Chisel Curved To Left With Side Probe Guard 180mm | 133 |
| Silver Nasal Chisel Curved To Right With Side Probe Guard 180mm | 133 |
| Silver Nasal Chisel Straight With Side Probe Guard 180mm | 133 |
| Sims Diathermy Vaginal Speculum Double Ended Large 80mm x 40mm | 200 |
| Sims Diathermy Vaginal Speculum Double Ended Medium 70mm x 35mm | 200 |
| Sims Diathermy Vaginal Speculum Double Ended Small 60mm x 30mm | 200 |
| Sims Scissors (Uterine) 200mm Curved | 18 |
| Sims Scissors (Uterine) 200mm Straight | 18 |
| Sims Scissors (Uterine) 230mm Curved | 18 |
| Sims Scissors (Uterine) 230mm Straight | 18 |

| Name | Page |
|---|------|
| Sims Scissors (Uterine) Angled To One Side 200mm | 18 |
| Sims Scissors (Uterine) Angled To One Side 230mm | 18 |
| Sims Uterine Curette 12mm Single Ended Sharp 255mm | 194 |
| Sims Uterine Curette 15mm Single Ended Sharp 255mm | 194 |
| Sims Uterine Curette 19mm Single Ended Sharp 255mm | 194 |
| Sims Uterine Curette 3mm Single Ended Sharp 255mm | 194 |
| Sims Uterine Curette 5mm Single Ended Sharp 255mm | 194 |
| Sims Uterine Curette 6mm Single Ended Sharp 255mm | 194 |
| Sims Uterine Curette 7mm Single Ended Sharp 255mm | 194 |
| Sims Uterine Curette 9mm Single Ended Sharp 255mm | 194 |
| Sims Uterine Curette Large Double Ended Sharp / Blunt 255mm | 194 |
| Sims Uterine Curette Large Double Ended Sharp / Sharp 255mm | 194 |
| Sims Uterine Curette Medium Double Ended Sharp / Blunt 255mm | 194 |
| Sims Uterine Curette Medium Double Ended Sharp / Sharp 255mm | 194 |
| Sims Uterine Curette Small Double Ended Sharp / Blunt 255mm | 194 |
| Sims Uterine Curette Small Double Ended Sharp / Sharp 255mm | 194 |
| Sims Uterine Sound Malleable Cm Graduated | 194 |
| Sims Vaginal Speculum Large Double Ended 75mm x 30mm / 85mm x 30mm 185mm | 198 |
| Sims Vaginal Speculum Large Single End 115mm x 43mm 180mm | 198 |
| Sims Vaginal Speculum Medium Double Ended 70mm x 30mm / 75mm x 30mm 155mm | 198 |
| Sims Vaginal Speculum Medium Single End 86mm x 29mm 180mm | 198 |
| Sims Vaginal Speculum Small Double Ended 60mm x 30mm / 70mm x 30mm 130mm | 198 |
| Sims Vaginal Speculum Small Single Ended 76mm x 21mm 180mm | 198 |

Product Name A-Z

| Name | Page |
|--|------|
| Small Pliers | 31 |
| Small Purse String Applicator | 71 |
| Small Teale Gorgets | 98 |
| Small Wire Cutters Tungsten Carbide Jaws 125mm Straight | 28 |
| Smith Peterson Bone Osteotomes 10mm Wide 200mm | 158 |
| Smith Peterson Bone Osteotomes 13mm Wide 200mm | 158 |
| Smith Peterson Bone Osteotomes 20mm Wide 200mm | 158 |
| Smith Peterson Bone Osteotomes 25mm Wide 200mm | 158 |
| Smith Peterson Bone Osteotomes 30mm Wide 200mm | 158 |
| Smith Peterson Bone Osteotomes 6mm Wide 200mm | 158 |
| Spencer Wells Artery Forceps With Box Joint 180mm Straight | 50 |
| Spencer Wells Artery Forceps With Box Joint 130mm Straight | 50 |
| Spencer Wells Artery Forceps With Box Joint 150mm Curved | 50 |
| Spencer Wells Artery Forceps With Box Joint 150mm Straight | 50 |
| Spencer Wells Artery Forceps With Box Joint 180mm Curved | 50 |
| Spencer Wells Artery Forceps With Box Joint 200mm Curved | 50 |
| Spencer Wells Artery Forceps With Box Joint 200mm Straight | 50 |
| Spencer Wells Artery Forceps With Box Joint 230mm Curved | 50 |
| Spencer Wells Artery Forceps With Box Joint 230mm Straight | 50 |
| Spencer Wells Artery Forceps With Box Joint 300mm Curved | 50 |
| Spencer Wells Artery Forceps With Box Joint 300mm Straight | 50 |
| Split Proctoscope | N/A |
| Square Jaw Tubing Clamp Angled 115mm | 76 |
| Square Jaw Tubing Clamp Angled 180mm | 76 |
| Square Jaw Tubing Clamp Straight 130mm | 76 |
| Square Jaw Tubing Clamp Straight 180mm | 76 |

| Name | Page |
|--|------|
| Square Jaw Tubing Clamp Straight 200mm | 76 |
| St Bartholomews Forceps For Tubing Clamp With Square Jaws And Box Joint 150mm Straight | 76 |
| St Bartholomews Wax Hook 180mm Angled | 93 |
| St Barts Blunt Hook Wool Carrier 170mm Angled | 125 |
| St Clair Thomson Adenoid Curette Size 0 With Cage And Metal Handle 8mm 220mm | 130 |
| St Clair Thomson Adenoid Curette Size 1 1/2 With Cage And Metal Handle 12mm 220mm | 130 |
| St Clair Thomson Adenoid Curette Size 1 With Cage And Metal Handle 10mm 220mm | 130 |
| St Clair Thomson Adenoid Curette Size 2 With Cage And Metal Handle 14mm 220mm | 130 |
| St Clair Thomson Adenoid Curette Size 3 With Cage And Metal Handle 16mm 220mm | 130 |
| St Clair Thomson Adenoid Curette Size 4 With Cage And Metal Handle 18mm 220mm | 130 |
| St Clair Thomson Nasal Speculum Blade 38mm Length. Small | 131 |
| St Clair Thomson Nasal Speculum Blade 51mm Length. Medium | 131 |
| St Clair Thomson Nasal Speculum Blade 64mm Length. Large | 131 |
| St Clair Thomson Nasal Speculum Blade 89mm Length. Extra Large | 131 |
| St Georges Pattern Vaginal Speculum Double Ended 100mm x 25mm And 120mm x 30mm | 200 |
| St Marks / Lloyd Davis Retractor With Blade Length 130mm x 60mm Width With Lip | 105 |
| St Marks / Lloyd Davis Retractor With Blade Length 130mm x 60mm Width Without Lip | 105 |
| St Marks / Lloyd Davis Retractor With Blade Length 170mm Into Taper 64mm | 105 |
| St Marks / Lloyd Davis Retractor With Blade Length 170mm Into Taper 64mm | 105 |
| St Marks Hospital Dissecting Forceps With 1 Into 2 Teeth 290mm Straight | 59 |
| St Marks Rectal Dilator Large Hollow Stainless Steel 75mm | 187 |
| St Marks Rectal Dilator Medium Hollow Stainless Steel 75mm | 187 |
| St Marks Rectal Dilator Small Hollow Stainless Steel 75mm | 187 |

Product Name A-Z

| Name | Page |
|---|------|
| Stamm Bone Cutting Forceps 150mm | 68 |
| Stanley Boyd Bone Curette Large Double End Oval Scoops 280mm | 153 |
| Stanley Boyd Bone Curette Medium Double End Oval Scoops 200mm | 153 |
| Stanley Boyd Bone Curette Small Double End Oval Scoops 150mm | 153 |
| Sternum Needle Holder / Wire Twister With Heavy Tungsten Carbide Jaws And Box Joint | 88 |
| Sternum Needle Holder / Wire Twister With Heavy Tungsten Carbide Jaws And Box Joint | 88 |
| Sternum Needleholder Tc Jaws 200mm | 88 |
| Stevens Supercut Tenotomy Scissors 115mm Curved | 23 |
| Stevens Supercut Tenotomy Scissors 115mm Straight | 23 |
| Stevens Supercut Tenotomy Scissors 135mm Curved | 23 |
| Stiles Tissue Forceps With 6 Into 7 Teeth And Box Joint 200mm | 40 |
| Stille Bone Chisel 10mm Wide x 200mm | 156 |
| Stille Bone Chisel 15mm Wide 200mm | 156 |
| Stille Bone Chisel 20mm Wide 200mm | 156 |
| Stille Bone Chisel 25mm Wide 200mm | 156 |
| Stille Bone Chisel 30mm Wide 200mm | 156 |
| Stille Bone Gouge 10mm Wide 200mm | 157 |
| Stille Bone Gouge 13mm Wide 200mm | 157 |
| Stille Bone Gouge 15mm Wide 200mm | 157 |
| Stille Bone Gouge 20mm Wide 200mm | 157 |
| Stille Bone Gouge 25mm Wide 200mm | 157 |
| Stille Bone Gouge 30mm Wide 200mm | 157 |
| Stille Bone Osteotomes 10mm Wide 200mm | 158 |
| Stille Bone Osteotomes 10mm Wide Curved 200mm | 158 |
| Stille Bone Osteotomes 11mm Wide 200mm | 158 |
| Stille Bone Osteotomes 13mm Wide 200mm | 158 |
| Stille Bone Osteotomes 15mm Wide 200mm | 158 |
| Stille Bone Osteotomes 15mm Wide Curved | 158 |

| Name | Page |
|--|------|
| Stille Bone Osteotomes 19mm Wide (3/4inch) Curved | 158 |
| Stille Bone Osteotomes 20mm Wide 200mm | 158 |
| Stille Bone Osteotomes 25mm Wide 200mm | 158 |
| Stille Bone Osteotomes 2mm Wide 200mm | 158 |
| Stille Bone Osteotomes 30mm Wide 200mm | 158 |
| Stille Bone Osteotomes 5mm Wide 200mm | 158 |
| Stille Bone Osteotomes 6mm Wide 200mm | 158 |
| Stille Bone Osteotomes 7mm Wide 200mm | 158 |
| Stille Bone Osteotomes 8mm Wide 200mm | 158 |
| Stille Bone Osteotomes 9mm Wide 200mm | 158 |
| Stille Bone Rongeur Angled To One Side With A Compound Double Action 230mm Angle | 147 |
| Stille Lorenz Plaster Shears 230mm | 167 |
| Stille Lorenz Plaster Shears 265mm | 167 |
| Stille Lorenz Plaster Shears 300mm | 167 |
| Stille Lorenz Plaster Shears 355mm | 167 |
| Stille Lorenz Plaster Shears 380mm | 167 |
| Stille Luer Bone Rongeur Curved With A Compound Action (Or Luer Stille) 220mm Curved | 147 |
| Stille Luer Bone Rongeur Straight With A Compound Action (Or Luer Stille) 220mm | 146 |
| Strabismus Scissors 115mm Curved | 15 |
| Strabismus Scissors 115mm Straight | 15 |
| Strabismus Scissors With Tungsten Carbide Jaws 115mm Curved | 21 |
| Strabismus Scissors With Tungsten Carbide Jaws 115mm Straight | 21 |
| Strabismus Supercut Scissors Curved | 25 |
| Strabismus Supercut Scissors Straight | 25 |
| Straight Nail Nipper Compound Action 150mm Straight | 34 |
| Stratte Needle Holder With Tungsten Carbide Jaws And Curved Shanks 230mm | 88 |
| Syme Double Ended Dissector One Blunt End And One Semi Sharp End 180mm | 136 |
| Symes Aneurysm Needle Flat Handle And Curved Laterally 160mm | 92 |

Product Name A-Z

| Name | Page |
|---|------|
| Teale Uterine Vulsellum Forceps Curved With 3 Into 4 Teeth And A Box Joint 180mm | 72 |
| Teale Uterine Vulsellum Forceps Curved With 3 Into 4 Teeth And A Box Joint 230mm | 72 |
| Thomson Walker Needle Holder Curved Shanks And Box Joint 200mm Straight | 84 |
| Thomson Walker Penile Clamp With Cross Action 115mm | 187 |
| Thudichum Nasal Speculum Size 1 Blade (Smallest) 70mm | 131 |
| Thudichum Nasal Speculum Size 2 Blade 70mm | 131 |
| Thudichum Nasal Speculum Size 3 Blade 70mm | 131 |
| Thudichum Nasal Speculum Size 4 Blade (Medium) 70mm | 131 |
| Thudichum Nasal Speculum Size 5 Blade 70mm | 131 |
| Thudichum Nasal Speculum Size 6 Blade 70mm | 131 |
| Thudichum Nasal Speculum Size 7 Blade (Biggest) 70mm | 131 |
| Thwaites Nail Nipper Pistol Grip Single Spring 140mm | 34 |
| Tilley Aural Dressing Forceps With Serrated Jaws 150mm | 63 |
| Tilley Gouge V Shaped 125mm | 157 |
| Tilley Henckel Nasal Turbinate Forceps 4.5mm Ethmoid Punches | 65 |
| Tilley Henckel Nasal Turbinate Forceps 5.5mm Ethmoid Punches | 65 |
| Tilley Henckel Nasal Turbinate Forceps 6.5mm Ethmoid Punches | 65 |
| Tilley Henckel Nasal Turbinate Forceps 8mm Ethmoid Punches | 65 |
| Tilley Nasal Dressing Forceps With Serrated Jaws And Angled Shanks 165mm | 64 |
| Tischler Biopsy Punch Forceps Crocodile Action 210mm Effective Length With Single | 201 |
| Top Cutters Tc Jaws 150mm | 28 |
| Towel Clip English Pattern With Cross Action 90mm | 78 |
| Towel Clip Non-Tearing With Ball And Socket. Box Joint 100mm Curved | 78 |
| Towel Clip Non-Tearing With Ball And Socket. Box Joint 115mm Curved | 78 |
| Towel Clip Non-Tearing With Ball And Socket. Box Joint 140mm Curved | 78 |

| Name | Page |
|---|------|
| Towel Clip Non-Tearing With Ball And Socket. Box Joint 90mm Curved | 78 |
| Towel Clip Non-Tearing With Serrated Jaws And Box Joint 100mm Curved | 78 |
| Towel Clip Non-Tearing With Serrated Jaws And Box Joint 115mm Curved | 78 |
| Towel Clip Non-Tearing With Serrated Jaws And Box Joint 140mm Curved | 78 |
| Toynbee Aural Speculum 4mm Dia Round 41mm | 127 |
| Toynbee Aural Speculum 5mm Dia Round 41mm | 127 |
| Toynbee Aural Speculum 6mm Dia Round 41mm | 127 |
| Toynbee Aural Speculum 7mm Dia Round 41mm | 127 |
| Travers Retractor With 4 Into 5 Teeth Blunt With Ratchet 210mm | 111 |
| Trethowen Angled Bone Lever Angled | 149 |
| Trethowen Bone Lever 215mm | 149 |
| Treves Dissecting Forceps With 1 Into 2 Teeth 130mm 130mm Straight | 59 |
| Treves Dissecting Forceps With 1 Into 2 Teeth 150mm Straight | 59 |
| Treves Dissecting Forceps With 1 Into 2 Teeth 180mm Straight | 59 |
| Treves Dissecting Forceps With 1 Into 2 Teeth 200mm Straight | 59 |
| Trinket Box Perforated 105mm x 70mm x 25mm | 210 |
| Trinket Box Perforated 120mm x 80mm x 25mm | 210 |
| Trinket Box Perforated 70mm x 30mm x 15mm | 210 |
| Trinket Box Perforated 90mm x 60mm x 25mm | 210 |
| Trinket Box Perforated Sides And Wire Mesh Lid 105mm x 70mm x 25mm | 210 |
| Trinket Box Perforated Sides And Wire Mesh Lid 120mm x 80mm x 25mm | 210 |
| Trinket Box Perforated Sides And Wire Mesh Lid 150mm x 100mm x 25mm | 210 |
| Trinket Box Perforated Sides And Wire Mesh Lid 180mm x 130mm x 25mm | 210 |
| Trinket Box Perforated Sides And Wire Mesh Lid 90mm x 60mm x 25mm | 210 |
| Trotter Bone Rongeur Angled To One Side With A 10mm Bite And Simple Action With | 147 |

Product Name A-Z

| Name | Page |
|--|------|
| Tumarkin Tapered Spec Slotted Size 1 = 6.5mm Black | 127 |
| Tudor Edwards Scapula Retractor With 45mm Wide Blade x 50mm Deep 160mm | 112 |
| Tufnol Handled Curved Osteotome 10mm Wide 265mm Curved | 160 |
| Tufnol Handled Curved Osteotome 12mm Wide 265mm Curved | 160 |
| Tufnol Handled Curved Osteotome 15mm Wide 265mm Curved | 160 |
| Tufnol Handled Curved Osteotome 20mm Wide 265mm Curved | 160 |
| Tufnol Handled Curved Osteotome 25mm Wide 265mm Curved | 160 |
| Tufnol Handled Curved Osteotome 30mm Wide 265mm Curved | 160 |
| Tufnol Handled Straight Gouge 10mm Wide 270mm | 157 |
| Tufnol Handled Straight Gouge 12mm Wide 270mm | 157 |
| Tufnol Handled Straight Gouge 15mm Wide 270mm | 157 |
| Tufnol Handled Straight Gouge 20mm Wide 270mm | 157 |
| Tufnol Handled Straight Gouge 25mm Wide 270mm | 157 |
| Tufnol Handled Straight Gouge 30mm Wide 270mm | 157 |
| Tufnol Handled Straight Osteotome 10mm Wide 265mm Straight | 160 |
| Tufnol Handled Straight Osteotome 12mm Wide 265mm Straight | 160 |
| Tufnol Handled Straight Osteotome 15mm Wide 265mm Straight | 160 |
| Tufnol Handled Straight Osteotome 20mm Wide 265mm Straight | 160 |
| Tufnol Handled Straight Osteotome 25mm Wide 265mm Straight | 160 |
| Tufnol Handled Straight Osteotome 30mm Wide 265mm Straight | 160 |
| Tumarkin Aural Speculum Tapered Slot Set Of 8 | 127 |
| Tumarkin Aural Speculum With Tapered Slot And A Black Finish. Individual Size 0 4.5mm 35mm | 127 |
| Tumarkin Nasal Speculum With Tapered Slot And A Black Finish. | 127 |
| Tumarkin Nasal Speculum With Tapered Slot And A Black Finish. | 127 |

| Name | Page |
|---|------|
| Tumarkin Nasal Speculum With Tapered Slot And A Black Finish. | 127 |
| Tumarkin Tapered Spec Slotted Size 2 = 7.0mm Black | 127 |
| Tumarkin Tapered Spec Slotted Size 3 = 7.5mm Black | 127 |
| Tumarkin Tapered Spec Slotted Size 4 = 8.0mm Black | 127 |
| Umbilical Dilator 95mm | 190 |
| Uni Eppendorf Biopsy Punch Forceps Crocodile Action 230mm Effective Length With Rotating Punch Complete Straight Cut 3 x 6mm 230mm | 202 |
| Uni Tischler Biopsy Punch Forceps Crocodile Action 230mm Effective Length With R | 202 |
| Uni Tischler Biopsy Punch Forceps Crocodile Action 230mm Effective Length With Rotating Punch Complete Upward Angle Cut 3 x 7mm 230mm | 202 |
| Universal Scissors For Ligature And Wire Cutting 115mm Angled | 15 |
| Universal Scissors With Tungsten Carbide Jaws For Ligature And Wire Cutting 115mm | 22 |
| Universal Scope Basket 290mm x 160mm x 52mm | 207 |
| Universal Scope Basket 310mm x 150mm x 77mm | 207 |
| Universal Scope Basket 460mm x 160mm x 52mm | 207 |
| Universal Scope Basket 580mm x 160mm x 52mm | 207 |
| Universal Tough Cut Scissors With Plastic Handles 140mm Angled (Tubing Scissors) | 19 |
| Universal Tough Cut Scissors With Plastic Handles 190mm Angled (Tubing Scissors) | 19 |
| Vehmehren Rib Shears With Sliding Blade Action | 161 |
| Volkman Aural Curette Single End Size 0 Hollow Handle 180mm Straight | 124 |
| Volkman Aural Curette Single End Size 00 Hollow Handle 180mm Straight | 124 |
| Volkman Aural Curette Single End Size 1 Hollow Handle 180mm Straight | 124 |
| Volkman Aural Curette Single End Size 2 Hollow Handle 180mm Straight | 124 |
| Volkman Aural Curette Single End Size 3 Hollow Handle 180mm Straight | 124 |

Product Name A-Z

| Name | Page |
|--|------|
| Volkman Aural Curette Single End Size 4 Hollow Handle 180mm Straight | 124 |
| Volkman Aural Curette Single End Size 5 Hollow Handle 180mm Straight | 124 |
| Volkman Aural Curette Single End Size 6 Hollow Handle 180mm Straight | 124 |
| Volkman Aural Large Curette Double Ended Scoops 215mm Straight | 124 |
| Volkman Aural Medium Curette Double Ended Scoops 215mm Straight | 124 |
| Volkman Aural Small Curette Double Ended Scoops 215mm Straight | 124 |
| Volkman Retractor With 1 Blunt Prong 215mm | 106 |
| Volkman Retractor With 1 Sharp Prong 215mm | 106 |
| Volkman Retractor With 2 Blunt Prongs 215mm | 106 |
| Volkman Retractor With 2 Sharp Prongs 215mm | 106 |
| Volkman Retractor With 3 Blunt Prongs 215mm | 106 |
| Volkman Retractor With 3 Sharp Prongs 215mm | 106 |
| Volkman Retractor With 4 Blunt Prongs 215mm | 106 |
| Volkman Retractor With 4 Sharp Prongs 215mm | 106 |
| Volkman Retractor With 5 Blunt Prongs 215mm | 106 |
| Volkman Retractor With 5 Sharp Prongs 215mm | 106 |
| Volkman Retractor With 6 Blunt Prongs 215mm | 106 |
| Volkman Retractor With 6 Sharp Prongs 215mm | 106 |
| Volkman Spoon Bone Curette Large Double Ended 215mm | 153 |
| Volkman Spoon Bone Curette Medium Double Ended 215mm | 153 |
| Volkman Spoon Bone Curette Small Double Ended 215mm | 153 |
| Walsham Nasal Septum Forceps For Redressing Left With Box Joint 220mm | 66 |
| Walsham Nasal Septum Forceps For Redressing Right With Box Joint 220mm | 66 |
| Walsham Nasal Septum Forceps For Straightening Box Joint 220mm | 66 |

| Name | Page |
|---|------|
| Ward Bone Rongeur Small Slightly Curved On Flat With A Compound Action 180mm Curved | 147 |
| Ward Double Ended Cheek Retractor Used In The Removal Of Wisdom Teeth 150mm | 102 |
| Ward Osteotome 3mm Wide 180mm | 159 |
| Ward Osteotome 5mm Wide 180mm | 159 |
| Warwick James Single End Elevator Curved To Left Light Model 170mm Curved | 137 |
| Warwick James Single End Elevator Curved To Right Light Model 170mm Curved | 137 |
| Warwick James Single End Elevator Straight Light Model 170mm Straight | 137 |
| Watersons Rib Spreader Infant Size With 40mm Wide x 20mm Deep Blades And 105mm Spread | 163 |
| Watson Cheyne Double Ended Elevator Fine With Probe 130mm | 137 |
| Watson Cheyne Double Ended Elevator Fine With Probe 170mm 170mm | 137 |
| Watson Cheyne Double Ended Elevator Fine With Probe 230mm | 137 |
| Waugh's Dissecting Forceps With 1 Into 2 Teeth 150mm Straight | 59 |
| Waugh's Dissecting Forceps With 1 Into 2 Teeth 180mm 180mm Straight | 59 |
| Waugh's Dissecting Forceps With 1 Into 2 Teeth 200mm Straight | 59 |
| Waugh's Dissecting Forceps With Serrated Jaws 150mm (Tonsil) Straight | 59 |
| Waugh's Dissecting Forceps With Serrated Jaws 180mm Straight | 59 |
| Waugh's Dissecting Forceps With Serrated Jaws 200mm Straight | 59 |
| Webster Needle Holder With Tungsten Carbide Jaws And Box Joint 130mm Straight | 88 |
| Weil Blakley Forcep 90deg Up | 65 |
| Weil Straight Blakley Forcep | 65 |
| Weil Upwards Blakley Forcep | 65 |
| Weitlander Retractor With 2 Into 3 Teeth Blunt 115mm | 111 |
| Weitlander Retractor With 2 Into 3 Teeth Sharp (Weislander) 115mm | 111 |
| Weitlander Retractor With 3 Into 4 Teeth Blunt 140mm | 111 |

Product Name A-Z

| Name | Page |
|---|------|
| Weitlander Retractor With 3 Into 4 Teeth Blunt 160mm | 111 |
| Weitlander Retractor With 3 Into 4 Teeth Sharp 140mm | 111 |
| Weitlander Retractor With 3 Into 4 Teeth Sharp 160mm | 111 |
| West Retractor With 2 Into 3 Teeth Blunt 115mm | 111 |
| West Retractor With 2 Into 3 Teeth Sharp 115mm | 111 |
| West Retractor With 3 Into 4 Teeth Blunt 140mm | 111 |
| West Retractor With 3 Into 4 Teeth Sharp 140mm | 111 |
| Whitfield Epilation Forceps With Angled (Oblique) End 90mm Angled | |
| Whyndham Powell Urethral Bougies (Sound) Size 10fg | 188 |
| Whyndham Powell Urethral Bougies (Sound) Size 8fg | 188 |
| Whyndham Powell Urethral Bougies (Sounds) | 188 |
| Whyndham Powell Urethral Bougies (Sounds) Individual Size 12 Fg | 188 |
| Whyndham Powell Urethral Bougies (Sounds) Individual Size 13 Fg | 188 |
| Whyndham Powell Urethral Bougies (Sounds) Individual Size 14 Fg | 188 |
| Whyndham Powell Urethral Bougies (Sounds) Individual Size 15 Fg | 188 |
| Whyndham Powell Urethral Bougies (Sounds) Individual Size 16 Fg | 188 |
| Whyndham Powell Urethral Bougies (Sounds) Individual Size 17 Fg | 188 |
| Whyndham Powell Urethral Bougies (Sounds) Individual Size 18 Fg | 188 |
| Whyndham Powell Urethral Bougies (Sounds) Individual Size 19 Fg | 188 |
| Whyndham Powell Urethral Bougies (Sounds) Individual Size 20 Fg | 188 |
| Whyndham Powell Urethral Bougies (Sounds) Individual Size 21 Fg | 188 |
| Whyndham Powell Urethral Bougies (Sounds) Individual Size 22 Fg | 188 |
| Whyndham Powell Urethral Bougies (Sounds) Individual Size 23 Fg | 188 |
| Whyndham Powell Urethral Bougies (Sounds) Individual Size 24 Fg | 188 |

| Name | Page |
|---|------|
| Whyndham Powell Urethral Bougies (Sounds) Individual Size 25 Fg | 188 |
| Whyndham Powell Urethral Bougies (Sounds) Individual Size 26 Fg | 188 |
| Whyndham Powell Urethral Bougies (Sounds) Individual Size 27 Fg | 188 |
| Whyndham Powell Urethral Bougies (Sounds) Individual Size 28 Fg | 188 |
| Whyndham Powell Urethral Bougies (Sounds) Individual Size 29 Fg | 188 |
| Whyndham Powell Urethral Bougies (Sounds) Individual Size 30 Fg | 188 |
| Whyndham Powell Urethral Bougies (Sounds) Individual Size 31 Fg | 188 |
| Whyndham Powell Urethral Bougies (Sounds) Individual Size 32 Fg | 188 |
| Whyndham Powell Urethral Bougies (Sounds) Individual Size 36 Fg | 188 |
| Whyndham Powell Urethral Bougies (Sounds) Individual Size 38 Fg | 188 |
| Whyndham Powell Urethral Bougies (Sounds) Individual Size 40 Fg | 188 |
| Wilde Aural Spring Pattern Dressing Forceps With Angular Serrated Jaws 130mm | 63 |
| Williger Bone Curette Large Double Ended 125mm | 153 |
| Williger Bone Curette Medium Double Ended 125mm | 153 |
| Williger Bone Curette Small Double Ended 125mm | 153 |
| Wilms Bone Rongeur Angled On Flat With Narrow Jaws And A 4mm Bite With A Simple Action 200mm Angled | 147 |
| Wilson Tonsil Artery Forceps With Angled Shanks Screw Joint 200mm Curved | 50 |
| Wilson Tonsil Scissors 200mm Curved | 16 |
| Wilson Tonsil Scissors 200mm Straight | 16 |
| Winterton Vaginal Speculum | 198 |
| Winterton Vaginal Speculum Diathermy Insulated With Smoke Tube | 199 |
| Winterton Vaginal Speculum With Smoke Tube And Bead Finish | 199 |
| Winterton Vaginal Speculum With Smoke Tube And Black Finish | 199 |
| Wire Cutter Tc 200mm | 29 |

Product Name A-Z

| Name | Page |
|---|------|
| Wire Mesh Base Basket With Perforated Sides 230mm x 130mm x 50mm | 208 |
| Wire Mesh Base Basket With Perforated Sides 230mm x 130mm x 60mm | 208 |
| Wire Mesh Base Basket With Perforated Sides 230mm x 130mm x 70mm | 208 |
| Wire Mesh Base Basket With Perforated Sides 240mm x 130mm x 50mm | 208 |
| Wire Mesh Base Basket With Perforated Sides 240mm x 130mm x 60mm | 208 |
| Wire Mesh Base Basket With Perforated Sides 240mm x 250mm x 100mm | 208 |
| Wire Mesh Base Basket With Perforated Sides 240mm x 250mm x 50mm Half Din | 208 |
| Wire Mesh Base Basket With Perforated Sides 240mm x 250mm x 60mm | 208 |
| Wire Mesh Base Basket With Perforated Sides 240mm x 250mm x 70mm | 208 |
| Wire Mesh Base Basket With Perforated Sides 280mm x 200mm x 60mm | 208 |
| Wire Mesh Base Basket With Perforated Sides 300mm x 300mm x 100mm | 208 |
| Wire Mesh Base Basket With Perforated Sides 300mm x 300mm x 50mm | 208 |
| Wire Mesh Base Basket With Perforated Sides 320mm x 320mm x 60mm | 208 |
| Wire Mesh Base Basket With Perforated Sides 350mm x 350mm x 60mm | 208 |
| Wire Mesh Base Basket With Perforated Sides 380mm x 240mm x 50mm | 208 |
| Wire Mesh Base Basket With Perforated Sides 450mm x 230mm x 50mm | 208 |
| Wire Mesh Base Basket With Perforated Sides 460mm x 320mm x 60mm | 208 |
| Wire Mesh Base Basket With Perforated Sides 470mm x 400mm x 60mm | 208 |
| Wire Mesh Base Basket With Perforated Sides 480mm x 120mm x 60mm | 208 |
| Wire Mesh Base Basket With Perforated Sides 480mm x 150mm x 50mm | 208 |
| Wire Mesh Base Basket With Perforated Sides 480mm x 250mm x 100mm | 208 |
| Wire Mesh Base Basket With Perforated Sides 480mm x 250mm x 120mm | 208 |
| Wire Mesh Base Basket With Perforated Sides 480mm x 250mm x 50mm One Din | 208 |
| Wire Mesh Base Basket With Perforated Sides 480mm x 250mm x 60mm | 208 |

| Name | Page |
|--|------|
| Wire Mesh Base Basket With Perforated Sides 480mm x 250mm x 70mm | 208 |
| Wire Mesh Base Basket With Perforated Sides 480mm x 350mm x 60mm | 208 |
| Wire Mesh Base Basket With Perforated Sides 480mm x 480mm x 100mm | 208 |
| Wire Mesh Base Basket With Perforated Sides 480mm x 480mm x 120mm | 208 |
| Wire Mesh Base Basket With Perforated Sides 480mm x 480mm x 50mm Two Din | 208 |
| Wire Mesh Base Basket With Perforated Sides 480mm x 480mm x 60mm | 208 |
| Wire Mesh Base Basket With Perforated Sides 480mm x 480mm x 70mm | 208 |
| Wire Mesh Base Basket With Perforated Sides 500mm x 500mm x 50mm | 208 |
| Wire Mesh Base Basket With Perforated Sides 520mm x 250mm x 50mm | 208 |
| Wire Mesh Base Basket With Perforated Sides 520mm x 250mm x 60mm | 208 |
| Wire Mesh Base Basket With Perforated Sides 520mm x 250mm x 75mm | 208 |
| Wire Mesh Base Basket With Perforated Sides 520mm x 250mm x 95mm | 208 |
| Wire Mesh Base Basket With Perforated Sides 520mm x 300mm x 60mm | 208 |
| Wire Mesh Base Basket With Perforated Sides 540mm x 250mm x 95mm | 208 |
| Wire Mesh Base Basket With Perforated Sides 550mm x 250mm x 60mm | 208 |
| Wire Mesh Base Basket With Perforated Sides 550mm x 280mm x 50mm | 208 |
| Wire Mesh Base Basket With Perforated Sides 550mm x 325mm x 50mm | 208 |
| Wire Mesh Base Basket With Perforated Sides 550mm x 430mm x 50mm | 208 |
| Wire Mesh Base Basket With Perforated Sides 600mm x 250mm x 150mm | 208 |
| Wire Mesh Base Basket With Perforated Sides 600mm x 300mm x 60mm | 208 |
| Wire Mesh Base Basket With Perforated Sides 600mm x 350mm x 60mm | 208 |
| Wire Mesh Base Basket With Perforated Sides 670mm x 300mm x 70mm | 208 |
| Wire Mesh Base Basket With Perforated Sides 700mm x 250mm x 60mm | 208 |
| Wire Mesh Base Basket With Perforated Sides 700mm x 300mm x 90mm | 208 |

Notes

| Name | Page |
|---|------|
| X-Ray Bag Clip 1inch | 77 |
| Zielke Bone Nibbler 250mm x 3mm Bite 250mm Curved | 148 |
| Zimmer Clip 2.5 inch (63.5mm) Wide With Serrated Jaws. Suitable For Holding Towelling | 77 |

Notes

[illegible]



STERIS IMS INSTRUMENT SALES AND REPAIRS

14 Pindar Road, Hoddesdon,
Hertfordshire EN11 0BZ
United Kingdom

+44 (0)3452 414818
instrument-sales@steris.com
www.healthcare.steris.com/en-gb/gb